



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 3 Stretto Order Form with TRU Balance® Flex Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked For: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and off board charger. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 3 Stretto 3S-SS..... \$8,870

Part: Q6 EDGE 3 STRETTO 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 6 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology \$1,360

Part: XRF000268. HCPCS: K0108

3. BASE COLOR**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.**3.1a. Standard Shrouds**


 Red Corvette Part: DGN1900199-PRD1	 Deep Purple Part: DGN1900199-PPR5	 Midnight Blue Part: DGN1900199-PBU10	 Green Machine Part: DGN1900199-PGN11
 Back in Black Part: DGN1900199-PBK1	 Raspberry Beret Part: DGN1900199-PPK4	 Electric Blue Part: DGN1900199-PBU9	 Military Green Part: DGN1900199-PGN8
 Night Train Part: DGN1900199-PBK3	 Pretty in Pink Part: DGN1900199-PPK3	 Blue Jean Part: DGN1900199-PBU11	
 Grey Street Part: DGN1900199-PSL3	 Evolution Orange Part: DGN1900199-POR7	 Lemon Crush Part: DGN1900199-PYL4	

3.1b. Hydrographic Shrouds

 Military Digital Camo\$263 Part: DGN1801101-WCM7	 Faded Glory\$263 Part: DGN1801101-WFL2
--	--

4. BASE MOUNTED LIGHTS**Optional.** Select one option.

Switch Fender Lights..... No Charge
Part: DGN1810041
Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

**5. QUANTUM BACKUP CAMERA****Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket..... \$21 **Backup Camera Track Mount Bracket** \$27
Part: MEC1910373 Part: MEC2104227
Works with .5" x .5" vertical supports only.

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket\$105
Part: MEC2106123

**6. DRIVE WHEELS****REQUIRED.** Select a type of drive wheel. Accent colors are available for 14" tires.**6.1. Drive Wheels**

3"x14" Pneumatic Tires No Charge Part: WHL1713925	2"x12.5" Pneumatic Tires No Charge Part: WHL1810666
3"x14" Flat Free Tires No Charge Part: WHL1704058	2"x12.5" Flat Free Tires No Charge Part: WHL1810665
4"x14" Knobby, Flat Free Black Tires \$580 Part: XRF000368	4"x14" Knobby, Pneumatic Black Tires \$580 Part: XRF000367

6.2. Drive Wheel Accent Color for 14" Tires

Not available with 12.5" tires.

Back in Black No Charge Part: DGN164245-PBK1 x2	Deep Purple No Charge Part: DGN164245-PPR5 x2	Midnight Blue No Charge Part: DGN164245-PBU10 x2	Green Machine No Charge Part: DGN164245-PGN11 x2
Red Corvette No Charge Part: DGN164245-PRD1 x2	Raspberry Beret No Charge Part: DGN164245-PPK4 x2	Electric Blue No Charge Part: DGN164245-PBU9 x2	Military Green No Charge Part: DGN164245-PGN8 x2
Night Train No Charge Part: DGN164245-PBK3 x2	Pretty in Pink No Charge Part: DGN164245-PPK3 x2	Blue Jean No Charge Part: DGN164245-PBU11 x2	Military Digital Camo \$58 ea. Part: DGN164245-WCM7 x2
Grey Street No Charge Part: DGN164245-PSL3 x2	Evolution Orange No Charge Part: DGN164245-POR7 x2	Lemon Crush No Charge Part: DGN164245-PYL4 x2	Faded Glory \$58 ea. Part: DGN164245-WFL2 x2

7. BATTERIES**Optional.** The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor 2 Battery \$355 ea.
 Part: BAT2004376. HCPCS: E2361
 Uses 8 AMP charger

NF-22, AGM, MK Flat Top Battery \$432 ea.
 Part: BAT2003982. HCPCS: E2361
 Uses 8 AMP charger

Gel Cell, 40AH, Interceptor 2 Battery.....\$191 ea.
 Part: BAT2100316. HCPCS: K0108
 Uses 5 AMP charger

8. PATIENT'S WEIGHT**REQUIRED.** The patient's weight is required to quote and order the power chair.

Up to 150 lb User Weight
 Part: PTOINDV3977

250-300 lb User Weight
 Part: PTOINDV3984

151-250 lb User Weight
 Part: PTOINDV3985

9. ELECTRONICS**REQUIRED.** Select one type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control found on the [Alternative Controls and Electronics Order Form](#).**9.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key JoystickStandard
 Part: CTL2205379
 Controls up to 2 actuators.

QLNES 10 Key Joystick\$420
 Part: CTL2205380
 Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

**9.1b. Expandable Controls**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module\$730
 Part: XRF000002. HCPCS: E2377

Q-Logic 3 Expandable Controller Harness.....\$458
 Part: XRF000264. HCPCS: E2313

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick includes XLR USB charger.

Q-Logic 3 EX Joystick\$840
 Part: CTL167724

Q-Logic 3 EX Joystick for Full Lights\$840
 Part: CTL167725

Q-Logic 3e Joystick\$735
 Part: CTL1910600
 Can operate full light package. Not available with alternative or specialty controls.

No Joystick, Setup for Q-Logic 3
 Part: XRF000010
 Requires selection of Enhanced Display

**9.2. Q-Logic Magnetic Charger**

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$60
 Part: ELE2103635

9.3. Q-Logic 3 Joystick Armor**REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft TouchStandard
 Part: DGN124307
 Also matches Night Train shrouds.

Raspberry Beret..... \$63
 Part: DGN124336-PPK4

Lemon Crush..... \$63
 Part: DGN124336-PYL4

Red Corvette..... \$63
 Part: DGN124336-PRD1

Pretty in Pink \$63
 Part: DGN124336-PPK3

Green Machine \$63
 Part: DGN124336-PGN11

Back in Black..... \$63
 Part: DGN124336-PBK1

Evolution Orange..... \$63
 Part: DGN124336-POR7

Military Green \$63
 Part: DGN124336-PGN8

Grey Street..... \$63
 Part: DGN124336-PSL3

Midnight Blue \$63
 Part: DGN124336-PBU10

Military Digital Camo..... \$63
 Part: DGN124336-WCM7

Deep Purple \$63
 Part: DGN124336-PPR5

Electric Blue \$63
 Part: DGN124336-PBU9

Faded Glory \$63
 Part: DGN124336-WFL2

Blue Jean \$63
 Part: DGN124336-PBU11

10. JOYSTICK HANDLES**Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.



- | | |
|---|---|
| ☐ Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323 | ☐ Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323 |
| ☐ Chin Cup (C) \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control) | ☐ Bodypoint 3" Goal Post Handle \$104
Part: PLSKNCP1057. HCPCS: E2323 |
| ☐ T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323 | ☐ Bodypoint 4" Goal Post Handle \$104
Part: PLSKNCP1058. HCPCS: E2323 |
| ☐ Medium Ball Handle (E) \$69
Part: ACCASMB1859. HCPCS: E2323 | ☐ Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308 |
| ☐ Soft Ball Handle (F) \$80
Part: ACCASMB1856. HCPCS: E2323 | |

**11. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.**11.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617**11.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1032**11.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1032**11.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$368
Part: XRF000003. HCPCS: E1032
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.Multi-Axis Mount \$126
Part: XRF000004
Allows height adjustment of joystick.**11.2. Joystick Mount Location**Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**12. TRU BALANCE FLEX SEATING SYSTEM**

Seat width and depth can be adjusted in 1" increments from 15"W-20"W and 16"D-20"D. Maximum seat size 20"W x 20"D.

12.1. TRU Balance Flex Frame Selection**REQUIRED.**TRU Balance Flex Seating System No Charge
Part: PTOINDV3853**12.2. Occupied Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.WC-19 Occupied Transit Option \$289
Part: KIT152047To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>**12.3. Seat Width****REQUIRED.** Maximum seat width is 20". Part number is SETWDTH plus number following # below.

10" W (#1000)	11" W (#1001)	12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.4. Seat Depth**REQUIRED.** Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

10" D (#1000)	12" D (#1002)	14" D (#1004)	16" D (#1006)	18" D (#1008)	20" D (#1010)
11" D (#1001)	13" D (#1003)	15" D (#1005)	17" D (#1007)	19" D (#1009)	


TRU Balance Flex seating section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**12.5. Seat to Floor Height****Lowest Seat to Floor**
Part: PTOINDV3558**Middle Seat to Floor**
Part: PTOINDV3559**Highest Seat to Floor**
Part: PTOINDV3560

High	18.5" at front with 4" posterior angle
Med	17.5" at front with 2" posterior angle
TRU Balance Flex Seat to Floor Heights	

12.6. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608**Essence SPP**Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608**Solution SPP**Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608**Spectrum Foam SPP**Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608**Spectrum Gel SPP**Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608**Premiere P**Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606**Glacial SP**Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604**Zen SP**Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604**Simplicity G**Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602**2. Essence Cushion Wedge Kit**

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller**Medium Wedge Positioning Kit.....\$72**
Part: POP1909006-M
14"W and smaller**Large Wedge Positioning Kit..... \$105**
Part: POP1909006-L
14"W and smaller**Small Wedge Positioning Kit\$39**
Part: POP1712340-S
15"W and larger**Medium Wedge Positioning Kit.....\$72**
Part: POP1712340-M
15"W and larger**Large Wedge Positioning Kit..... \$105**
Part: POP1712340-L
15"W and larger**3. Cushion Top Foam****Visco Top Foam**
Reference part: XRF000380**FRG Top Foam**
Reference part: XRF000381**No Top Foam**
Reference part: XRF000382**4. Cushion Liner****Stealth Urethane Liner**
Reference part: XRF000383**Stealth Polyurethane Liner**
Reference part: XRF000384**No Liner**
Reference part: XRF000385**5. Cushion Cover****CoolCore Cover**
Reference part: XRF000386**Mesh Cover**
Reference part: XRF000387**Stealth Tek Cover**
Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

TRU Balance Flex seating section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**6. Optional Base Insert**

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer

Reference part: XRF000389

MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner

Reference Part: XRF000390

QTY

Additional Stealth Polyurethane Liner

Reference Part: XRF000391

QTY

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover

Reference Part: XRF000392

QTY

Additional Stealth Tek Cover

Reference Part: XRF000394

QTY

Additional Mesh Cover

Reference Part: XRF000393

QTY

12.7. TRU Balance Flex Back Cane Selection**REQUIRED.**

Angle Adjustable Back Cane..... No Charge

Part: MEC140266

2" Endomorph (1" per side) Back Cane..... \$205

Part: XRF000270

Flip Forward Back Cane \$184

Part: XRF000274

Not available with oxygen holder, or the flip back, cane mount armrest.

4" Endomorph (2" per side) Back Cane..... \$205

Part: XRF000272

2" Mesomorph (1" per side) Back Cane..... \$205

Part: MEC1804736

Manual Recline Back \$788

Part: FRMASMB12910. HCPCS: E1226

Select Synergy Back and an overall back height, or to omit back. Back canes will be 22" H. Not available with flip back armrests, TRU-Comfort Back, Sport Back, rear accessories, or articulating foot platform. Pivot point is 2.5" above the seat pan.

12.8. TRU Balance Flex Back Cane Height Selection**REQUIRED.** 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests.

14" High Canes

Part: FRMASMB15832

17" High Canes

Part: FRMASMB15835

20" High Canes

Part: FRMASMB15838

23" High Canes

Part: FRMASMB15841

15" High Canes

Part: FRMASMB15833

18" High Canes

Part: FRMASMB15836

21" High Canes

Part: FRMASMB15839

24" High Canes

Part: FRMASMB15842

16" High Canes

Part: FRMASMB15834

19" High Canes

Part: FRMASMB15837

22" High Canes

Part: FRMASMB15840

25" High Canes

Part: FRMASMB15843

12.9. TRU Balance Flex Back Cushion Options

Select one type of back. See the seat back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

1a. Synergy Back Selection

Synergy Back..... \$289

Part: XRF000363

1b. Sport Back Selection

Sport Back Cushion, 21" Wide or Less..... \$394

Part: XRF000361

Sport Back Cushion, 22" Wide or More..... \$604

Part: XRF000362

		Back Width											
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
Overall Height	13"	-	-	S	-	-	-	-	-	-	-	-	
	14"	P	T	-	S	P	-	-	-	-	-	-	
	15"	S	P	T	-	S	P	-	-	-	-	-	
	16"	S	P	T	-	S	P	T	P	T	P	T	P
	17"	P	T	S	S	P	T	S	P	T	S	P	T
	18"	-	-	S	S	S	P	T	S	P	T	S	P
	19"	-	-	-	-	S	S	P	T	S	P	T	S
	20"	-	-	-	-	S	S	P	T	S	P	T	S
	21"	-	-	-	-	S	S	P	T	S	P	T	S
	22"	-	-	-	-	S	S	P	T	S	P	T	S
	23"	-	-	-	-	S	S	P	T	S	P	T	S
	24"	-	-	-	-	S	S	P	T	S	P	T	S
	25"	-	-	-	-	S	S	P	T	S	P	T	S

S

Synergy Back

P

Sport Back

T

TRU-Comfort Back

TRU Balance Flex Seat Back Availability Chart

TRU Balance Flex seating section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**1c. TRU-Comfort Back Selection**

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width..... \$744
Part: TC4V10W21WBC. HCPCS: E2620TRU-Comfort Back, Back Cane Option, 22" Width \$855
Part: TC4V22W32WBC. HCPCS: E2621**1d. ADI Back, AL Series, 4-Point, Active Contour**

Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 16" Wide x 19" or 20" Overall Height \$622
Part: ALBR16SA4. HCPCS: E2613Active Contour, 18" Wide x 23" or 24" Overall Height \$622
Part: ALBR18XA4. HCPCS: E2613Active Contour, 16" Wide x 21" or 22" Overall Height \$622
Part: ALBR16TA4. HCPCS: E2613Active Contour, 20" Wide x 19" or 20" Overall Height..... \$622
Part: ALBR20SA4. HCPCS: E2613Active Contour, 16" Wide x 23" or 24" Overall Height \$622
Part: ALBR16XA4. HCPCS: E2613Active Contour, 20" Wide x 21" or 22" Overall Height \$622
Part: ALBR20TA4. HCPCS: E2613Active Contour, 18" Wide x 19" or 20" Overall Height \$622
Part: ALBR18SA4. HCPCS: E2613Active Contour, 20" Wide x 23" or 24" Overall Height..... \$622
Part: ALBR20XA4. HCPCS: E2613Active Contour, 18" Wide x 21" or 22" Overall Height \$622
Part: ALBR18TA4. HCPCS: E2613**1e. ADI Back, AL Series, 4-Point, Deep Contour**

Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 16" Wide x 19" or 20" Overall Height..... \$735
Part: ALBR16SD4. HCPCS: E2615Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$735
Part: ALBR18TD4. HCPCS: E2615Deep Contour, 16" Wide x 21" or 22" Overall Height..... \$735
Part: ALBR16TD4. HCPCS: E2615Deep Contour, 20" Wide x 19" or 20" Overall Height \$735
Part: ALBR20SD4. HCPCS: E2615Deep Contour, 16" Wide x 23" or 24" Overall Height..... \$735
Part: ALBR16XD4. HCPCS: E2615Deep Contour, 20" Wide x 23" or 24" Overall Height \$735
Part: ALBR20XD4. HCPCS: E2615Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$735
Part: ALBR18SD4. HCPCS: E2615**1f. ADI Back, AL Series, 4-Point, Extra Deep Contour**

Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$822
Part: ALBR18S64. HCPCS: E2615Extra Deep Contour, 18" Wide x 23" or 24" Overall Height..... \$822
Part: ALBR18X64. HCPCS: E2615Extra Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$822
Part: ALBR18T64. HCPCS: E2615**1h. Omit Back Selection**Omit Back..... No Charge
Part: PTOINDV1287**12.10. Back Width Package****REQUIRED with Sport and TRU-Comfort back selections.**16"-20" Back Width Package
Part: XRF00003818"-22" Back Width Package
Part: XRF000039**12.11. Overall Back Height****REQUIRED with back selection.** Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges.14" Height
Part: SETHGHT102717" Height
Part: SETHGHT103020" Height
Part: SETHGHT103323" Height
Part: SETHGHT103626" Height
Part: SETHGHT103915" Height
Part: SETHGHT102818" Height
Part: SETHGHT103121" Height
Part: SETHGHT103424" Height
Part: SETHGHT103727" Height
Part: SETHGHT104016" Height
Part: SETHGHT102919" Height
Part: SETHGHT103222" Height
Part: SETHGHT103525" Height
Part: SETHGHT103828" Height
Part: SETHGHT1041**12.12. Armrests****REQUIRED.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. TRU Balance Flex Seat, Drop-In Armrests**

Includes full length, flat arm pads. Available with 14" deep seats and greater.

☐ ☐ Full Length, Drop-In Armrests..... No Charge
L R Left part: ARM1803253. Right part: ARM156742☐ ☐ Desk Length, Drop-In Armrests No Charge
L R Left part: ARM1803254. Right part: ARM156743**1b. Pediatric, Height Adjustable, Removable Armrests**

Available with straight arm pads only.

☐ ☐ Desk Length, Pediatric Armrests \$79
L R Left part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973**1c. Flip Back, Cane Mounted Armrests**

May interfere with swing-away laterals.

☐ ☐ Full Length, Flip Back Armrests..... \$63
L R Left part: XRF000371. Right part: XRF000372☐ ☐ Desk Length, Flip Back Armrests \$63
L R Left part: XRF000373. Right part: XRF000374

TRU Balance Flex seating section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**1d. Two-Post, Flip Back, Height Adjustable, Removable Armrests**

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐ Full Length, Two-Post Armrests.....\$158 ea.
☐ ☐ Left part: XRF000051. Right part: XRF000052. HCPCS: E0973

☐ ☐ Desk Length, Two-Post Armrests.....\$158 ea.
☐ ☐ Left part: XRF000053. Right part: XRF000054. HCPCS: E0973
1e. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests\$79 ea.
☐ ☐ Left part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973

☐ ☐ Desk Length, Single Post Armrests.....\$79 ea.
☐ ☐ Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973
12.13. Armpads**REQUIRED with armrest selection except Drop-In Armrests.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.**1a. Straight Armpads**Straight, Full Length Left and Right No Charge
Part: ARM134275Straight, Desk Length Left, Full Length Right No Charge
Part: ARM134280Straight, Full Length Left, Desk Length Right No Charge
Part: ARM134281Straight, Desk Length Left and Right..... No Charge
Part: ARM134254**1b. Waterfall Armpads**Waterfall, Full Length Left and Right..... No Charge
Part: ARM134285Waterfall, Desk Length Left, Full Length Right No Charge
Part: ARM134288Waterfall, Full Length Left, Desk Length Right No Charge
Part: ARM134289Waterfall, Desk Length Left and Right No Charge
Part: ARM134282**1c. Falcon Gel Armpads**Falcon Gel, Full Length Left and Right.....\$311.22
Part: ARM136438Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
Part: ARM136439Falcon Gel, Full Length Left, Desk Length Right.....\$287.93
Part: ARM136442Falcon Gel, Desk Length Left and Right \$265.02
Part: ARM136437**13. HEADRESTS****Optional.** Select a headrest pad and bracket.**13.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q. HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q. HCPCS: E0955**2. Headrest Brackets**Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033**13.2. Q Cam Lock****Optional.**Q Cam Lock.....\$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

14. LAP BELTS**REQUIRED.** Select one lap belt below.
☐ 38" Lap Belt Part: ACCBELT1006
 ☐ 50" Lap Belt Part: ACCBELT1000
 ☐ 60" Lap Belt Part: ACCBELT1002
 ☐ 70" Lap Belt Part: ACCBELT1003
 ☐ 80" Lap Belt Part: ACCBELT1004
 ☐ 90" Lap Belt Part: ACCBELT1001
15. LEGRESTS**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**15.1a. Basic Center Mount Foot Platform**

Select legrest length in following section.

1. Center Mount Foot Platform TypeBasic Center Mount..... No Charge
Part: XRF000323**2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates**Short Extension, 5"x9" Footplates No Charge
Part: MEC1906675Medium Extension, 5"x9" Footplates No Charge
Part: MEC1906680**3. Basic Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount.....\$158
Part: RIG2104369

Short	Medium
8"-10"	10"-14"
Basic Center Mount Extension Range	

Legrests seating section continued on following pages.

15. LEGRESTS CONTINUED**4. Basic Center Mount Additional Accessories**

Basic Center Mount Footplate Connector \$16
 Part: RIG2104240

15.1b. Center Mount Foot Platform

Select legrest length in following section.

1. Foot Platform Assembly

Quick Adjustable Center Mount Foot Platform \$263 Power Articulating Foot Platform \$3,145
 Part: XRF000311. HCPCS: K0040 Part: PTOINDV3754. HCPCS: E1012+K0040

2. Extension and Angle Adjustable Footplates

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates
 Part: RIG152577

Medium Extension, 5"x11" Footplates
 Part: RIG2404018

Pediatric Extension, 5"x8" Footplates
 Part: RIG2404015

Long Extension, 5"x6" Footplates
 Part: RIG151505

Medium Extension, 5"x6" Footplates
 Part: RIG150206

Long Extension, 5"x8" Footplates
 Part: RIG2404017

Medium Extension, 5"x8" Footplates
 Part: RIG2404016

Long Extension, 5"x11" Footplates
 Part: RIG2404019

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
 Part: RIGASMB7120018

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
 Part: RIGASMB7120019

Small Calf Pads, Quick Adjustable Center Mount Only, \$158
 6"Wx6"H (pair)
 Part: XRF000156

Large Calf Pads, Quick Adjustable Center Mount Only, \$158
 8" taper to 6"Wx8"H (pair)
 Part: XRF000157

3b. Center Mount Calf Panel

Small Calf Panel, 18"Wx8"H \$184
 Part: FRMASMB7110012

Large Calf Panel, 24"Wx8"H \$237
 Part: RIG123825
 Only available on widths of 20" or greater.

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
 Part: RIG152519

Only available with the 5x6 footplates.

Width Extension for 8" \$69
 Part: RIG2404919

Width Extension for 11" \$69
 Part: RIG2404920

Outside Edge \$69
 Part: RIG2404921

Available with 5"x8" footplates.

Front and Rear Cap \$16
 Part: RIG2404923

Available with 5"x8" footplates.

5. Center Mount Additional Accessories

Small Foam Calf Pad Cover (pair) \$67
 Part: ST-CSCS0775-Q

Small Foam Panel Cover (pair) \$176
 Part: ST-CSCP0920-Q

Footplate Connector \$16
 Part: RIGASMB7120023

Footplate Riser \$75
 Part: XRF000348

Select legrest length between 2" to 9" at 1" increments.

Footplate Riser Calf Panel \$50
 Part: FRMASMB7110011

Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222
 Part: ST-FPG5008T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers \$270
 Part: ST-AFPG5011T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension \$247
 Part: ST-AFPG5008E-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Cover for Width Extension \$270
 Part: ST-AFPG5011E-Q. HCPCS: K0108

15.1c. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

☐ P ☐ L ☐ R 60°-90° Swing-Away Legrest, Static Seat \$86 ea.
 Pair part: XRF000082. Left: XRF000083. Right: XRF000084
 Not available with 11"-13" extensions.

☐ P ☐ L ☐ R 60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea.
 Pair part: XRF000085. Left: XRF000086. Right: XRF000087
 For use only with 11"-13" lower extensions.

2. Lower Extension

☐ P ☐ L ☐ R 11"-13" Length
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.

☐ P ☐ L ☐ R 13"-16.5" Length
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.

☐ P ☐ L ☐ R 15.5"-19.5" Length
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

Legrests section continued on following pages.

15. LEGRESTS CONTINUED**3. Angle Adjustable Footplate(s)**

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐	5"Wx8"D with Heel Loops..... \$90 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951	☐ ☐ ☐	6"Wx11"D¹ with Heel Loops..... \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951
☐ ☐ ☐	6"Wx8"D¹ with Heel Loops..... \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951	☐ ☐ ☐	8"Wx11"D² with Heel Loops..... \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea. Part: INDPART2822. HCPCS: E1020	Therafin 32387 Amputee Support Gel Cover \$155 ea. Part: INDPART2823. HCPCS: K0108
--	---

15.1d. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐	Elevating Legrest Assembly \$171 ea. Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990 Select lower extension and footplates
--	---

2a. Lower Extension

☐ ☐ ☐	11"-14" Length Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030 For reinforced only	☐ ☐ ☐	14"-17" Length Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035 For reinforced only
☐ ☐ ☐	8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103 For angle adjustable only	☐ ☐ ☐	13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106 For angle adjustable only

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad..... \$342 Part: FRMASMB5301. HCPCS: E1020	Right Residual Limb Support with Calf Pad \$342 Part: FRMASMB5300. HCPCS: E1020
--	--

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐	7.25"Wx6"D Reinforced Footplate..... Standard Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056	☐ ☐ ☐	8.5"Wx6"D Reinforced Footplate Standard Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062 May cause interference on seat widths less than 18."
--	---	--	---

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐	5"Wx8"D with Heel Loops..... \$90 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951	☐ ☐ ☐	6"Wx11"D¹ with Heel Loops..... \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951
☐ ☐ ☐	6"Wx8"D¹ with Heel Loops..... \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951	☐ ☐ ☐	8"Wx11"D² with Heel Loops..... \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25 Part: INDPART2242. HCPCS: E0953
--

15.2. Legrest Accessories

Heel Loops (pair)..... \$53 Part: FRMASMB7873. HCPCS: E0951 Standard on angle adjustable and multi axis footplates when noted.	Wheel Bumper..... \$19 ea. Part: XRF000100
--	---

15.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	4" L (#3991)	6" L (#3993)	8" L (#3680)	9" L (#3682)	10" L (#3684)	11" L (#3686)	12" L (#3688)	13" L (#3690)	14" L (#3692)	15" L (#3694)
3" L (#3990)	5" L (#3992)	7" L (#3994)	8.5" L (#3681)	9.5" L (#3683)	10.5" L (#3685)	11.5" L (#3687)	12.5" L (#3689)	13.5" L (#3691)	14.5" L (#3693)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)
7.5" L (#3471)	8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	

16. POWER OPTION ELECTRONICS**REQUIRED** for articulating foot platform.Single Actuator Through Controller Interface \$2,163
Part: XRF000335. HCPCS: E2310Single Actuator Through Controller and Interface (AAM).....\$3,040
Part: XRF000068. HCPCS: E2310**17. SEAT MOUNTED LED LIGHTS****Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045

Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

18. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.Flex Seat Cup Holder No Charge
Part: ACCASMB2634

Not available on TRU Balance Flex seating without two-post armrests.

Oxygen Holder for Flex Seat \$210
Part: ACC125006. HCPCS: E2208Cane/Crutch Holder \$105
Flex part: ACC152515. HCPCS: E2207Cell Phone Holder for 12" Lift Captain Seat.....\$70
Part: ACC1611300Flag Holder \$55
Part: ACCASMB2415Saddle Bag \$40
Part: DWR1020D010Stealth Cell Phone/Media Device Holder \$215
Part: ACC7120006Quantum Backpack \$168
Part: ACC1905853Grip Lap Board (20005).....\$38.95
Part: ACC152992Grip Drink Holder (20085) \$15.50
Part: ACC152993Accessory Bar \$132
Part: XRF000128Large, Black, Weather Cover \$202
Part: ACCCOVR1021TRU Balance 3 Transfer Handle \$105 ea.
Part: MEC2205680. HCPCS: E0705
Not available with Heavy Duty armrests.Q-Logic EX/NE Programmer \$573
Part: ELEASMB6976Q-Logic PC Programming Disk with Cable \$445
Part: ELEASMB5216Privacy Bag \$42
Part: ACC2005723Personal Item Pouch \$84
Part: ACC2101017Remote Stop Switch \$1,295
Part: XRF000379**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
QTY Part: ACC2407046Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
QTY Part: ACC2405496**19. BRANCH PATCHES****Optional.**Air Force Patch\$21
Part: ACC140240Coast Guard Patch\$21
Part: ACC140239Navy Patch\$21
Part: ACC140244Army Patch\$21
Part: ACC140241Marines Patch.....\$21
Part: ACC140243Purple Heart Patch\$21
Part: ACC140238**20. POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#).**20.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with 2-post, flip up armrests. These laterals are not available with the back cane option.

1. BracketSwing-away, Rail Mounts, pair (TWBSTL) \$562
Part: INDPART2896. HCPCS: E1034
For use with Synergy Back.Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034
Only available with back cane and omit back configurations.**2. Pads**Curved Pads, 3"x5", pair (TWBL-3X5) \$250
Part: INDPART2663. HCPCS: E0956Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250
Part: INDPART2664. HCPCS: E0956**3. Installation**Install Laterals
Part: PTOINDV3499Do Not Install
Part: PTOINDV3500**4. Lateral Extensions****REQUIRED** with Synergy Back and thoracic lateral selections.Lateral Extension for Synergy Back (pair)..... \$28 ea.
Part: XRF000150

Positioning components section continued on following page.

20. POSITIONING COMPONENTS CONTINUED**20.2. Stealth UniLink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea. Removable Mounting Hardware QTY Part: ST-UTTP5-Q. HCPCS: E1028/E1034^	Flip Down with LEFT Full Surface Contact \$281 ea. Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028/E1034^
UniLink, 7" Height Adjustable, \$276 ea. Removable Mounting Hardware QTY Part: ST-UTTP7-Q. HCPCS: E1028/E1034^	Flip Down with RIGHT Full Surface Contact..... \$281 ea. Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028/E1034^
UniLink Multi-Angle Thigh Guide, Left \$358.87 ea. Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^	UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea. Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea. QTY Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^	UniLink, 4" x 12" Pad \$106 ea. QTY Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^
UniLink, 4" x 8" Pad..... \$106 ea. QTY Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^	

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea. QTY Part: ST-UIL00150A000-Q	UniLink, 2" Link..... \$81 ea. QTY Part: ST-UIL001200A000-Q
---	---

4. UniLink Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
---	--	--

20.3. Stealth Position Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea. Part: ST-IFPS1103-Q. HCPCS: E0960	Upper Contour, Medium (IFPS2103)..... \$154 ea. Part: ST-IFPS2103-Q. HCPCS: E0960
Upper Classic, Large (IFPS1104) \$154 ea. Part: ST-IFPS1104-Q. HCPCS: E0960	Upper Contour, Large (IFPS2104) \$154 ea. Part: ST-IFPS2104-Q. HCPCS: E0960

20.4. Stealth Position Chest Straps

Chest Strap, Medium.....\$97 ea. w/ Belt Anchor, (IFCS4113) Part: ST-IFCS4113-Q. HCPCS: E0960	Chest Strap, Large.....\$97 ea. w/ Belt Anchor, (IFCS4114) Part: ST-IFCS4114-Q. HCPCS: E0960
--	---

20.5. Stealth Swing-Away Headrest Mounts


These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right..... \$318 ea. Part: ST-SAR2A-R-Q. HCPCS: E1028	UniLink Facial Lateral, Dual Swing Away \$492 ea. 2 Upper & 1 Lower Left Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028
UniLink Facial Lateral, Swing Away, Left \$318 ea. Part: ST-SAR2A-L-Q. HCPCS: E1028	UniLink Facial Lateral, Dual Swing Away \$434 ea. 2 Upper Part: ST-SAR4A-DBU-Q. HCPCS: E1028
UniLink Facial Lateral, Dual Swing Away \$492 ea. 2 Upper & 1 Lower Right Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028	

2. Facial Lateral Pads and Accessories

Egg Switch Mount (ESM) (A)\$57 ea. QTY Part: INDPART2815	Swing-Away Joystick Mount, 18" Rod (C) \$276 ea. QTY Part: ST-SUS18-Q Allows mounting of facial pads and attaches to the headrest
Facial Pad, 4".5x2", Cool Core (900M-Q)\$61 ea. QTY Part: ST-900M-Q. HCPCS: K0108	
Facial Pad, 6"x2", Cool Core (920M-Q).....\$61 ea. QTY Part: ST-920M-Q. HCPCS: K0108	

**20.6. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket, Tube\$199 ea. Style Arms (ARMS-ES-L) Part: POS124064. HCPCS: E1028	Right, Swing-Away, Elbow Support Bracket, Tube.....\$199 ea. Style Arms (ARMS-ES-R) Part: POS124066. HCPCS: E1028
--	--

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea. QTY Part: POS124065. HCPCS: K0108	4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea. QTY Part: POS124067. HCPCS: K0108
--	---

Positioning components section continued on following page.

20. POSITIONING COMPONENTS CONTINUED

20.7. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

20.8. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

21. LAP TRAY

21.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out..... \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware..... \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	ML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.