



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 3 Stretto Order Form TRU Balance® 4 Positioning

Account Number: _____ Date: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
PO Number: _____
Marked for: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and off board charger. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 3 Stretto 3SP-SS\$8,559

Part: Q6 EDGE 3 STRETTO 3SP-SS. HCPCS: K0856

Single power option.

Quantum Q6 Edge 3 Stretto 3MP-SS\$8,373

Part: Q6 EDGE 3 STRETTO 3MP-SS. HCPCS: K0861

Multiple actuator power option.

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 6 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology\$1,360

Part: XRF000268. HCPCS: K0108

Account Number: _____ Date: _____ Marked for: _____

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette Part: DGN1900199-PRD1	 Deep Purple Part: DGN1900199-PPR5	 Midnight Blue Part: DGN1900199-PBU10	 Green Machine Part: DGN1900199-PGN11
 Back in Black Part: DGN1900199-PBK1	 Raspberry Beret Part: DGN1900199-PPK4	 Electric Blue Part: DGN1900199-PBU9	 Military Green Part: DGN1900199-PGN8
 Night Train Part: DGN1900199-PBK3	 Pretty in Pink Part: DGN1900199-PPK3	 Blue Jean Part: DGN1900199-PBU11	
 Grey Street Part: DGN1900199-PSL3	 Evolution Orange Part: DGN1900199-POR7	 Lemon Crush Part: DGN1900199-PYL4	

3.1b. Hydrographic Shrouds



Military Digital Camo \$263
Part: DGN1801101-WCM7



Faded Glory..... \$263
Part: DGN1801101-WFL2


4. BASE MOUNTED LIGHTS

Optional. Select one option.

Switch Fender Lights..... No Charge

Part: DGN1810041

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.



5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Mounting Bracket

Backup Camera Headrest Bracket..... \$21

Part: MEC1910373

Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$27

Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



6. DRIVE WHEELS

REQUIRED. Select a type of drive wheel. Accent colors are available for 14" tires.

6.1. Drive Wheels

3"x14" Pneumatic Tires..... No Charge
Part: WHL1713925

2"x12.5" Pneumatic Tires No Charge
Part: WHL1810666

3"x14" Flat Free Tires..... No Charge
Part: WHL1704058

2"x12.5" Flat Free Tires No Charge
Part: WHL1810665

6.2. Drive Wheel Accent Color for 14" Tires

Not available with 12.5" tires.

Back in Black..... No Charge
Part: DGN164245-PBK1 x2

Pretty in Pink No Charge
Part: DGN164245-PPK3 x2

Green Machine No Charge
Part: DGN164245-PGN11 x2

Red Corvette..... No Charge
Part: DGN164245-PRD1 x2

Evolution Orange..... No Charge
Part: DGN164245-POR7 x2

Military Green No Charge
Part: DGN164245-PGN8 x2

Night Train..... No Charge
Part: DGN164245-PBK3 x2

Midnight Blue No Charge
Part: DGN164245-PBU10 x2

Military Digital Camo..... \$58 ea.
Part: DGN164245-WCM7 x2

Grey Street..... No Charge
Part: DGN164245-PSL3 x2

Electric Blue No Charge
Part: DGN164245-PBU9 x2

Faded Glory \$58 ea.
Part: DGN164245-WFL2 x2

Deep Purple..... No Charge
Part: DGN164245-PPR5 x2

Blue Jean No Charge
Part: DGN164245-PBU11 x2

Raspberry Beret..... No Charge
Part: DGN164245-PPK4 x2

Lemon Crush..... No Charge
Part: DGN164245-PYL4 x2

7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor 2 Battery \$355 ea.
Part: BAT2004376 x2. HCPCS: E2361
Uses 8 AMP charger

NF-22, AGM, MK Flat Top Battery \$432 ea.
Part: BAT2003982 x2. HCPCS: E2361
Uses 8 AMP charger

8. PATIENT'S WEIGHT**REQUIRED.** The patient's weight is required to quote and order the power chair.Up to 150 lb User Weight
Part: PTOINDV3977151 lb-250 lb User Weight
Part: PTOINDV3985251 lb-300 lb User Weight
Part: PTOINDV3984**9. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**9.1. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module\$730
Part: XRF000002. HCPCS: E2377Q-Logic 3 Expandable Controller Harness.....\$458
Part: XRF000264. HCPCS: E2313**2. Q-Logic 3 Joystick**

Q-Logic 3 EX Joystick includes XLR USB charger.

Q-Logic 3 EX Joystick\$840
Part: CTL167724No Joystick, Setup for Q-Logic 3
Part: XRF000010Q-Logic 3e Joystick\$735
Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.




Q-Logic 3 EX Joystick

9.2. Q-Logic 3 Joystick Armor**REQUIRED with Q-Logic 3 Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2.Black, Soft TouchStandard
Part: DGN124307Raspberry Beret.....\$63
Part: DGN124336-PPK4Blue Jean\$63
Part: DGN124336-PBU11Red Corvette.....\$63
Part: DGN124336-PRD1Pretty in Pink\$63
Part: DGN124336-PPK3Lemon Crush.....\$63
Part: DGN124336-PYL4Back in Black.....\$63
Part: DGN124336-PBK1Evolution Orange.....\$63
Part: DGN124336-POR7Green Machine\$63
Part: DGN124336-PGN11Grey Street.....\$63
Part: DGN124336-PSL3Midnight Blue\$63
Part: DGN124336-PBU10Military Green\$63
Part: DGN124336-PGN8Deep Purple.....\$63
Part: DGN124336-PPR5Electric Blue\$63
Part: DGN124336-PBU9**10. JOYSTICK HANDLES****Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.Mushroom Handle (B)\$90
Part: ACCASMB1860. HCPCS: E2323Extended I-Handle (G)\$90
Part: ACCASMB1857. HCPCS: E2323Chin Cup (C)\$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)Bodypoint 3" Goal Post Handle\$104
Part: PLSKNCP1057. HCPCS: E2323T-Bar (D)\$126
Part: ACCASMB1858. HCPCS: E2323Bodypoint 4" Goal Post Handle\$104
Part: PLSKNCP1058. HCPCS: E2323Medium Ball Handle (E)\$69
Part: ACCASMB1859. HCPCS: E2323Q-Logic 3 Carrot Handle (A)No Charge
Part: DGN124308Soft Ball Handle (F)\$80
Part: ACCASMB1856. HCPCS: E2323


Joystick Handles

11. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**11.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for Q-Logic 3Standard
Part: FRMASMB10617**11.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3\$319
Part: ACC1908434. HCPCS: E1032**11.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side\$245
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side\$245
Part: ACC2201042. HCPCS: E1032

Swing Away Joystick



Retract4

Joystick mounting brackets section continued on following page.

11. JOYSTICK MOUNTING BRACKETS CONTINUED**11.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$368
 Part: XRF000003. HCPCS: E1032
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount \$126
 Part: XRF000004
 Allows height adjustment of joystick.

**11.2. Joystick Mount Location**

Left Side Mount
 Part: PTOINDV1225

Right Side Mount
 Part: PTOINDV1226

12. TRU BALANCE® 4 SEATING

Maximum seat size is 22"W x 22"D.

12.1. TRU Balance 4 Configuration**REQUIRED.**

TRU-Balance 4 Power Tilt Only \$8,595
 Part: PTOINDV4014. HCPCS: E1002
 50° Center of Gravity (CG) power tilt seating system.

TRU-Balance 4 Power Tilt and Recline \$12,790
 Part: PTOINDV3954. HCPCS: E1007
 50° of posterior tilt and 165° of recline.

12.2. 12" Power Lift

Requires combined selection of TRU-Balance 4 Power Tilt and Recline and Articulating Foot Platform. Requires selection in 12.3 Anterior Tilt. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

Power Lift with iLevel \$3,725
 Part: PTOINDV3956. HCPCS: E2298/E2300^
 12" scissor lift with positional feedback.

12.3. Anterior Tilt

Chest support is required with anterior tilt 20 and 30 degree functions. See chest support section. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.

1. Anterior Tilt Angle

Only available with Power Lift.

10° Anterior Tilt Angle No Charge
 Part: MEC2009356

30° Anterior Tilt Angle \$2,595
 Part: MEC2009400. HCPCS: K0108
 Requires chest support and knee block.

20° Anterior Tilt Angle \$1,995
 Part: MEC2009399. HCPCS: K0108
 Requires chest support and knee block.

0° Anterior Tilt No Charge
 Part: PTOINDV3980

12.4. Knee Block Assembly**REQUIRED with 20° or 30° anterior tilt. Optional with 10° and 0° anterior tilt.**

Anterior Knee Block Assembly (ST-KBATB4ATC-1-Q) \$314
 Part: ST-KBATB4ATC-1-Q. HCPCS: K0108

Anterior Knee Block Assembly, Curved (ST-KBATB4ATC-U-Q) .. \$354
 Part: ST-KBATB4ATC-U-Q. HCPCS: K0108

12.5. Occupied Transit Kit**Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option \$289
 Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

12.6. Seat Width**REQUIRED with seat selection. 12"-15" widths not available with Power Lift. Maximum seat width is 22". Part number is SETWDTH plus number following # below.**

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.7. Seat Width Package**REQUIRED with seat selection.**

12"-16" Width Package
 Part: XRF000006
 Available with Power Tilt Only.

16"-20" Width Package
 Part: XRF000008

20"-24" Width Package
 Part: XRF000158
 Maximum width is 22".

12.8. Seat Depth**REQUIRED with seat selection. 12"-13" depth not available with Power Lift. Maximum seat depth is 22". Part number is SETDPTH plus number following # below.**

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------


TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**12.9. Seat Depth Package****REQUIRED with seat selection.****1a. Power Tilt Depth Package**12"-16" Depth Package
Part: XRF00001116"-20" Depth Package
Part: XRF00001320"-24" Depth Package
Part: XRF000159
Maximum depth is 22".**1b. Power Lift Depth Package**16"-22" Depth Package
Part: XRF000309**12.10. Seat to Floor Height****REQUIRED.**Low Seat to Floor
Part: PTOINDV3976Middle Seat to Floor
Part: PTOINDV3559Highest Seat to Floor
Part: PTOINDV3560

Low	Med	High
Q6 Edge 3 Stretto		
18.0"	18.5"	19.0"
TRU Balance 4 Seat to Floor Heights		

12.11. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type**

TRU-Comfort 2 SPP
Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-####XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPP
Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-####XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPP
Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-####XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPP
Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-####XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPP
Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-####XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere P
Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-####XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SP
Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-####XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SP
Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-####XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity G
Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-####XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller

Medium Wedge Positioning Kit.....\$72
Part: POP1909006-M
14"W and smaller

Large Wedge Positioning Kit..... \$105
Part: POP1909006-L
14"W and smaller

Small Wedge Positioning Kit\$39
Part: POP1712340-S
15"W and larger

Medium Wedge Positioning Kit.....\$72
Part: POP1712340-M
15"W and larger

Large Wedge Positioning Kit..... \$105
Part: POP1712340-L
15"W and larger

TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**3. Cushion Top Foam**Visco Top Foam
Reference part: XRF000380No Top Foam
Reference part: XRF000382FRG Top Foam
Reference part: XRF000381**4. Cushion Liner**Stealth Urethane Liner
Reference part: XRF000383No Liner
Reference part: XRF000385Stealth Polyurethane Liner
Reference part: XRF000384**5. Cushion Cover**CoolCore Cover
Reference part: XRF000386Stealth Tek Cover
Reference part: XRF000388Mesh Cover
Reference part: XRF000387

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	—	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389
MSRP: \$48-\$69**7. Additional Cushion Liner**

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390Additional Stealth Polyurethane Liner
Reference Part: XRF000391**8. Additional Cushion Cover**

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392Additional Stealth Tek Cover
Reference Part: XRF000394Additional Mesh Cover
Reference Part: XRF000393**12.12. Back Type Selection****REQUIRED.** Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back

Part: XRF000014

Select a back cushion configuration in the subsequent section.

Back Cane Option

Part: XRF000015

Available only with TRU-Balance 3 Power Tilt (1 actuator) selections.
Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back..... \$184

Part: XRF000016

Not available with power recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

12.13. Back Cane Selection**REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type**Angle Adjustable Back Cane..... No Charge
Part: MEC1402662" Endomorph (1" per side) Back Cane..... \$74
Part: XRF000269Flip Forward Back Cane \$53
Part: XRF0002734" Endomorph (2" per side) Back Cane..... \$74
Part: XRF000271**2. Back Cane Height**16"H Cane
Part: FRMASMB1583420"H Cane
Part: FRMASMB1583824"H Cane
Part: FRMASMB1584217"H Cane
Part: FRMASMB1583521"H Cane
Part: FRMASMB1583925"H Cane
Part: FRMASMB15843**12.14. Back Cushion Configuration****REQUIRED.** Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1a. Sport Back/Omit/After-Market Interface**Sport Back..... No Charge
Part: XRF000017
Complete following width package and size.

After-market Back Interface \$535

Part: ST-SSAPT84-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.

Omit Back Pan/Cushion..... No Charge
Part: PTOINDV1287
Select a following width package.

TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**1b. TRU-Comfort Back**

TRU-Comfort Back (Static/Flip Forward), 12" to 20" Width..... \$744
 Part: TC4V10W21W. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$855
 Part: TC4V22W32W. HCPCS: E2621
 Complete following width package and size.

1c. Q-Back 2.0, AL Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568
 Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$568
 Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$568
 Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL20AXX. HCPCS: E2613

1d. Q-Back 2.0, AL Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643
 Part: QB2AL14DS. HCPCS: E2620

Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643
 Part: QB2AL14DT. HCPCS: E2620

Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643
 Part: QB2AL16DS. HCPCS: E2620

Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643
 Part: QB2AL16DT. HCPCS: E2620

Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643
 Part: QB2AL16DX. HCPCS: E2620

Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL16DXX. HCPCS: E2620

Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
 Part: QB2AL18DS. HCPCS: E2620

Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643
 Part: QB2AL18DT. HCPCS: E2620

Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643
 Part: QB2AL18DX. HCPCS: E2620

Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL18DXX. HCPCS: E2620

Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643
 Part: QB2AL20DS. HCPCS: E2620

Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
 Part: QB2AL20DT. HCPCS: E2620

Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643
 Part: QB2AL20DX. HCPCS: E2620

Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL20DXX. HCPCS: E2620

1e. Q-Back 2.0, AL Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703
 Part: QB2AL14ES. HCPCS: E2620

Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$703
 Part: QB2AL16ES. HCPCS: E2620

Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$703
 Part: QB2AL16ET. HCPCS: E2620

Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$703
 Part: QB2AL16EX. HCPCS: E2620

Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703
 Part: QB2AL18ES. HCPCS: E2620

Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$703
 Part: QB2AL18ET. HCPCS: E2620

Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703
 Part: QB2AL18EX. HCPCS: E2620

Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL18EXX. HCPCS: E2620

Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703
 Part: QB2AL20ES. HCPCS: E2620

Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703
 Part: QB2AL20ET. HCPCS: E2620

Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703
 Part: QB2AL20EX. HCPCS: E2620

Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL20EXX. HCPCS: E2620

1f. Q-Back 2.0, AFT Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AFT16AS. HCPCS: E2613

AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AFT18AS. HCPCS: E2613

1g. Q-Back 2.0, AFT Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643
 Part: QB2AFT16DS. HCPCS: E2620

AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643
 Part: QB2AFT16DT. HCPCS: E2620

AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
 Part: QB2AFT18DS. HCPCS: E2620

AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643
 Part: QB2AFT18DT. HCPCS: E2620

AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
 Part: QB2AFT20DT. HCPCS: E2620

TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**1h. Q-Back 2.0, AFT Series, Extra Deep Contour**

Select a back type and the specific overall back height in subsequent section.

AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline\$703
Part: QB2AFT16EXX. HCPCS: E2620**12.15. TRU Balance 4 Back Shroud****REQUIRED.** The back shroud is a feature of the TRU Balance 4 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette Part: XRF000019	Grey Street Part: XRF000023	Pretty in Pink Part: XRF000027	Electric Blue Part: XRF000030	Green Machine Part: XRF000033
Back in Black Part: XRF000021	Deep Purple Part: XRF000025	Evolution Orange Part: XRF000028	Blue Jean Part: XRF000031	Military Green Part: XRF000034
Night Train Part: XRF000022	Raspberry Beret Part: XRF000026	Midnight Blue Part: XRF000029	Lemon Crush Part: XRF000032	

12.16. Lateral Positioning Wedge**Optional.** Available with Sport and TRU-Comfort backs. Not available on power recline back 19"-20"H with limited adjustability on 21"H.Lateral Positioning Wedge \$16 ea.
Part: SET123465 x2**12.17. Back Width Package****REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections Part: PTOINDV3468	14"W-16"W Package Part: XRF000037	16"W-20"W Package Part: XRF000038	18"W-22"W Package Part: XRF000039
--	--------------------------------------	--------------------------------------	--------------------------------------

12.18. Back Width**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

12" Back Width Part: PTOINDV3562	16" Back Width Part: PTOINDV3564	18" Back Width Part: PTOINDV3566	20" Back Width Part: PTOINDV3568	22" Back Width Part: PTOINDV3570
14" Back Width Part: PTOINDV3563	17" Back Width Part: PTOINDV3565	19" Back Width Part: PTOINDV3567	21" Back Width Part: PTOINDV3569	

Overall Height	Tilt Only Width												Recline Width																			
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"		
	14"	S	T	S	-	-	-	-	-	-	-	12"	-	-	-	-	-	-	-	-	-	12"	-	-	-	-	-	-	-	-	14"	
15"	S	T	S	-	-	-	-	-	-	-	-	15"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"	
16"	S	S	T	S	T	S	T	S	T	S	T	S	T	-	-	-	-	-	-	-	16"	-	-	-	-	-	-	-	-	-	-	16"
17"	S	S	T	S	T	S	T	S	T	S	T	S	T	T	-	-	-	-	-	-	17"	-	-	-	-	-	-	-	-	-	-	17"
18"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	-	-	-	-	-	-	18"	-	-	-	-	-	-	-	-	-	-	18"
19"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	S	S	T	S	T	S	19"	S	S	T	S	T	S	T	S	T	S	19"
20"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	S	S	T	S	T	S	20"	S	S	T	S	T	S	T	S	T	S	20"
21"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	-	S	T	S	T	S	21"	S	S	T	S	T	S	T	S	T	S	21"
22"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	-	S	T	S	T	S	22"	S	S	T	S	T	S	T	S	T	S	22"
23"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	-	S	T	S	T	S	23"	S	S	T	S	T	S	T	S	T	S	23"
24"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	-	S	T	S	T	S	24"	S	S	T	S	T	S	T	S	T	S	24"
25"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	-	S	T	S	T	S	25"	S	S	T	S	T	S	T	S	T	S	25"
26"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	-	S	T	S	T	S	26"	S	S	T	S	T	S	T	S	T	S	26"
27"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	-	S	T	S	T	S	27"	S	S	T	S	T	S	T	S	T	S	27"
28"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	-	S	T	S	T	S	28"	S	S	T	S	T	S	T	S	T	S	28"
29"	-	S	T	S	T	S	T	S	T	S	T	S	T	T	-	S	T	S	T	S	29"	S	S	T	S	T	S	T	S	T	S	29"
30"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	30"	S	T	S	T	S	T	S	T	S	T	30"
31"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	31"	S	T	S	T	S	T	S	T	S	T	31"
32"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	32"	S	T	S	T	S	T	S	T	S	T	32"

S

 Sport Back

Q

 Q-Back

T

 TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

12.19. Back Cushion Overall Back Height**Power Recline available height range: 19"-32"H. Odd heights include 6" gap and even heights include 7" gap between seat pan and bottom of back cushion. See the availability chart for stock order back size ranges.**

19" Height Part: SETHGHT1032	22" Height Part: SETHGHT1035	25" Height Part: SETHGHT1038	28" Height Part: SETHGHT1041	31" Height Part: SETHGHT1102
20" Height Part: SETHGHT1033	23" Height Part: SETHGHT1036	26" Height Part: SETHGHT1039	29" Height Part: SETHGHT1042	32" Height Part: SETHGHT1103
21" Height Part: SETHGHT1034	24" Height Part: SETHGHT1037	27" Height Part: SETHGHT1040	30" Height Part: SETHGHT1043	

12.20. Optional Back Height Extension**Optional.** Adds 2" of seat back pan height. 25"H selection includes a back height extension.Back Height Extension \$48
Part: XRF000041

13. ARMRESTS**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.**13.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

☐ L	☐ R	Flip-Back, 8"-16"H Armrests Standard
Left part: ARM1806662. Right part: ARM1806663		

Armrest Support**Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

☐ L	☐ R	Flip-Back Armrest Downpost \$92 ea.
Left part: XRF000265. Right part: XRF000266		

13.1b. Two-Post, Height Adjustable, Removable Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

☐ L	☐ R	Full Length, Two-Post Armrest for Power Seating No Charge	☐ L	☐ R	Desk Length, Two-Post Armrest for Power Seating.... No Charge
Left part: XRF000047. Right part: XRF000048			Left part: XRF000049. Right part: XRF000050		

13.1c. Pediatric, Removable, Height Adjustable ArmrestsHeight adjustment from 7.75" to 9.75" in .5" increments. **Not available with full length armpads.**

☐ L	☐ R	Pediatric, Removable Armrests \$79 ea.
Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973		

13.2. Armrest Height**REQUIRED.** Set left and right armrest height. Note the available height range in the option description.**1a. Flip-Back and Two-Post Left Armrest Height**

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

1b. Pediatric, Left Armrest Height

7.75" H PTOINDV3722	8.25" H PTOINDV3724	8.75" H PTOINDV3726	9.25" H PTOINDV3728	9.75" H PTOINDV3730
------------------------	------------------------	------------------------	------------------------	------------------------

2a. Flip-Back and Two-Post Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

2b. Pediatric, Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	8.75" H PTOINDV3725	9.25" H PTOINDV3727	9.75" H PTOINDV3729
------------------------	------------------------	------------------------	------------------------	------------------------

13.3. Armrest Size**REQUIRED.** Select a left (L) and right (R) armrest size. This also affects armpad selection.

☐ L	☐ R	10" Length L part: XRF000057. R part: XRF000123	☐ L	☐ R	14" Length L part: XRF000058. R part: XRF000124	☐ L	☐ R	17" Length L part: XRF000059. R part: XRF000125 Only available with Flip-Back, 8"-16" H Armrests and straight arm pads.
----------------------------	----------------------------	--	----------------------------	----------------------------	--	----------------------------	----------------------------	---

14. ARMPADS**REQUIRED.****14.1a. Omit Armpad**Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection..... No Charge	Interface for Aftermarket Pad No Charge
Part: XRF000060	Part: XRF000308
Not available with Flip-back, 8"-16" H armrests.	Only available with Flip-back, 8"-16" armrests. Not available with chest bar.

14.1b. Foam Armpads

☐ L	☐ R	Straight, Armpad Standard	☐ L	☐ R	Waterfall, Armpad No Charge
Left part: XRF000061. Right part: XRF000126			Left part: XRF000062. Right part: XRF000127		
			Not available with Flip-back, 8"-16" armrests.		

14.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right\$195	Falcon Gel, Desk Length Left, Full Length Right.....\$171
Part: ARM1806807	Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right.....\$171	Falcon Gel, Desk Length Left and Right\$147
Part: ARM1806813	Part: ARM1806806

Armpads section continued on following page.

14. ARMPADS CONTINUED**b. Falcon Gel Armpads for Tube Style Armrests**Falcon Gel, Full Length Left and Right\$311.22
Part: ARM136438Falcon Gel, Full Length Left, Desk Length Right.....\$287.93
Part: ARM136442Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
Part: ARM136439Falcon Gel, Desk Length Left and Right \$265.02
Part: ARM136437**15. ARMREST BRACKETS**Standard Fixed Receiver Standard
Part: XRF000063Stealth Height Adjustable Receiver \$263
Part: XRF000064. HCPCS: K0108**16. HEADRESTS****Optional.** Select a headrest pad and bracket.**16.1. Headrest Pads, Redesign Back Assembly**

Only available with 14" to 22" back widths of Static Type or Flip Forward back, without occupied transit.

1. Headrest PadsStealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q. HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q. HCPCS: E0955**2. Headrest Brackets**Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033**16.2. Q Cam Lock****Optional.**Q Cam Lock.....\$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

17. LAP BELTS**REQUIRED.** Select a lap belt below.**17.1a. Standard Lap Belts**
☐ 38" Lap Belt Part: ACCBELT1006

 ☐ 50" Lap Belt Part: ACCBELT1000

 ☐ 60" Lap Belt Part: ACCBELT1002

 ☐ 70" Lap Belt Part: ACCBELT1003

 ☐ 80" Lap Belt Part: ACCBELT1004

 ☐ 90" Lap Belt Part: ACCBELT1001
17.1b. Optional Lap BeltsStealth Positioning, 60" Single Pull, Push Button, Medium.....\$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-QStealth Positioning, 60" Single Pull, Push Button, Large\$50
(9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-QStealth Positioning, 60" Single Pull, Push Button, Extra Large.....\$50
(11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q60" Sure-Lok Lap Belt (Non-retractable)..... \$120
Part: ACCBELT1016
Required with ISO 7176-19 Compliant Transit Option.**18. LEGRESTS**

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

18.1a. Basic Center Mount Foot Platform**1. Basic Center Mount Foot Platform Assembly**Basic Center Mount..... No Charge
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Extension and Angle Adjustable FootplatesShort Extension, 5"x9" Footplates \$180
Part: XRF000312. HCPCS: K0040Medium Extension, 5"x9" Footplates \$180
Part: XRF000313. HCPCS: K0040Long Extension, 5"x9" Footplates \$180
Part: XRF000314. HCPCS: K0040**3. Basic Center Mount Calf Pads**Small Calf Pads, Basic Center Mount..... \$158
Part: RIG2104369Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997
Only available with Long Extension.**4. Basic Center Mount Additional Accessory**Basic Center Mount Footplate Connector \$16
Part: RIG2104240

Legrests section continued on following page.

18. LEGRESTS CONTINUED**18.1b. Center Mount Foot Platform****1. Foot Platform Assembly**

Dual Actuator, Power Articulating Foot Platform..... \$3,595

Part: RIG2009246. HCPCS: E1012+K0040

Positional feedback with independent control of articulation.

Quick Adjustable Center Mount Foot Platform..... \$263

Part: XRF000311. HCPCS: K0040

Not available with power lift.

2. Extension and Angle Adjustable Footplates

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x8" Footplates

Part: RIG2404015

Medium Extension, 5"x11" Footplates

Part: RIG2404018

Long Extension, 5"x11" Footplates

Part: RIG2404019

Medium Extension, 5"x6" Footplates

Part: RIG150206

Long Extension, 5"x6" Footplates

Part: RIG151505

Medium Extension, 5"x8" Footplates

Part: RIG2404016

Long Extension, 5"x8" Footplates

Part: RIG2404017

Legrest Compatibility Articulating Foot Platform TRU Balance 4	Seat to Floor Option		
	Low	Middle	High
Medium Length Extension	11.5"-14"	11.5"-15"	11.5"-16"
Long Length Extension	N/A	14.5"-15"	14.5"-16"

Legrest Compatibility Quick Adjustable CM TRU Balance 4	Seat to Floor Option		
	Low	Middle	High
Pediatric Length Extension	8.5"-12.5"	8.5"-12.5"	8.5"-12.5"
Medium Length Extension	11"-14.5"	11"-15.5"	11"-16.5"
Long Length Extension	N/A	14"-15.5"	14"-16.5"

3. Center Mount Calf Support

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard

Part: RIGASMB7120018

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69

Part: RIG152519

Only available with the 5x6 footplates.

Outside Edge \$69

Part: RIG2404921

Available with 5"x8" footplates.

Width Extension for 8" \$69

Part: RIG2404919

Front and Rear Cap \$16

Part: RIG2404923

Available with 5"x8" footplates.

Width Extension for 11" \$69

Part: RIG2404920

5. Center Mount Additional Accessories

Small Reversed Dartex Calf Pad Cover (pair) \$158

Part: RIG133717

Footplate Connector \$16

Part: RIGASMB7120023

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222

Part: ST-FPG5008T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension \$247

Part: ST-AFPG5008E-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers \$270

Part: ST-AFPG5011T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Cover for Width Extension \$270

Part: ST-AFPG5011E-Q. HCPCS: K0108

18.2. Legrest Accessories

Wheel Bumper \$37

Part: XRF000100

18.3. Legrest Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)
11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	

19. POWER POSITIONING ELECTRONICS**REQUIRED.****19.1a. Q-Logic 3 Controls through Joystick**Single Actuator Through Controller and Interface \$2,163
Part: XRF000335. HCPCS: E2310Single Actuator Through Controller and Interface (AAM) \$3,040
Part: XRF000068. HCPCS: E2310Multiple Actuator Through Controller and Interface \$3,040
Part: XRF000336. HCPCS: E2311**19.1b. Q-Logic 3 Controls with iAccess Switch Control**Single Actuator Control Through Flat Button \$2,547
iAccess and Controller and Interface
Part: XRF000337. HCPCS: E2310Single Actuator Control Through Toggle Button \$2,547
iAccess and Controller and Interface
Part: XRF000338. HCPCS: E2310Single Actuator Control Through Flat Button \$3,460
iAccess and Controller and Interface (AAM)
Part: XRF000070. HCPCS: E2310Single Actuator Control Through Toggle Button \$3,460
iAccess and Controller and Interface (AAM)
Part: XRF000074. HCPCS: E2310Multiple Actuator Control Through Flat Button \$3,460
iAccess and Controller and Interface
Part: XRF000339. HCPCS: E2311Multiple Actuator Control Through Toggle Button \$3,460
iAccess and Controller and Interface
Part: XRF000340. HCPCS: E2311**19.2. iAccess Mount****REQUIRED with iAccess options.****1. iAccess Mount Selection**Joystick Mount No Charge
Part: MEC1907500Under Arm Mount, Flip-Back Armrest No Charge
Part: MEC1906531Under Arm Mount, Tube Style Armrest No Charge
Part: MEC1906530Attendant Mount \$26
Part: MEC1906564Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$26
Part: MEC2001392Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$26
Part: MEC2001391**2. Under Arm Mount Side**Left Side Under Arm Mount
Part: PTOINDV3678Right Side Under Arm Mount
Part: PTOINDV3679**20. ATTENDANT CONTROL****Optional.** Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.Attendant Control \$525
Part: XRF000078. HCPCS: E2331

Select this option if no charge option criteria is not met.

Attendant Control with Tilt, Recline, Alternative No Charge
Controls and Q-Logic Display or SCIM Module
Part: XRF000079**21. ACCESSORIES****Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.Accessory Bar \$132
Part: XRF000128


Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.

Backpack Holder (Pair) (A) \$69
Part: ACC147774Cane/Crutch Holder \$105
Part: ACC152515. HCPCS: E2207Clothing Guard (Pair) \$237
Part: ACC125007Compact Vent Tray \$1,465
Part: ACC133889. HCPCS: E1029Cup Holder No Charge
Part: ACC125003Flag Kit \$69
Part: ACC1700443Flip-Back Armrest, Side Track Accessory Bracket ... \$11
Part: MEC1809178Front Mounted, Cell Phone/ Media Holder, \$69
X-Grip (E)
Part: ACC155835Glove Box (C) \$158
Part: ACC125011Grip Drink Holder (20085) \$15.50
Part: ACC152993Grip Lap Board (20005) \$38.95
Part: ACC152992Hydration System (B) \$237
Part: ACC125010Large, Black, Weather Cover \$202
Part: ACCCOVR1021Oxygen Holder \$210
Part: ACC125006. HCPCS: E2208

Includes accessory bar, a credit will be applied PTOINDV3769.

Personal Item Hook \$69
Part: ACC125008Personal Item Pouch \$84
Part: ACC2101017Privacy Bag \$42
Part: ACC2005723Push Handles (pair) \$195
Part: ACC2005027

Not available with oxygen holder or vent tray.

Q-Logic EX/NE Programmer \$573
Part: ELEASMB6976Quantum Backpack \$168
Part: ACC1905853Rear View Mirror \$69
Part: ACC155833

Accessories section continued on following page.

21. ACCESSORIES CONTINUED

Reflector Kit \$53
Part: ACC124670

Stealth Large Tablet Holder \$625.28
Part: ACC141818

Transfer Handle (D).....\$105 ea.
Part: MEC2205680 x2. HCPCS: E0705

Depth Adjustable Vent Tray \$595
Part: ACC2302500. HCPCS: E1029
Includes accessory bar, a credit will be applied PTOINDV3769.

Remote Stop Switch \$1,295
Part: XRF000379

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
Part: ACC2407046

Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
Part: ACC2405496

22. BRANCH PATCHES

Optional.

Air Force Patch\$21
Part: ACC140240

Army Patch\$21
Part: ACC140241

Coast Guard Patch\$21
Part: ACC140239

Marines Patch\$21
Part: ACC140243

Navy Patch\$21
Part: ACC140244

Purple Heart Patch\$21
Part: ACC140238

**23. POSITIONING COMPONENTS**Optional. More options are available on the [Positioning Components Order Form](#).**23.1. Stealth Unilink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable & Removable Hardware \$265 ea.
Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

UniLink, 7" Height Adjustable & Removable Hardware \$276 ea.
Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$281 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact..... \$281 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Left..... \$332 ea.
Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Right \$332 ea.
Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
Part: ST-UTT8A116O2-L-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea.
Part: ST-UTT8A116O2-R-Q. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.
Part: ST-UTP544056100-Q. HCPCS: E0956/E0953^

UniLink, 4" x 8" Pad..... \$106 ea.
Part: ST-UTP548056100-Q. HCPCS: E0956/E0953^

UniLink, 4" x 12" Pad \$106 ea.
Part: ST-UTP54D056100-Q. HCPCS: E0956/E0953^

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea.
Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$81 ea.
Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
Part: PTOINDV3585

Install on Right
Part: PTOINDV3586

Do Not Install
Part: PTOINDV3719

23.2. Stealth Thoracic Laterals

Swing-away laterals may interfere with 2-post, flip up armrests.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$664
Part: ST-ULS4042000-Q. HCPCS: E1034
Only available with Redesigned back assemblies.

Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664
Part: ST-ULPK404108B00-Q. HCPCS: E1034
Only available with Sport or TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034
Only available with back cane and omit back configurations.

Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$332
Part: ST-ULPK404103L00-Q. HCPCS: E1034
Only available with Active Contour Q-Backs.

Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$332
Part: ST-ULPK404103R00-Q. HCPCS: E1034
Only available with Active Contour Q-Backs.

Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**2. Pads**

Curved Pads, 3"x5", pair (TWBL3X5)..... \$250
 Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL4X6)..... \$250
 Part: INDPART2664. HCPCS: E0956

3. Installation

Install Laterals
 Part: PTOINDV3499

Do Not Install
 Part: PTOINDV3500

23.4. Chest Bar Assembly

Chest support is required with anterior tilt 20 and 30 degree functions. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.

Small Chest Bar, 14" Foam Pad (ST-CBATB4ATT-S-Q)..... \$368
 Part: ST-CBA-TB4-S. HCPCS: K0108
 Available with 14" to 20" back widths.

Medium Chest Bar, 16" Foam Pad (ST-CBATB4ATT-M-Q)..... \$368
 Part: ST-CBA-TB4-M. HCPCS: K0108
 Available with 17" to 22" back widths.

23.5. Stealth Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea.
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$154 ea.
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$154 ea.
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$154 ea.
 Part: ST-IFPS2104-Q. HCPCS: E0960

23.6. Stealth Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor, (IFCS4113)..... \$97 ea.
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor, (IFCS4114)..... \$97 ea.
 Part: ST-IFCS4114-Q. HCPCS: E0960

23.7. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right..... \$318
 Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left \$318
 Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492
 2 Upper & 1 Lower Right
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492
 2 Upper & 1 Lower Left
 Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$434
 2 Upper
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028

2. Facial Lateral Pads and Accessories

Egg Switch Mount (ESM) (A) \$57 ea.
 Part: INDPART2815

Facial Pad, 4".5x2", Cool Core (900M-Q) \$61 ea.
 Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2", Cool Core (920M-Q)..... \$61 ea
 Part: ST-920M-Q. HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C)..... \$276 ea.
 Part: ST-SUS18-Q
 Allows mounting of facial pads and attaches to the headrest




Contoured and Classic Chest Harness



Chest Strap



Left and Right Facial Lateral with Pads



Headrest Accessories



Elbow Supports

23.8. Stealth Elbow Support**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support \$199 ea.
 Bracket Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket..... \$211
 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
 Part: ST-UETB4021FL-Q

Right, Swing-Away, Elbow Support..... \$199 ea.
 Bracket Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

Right, UniLink, Fixed, Elbow Support Bracket..... \$211
 for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
 Part: ST-UETB4021FR-Q

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
 Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.
 Part: POS124067. HCPCS: K0108

23.9. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

23.10. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
 Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
 Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
 Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Positioning components section continued on following page.

24. LAP TRAY**24.1. Quantum Lap Tray**

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

Black Lap Tray, 11" W \$200
 Part: MEC1905328-11B. HCPCS: E0950
 For use with 12"-13" back width.

Black Lap Tray, 13" W \$200
 Part: MEC1905328-13B. HCPCS: E0950
 For use with 14"-15" back width.

Black Lap Tray, 15" W \$200
 Part: MEC1905328-15B. HCPCS: E0950
 For use with 16"-17" back width.

Black Lap Tray, 17" W \$200
 Part: MEC1905328-17B. HCPCS: E0950
 For use with 18"-19" back width.

Black Lap Tray, 19" W \$200
 Part: MEC1905328-19B. HCPCS: E0950
 For use with 20"-22" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$200
 Part: MEC1905328-11C. HCPCS: E0950
 For use with 12"-13" back width.

Clear Lap Tray, 13" W \$200
 Part: MEC1905328-13C. HCPCS: E0950
 For use with 14"-15" back width.

Clear Lap Tray, 15" W \$200
 Part: MEC1905328-15C. HCPCS: E0950
 For use with 16"-17" back width.

Clear Lap Tray, 17" W \$200
 Part: MEC1905328-17C. HCPCS: E0950
 For use with 18"-19" back width.

Clear Lap Tray, 19" W \$200
 Part: MEC1905328-19C. HCPCS: E0950
 For use with 20"-22" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray \$75
 Part: MEC2104886

24.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
 Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
 Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	S	M	M	M	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	S	S	M	M	L	L	
Drop-In, T-Arm	S	S	S	S	S	S	S	M	M	L	L	
Two-Post, Flip-Back	S	S	S	S	S	S	S	M	M	L	L	
Cantilever	S	S	S	S	S	M	M	M	M	L	L	
Q-Captain, Flip-Up						S	M	M	L	L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	S	S	M	M	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

S:Compatible

M:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

25. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**25.1. Enhanced Display**Q-Logic 3 Enhanced Display Kit.....\$1,980
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard Part: DGN137728	Deep Purple \$63 Part: DGN137729-PPR5	Electric Blue \$63 Part: DGN137729-PBU9	Military Digital Camo \$63 Part: DGN137729-WCM7
Red Corvette \$63 Part: DGN137729-PRD1	Raspberry Beret \$63 Part: DGN137729-PPK4	Blue Jean \$63 Part: DGN137729-PBU11	Faded Glory \$63 Part: DGN137729-WFL2
Back in Black \$63 Part: DGN137729-PBK1	Pretty in Pink \$63 Part: DGN137729-PPK3	Lemon Crush \$63 Part: DGN137729-PYL4	Carbon Fiber \$63 Part: DGN137729-WCB1
Night Train \$63 Part: DGN137729-PBK3	Evolution Orange \$63 Part: DGN137729-POR7	Green Machine \$63 Part: DGN137729-PGN11	
Grey Street \$63 Part: DGN137729-PSL3	Midnight Blue \$63 Part: DGN137729-PBU10	Military Green \$63 Part: DGN137729-PGN8	



Enhanced Display

25.2. SCIM Module**Optional.**SCIM Module for Q-Logic 3 \$525
Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

25.3. Q-Logic 3 SAJQ-Logic 3 EX SAJ \$840
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

25.4. Head ArrayAdult, i-Drive Tri-Array Series, 3-Switch \$4,774
Head Array (IDHBT200-1)

Part: ST-IDHBT200-1-Q. HCPCS: E2330

Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

25.5. Sip and Puff Drive**1. Sip and Puff Module**Sip and Puff Module for Q-Logic 3 \$1,990
Part: CTL168056. HCPCS: E2325

Includes module and mounting hardware.

2. Sip and Puff Straw and Mount**The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.**Stealth Sip and Puff for Comfort Series, \$318
Left Mount (SP625-CL)

Part: ST-SP625-CL-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, \$318
Head Array (SP625-U)

Part: ST-SP625-U-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series, \$318
Right Mount (SP622-CR)

Part: ST-SP625-CR-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

3. Sip & Puff AccessoriesSip & Puff Straws, 12 pack (SP625-S12) \$35
Part: ST-SP625-S12-Q

Sip & Puff Module



Stealth Sip & Puff

25.6. Egg SwitchesRed Egg Switch (ESRED) \$87
Part: SWTMCRO1077Blue Egg Switch (ESBLU) \$87
Part: SWTMCRO1079Black Egg Switch (ESBLK) \$87
Part: SWTMCRO1081Green Egg Switch (ESGRN) \$87
Part: SWTMCRO1078Yellow Egg Switch (ESYEL) \$87
Part: SWTMCRO1080

Egg Switch

Alternative controls and electronics section continued on following page.

25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

25.7. Switch Adapter

2 to 1 Switch Adapter \$15

Part: DWR1361H037

Connects 2 single switches to a stereo jack to act as a dual switch.

25.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520

Part: CTL168057

25.9. USB Mobile Device Charger

USB Mobile Device Charger \$184

Part: ELE168069

Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

