



Q6 Edge® 4 Order Form with Contoured Seating

Account Number: _____ Date: _____
 Provider Name: _____
 Contact: _____
 Phone: _____ Fax: _____
 Email: _____
 PO Number: _____
 Quote Name/Marked For: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
 B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
 C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and off board charger. XLR USB charger is included. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 4 3S-C \$8,870
 Part: Q6 EDGE 4 3S-C. HCPCS: K0849, Pending PDAC Code Verification
 Contoured seat (captain) standard.

Quantum Q6 Edge 4 3S-SS \$8,870
 Part: Q6 EDGE 4 3S-SS. HCPCS: K0848, Pending PDAC Code Verification
 Solid seat pan standard with no power option.

Quantum Q6 Edge 4 3SP-C \$8,559
 Part: Q6 EDGE 4 3SP-C. HCPCS: K0857, Pending PDAC Code Verification
 Contoured seat (captain), single power option for 12" lift.

Quantum Q6 Edge 4 3SP-SS \$8,559
 Part: Q6 EDGE 4 3SP-SS. HCPCS: K0856
 Solid seat pan standard, single power option for 12" lift.

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 7 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology \$1,360
Coming soon. Part: XRF000268. HCPCS: K0108

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

- | | | | |
|---|--|---|---|
|  Inferno Red
Part: DGN2401066-PRD5 |  Cosmic Purple
Part: DGN2401066-PPR8 |  Steel Blue
Part: DGN2401066-PBU13 |  Lemon Crush
Part: DGN2401066-PYL4 |
|  Back in Black
Part: DGN2401066-PBK1 |  Ice Pink
Part: DGN2401066-PPK5 |  Royale Blue
Part: DGN2401066-PBU14 |  Forest Green
Part: DGN2401066-PGN13 |
|  Night Train
Part: DGN2401066-PBK3 |  Pretty in Pink
Part: DGN2401066-PPK3 |  Bermuda Blue
Part: DGN2401066-PBU15 | |
|  Ghost Grey
Part: DGN2401066-PGY6 |  Evolution Orange
Part: DGN2401066-POR7 |  Blue Jean
Part: DGN2401066-PBU11 | |

3.1b. Hydrographic Shrouds

- | | | |
|---|---|--|
|  Military Digital Camo \$263
Part: DGN2401066-WCM7 |  Faded Glory \$263
Part: DGN2401066-WFL2 |  Black Camouflage \$263
Part: DGN2401066-WCM9-PCL2
<i>Coming soon.</i> |
|---|---|--|

4. POWER ADJUSTABLE SEAT HEIGHT

Optional. The seat to floor height will be increased. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

12" Power Lift for Captain Seat with iLevel..... \$3,725
 Part: PTOINDV3857. HCPCS: E2298/E2300^
 Available with swing away legrests or manually adjustable center mount foot platform.

5. BASE MOUNTED LIGHTS

Optional. Select one option.

Accent Light Package Standard Part: ELE2404196	Full Light Package..... \$1,045 Part: ELE2404199 Requires QLNES 10 Key, Q-Logic 3 EX Joystick for Full Lights, or Q-Logic 3e Joystick.
---	--

6. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
 Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Captain's Headrest Bracket \$21 Part: MEC1901231	Backup Camera Track Mount Bracket \$27 Part: MEC2104227
---	--

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
 Part: MEC2106123



7. DRIVE WHEELS

REQUIRED. Select a type of drive wheel.

3"x14" Flat Free Tires Standard Part: PTOINDV3525	4"x14" Knobby, Flat Free Black Tires \$580 Part: XRF000368
3"x14" Pneumatic Tires No Charge Part: WHL1713925	

8. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

Group 24, Gel Cell, Interceptor 2 Battery..... \$399 ea. Part: BAT2404214. HCPCS: E2361 Battery tray configuration.	NF-22, AGM, MK Flat Top Battery \$432 ea. <i>Coming soon.</i> Part: BAT2501912. HCPCS: E2361 Without battery tray configuration
NF-22, Gel Cell, Interceptor 2 Battery \$355 ea. <i>Coming soon.</i> Part: BAT2500310. HCPCS: E2361 Without battery tray configuration.	

9. ELECTRONICS

REQUIRED. Select one type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control found on the [Alternative Controls and Electronics Order Form](#).

9.1a. Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

- | | |
|---|--|
| QLNES 6 Key JoystickStandard
Part: CTL2205379
Controls up to 2 actuators. | QLNES 10 Key Joystick\$420
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package. |
|---|--|



9.1b. Expandable Controls

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

- | | |
|--|--|
| Q-Logic 3 EX Controller Module\$730
Part: XRF000002. HCPCS: E2377 | Q-Logic 3 Expandable Controller Harness.....\$458
Part: XRF000264. HCPCS: E2313 |
|--|--|

2. Q-Logic 3 Joystick

- | | |
|---|---|
| Q-Logic 3 EX Joystick\$840
Part: CTL167724 | Q-Logic 3 EX Joystick for Full Lights\$840
Part: CTL167725 |
| Q-Logic 3e Joystick\$735
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls. | No Joystick, Setup for Q-Logic 3
Part: XRF000010
Requires selection of Enhanced Display |



9.2. Q-Logic Magnetic Charger

Available with QLNES. Included with Q-Logic 3e Joystick.

- | |
|--|
| Magnetic to XLR Charger Adapter with Wall Holder\$60
Part: ELE2103635 |
|--|

9.3. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

- | | | |
|---|---|---|
| Black, Soft TouchStandard
Part: DGN124307
Also matches Night Train shrouds. | Pretty in Pink\$63
Part: DGN124336-PPK3 | Lemon Crush.....\$63
Part: DGN124336-PYL4 |
| Inferno Red.....\$63
Part: DGN124336-PRD5 | Evolution Orange.....\$63
Part: DGN124336-POR7 | Forest Green\$63
Part: DGN124336-PGN13 |
| Back in Black.....\$63
Part: DGN124336-PBK1 | Steel Blue\$63
Part: DGN124336-PBU3 | Military Digital Camo\$63
Part: DGN124336-WCM7 |
| Ghost Grey\$63
Part: DGN124336-PGY6 | Royale Blue.....\$63
Part: DGN124336-PBU14 | Faded Glory\$63
Part: DGN124336-WFL2 |
| Cosmic Purple\$63
Part: DGN124336-PPR8 | Bermuda Blue.....\$63
Part: DGN124336-PBU15 | Black Camouflage.....\$63
Part: DGN124336-WCM9-PCL2 |
| Ice Pink.....\$63
Part: DGN124336-PPK5 | Blue Jean\$63
Part: DGN124336-PBU11 | |

10. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

- | | |
|--|--|
| <input type="checkbox"/> Mushroom Handle (B)\$90
Part: ACCASMB1860. HCPCS: E2323 | <input type="checkbox"/> Extended I-Handle (G)\$90
Part: ACCASMB1857. HCPCS: E2323 |
| <input type="checkbox"/> Chin Cup (C)\$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control) | <input type="checkbox"/> Bodypoint 3" Goal Post Handle\$104
Part: PLSKNCP1057. HCPCS: E2323 |
| <input type="checkbox"/> T-Bar (D)\$126
Part: ACCASMB1858. HCPCS: E2323 | <input type="checkbox"/> Bodypoint 4" Goal Post Handle\$104
Part: PLSKNCP1058. HCPCS: E2323 |
| <input type="checkbox"/> Medium Ball Handle (E)\$69
Part: ACCASMB1859. HCPCS: E2323 | <input type="checkbox"/> Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308 |
| <input type="checkbox"/> Soft Ball Handle (F)\$80
Part: ACCASMB1856. HCPCS: E2323 | |



11. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

11.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket for Q-Logic 3 Standard
 Part: FRMASMB10617

11.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
 Part: ACC1908434. HCPCS: E1032



11.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side \$245
 Part: ACC2201041. HCPCS: E1032


Basic Swing Away, Q-Logic, Right Side \$245
 Part: ACC2201042. HCPCS: E1032

11.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount \$368
 Part: XRF000003. HCPCS: E1032
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount \$126
 Part: XRF000004
 Allows height adjustment of joystick.



11.2. Joystick Mount Location

Left Side Mount
 Part: PTOINDV1225

Right Side Mount
 Part: PTOINDV1226

12. CONTOURED SEATING

Optional. Select only one seating type.

12.1. Occupied Transit

Optional.

Occupied Transit Option for Contoured Seat \$368
 Part: KIT122609
 ISO standard occupied transit.

12.2a. Q-Captain Seat

16"W x 16"D No Charge Part: SET2005145	18"W x 16"D No Charge Part: SET2005149	20"W x 16"D No Charge Part: SET2005152
16"W x 18"D No Charge Part: SET2005146	18"W x 18"D No Charge Part: SET2005150	20"W x 18"D No Charge Part: SET2005153
16"W x 20"D No Charge Part: SET2005147	18"W x 20"D No Charge Part: SET2005151	20"W x 20"D No Charge Part: SET2005154

12.2b. Q-Captain Seat with Solid Seat Pan

16"W x 16"D No Charge Part: SET2005394	18"W x 16"D No Charge Part: SET2005397	20"W x 16"D No Charge Part: SET2005400
16"W x 18"D No Charge Part: SET2005395	18"W x 18"D No Charge Part: SET2005398	20"W x 18"D No Charge Part: SET2005401
16"W x 20"D No Charge Part: SET2005396	18"W x 20"D No Charge Part: SET2005399	20"W x 20"D No Charge Part: SET2005402

12.3. Seat to Floor Height

REQUIRED.

Lowest Seat to Floor
 Part: PTOINDV3558

Middle Seat to Floor
 Part: PTOINDV3559

Highest Seat to Floor
 Part: PTOINDV3560

	Cushion	SSP
High	22.5"	18.5"
Mid	21.5"	17.5"
Low	20.5"	16.5"


Q-Captain
 Seat to Floor Height

Contoured seating section continued on following pages.

12. CONTOURED SEATING CONTINUED

12.4. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.



1. Cushion Type

- TRU-Comfort 2 SPP**
Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608
- Essence SPP**
Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608
- Solution SPP**
Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608
- Spectrum Foam SPP**
Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608
- Spectrum Gel SPP**
Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608
- Premiere P**
Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606
- Glacial SP**
Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604
- Zen SP**
Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604
- Simplicity G**
Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

- Small Wedge Positioning Kit.....\$39**
Part: POP1909006-S
14"W and smaller
- Medium Wedge Positioning Kit.....\$72**
Part: POP1909006-M
14"W and smaller
- Large Wedge Positioning Kit..... \$105**
Part: POP1909006-L
14"W and smaller
- Small Wedge Positioning Kit.....\$39**
Part: POP1712340-S
15"W and larger
- Medium Wedge Positioning Kit.....\$72**
Part: POP1712340-M
15"W and larger
- Large Wedge Positioning Kit..... \$105**
Part: POP1712340-L
15"W and larger

3. Cushion Top Foam

- Visco Top Foam**
Reference part: XRF000380
- FRG Top Foam**
Reference part: XRF000381
- No Top Foam**
Reference part: XRF000382

4. Cushion Liner

- Stealth Urethane Liner**
Reference part: XRF000383
- Stealth Polyurethane Liner**
Reference part: XRF000384
- No Liner**
Reference part: XRF000385

5. Cushion Cover

- CoolCore Cover**
Reference part: XRF000386
- Mesh Cover**
Reference part: XRF000387
- Stealth Tek Cover**
Reference part: XRF000388

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

- ABS Flat Rigidizer**
Reference part: XRF000389
MSRP: \$48-\$69

	None	Visco/V/Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam	Liner	Cover						
TRU-Comfort 2	-	✓	-	✓	-	✓	✓	-	-
Essence	-	✓	-	✓	-	✓	✓	✓	-
Solution	-	✓	-	✓	-	✓	✓	✓	✓
Spectrum Foam	-	✓	✓	-	✓	✓	✓	✓	✓
Spectrum Gel	-	✓	✓	-	✓	✓	✓	✓	✓
Premiere	-	✓	✓	✓	✓	✓	✓	-	-
Glacial	✓	-	-	✓	✓	✓	✓	✓	✓
Zen	-	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	-	-	✓	-	-	-	-	✓

Contoured seating section continued on following pages.

14. CONTOURED SEATING CONTINUED

12.5. Armrests

Optional. Configured in pairs. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Q-Captain, Flip-Up, Height Adjustable Armrests

Includes full length, flat armpads.

Full Length, Left Armrest \$90 ea. Full Length, Right Armrest..... \$90 ea.
 Part: MEC2201304. HCPCS: K0020 Part: MEC2201303. HCPCS: K0020

1b. Drop-In, Height Adjustable, Armrests (T-Arm)

Includes flat armpads.

Full Length, Drop-In Armrests..... \$79 ea. Desk Length, Drop-In Armrests \$79 ea.
 L R Left part: XRF000298. Right part: XRF000299. HCPCS: E0973 L R Left part: XRF000300. Right part: XRF000301. HCPCS: E0973

1c. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

Full Length, Single Post Armrests \$79 ea. Desk Length, Single Post Armrests \$79 ea.
 L R Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973 L R Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973

1d. Bariatric, Height Adjustable, Removable Armrests

Select armrests and armpads. Armrests are sold in same-size pairs.

Full Length, Bariatric Armrests, 8"-12" H \$90 ea. Full Length, Bariatric Armrests, 10"-14"H \$90 ea.
 L R Left part: FRMASMB6728. Right part: FRMASMB6729. HCPCS: E0973 L R Left part: FRMASMB6730. Right part: FRMASMB6731. HCPCS: E0973

1e. Heavy Duty, Height Adjustable, Removable Armrests

Select armrests and armpads. Armrests are sold in same-size pairs.

Full Length, Heavy Duty Armrests, 8"-12" H \$90 ea. Full Length, Heavy Duty Armrests, 10"-14"H..... \$90 ea.
 L R Left part: XRF000253. Right part: XRF000254. HCPCS: E0973 L R Left Part: XRF000257. Right part: XRF000258. HCPCS: E0973

1f. Pediatric, Removable, Height Adjustable Armrests

Height adjustment from 7.75" to 8.25". Includes desk length, straight armpads.

Pediatric, Removable Armrests \$79 ea.
 L R Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973

12.6. Q-Captain Armrest Spacer

Optional.

Armrest Spacer \$35
 Part: XRF000370
 1" spacer for Q-Captain, Flip-Up, Height Adjustable Armrests only.

12.7. Armrest Height

REQUIRED with Bariatric, Heavy Duty, and Pediatric Height Adjustable armrests.

1. Left Armrest Height

7.75" H PTOINDV3722	8.25" H PTOINDV3724	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
8" H PTOINDV3535	9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

2. Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
8" H PTOINDV3450	9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

12.8. Armpads

REQUIRED except with Flip-Up, Drop-In, and Pediatric armrests.

1a. Straight Armpads for Single Post Armrests

Straight, Full Length Left and Right No Charge Straight, Desk Length Left, Full Length Right No Charge
 Part: ARM134275 Part: ARM134280

Straight, Full Length Left, Desk Length Right No Charge Straight, Desk Length Left and Right..... No Charge
 Part: ARM134281 Part: ARM134254

1b. Waterfall Armpads for Single Post Armrests

Waterfall, Full Length Left and Right..... No Charge Waterfall, Desk Length Left, Full Length Right No Charge
 Part: ARM134285 Part: ARM134288

Waterfall, Full Length Left, Desk Length Right No Charge Waterfall, Desk Length Left and Right No Charge
 Part: ARM134289 Part: ARM134282

1c. Falcon Gel Armpads for Single Post Armrests

Falcon Gel, Full Length Left and Right\$311.22 Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
 Part: ARM136438 Part: ARM136439

Falcon Gel, Full Length Left, Desk Length Right.....\$287.93 Falcon Gel, Desk Length Left and Right \$265.02
 Part: ARM136442 Part: ARM136437

Contoured seating section continued on following page.

12. CONTOURED SEATING CONTINUED

1d. Armpads for Bariatric or Heavy Duty Armrests

Bariatric Armpads No Charge
 Part: UPSASMB1027
 Only for use with Bariatric Armrests.

Straight Armpads..... No Charge
 Part: UPSASMB1021
 Only for use with Heavy Duty Armrests

Waterfall Armpads..... No Charge
 Part: UPSASMB1024
 Only for use with Heavy Duty Armrests

12.9. Quad Release Lever

Optional. Available with Single Post and Pediatric armrests.

Quad Release Lever, Pair.....No Charge
 Part: FRMASMB10517

15. LAP BELTS

REQUIRED. Select one lap belt below.

15.1a. Standard Lap Belts

38" Lap Belt Part: ACCBELT1006	50" Lap Belt Part: ACCBELT1000	60" Lap Belt Part: ACCBELT1002	70" Lap Belt Part: ACCBELT1003	80" Lap Belt Part: ACCBELT1004	90" Lap Belt Part: ACCBELT1001
--	--	--	--	--	--

15.1b. Optional Lap Belts

Stealth Positioning, 60" Single Pull, Push Button, Medium \$50 (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q	Stealth Positioning, 60" Single Pull, Push Button, Extra Large \$50 Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q
Stealth Positioning, 60" Single Pull, Push Button, Large \$50 Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q	Stealth Positioning, 60" Single Pull, Push Button, Non-padded ... \$50 2" Belt (IFPB5310) Part: ST-IFPB5310-Q

16. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

16.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Type

Basic Center Mount.....No Charge
 Part: XRF000323

2. Basic Center Mount Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates \$180 Part: XRF000312. HCPCS: K0040	Long Extension, 5"x9" Footplates \$180 Part: XRF000314. HCPCS: K0040
Medium Extension, 5"x9" Footplates \$180 Part: XRF000313. HCPCS: K0040	

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount.....\$158
 Part: RIG2104369

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16
 Part: RIG2104240

16.1b. Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Assembly

Power Articulating Foot Platform \$3,145 Part: PTOINDV3754. HCPCS: E1012+K0040 Only available with power lift option.	Quick Adjustable Center Mount Foot Platform \$263 Part: XRF000311. HCPCS: K0040
--	---

2. Center Mount Lower Extension and Angle Adjustable Footplates

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates Part: RIG152577	Medium Extension, 5"x11" Footplates Part: RIG2404018	Long Extension, 5"x11" Footplates Part: RIG2404019
Pediatric Extension, 5"x8" Footplates Part: RIG2404015	Long Extension, 5"x6" Footplates Part: RIG151505	
Medium Extension, 5"x8" Footplates Part: RIG2404016	Long Extension, 5"x8" Footplates Part: RIG2404017	

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard Part: RIGASMB7120018	Small Calf Pads, Quick Adjustable Center Mount Only, \$158 6"Wx6"H (pair) Part: XRF000156
---	--

Legrests seating section continued on following page.

16. LEGRESTS CONTINUED

3b. Center Mount Calf Panel

Small Calf Panel, 18"Wx8"H.....\$184
 Part: FRMASMB7110012

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69 Part: RIG152519 Only available with the 5x6 footplates.	Outside Edge \$69 Part: RIG2404921 Available with 5"x8" and 5"x11" footplates.
Width Extension for 8" \$69 Part: RIG2404919	Front and Rear Cap \$16 Part: RIG2404923 Available with 5"x8" and 5"x11" footplates.
Width Extension for 11" \$69 Part: RIG2404920	

5. Center Mount Additional Accessories

Small Foam Calf Pad Cover (pair) \$67 Part: ST-CSCS0775-Q	Footplate Riser \$75 Part: XRF000348 Not available with long extension. Select legrest length between 2" to 9" at 1" increments.
Small Foam Panel Cover (pair) \$176 Part: ST-CSCP0920-Q	Footplate Riser Calf Panel \$50 Part: FRMASMB7110011 Only available and required with Footplate Riser selection and extension lengths of 6" to 9".
Footplate Connector \$16 Part: RIGASMB7120023	

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222 Part: ST-FPG5008T-Q. HCPCS: K0108	5"x8" Tapered Gel Footplate Covers for Width Extension \$247 Part: ST-AFPG5008E-Q. HCPCS: K0108
5"x11" Tapered Gel Footplate Covers \$270 Part: ST-AFPG5011T-Q. HCPCS: K0108	5"x11" Tapered Gel Footplate Cover for Width Extension \$270 Part: ST-AFPG5011E-Q. HCPCS: K0108

16.1c. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 60°-90° Swing-Away Legrest, Static Seat \$86 ea. Pair part: XRF000082. Left: XRF000083. Right: XRF000084 Not available with 11"-13" extensions.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea. Pair part: XRF000085. Left: XRF000086. Right: XRF000087 For use only with 11"-13" lower extensions.
--	--

2. Lower Extension

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 11"-13" Length Pair part: XRF000088. Left: XRF000089. Right: XRF000090.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 13"-16.5" Length Pair part: XRF000094. Left: XRF000095. Right: XRF000096.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 15.5"-19.5" Length Pair part: XRF000097. Left: XRF000098. Right: XRF000099.
---	---	---

3a. Extruded Footplate(s)

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 7.25"Wx6"D Extruded Footplate..... Standard Part pair: FRMASMB7847. Left: FRMASMB7861 Right: FRMASMB7862	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 8.5"Wx6"D Extruded Footplate Standard Part pair: FRMASMB7848. Left: FRMASMB7863 Right: FRMASMB7864 May cause interference on seat widths less than 18."
---	---

3b. Swing-Away Legrest, Angle Adjustable Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 4"Wx6"D \$90 ea. Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560 HCPCS: K0040+E0951	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 6"Wx11"D ¹ \$90 ea. Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564 HCPCS: K0040+E0951
<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 5"Wx6"D..... \$90 ea. Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561 HCPCS: K0040+E0951	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 8"Wx11"D ² \$90 ea. Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565 HCPCS: K0040+E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

<input type="checkbox"/> QTY Therafin 31515 Multi-Axis Amputee Support \$230 ea. Part: INDPART2822. HCPCS: E1020	<input type="checkbox"/> QTY Therafin 32387 Amputee Support Gel Cover \$155 ea. Part: INDPART2823. HCPCS: K0108
---	--

Legrests seating section continued on following page.

16. LEGRESTS CONTINUED

16.1d. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

P L R Elevating Legrest Assembly \$171 ea.
 Pair part: FRMASMB483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

P L R **11"-14" Length**
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

P L R **8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length**
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

P L R **14"-17" Length**
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

P L R **13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length**
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

P L R **7.25"Wx6"D Reinforced Footplate**..... Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P L R **8.5"Wx6"D Reinforced Footplate**..... Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

P L R **4"Wx6"D** \$90 ea.
 Pair parts: FRMASMB8812 x2.
 HCPCS: K0040

P L R **6"Wx8"D¹ with Heel Loops**..... \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

P L R **5"Wx6"D**..... \$90 ea.
 Pair parts: FRMASMB8811 x2.
 HCPCS: K0040

P L R **6"Wx11"D¹ with Heel Loops**..... \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

P L R **5"Wx8"D with Heel Loops**..... \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

P L R **8"Wx11"D² with Heel Loops**..... \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

16.2. Legrest Accessories

Heel Loops (pair)..... \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$19 ea.
 Part: XRF000100

16.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	4" L (#3991)	6" L (#3993)	8" L (#3680)	9" L (#3682)	10" L (#3684)	11" L (#3686)	12" L (#3688)	13" L (#3690)	14" L (#3692)	15" L (#3694)
3" L (#3990)	5" L (#3992)	7" L (#3994)	8.5" L (#3681)	9.5" L (#3683)	10.5" L (#3685)	11.5" L (#3687)	12.5" L (#3689)	13.5" L (#3691)	14.5" L (#3693)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)
7.5" L (#3471)	8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	

17. POWER ELECTRONICS

REQUIRED with actuator seating.

Single Actuator Through Controller and Interface \$2,163
 Part: XRF000335. HCPCS: E2310

Multiple Actuator Through Controller and Interface.....\$3,040
 Part: XRF000336. HCPCS: E2311

Single Actuator Through Controller and Interface (AAM).....\$3,040
 Part: XRF000068. HCPCS: E2310

19. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Cup Holder, Q-Captain No Charge Part: ACC2006355	Grip Lap Board (20005)..... \$38.95 Part: ACC152992
Front Mounted Cell Phone/ Media Holder, X-Grip, Q-Captain..... \$74 Part: ACC2103691	Grip Drink Holder (20085) \$15.50 Part: ACC152993
Rear View Mirror, Q-Captain \$69 Part: ACC2103692	Large, Black, Weather Cover \$202 Part: ACCCOVR1021
Oxygen Holder..... \$118 Part: ACC2006353. HCPCS: E2208 Not available with Power Lift option.	Reflector Kit \$53 Part: ACC124670
Oxygen Holder..... \$180 Part: ACC2400456. HCPCS: E2208 Select with Power Lift option.	Clothing Guard (Pair)..... \$237 Part: ACC125007 Not available with Heavy Duty armrests.
Cane/Crutch Holder \$110 Part: ACC2006352. HCPCS: E2207 Not available with Power Lift option.	Transfer Handles \$105 ea. Part: MEC2205680. HCPCS: E0705 Not available with Heavy Duty armrests.
Cane/Crutch Holder \$250 Part: ACC2400455. HCPCS: E2207 Select with Power Lift option.	Privacy Bag \$84 Part: ACC2005723
Rear Basket \$118 Part: ACC2006351 Not available with Power Lift option.	Personal Item Pouch, Rail Mounted..... \$84 Part: ACC2101017
Rear Basket \$124 Part: ACC2400454 Select with Power Lift option.	Personal Item Pouch, Arm Mounted \$84 Part: ACC2103694
Saddle Bag..... \$40 Part: DWR1020D010	Glove Box, Rail Mounted \$158 Part: ACC125011
Quantum Backpack \$168 Part: ACC1905853	Glove Box, Arm Mounted..... \$158 Part: ACC2103693
	Remote Stop Switch \$1,295 Part: XRF000379

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea. QTY Part: ACC2407046	Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea. QTY Part: ACC2405496
--	---

20. BRANCH PATCHES

Optional.

Air Force Patch \$21 Part: ACC140240	Coast Guard Patch \$21 Part: ACC140239	Navy Patch \$21 Part: ACC140244
Army Patch \$21 Part: ACC140241	Marines Patch \$21 Part: ACC140243	Purple Heart Patch \$21 Part: ACC140238



21. POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.

21.1. Stealth Unilink Pelvic/Thigh Guide

Available with Flip-Up armrests only. Not available with Elevating Legrests.

1. UniLink Pelvic/Thigh Hardware

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea. Removable Mounting Hardware QTY Part: ST-UTTP5-Q. HCPCS: E1028/E1034^	Flip Down with LEFT Full Surface Contact \$281 ea. Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028/E1034^
UniLink, 7" Height Adjustable, \$276 ea. Removable Mounting Hardware QTY Part: ST-UTTP7-Q. HCPCS: E1028/E1034^	Flip Down with RIGHT Full Surface Contact..... \$281 ea. Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea. QTY Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^	UniLink, 4" x 12" Pad \$106 ea. QTY Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^
UniLink, 4" x 8" Pad..... \$106 ea. QTY Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^	

Positioning components section continued on following page.

21. POSITIONING COMPONENTS CONTINUED

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea. UniLink, 2" Link..... \$81 ea.
QTY Part: ST-UIL001150A000-Q QTY Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left Install on Right Do Not Install
 Part: PTOINDV3585 Part: PTOINDV3586 Part: PTOINDV3719

21.2. Stealth Elbow Support

1. Elbow Support Bracket

Left, Swing-Away, Elbow Support Bracket, Tube \$199 ea. Right, Swing-Away, Elbow Support Bracket, Tube..... \$199 ea.
 Style Arms (ARMS-ES-L) Style Arms (ARMS-ES-R)
 Part: POS124064. HCPCS: E1028 Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea. 4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.
QTY Part: POS124065. HCPCS: K0108 QTY Part: POS124067. HCPCS: K0108

21.3. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539 Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
 (SPLE107) (pair) (SPLE109) (pair)
 Part: INDPART2786. HCPCS: E0954 Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

21.4. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea. Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
QTY Part: ST-IFSH4211-Q. HCPCS: E0951+E0952 QTY Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
QTY Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

22. LAP TRAY

22.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309 Simple Solutions Large Tray with Cut Out..... \$309
 Part: ST-SSLXS010001171-Q. HCPCS: E0950 Part: ST-SSLXL010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out..... \$309
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware..... \$302
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width												
	10	11	12	13	14	16	17	18	19	20	21	22	
Flip-Back	S	S	S	S	S	S	S	M	M	M	L	L	
Flip-Back, Cane Mounted	S	S	S	S	S	S	M	M	M	L	L		
Drop-In, T-Arm	S	S	S	S	S	S	M	M	M	L	L		
Two-Post, Flip-Back	S	S	S	S	S	S	M	M	M	L	L		
Cantilever	S	S	S	S	S	M	M	M	M	L	L		
Q-Captain, Flip-Up						S	M			L			
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	M	S	M	L	M	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q
■:Compatible ■:Armrest configuration may not hold tray position
 *Some armrests listed may not be available on this order form.