



Quantum Rehab®
 401 York Ave., Duryea, PA 18642
 Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 4 Order Form TRU Balance® 3 Positioning

Account Number: _____ Date: _____
 Provider Name: _____
 Contact: _____
 Phone: _____ Fax: _____
 Email: _____
 PO Number: _____
 Quote Name/Marked For: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
 B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
 C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 7mph programmable motors and off board charger. XLR USB charger is included. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 4 3SP-SS\$8,559
 Part: Q6 EDGE 4 3SP-SS. HCPCS: K0856
 Single power option.

Quantum Q6 Edge 4 3MP-SS\$8,373
 Part: Q6 EDGE 4 3MP-SS. HCPCS: K0861
 Multiple actuator power option.

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 7 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology \$1,360
Coming soon. Part: XRF000268. HCPCS: K0108

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

- | | | | |
|---|--|---|---|
|  Inferno Red
Part: DGN2401066-PRD5 |  Cosmic Purple
Part: DGN2401066-PPR8 |  Steel Blue
Part: DGN2401066-PBU13 |  Lemon Crush
Part: DGN2401066-PYL4 |
|  Back in Black
Part: DGN2401066-PBK1 |  Ice Pink
Part: DGN2401066-PPK5 |  Royale Blue
Part: DGN2401066-PBU14 |  Forest Green
Part: DGN2401066-PGN13 |
|  Night Train
Part: DGN2401066-PBK3 |  Pretty in Pink
Part: DGN2401066-PPK3 |  Bermuda Blue
Part: DGN2401066-PBU15 | |
|  Ghost Grey
Part: DGN2401066-PGY6 |  Evolution Orange
Part: DGN2401066-POR7 |  Blue Jean
Part: DGN2401066-PBU11 | |

3.1b. Hydrographic Shrouds

- | | | |
|---|---|--|
|  Military Digital Camo \$263
Part: DGN2401066-WCM7 |  Faded Glory \$263
Part: DGN2401066-WFL2 |  Black Camouflage \$263
Part: DGN2401066-WCM9-PCL2
<i>Coming soon.</i> |
|---|---|--|

4. BASE MOUNTED LIGHTS

Optional. Select one option.

Accent Light Package Standard
Part: ELE2404196

Full Light Package..... \$1,045
Part: ELE2404199
Requires QLNES 10 Key, Q-Logic 3 EX Joystick for Full Lights, or Q-Logic 3e Joystick.

5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket..... \$21
Part: MEC1910373
Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$27
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



6. DRIVE WHEELS

REQUIRED. Select a type of drive wheel.

3"x14" Flat Free Tires Standard
Part: PTOINDV3525

4"x14" Knobby, Flat Free Black Tires \$580
Part: XRF000368

3"x14" Pneumatic Tires No Charge
Part: WHL1713925

7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered. Uses 8 AMP charger.

Group 24, Gel Cell, Interceptor 2 Battery..... \$399 ea.
Part: BAT2404214. HCPCS: E2361
Battery tray configuration.

NF-22, AGM, MK Flat Top Battery \$432 ea.
Coming soon. Part: BAT2501912. HCPCS: E2361
Without battery tray configuration

NF-22, Gel Cell, Interceptor 2 Battery \$355 ea.
Coming soon. Part: BAT2500310. HCPCS: E2361
Without battery tray configuration.

10. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

10.1a. Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

- | | |
|---|---|
| QLNES 6 Key JoystickStandard
Part: CTL2205379
Controls up to 2 actuators. | QLNES 10 Key Joystick\$420
Part: CTL2205380
Controls up to 2 actuators through joystick and full light package. |
|---|---|



10.1b. Expandable Controls

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

- | | |
|--|--|
| Q-Logic 3 EX Controller Module\$730
Part: XRF000002. HCPCS: E2377 | Q-Logic 3 Expandable Controller Harness.....\$458
Part: XRF000264. HCPCS: E2313 |
|--|--|

2. Q-Logic 3 Joystick

- | | |
|---|---|
| Q-Logic 3 EX Joystick\$840
Part: CTL167724 | Q-Logic 3 EX Joystick for Full Lights\$840
Part: CTL167725 |
| Q-Logic 3e Joystick\$735
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls. | No Joystick, Setup for Q-Logic 3
Part: XRF000010 |



10.2. Q-Logic Magnetic Charger

Available with QLNES. Included with Q-Logic 3e Joystick.

- | |
|--|
| Magnetic to XLR Charger Adapter with Wall Holder\$60
Part: ELE2103635 |
|--|

10.3. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

- | | | |
|---|--|--|
| Black, Soft TouchStandard
Part: DGN124307
Also matches Night Train shrouds. | Pretty in Pink \$63
Part: DGN124336-PPK3 | Lemon Crush \$63
Part: DGN124336-PYL4 |
| Inferno Red..... \$63
Part: DGN124336-PRD5 | Evolution Orange..... \$63
Part: DGN124336-POR7 | Forest Green \$63
Part: DGN124336-PGN13 |
| Back in Black..... \$63
Part: DGN124336-PBK1 | Steel Blue \$63
Part: DGN124336-PBU13 | Military Digital Camo \$63
Part: DGN124336-WCM7 |
| Ghost Grey \$63
Part: DGN124336-PGY6 | Royale Blue..... \$63
Part: DGN124336-PBU14 | Faded Glory \$63
Part: DGN124336-WFL2 |
| Cosmic Purple \$63
Part: DGN124336-PPR8 | Bermuda Blue..... \$63
Part: DGN124336-PBU15 | Black Camouflage \$63
Part: DGN124336-WCM9-PCL2 |
| Ice Pink..... \$63
Part: DGN124336-PPK5 | Blue Jean \$63
Part: DGN124336-PBU11 | |

11. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

- | | |
|--|--|
| <input type="checkbox"/> Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323 | <input type="checkbox"/> Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323 |
| <input type="checkbox"/> Chin Cup (C)..... \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control) | <input type="checkbox"/> Bodypoint 3" Goal Post Handle\$104
Part: PLSKNCP1057. HCPCS: E2323 |
| <input type="checkbox"/> T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323 | <input type="checkbox"/> Bodypoint 4" Goal Post Handle\$104
Part: PLSKNCP1058. HCPCS: E2323 |
| <input type="checkbox"/> Medium Ball Handle (E)..... \$69
Part: ACCASMB1859. HCPCS: E2323 | <input type="checkbox"/> Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308 |
| <input type="checkbox"/> Soft Ball Handle (F)..... \$80
Part: ACCASMB1856. HCPCS: E2323 | |



12. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

12.1a. Fixed, In-line Joystick Bracket


Fixed, In-line, Bracket for Q-Logic 3Standard
Part: FRMASMB10617



12.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1032



12.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side\$245
Part: ACC2201041. HCPCS: E1032

Basic Swing Away, Q-Logic, Right Side\$245
Part: ACC2201042. HCPCS: E1032

12.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount\$368
Part: XRF000003. HCPCS: E1032
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount\$126
Part: XRF000004
Allows height adjustment of joystick.



12.2. Joystick Mount Location

Left Side Mount
Part: PTOINDV1225

Right Side Mount
Part: PTOINDV1226

13. TRU BALANCE® 3 SEATING

REQUIRED. Maximum seat size is 22"W x 22"D.

13.1. TRU Balance 3 Configuration

TRU Balance 3 Static Seat.....\$2,410
Part: PTOINDV3442
Only available with Power Lift. 20° of limited recline adjustment (-5/+15) via gas strut within back cover.

TRU Balance 3 Power Tilt Only\$8,202
Part: PTOINDV3443. HCPCS: E1002
50° Center of Gravity (CG) power tilt seating system.

TRU Balance 3 Power Recline Only\$8,234
Part: PTOINDV3444. HCPCS: E1004
Only available with Power Lift. 175° of power recline.

TRU Balance 3 Power Tilt & Recline\$12,615
Part: PTOINDV3446. HCPCS: E1007
50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

13.2. 12" Power Lift

Optional. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

TRU Balance 3 12" Power Lift with iLevel® \$3,725
Part: PTOINDV3773. HCPCS: E2298/E2300^
12" of scissor lift with iLevel System.



13.3. Occupied Transit Kit

Optional. The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

Occupied Transit Option \$289
Part: KIT152047
ISO standard occupied transit.

13.4. Seat Width

REQUIRED with seat selection. Maximum seat width is 20". Part number is SETWDTH plus number following # below. 12" Width is not available with 16" and 17" Depth. 13" and 15" widths are only available with non-matching back width.

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

13.5. Seat Width Package

REQUIRED with seat selection.

12"-16" Width Package Part: XRF000006	14"-16" Width Package Part: XRF000007	16"-20" Width Package Part: XRF000008	20"-24" Width Package Part: XRF000158 Maximum width is 22".
--	--	--	---

13.6. Seat Depth

REQUIRED with seat selection. Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

12" D (#1002)	14" D (#1004)	16" D (#1006)	18" D (#1008)	20" D (#1010)	22" D (#1012)
13" D (#1003)	15" D (#1005)	17" D (#1007)	19" D (#1009)	21" D (#1011)	

Low	Med	High
Static Seat/ Recline Only		
16"	17"	18"
Power Tilt/Tilt & Recline		
17.5"	18.5"	19.5"
Any Power Lift		
18"	19"	20"
TRU Balance 3 Seat to Floor Heights		

13.7. Seat to Floor Height

REQUIRED.


Lowest Seat to Floor Part: PTOINDV3558	Middle Seat to Floor Part: PTOINDV3559	Highest Seat to Floor Part: PTOINDV3560
---	---	--

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED

13.8. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.



1. Cushion Type

- | | |
|--|--|
| <p>TRU-Comfort 2 SPP
Reference part: XRF000116. See cushion order form.
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608</p> <p>Essence SPP
Reference part: XRF000118. See cushion order form.
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608</p> <p>Solution SPP
Reference part: XRF000113. See cushion order form.
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608</p> <p>Spectrum Foam SPP
Reference part: XRF000112. See cushion order form.
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608</p> <p>Spectrum Gel SPP
Reference part: XRF000117. See cushion order form.
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608</p> | <p>Premiere P
Reference part: XRF000122. See cushion order form.
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606</p> <p>Glacial SP
Reference part: XRF000114. See cushion order form.
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604</p> <p>Zen SP
Reference part: XRF000115. See cushion order form.
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604</p> <p>Simplicity G
Reference part: XRF000111. See cushion order form.
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602</p> |
|--|--|

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

- | | | |
|--|--|--|
| <p>Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller</p> <p>Small Wedge Positioning Kit\$39
Part: POP1712340-S
15"W and larger</p> | <p>Medium Wedge Positioning Kit.....\$72
Part: POP1909006-M
14"W and smaller</p> <p>Medium Wedge Positioning Kit.....\$72
Part: POP1712340-M
15"W and larger</p> | <p>Large Wedge Positioning Kit..... \$105
Part: POP1909006-L
14"W and smaller</p> <p>Large Wedge Positioning Kit..... \$105
Part: POP1712340-L
15"W and larger</p> |
|--|--|--|

3. Cushion Top Foam

- | | |
|---|---|
| <p>Visco Top Foam
Reference part: XRF000380</p> <p>FRG Top Foam
Reference part: XRF000381</p> | <p>No Top Foam
Reference part: XRF000382</p> |
|---|---|

4. Cushion Liner

- | | |
|---|--|
| <p>Stealth Urethane Liner
Reference part: XRF000383</p> <p>Stealth Polyurethane Liner
Reference part: XRF000384</p> | <p>No Liner
Reference part: XRF000385</p> |
|---|--|

5. Cushion Cover

- | | |
|---|---|
| <p>CoolCore Cover
Reference part: XRF000386</p> <p>Mesh Cover
Reference part: XRF000387</p> | <p>Stealth Tek Cover
Reference part: XRF000388</p> |
|---|---|

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam	Liner		Cover					
TRU-Comfort 2	-	✓	-	✓	-	✓	✓	-	-
Essence	-	✓	-	✓	-	✓	✓	✓	-
Solution	-	✓	-	-	✓	✓	✓	✓	✓
Spectrum Foam	-	✓	✓	-	✓	✓	✓	✓	✓
Spectrum Gel	-	✓	✓	-	✓	✓	✓	✓	✓
Premiere	-	✓	✓	✓	✓	✓	✓	-	-
Glacial	✓	-	-	✓	✓	✓	✓	✓	✓
Zen	-	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	-	-	✓	-	-	-	-	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

- ABS Flat Rigidizer**
Reference part: XRF000389
MSRP: \$48-\$69

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390
QTY _____

Additional Stealth Polyurethane Liner
Reference Part: XRF000391
QTY _____

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392
QTY _____

Additional Stealth Tek Cover
Reference Part: XRF000394
QTY _____

Additional Mesh Cover
Reference Part: XRF000393
QTY _____

13.9. Back Type Selection

REQUIRED. Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back
Part: XRF000014
12" and 24" back width will use the TRU Balance 3 Back assembly. Select a back cushion configuration in the subsequent section.

Redesigned, Flip Forward Back \$184
Part: XRF000016
12" and 24" back width will use the TRU Balance 3 Back assembly. Available with TRU Balance 3 Static Seat and Power Tilt (1 actuator). Only available with 2-post or pediatric armrests. Not available with recline configurations. Select a back cushion type in the subsequent section.

Back Cane Option
Part: XRF000015
Available only with TRU Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below. Widths less than 15" must select Omit Back Pan/Cushion.

13.10. Back Cane Selection (Back Cane Option Only)

REQUIRED when Back Cane Option is selected above. Select type and height of back canes.

1. Back Cane Type

Angle Adjustable Back Cane No Charge
Part: MEC140266
Flip Forward Back Cane \$53
Part: XRF000273

2" Endomorph (1" per side) Back Cane \$74
Part: XRF000269
4" Endomorph (2" per side) Back Cane \$74
Part: XRF000271

2. Back Cane Height

16"H Cane
Part: FRMASMB15834
17"H Cane
Part: FRMASMB15835

20"H Cane
Part: FRMASMB15838
21"H Cane
Part: FRMASMB15839

24"H Cane
Part: FRMASMB15842
25"H Cane
Part: FRMASMB15843

13.11. Back Cushion Configuration

REQUIRED. Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

1a. Sport Back/Omit/After-Market Interface

Sport Back No Charge
Part: XRF000017
Complete following width package and size.
Omit Back Pan/Cushion No Charge
Part: PTOINDV1287
Select a following width package.

After-market Back Interface \$535
Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.

1b. TRU-Comfort Back

TRU-Comfort Back (Static/Flip Forward), 12" Width \$744
Part: SET123622. HCPCS: E2620
Complete following width package and size.
TRU-Comfort Back (Static/Flip Forward), 14" to 20" Width \$744
Part: TC4V10W21W. HCPCS: E2620
Complete following width package and size.
TRU-Comfort Back, Back Cane Option, 16" to 20" Width \$744
Part: TC4V10W21WBC. HCPCS: E2620
Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$855
Part: TC4V22W32W. HCPCS: E2621
Complete following width package and size.
TRU-Comfort Back, Back Cane Option, 22" Width \$855
Part: TC4V22W32WBC. HCPCS: E2621
Complete following width package and size.

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED

Overall Height	Tilt or Static Width												Recline Width											
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"		
	14"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
15"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
16"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-		
17"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-		
18"	-	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-		
19"	-	S	T	S	T	S	T	S	T	S	T	S	S	S	T	S	T	S	T	S	T	S		
20"	-	S	T	S	T	S	T	S	T	S	T	S	S	S	T	S	T	S	T	S	T	S		
21"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S		
22"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S		
23"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S		
24"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S		
25"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S		
26"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S		
27"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S		
28"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S		
29"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S		
30"	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T		
31"	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T		
32"	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T		

S Sport Back
 Q Q-Back
 T TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

1c. ADI Back, AL Series, 4-Point, Active Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

- | | |
|---|---|
| <p>Active Contour, 16" Wide x 19" or 20" Overall Height \$622
Part: ALBR16SA4. HCPCS: E2613</p> <p>Active Contour, 16" Wide x 21" or 22" Overall Height \$622
Part: ALBR16TA4. HCPCS: E2613</p> <p>Active Contour, 16" Wide x 23" or 24" Overall Height \$622
Part: ALBR16XA4. HCPCS: E2613</p> <p>Active Contour, 18" Wide x 19" or 20" Overall Height \$622
Part: ALBR18SA4. HCPCS: E2613</p> <p>Active Contour, 18" Wide x 21" or 22" Overall Height \$622
Part: ALBR18TA4. HCPCS: E2613</p> | <p>Active Contour, 18" Wide x 23" or 24" Overall Height \$622
Part: ALBR18XA4. HCPCS: E2613</p> <p>Active Contour, 20" Wide x 19" or 20" Overall Height \$622
Part: ALBR20SA4. HCPCS: E2613</p> <p>Active Contour, 20" Wide x 21" or 22" Overall Height \$622
Part: ALBR20TA4. HCPCS: E2613</p> <p>Active Contour, 20" Wide x 23" or 24" Overall Height \$622
Part: ALBR20XA4. HCPCS: E2613</p> |
|---|---|

1d. ADI Back, AL Series, 4-Point, Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

- | | |
|---|---|
| <p>Deep Contour, 16" Wide x 19" or 20" Overall Height \$735
Part: ALBR16SD4. HCPCS: E2615</p> <p>Deep Contour, 16" Wide x 21" or 22" Overall Height \$735
Part: ALBR16TD4. HCPCS: E2615</p> <p>Deep Contour, 16" Wide x 23" or 24" Overall Height \$735
Part: ALBR16XD4. HCPCS: E2615</p> <p>Deep Contour, 18" Wide x 19" or 20" Overall Height \$735
Part: ALBR18SD4. HCPCS: E2615</p> | <p>Deep Contour, 18" Wide x 21" or 22" Overall Height \$735
Part: ALBR18TD4. HCPCS: E2615</p> <p>Deep Contour, 20" Wide x 19" or 20" Overall Height \$735
Part: ALBR20SD4. HCPCS: E2615</p> <p>Deep Contour, 20" Wide x 23" or 24" Overall Height \$735
Part: ALBR20XD4. HCPCS: E2615</p> |
|---|---|

1e. ADI Back, AL Series, 4-Point, Extra Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

- | | |
|---|---|
| <p>Extra Deep Contour, 18" Wide x 19" or 20" Overall Height \$822
Part: ALBR18S64. HCPCS: E2615</p> <p>Extra Deep Contour, 18" Wide x 21" or 22" Overall Height \$822
Part: ALBR18T64. HCPCS: E2615</p> | <p>Extra Deep Contour, 18" Wide x 23" or 24" Overall Height \$822
Part: ALBR18X64. HCPCS: E2615</p> |
|---|---|

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED

1f. Q-Back 2.0, AL Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AL14AS. HCPCS: E2613	Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$568 Part: QB2AL18AT. HCPCS: E2613
Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$568 Part: QB2AL14AT. HCPCS: E2613	Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$568 Part: QB2AL18AX. HCPCS: E2613
Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AL16AS. HCPCS: E2613	Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18AXX. HCPCS: E2613
Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568 Part: QB2AL16AT. HCPCS: E2613	Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AL20AS. HCPCS: E2613
Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568 Part: QB2AL16AX. HCPCS: E2613	Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$568 Part: QB2AL20AT. HCPCS: E2613
Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL16AXX. HCPCS: E2613	Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$568 Part: QB2AL20AX. HCPCS: E2613
Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AL18AS. HCPCS: E2613	Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20AXX. HCPCS: E2613

1g. Q-Back 2.0, AL Series, Deep Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL14DS. HCPCS: E2620	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL18DT. HCPCS: E2620
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL14DT. HCPCS: E2620	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL18DX. HCPCS: E2620
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL16DS. HCPCS: E2620	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18DXX. HCPCS: E2620
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL16DT. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL16DX. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL20DT. HCPCS: E2620
Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL16DXX. HCPCS: E2620	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL20DX. HCPCS: E2620
Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL18DS. HCPCS: E2620	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20DXX. HCPCS: E2620

1h. Q-Back 2.0, AL Series, Extra Deep Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL14ES. HCPCS: E2620	Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18EXX. HCPCS: E2620
Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL16ES. HCPCS: E2620	Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL20ES. HCPCS: E2620
Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$703 Part: QB2AL16ET. HCPCS: E2620	Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703 Part: QB2AL20ET. HCPCS: E2620
Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL16EX. HCPCS: E2620	Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL20EX. HCPCS: E2620
Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL18ES. HCPCS: E2620	Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20EXX. HCPCS: E2620
Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$703 Part: QB2AL18ET. HCPCS: E2620	
Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL18EX. HCPCS: E2620	

1i. Q-Back 2.0, AFT Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT16AS. HCPCS: E2613	AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT18AS. HCPCS: E2613
--	--

1j. Q-Back 2.0, AFT Series, Deep Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AFT16DS. HCPCS: E2620	AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AFT18DT. HCPCS: E2620
AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AFT16DT. HCPCS: E2620	AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AFT20DT. HCPCS: E2620
AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AFT18DS. HCPCS: E2620	

1k. Q-Back 2.0, AFT Series, Extra Deep Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AFT16EXX. HCPCS: E2620	
---	--

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED

13.12. TRU Balance 3 Back Shroud

REQUIRED. The back shroud is a feature of the TRU Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Inferno Red Part: DGN1714085-PRD5	Ice Pink Part: DGN1714085-PPK5	Bermuda Blue Part: DGN1714085-PBU15	Faded Glory Part: DGN1714085-WFL2
Back in Black Part: XRF000021	Pretty in Pink Part: XRF000027	Blue Jean Part: XRF000031	Black Camouflage Part: DGN1714085-WCM9-PCL2
Night Train Part: XRF000022	Evolution Orange Part: XRF000028	Lemon Crush Part: XRF000032	
Ghost Grey Part: DGN1714085-PGY6	Steel Blue Part: DGN1714085-PBU13	Forest Green Part: DGN1714085-PGN13	
Cosmic Purple Part: DGN1714085-PPR8	Royale Blue Part: DGN1714085-PBU14	Military Digital Camo Part: DGN1714085-WCM7	

13.13. Lateral Positioning Wedge

Optional.

Lateral Positioning Wedge, pair \$32
 Part: SET123465
 Available with Sport and TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

13.14. Back Width Package

REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections Part: PTOINDV3468	14"W-16"W Package Part: XRF000037	18"W-22"W Package Part: XRF000039
10"W-14"W Package Part: XRF000036	16"W-20"W Package Part: XRF000038	20"W-24"W Package Part: XRF000040

13.15. Back Width

REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

12" Back Width Part: PTOINDV3562	16" Back Width Part: PTOINDV3564	19" Back Width Part: PTOINDV3567	22" Back Width Part: PTOINDV3570
14" Back Width Part: PTOINDV3563	17" Back Width Part: PTOINDV3565	20" Back Width Part: PTOINDV3568	
15" Back Width Part: PTOINDV3856	18" Back Width Part: PTOINDV3566	21" Back Width Part: PTOINDV3569	

13.16. Back Cushion Configuration - Overall Back Height

Static Seat and Power Tilt available height range: 14"-29"H. Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion. Power Recline available height range: 19"-32"H. Odd heights include 6" gap and even heights include 7" gap. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

14" Height Part: SETHGHT1027	18" Height Part: SETHGHT1031	22" Height Part: SETHGHT1035	26" Height Part: SETHGHT1039	30" Height Part: SETHGHT1043
15" Height Part: SETHGHT1028	19" Height Part: SETHGHT1032	23" Height Part: SETHGHT1036	27" Height Part: SETHGHT1040	
16" Height Part: SETHGHT1029	20" Height Part: SETHGHT1033	24" Height Part: SETHGHT1037	28" Height Part: SETHGHT1041	
17" Height Part: SETHGHT1030	21" Height Part: SETHGHT1034	25" Height Part: SETHGHT1038	29" Height Part: SETHGHT1042	

13.17. Power Tilt Privacy Flap

Optional. A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap \$48
 Part: XRF000351

13.18. Back Height Extension

Optional. Not available with back cane option, 12" back width, or 24" back width. Adds 2" of seat back pan height. 25"H selection includes a back height extension.

Back Height Extension \$48
 Part: XRF000041

14. TRU BALANCE 3 ARMRESTS

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

14.1a. Flip-Back, 8"-16" H Armrests

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

- | | | | |
|---|--|---|--|
| <input type="checkbox"/> L <input type="checkbox"/> R | Flip-Back, 8"-16"H Armrests Standard
Left part: ARM1806662. Right part: ARM1806663
Not available with TRU Balance 3 Static Seat. | <input type="checkbox"/> L <input type="checkbox"/> R | Flip-Back, 8"-16"H Armrests for Static Seat \$289
Left part: XRF000320. Right part: XRF000321 |
|---|--|---|--|

Armrest Support

Optional. Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

- | | |
|---|--|
| <input type="checkbox"/> L <input type="checkbox"/> R | Flip-Back Armrest Downpost \$92 ea.
Left part: XRF000265. Right part: XRF000266 |
|---|--|

14.1b. Drop-In, Height Adjustable, Removable, T-Arm Armrests

Only available with TRU Balance 3 Static Seat. Includes flat armpads.

- | | | | |
|---|---|---|---|
| <input type="checkbox"/> L <input type="checkbox"/> R | Full Length, Drop-In Armrests No Charge
Left part: ARM1803253. Right part: ARM156742 | <input type="checkbox"/> L <input type="checkbox"/> R | Desk Length, Drop-In Armrests No Charge
Left part: ARM1803254. Right part: ARM156743 |
|---|---|---|---|

14.1c. Two-Post, Height Adjustable, Removable Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

- | | | | |
|---|--|---|--|
| <input type="checkbox"/> L <input type="checkbox"/> R | Full Length, Two-Post Armrest for Power Seating No Charge
Left part: XRF000047. Right part: XRF000048 | <input type="checkbox"/> L <input type="checkbox"/> R | Desk Length, Two-Post Armrest for Power Seating No Charge
Left part: XRF000049. Right part: XRF000050 |
|---|--|---|--|

1b. Two-Post Armrest for Static Seating

- | | | | |
|---|---|---|---|
| <input type="checkbox"/> L <input type="checkbox"/> R | Full Length, Two-Post Armrest for Static Seating \$158 ea.
Left part: XRF000051. Right part: XRF000052 | <input type="checkbox"/> L <input type="checkbox"/> R | Desk Length, Two-Post Armrest for Static Seating \$158 ea.
Left part: XRF000053. Right part: XRF000054 |
|---|---|---|---|

14.1d. Height Adjustable, Removable, Heavy Duty Armrests

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

- | | |
|--|---|
| 8"-12"H, Single Post, Height Adjustable, Removable, \$79 ea.
Heavy Duty (Sold in pairs)
Part: XRF000055. HCPCS w/ Static Seat: E0973 | 10"-14"H, Single Post, Height Adjustable, Removable, \$79 ea.
Heavy Duty (Sold in pairs)
Part: XRF000056. HCPCS w/ Static Seat: E0973 |
|--|---|

14.2. Armrest Height

REQUIRED. Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

2. Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

14.3. Armrest Size

REQUIRED. Select a left (L) and right (R) armrest size. This also affects armpad selection.

- | | | | | | |
|---|--|---|--|---|--|
| <input type="checkbox"/> L <input type="checkbox"/> R | 10" Length
L part: XRF000057. R part: XRF000123 | <input type="checkbox"/> L <input type="checkbox"/> R | 14" Length
L part: XRF000058. R part: XRF000124 | <input type="checkbox"/> L <input type="checkbox"/> R | 17" Length
L part: XRF000059. R part: XRF000125
Only available with Flip-Back, 8"-16" H
Armrests and straight arm pads. |
|---|--|---|--|---|--|

15. TRU BALANCE 3 ARPAPDS

REQUIRED. Select to omit an armpad selection or select a left and right armpad type. Please note that waterfall and gel armpads may cause interference with the seat back.

- | | |
|--|---|
| No Armpad Selection No Charge
Part: XRF000060
Not available with Flip-back, 8"-16" H armrests. | Interface for Aftermarket Pad No Charge
Part: XRF000308
Only available with Flip-back, 8"-16" armrests. |
|--|---|

15.1a. Foam Armpads

- | | | | |
|---|--|---|--|
| <input type="checkbox"/> L <input type="checkbox"/> R | Straight, Armpad Standard
Left part: XRF000061. Right part: XRF000126 | <input type="checkbox"/> L <input type="checkbox"/> R | Waterfall, Armpad No Charge
Left part: XRF000062. Right part: XRF000127
Not available with Flip-back, 8"-16" armrests. |
|---|--|---|--|

TRU Balance 3 armpads section continued on following page.

15. TRU BALANCE 3 ARMPADS CONTINUED

15.1b. Gel Armpads
 Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right\$195 Part: ARM1806807	Falcon Gel, Desk Length Left, Full Length Right.....\$171 Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right.....\$171 Part: ARM1806813	Falcon Gel, Desk Length Left and Right\$147 Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right\$311.22 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right.....\$287.93 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right.....\$287.93 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$265.02 Part: ARM136437

16. TRU BALANCE 3 ARMREST BRACKETS

REQUIRED.

Standard Fixed Receiver Standard Part: XRF000063	Stealth Height Adjustable Receiver \$236 Part: XRF000064. HCPCS: K0108
---	---

17. TRU BALANCE 3 HEADRESTS

Optional. Select a headrest pad and bracket.

17.1. Headrests, Comfort Plus with Cool Core
 Only available with 14" to 22" back widths of Static Type or Flip Forward back, without occupied transit.

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179 Part: ST-CP08CCFPM-Q. HCPCS: E0955	Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179 Part: ST-CP14CCFPM-Q. HCPCS: E0955
Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)\$179 Part: ST-CP10CCFPM-Q. HCPCS: E0955	Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179 Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink.....\$311 Part: ST-UHB2D33200-Q. HCPCS: E1033	Multi-Axis Removable Headrest Bracket..... \$174 Part: ST-HMO475-17M-Q. HCPCS: E1033
--	---

17.2. Q Cam Lock

Optional.

Q Cam Lock..... \$50 Part: XRF000376	<table border="1"> <thead> <tr> <th>Q Cam Lock Availability (Overall height)</th> <th colspan="2">TRU-Comfort/Sport</th> <th colspan="2">Q-Back 2</th> </tr> <tr> <td></td> <th>Standard</th> <th>\$50</th> <th>Standard</th> <th>\$50</th> </tr> </thead> <tbody> <tr> <td>Tilt/Static/Back Canes</td> <td>14"-21"H</td> <td>22"-29"H</td> <td>19"-24"H</td> <td>25"-26"H</td> </tr> <tr> <td>Power Recline Back</td> <td>19"-21"H</td> <td>22"-32"H</td> <td>22"-25"H</td> <td>26"-29"H</td> </tr> </tbody> </table>	Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2			Standard	\$50	Standard	\$50	Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H	Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H
Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2																		
	Standard	\$50	Standard	\$50																	
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H																	
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H																	

18. LAP BELTS

REQUIRED. Select a lap belt below.

18.1a. Standard Lap Belts

38" Lap Belt Part: ACCBELT1006
 50" Lap Belt Part: ACCBELT1000
 60" Lap Belt Part: ACCBELT1002
 70" Lap Belt Part: ACCBELT1003
 80" Lap Belt Part: ACCBELT1004
 90" Lap Belt Part: ACCBELT1001

18.1b. Optional Lap Belts

Stealth Retractable Lap Belt (IFRB5350).....\$181 Part: ST-IFRB5350	Stealth Positioning, 60" Single Pull, Push Button, Extra Large \$50 Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q
Stealth Positioning, 60" Single Pull, Push Button, Medium..... \$50 (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q	60" Sure-Lok Lap Belt (Non-retractable)..... \$120 Part: ACCBELT1016 Required with ISO 7176-19 Compliant Transit Option.
Stealth Positioning, 60" Single Pull, Push Button, Large \$50 Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q	

19. TRU BALANCE 3 LEGREST CONTROLS

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (adds 2 actuator control requirement) Part: XRF000108
 Articulating Foot Platform (adds 1 actuator control requirement) Part: XRF000110

Combined Legrests (adds 2 actuator control requirement) Part: XRF000109

20. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

20.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

1. Center Mount Foot Platform Type

Basic Center Mount..... No Charge
 Part: XRF000323

2. Basic Center Mount Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates \$180 Long Extension, 5"x9" Footplates \$180
 Part: XRF000312. HCPCS: K0040 Part: XRF000314. HCPCS: K0040

Medium Extension, 5"x9" Footplates \$180
 Part: XRF000313. HCPCS: K0040

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount..... \$158
 Part: RIG2104369

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16
 Part: RIG2104240

20.1b. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Center Mount Foot Platform Type

Power Articulating Foot Platform \$3,145 Quick Adjustable Center Mount Foot Platform \$263
 Part: PTOINDV3754. HCPCS: E1012+K0040 Part: XRF000311. HCPCS: K0040

2. Extension and Angle Adjustable Footplates

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates Long Extension, 5"x6" Footplates
 Part: RIG152577 Part: RIG151505

Pediatric Extension, 5"x8" Footplates Long Extension, 5"x8" Footplates
 Part: RIG2404015 Part: RIG2404017

Medium Extension, 5"x8" Footplates Long Extension, 5"x11" Footplates
 Part: RIG2404016 Part: RIG2404019

Medium Extension, 5"x11" Footplates
 Part: RIG2404018

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
 Part: RIGASMB7120018 Part: RIGASMB7120019

Small Calf Pads, Manually Adjustable Center Mount Only, \$158 Large Calf Pads, Manually Adjustable Center Mount Only, \$158
 6"Wx6"H (pair) 8" taper to 6"Wx8"H (pair)
 Part: XRF000156 Part: XRF000157

3b. Center Mount Calf Panel

Small Calf Panel, 18"Wx8"H..... \$184 Large Calf Panel, 24"Wx8"H \$237
 Part: FRMASMB7110012 Part: RIG123825
 Only available on widths of 20" or greater.

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69 Outside Edge \$69
 Part: RIG152519 Part: RIG2404921
 Only available with the 5x6 footplates. Available with 5"x8" and 5"x11" footplates.

Width Extension for 8" \$69 Front and Rear Cap \$16
 Part: RIG2404919 Part: RIG2404923
 Available with 5"x8" and 5"x11" footplates.

Width Extension for 11" \$69
 Part: RIG2404920

5. Center Mount Additional Accessories

Small Foam Calf Pad Cover (pair) \$67 Footplate Riser \$75
 Part: ST-CSCS0775-Q Part: XRF000348
 Select legrest length between 2" to 9" at 1" increments. Not available with long extension.

Small Foam Panel Cover (pair) \$176 Footplate Riser Calf Panel \$50
 Part: ST-CSCP0920-Q Part: FRMASMB7110011
 Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

Footplate Connector \$16
 Part: RIGASMB7120023

Legrests section continued on following pages.

20. LEGRESTS CONTINUED

6. Stealth Gel Footplate Covers

- 5"x8" Tapered Gel Footplate Covers..... \$222
Part: ST-FPG5008T-Q. HCPCS: K0108
- 5"x11" Tapered Gel Footplate Covers..... \$270
Part: ST-AFPG5011T-Q. HCPCS: K0108

- 5"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K0108
- 5"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-AFPG5011E-Q. HCPCS: K0108

20.1c. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

- P L R 60°-90° Swing-Away Legrest, Power Seating..... No Charge
Pair part: RIG120062. Left: FRMASMB17558.
Right: FRMASMB17559
Available with Power Positioning only. Not available with 11"-13" extensions.
- P L R 60°-90° Swing-Away Legrest, Short, Power Seating.... No Charge
Pair part: RIG136481. Left: RIG132607.
Right: RIG132609.
Available with Power Positioning only. For use only with 11"-13" lower extensions.

- P L R 60°-90° Swing-Away Legrest, Static Seat \$86 ea.
Pair part: XRF000082. Left: XRF000083.
Right: XRF000084
Available with Static Seat only. Not available with 11"-13" extensions.
- P L R 60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea.
Pair part: XRF000085. Left: XRF000086.
Right: XRF000087.
Available with Static Seat only. For use only with 11"-13" lower extensions.

2. Lower Extension

- P L R 11"-13" Length
Pair part: XRF000088.
Left: XRF000089.
Right: XRF000090.
- P L R 13"-16.5" Length
Pair part: XRF000094.
Left: XRF000095.
Right: XRF000096.

- P L R 15.5"-19.5" Length
Pair part: XRF000097.
Left: XRF000098.
Right: XRF000099.

3a. Extruded Footplate(s)

- P L R 7.25"Wx6"D Extruded Footplate..... Standard
Part pair: FRMASMB7847. Left: FRMASMB7861
Right: FRMASMB7862

- P L R 8.5"Wx6"D Extruded Footplate Standard
Part pair: FRMASMB7848. Left: FRMASMB7863
Right: FRMASMB7864
May cause interference on seat widths less than 18."

3b. Swing-Away Legrest, Angle Adjustable Footplate(s) with Heel Loops

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

- P L R 4"Wx6"D \$90 ea.
Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560
HCPCS: K0040+E0951
- P L R 5"Wx6"D..... \$90 ea.
Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561
HCPCS: K0040+E0951
- P L R 6"Wx11"D¹ \$90 ea.
Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564
HCPCS: K0040+E0951
- P L R 8"Wx11"D² \$90 ea.
Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565
HCPCS: K0040+E0951

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-13"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓

11"-13" extension causes no interference at all seat-to-floor heights.

This rigging configuration causes no interference issues at this seat-to-floor height.

This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

- Therafin 31515 Multi-Axis Amputee Support \$230 ea.
QTY Part: INDPART2822. HCPCS: E1020

- Therafin 32387 Amputee Support Gel Cover \$155 ea.
QTY Part: INDPART2823. HCPCS: K0108

Legrests section continued on following pages.

20. LEGRESTS CONTINUED

20.1d. Power Elevating Legrests

Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Power Elevating Legrest Upper Assembly

Includes calf pads.

Power Elevating Legrests, Pair \$2,516
 Part: XRF000341. HCPCS: E1010+K0040+E0951

Legrest Length Matrix	Seat to Floor Option		
	Low	Middle	High
Power ELR: Static Seat, Recline Only			
8.5"-14.5" Lower Length	8.5"- 12.5"	8.5"- 13.5"	8.5"- 14.5"

Legrest Length Matrix	Seat to Floor Option		
	Low	Middle	High
Power ELR: Tilt Only, Tilt and Recline, Any Power Lift			
8.5"-14.5" Lower Length	8.5"- 14"	8.5"- 14.5"	8.5"- 14.5"

2. Lower Extension

8.5" to 14.5" Lower Length
 Part: XRF000344

3. Angle Adjustable Footplate with Heel Loop and Wheel Bumper

¹The 6"W footplates may cause interference with seat widths less than 16". ²The 8"W footplates may cause interference with seat widths less than 19".

4"Wx6"D, Pair Part: FRMASMB7448	5"Wx8"D, Pair Part: FRMASMB7450	6"Wx11"D ¹ , Pair Part: FRMASMB7452
5"Wx6"D, Pair Part: FRMASMB7449	6"Wx8"D ¹ , Pair Part: FRMASMB7451	8"Wx11"D ² , Pair Part: FRMASMB7453

20.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

P L R Elevating Legrest Assembly \$171 ea.
 Pair part: FRMASMB6483. Left: RIG138012
 Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 11"-14" Length Pair: FRMASMB4755 Left: RIG138026 Right: RIG138030 For reinforced only	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 14"-17" Length Pair: FRMASMB4756 Left: RIG138033 Right: RIG138035 For reinforced only
<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 8.5"-12" HD Drop-Ins/ 12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103. For angle adjustable only	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 13"-18" HD Drop-Ins/ 15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106 For angle adjustable only

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓

This rigging configuration causes no interference issues at this seat-to-floor height.
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 7.25"Wx6"D Reinforced Footplate Standard Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 8.5"Wx6"D Reinforced Footplate Standard Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062 May cause interference on seat widths less than 18."
--	---

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 4"Wx6"D \$90 ea. Pair parts: FRMASMB8812 x2. HCPCS: K0040	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 6"Wx8"D ¹ with Heel Loops \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951
<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 5"Wx6"D \$90 ea. Pair parts: FRMASMB8811 x2. HCPCS: K0040	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 6"Wx11"D ¹ with Heel Loops \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951
<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 5"Wx8"D with Heel Loops \$90 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 8"Wx11"D ² with Heel Loops \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

Legrests section continued on following pages.

20. LEGRESTS CONTINUED

20.2. Legrest Accessories

Heel Loops (pair)..... \$53 Wheel Bumper..... \$19 ea.
 Part: FRMASMB7873. HCPCS: E0951 Part: XRF000100
 Standard on angle adjustable and multi axis footplates when noted.

20.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	4" L (#3991)	6" L (#3993)	8" L (#3680)	9" L (#3682)	10" L (#3684)	11" L (#3686)	12" L (#3688)	13" L (#3690)	14" L (#3692)	15" L (#3694)
3" L (#3990)	5" L (#3992)	7" L (#3994)	8.5" L (#3681)	9.5" L (#3683)	10.5" L (#3685)	11.5" L (#3687)	12.5" L (#3689)	13.5" L (#3691)	14.5" L (#3693)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)
7.5" L (#3471)	8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	

21. TRU BALANCE 3 POWER POSITIONING ELECTRONICS

REQUIRED with power seating or legrest options. Select one option from section 21.1 or 21.3.

21.1a. Tilt Only through Toggle

Available with QLNEs only.

Single Actuator Control, Power Tilt Only, Through Toggle No Charge
 Part: ELE1905907
 Select toggle side.

Left Side Toggle Right Side Toggle
 Part: PTOINDV1839 Part: PTOINDV1838

21.1b. Actuator Control through Joystick

1. Actuator Control Selection

Single Actuator Through Controller and \$2,163 Multiple Actuator Through Controller and \$3,040
 Interface Interface
 Part: XRF000335. HCPCS: E2310 Part: XRF000336. HCPCS: E2311

Single Actuator Through Controller and \$3,040
 Interface (AAM)
 Part: XRF000068. HCPCS: E2310

2. Q-Logic 3 Optional Switch Control

Only available with Q-Logic 3, multiple actuator configurations.

4-Function Push Button 4-Way Toggle
 Part: CTL123218 Part: ELEASMB5925
 See section 21.2 for switch mounting.

4-Way Button Toggle
 Part: ELEASMB5926

3. User Controls Side Selection

Required with 4-Way Toggle and 4-Way Button Toggle.

Mount User Controls on Left Mount User Controls on Right
 Part: PTOINDV3996 Part: PTOINDV3997

21.1c. Actuator Control with iAccess Switch

Single Actuator Control Through Flat Button..... \$2,547 iAccess and Controller and Interface Part: XRF000337. HCPCS: E2310	Single Actuator Control Through Toggle Button \$3,460 iAccess and Controller and Interface (AAM) Part: XRF000074. HCPCS: E2310
Single Actuator Control Through Toggle Button \$2,547 iAccess and Controller and Interface Part: XRF000338. HCPCS: E2310	Multiple Actuator Control Through Flat Button \$3,460 iAccess and Controller and Interface Part: XRF000339. HCPCS: E2311
Single Actuator Control Through Flat Button..... \$3,460 iAccess and Controller and Interface (AAM) Part: XRF000070. HCPCS: E2310	Multiple Actuator Control Through Toggle Button..... \$3,460 iAccess and Controller and Interface Part: XRF000340. HCPCS: E2311

21.2. 4-Function Push Button/iAccess Mount

REQUIRED with 4-Function Push Button or iAccess options. Under arm mount is not available with TRU Balance 3 HD armrests.

Joystick Mount No Charge Part: MEC1907500	Attendant Mount \$26 Part: MEC1906564
Under Arm Mount, Flip-Back Armrest..... No Charge Part: MEC1906531	Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$26 Part: MEC2001392
Under Arm Mount, Tube Style Armrest No Charge Part: MEC1906530	Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$26 Part: MEC2001391



TRU Balance 3 power positioning electronics section continued on following page.

21. TRU BALANCE 3 POWER POSITIONING ELECTRONICS CONTINUED

1. Under Arm Mount Side Selection

Left Side Under Arm Mount
Part: PTOINDV3678

Right Side Under Arm Mount
Part: PTOINDV3679

21.3. Power Seating through Switches

Power seating through switches (mounting bracket included) is special order\$495
Part: PTOINDV3932

22. TRU BALANCE 3 ATTENDANT CONTROL

Optional. *Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control \$525
Part: XRF000078. HCPCS: E2331
Select this option if no charge option criteria is not met.

Attendant Control with Tilt, Recline, Alternative No Charge
Controls and Q-Logic Display or SCIM Module
Part: XRF000079

23. PROVIDER INSTALLED ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.

TRU Balance 3 Cane/Crutch Holder \$105
Part: ACC152515. HCPCS: E2207

Quantum Backpack \$168
Part: ACC1905853

Grip Lap Board (20005)..... \$38.95
Part: ACC152992

Grip Drink Holder (20085) \$15.50
Part: ACC152993

TRU Balance 3 Hydration System (B) \$237
Part: ACC125010

TRU Balance 3 Transfer Handle (D) \$105 ea.
Part: MEC2205680. HCPCS: E0705
Not available with Heavy Duty armrests.

TRU Balance 3 Cup Holder No Charge
Part: ACC125003
For use with TRU Balance 3 (track style) armrests.

Cup Holder (Tube-style Armrest)..... No Charge
Part: ACCASMB2634
For use with tube-style armrests.

TRU Balance 3 Personal Item Hook \$69
Part: ACC125008
Not available with back cane option.

TRU Balance 3 Reflector Kit..... \$53
Part: ACC124670

TRU Balance 3 Backpack Holder (Pair) (A) \$69
Part: ACC147774
Not available with back cane option.

TRU Balance 3 Front Mounted \$69
Cell Phone/ Media Holder, X-Grip (E)
Part: ACC155835
Not available with Heavy Duty armrests.

TRU Balance 3 Rear View Mirror..... \$69
Part: ACC155833
Not available with Heavy Duty armrests.

TRU Balance 3 Oxygen Holder \$210
Part: ACC125006. HCPCS: E2208
Includes accessory bar, a credit will be applied PTOINDV3769.

TRU Balance 3 Clothing Guard (Pair) \$237
Part: ACC125007
Not available with Heavy Duty armrests.

TRU Balance 3 Glove Box (C) \$158
Part: ACC125011

Stealth Large Tablet Holder\$625.28
Part: ACC141818
Not available with Heavy Duty armrests.

Push Handles (pair)..... \$195
Part: ACC2005027
Not available with oxygen holder, back cane option, or vent tray. Required with attendant manual recline.

Accessory Bar \$132
Part: XRF000128
Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder. Not available with back cane option.

Large, Black, Weather Cover \$202
Part: ACCCOVR1021


Flip-Back Armrest, Side Track Accessory Bracket ... \$11
Part: MEC1809178
Only available with flip-back, 8"-16" H armrests. This adds a 47mm channel to the side of the armrest to allow mounting of accessories with t-nuts.

Privacy Bag..... \$42
Part: ACC2005723

Personal Item Pouch \$84
Part: ACC2101017

Depth Adjustable Vent Tray \$595
Part: ACC2302500. HCPCS: E1029
Includes accessory bar, a credit will be applied PTOINDV3769.

Remote Stop Switch\$1,295
Part: XRF000379



Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
Part: ACC2407046

Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
Part: ACC2405496

24. BRANCH PATCHES

Optional.

Air Force Patch\$21 Part: ACC140240	Coast Guard Patch\$21 Part: ACC140239	Navy Patch\$21 Part: ACC140244
Army Patch\$21 Part: ACC140241	Marines Patch\$21 Part: ACC140243	Purple Heart Patch\$21 Part: ACC140238



25. POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

25.1. Stealth Thoracic Laterals

Swing-away laterals may interfere with TRU Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$664 Part: ST-ULS4042000-Q. HCPCS: E1034 Only available with Redesigned back assemblies.	Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$332 Part: ST-ULPK404103L00-Q. HCPCS: E1034 Only available with Active Contour Q-Backs.
Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664 Part: ST-ULPK404108B00-Q. HCPCS: E1034 Only available with Sport or TRU Comfort backs.	Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$332 Part: ST-ULPK404103R00-Q. HCPCS: E1034 Only available with Active Contour Q-Backs.
Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664 Part: ST-ULPK404104B00-Q. HCPCS: E1034 Only available with back cane and omit back configurations.	

2. Pads

Curved Pads, 3"x5", pair (TWBL3X5)..... \$250 Part: INDPART2663. HCPCS: E0956	Flat Pads, 3"x5", pair (ULP846050320)..... \$250 Part: ST-ULP846050320-Q. HCPCS: E0956
Curved Pads, 4"x6", pair (TWBL4X6)..... \$250 Part: INDPART2664. HCPCS: E0956	Flat Pads, 4"x6", pair (ULP835050320)..... \$250 Part: ST-ULP835050320-Q. HCPCS: E0956

3. Installation

Install Laterals Part: PTOINDV3499	Do Not Install Part: PTOINDV3500
---------------------------------------	-------------------------------------

25.2. Stealth UniLink Pelvic/Thigh Guide

1. UniLink Pelvic/Thigh Hardware

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea. Removable Mounting Hardware QTY Part: ST-UTTP5-Q. HCPCS: E1028/E1034^	UniLink Multi-Angle Thigh Guide, Left \$358.87 ea. Part: ST-UTT8A116O2-L-Q.
UniLink, 7" Height Adjustable, \$276 ea. Removable Mounting Hardware QTY Part: ST-UTTP7-Q. HCPCS: E1028/E1034^	UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea. Part: ST-UTT8A116O2-R-Q.
UniLink Flip Down Hip/Thigh Support, Left..... \$332 ea. Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^	Flip Down with LEFT Full Surface Contact \$281 ea. Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028/E1034^
UniLink Flip Down Hip/Thigh Support, Right \$332 ea. Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^	Flip Down with RIGHT Full Surface Contact..... \$281 ea. Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea. Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^	UniLink, 4" x 12" Pad \$106 ea. Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^
UniLink, 4" x 8" Pad..... \$106 ea. Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^	

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea. Part: ST-UIL001150A000-Q	UniLink, 2" Link..... \$81 ea. Part: ST-UIL001200A000-Q
---	--

4. UniLink Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

25.3. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea. Part: ST-IFPS1103-Q. HCPCS: E0960	Upper Contour, Medium (IFPS2103)..... \$154 ea. Part: ST-IFPS2103-Q. HCPCS: E0960
Upper Classic, Large (IFPS1104) \$154 ea. Part: ST-IFPS1104-Q. HCPCS: E0960	Upper Contour, Large (IFPS2104) \$154 ea. Part: ST-IFPS2104-Q. HCPCS: E0960



Positioning components section continued on following page.

25. POSITIONING COMPONENTS CONTINUED

25.4. Stealth Positioning Chest Straps

Chest Strap, Medium \$97 ea.
w/ Belt Anchor, (IFCS4113)
Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large \$97 ea.
w/ Belt Anchor, (IFCS4114)
Part: ST-IFCS4114-Q. HCPCS: E0960



25.5. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right..... \$318 ea.
Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left \$318 ea.
Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$434 ea.
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028




2. Facial Lateral Pads and Accessories

Egg Switch Mount (ESM) (A) \$57 ea.
Part: INDPART2815

Swing-Away Joystick Mount, 18" Rod (C) \$276 ea.
Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest

Facial Pad, 4".5x2", Cool Core (900M-Q) \$61 ea.
Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2", Cool Core (920M-Q) \$61 ea.
Part: ST-920M-Q. HCPCS:K0108



25.6. Stealth Elbow Support

1. Elbow Support Bracket

Left, Swing-Away, Elbow Support..... \$199 ea.
Bracket Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support..... \$199 ea.
Bracket Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket..... \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
Part: ST-UETB4021FL-Q

Right, UniLink, Fixed, Elbow Support Bracket..... \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
Part: ST-UETB4021FR-Q

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.
Part: POS124067. HCPCS: K0108



25.7. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

25.8. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

26. LAP TRAY

26.1. Quantum Lap Tray

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

Black Lap Tray, 11" W \$200
Part: MEC1905328-11B. HCPCS: E0950
For use with 12"-13" back width.

Black Lap Tray, 17" W \$200
Part: MEC1905328-17B. HCPCS: E0950
For use with 18"-19" back width.

Black Lap Tray, 13" W \$200
Part: MEC1905328-13B. HCPCS: E0950
For use with 14"-15" back width.

Black Lap Tray, 19" W \$200
Part: MEC1905328-19B. HCPCS: E0950
For use with 20"-22" back width.

Black Lap Tray, 15" W \$200
Part: MEC1905328-15B. HCPCS: E0950
For use with 16"-17" back width.

Lap tray section continued on following page.

26. LAP TRAY

1b. Clear Lap Tray

- Clear Lap Tray, 11" W \$200
Part: MEC1905328-11C. HCPCS: E0950
For use with 12"-13" back width.
- Clear Lap Tray, 13" W \$200
Part: MEC1905328-13C. HCPCS: E0950
For use with 14"-15" back width.
- Clear Lap Tray, 15" W \$200
Part: MEC1905328-15C. HCPCS: E0950
For use with 16"-17" back width.

- Clear Lap Tray, 17" W \$200
Part: MEC1905328-17C. HCPCS: E0950
For use with 18"-19" back width.
- Clear Lap Tray, 19" W \$200
Part: MEC1905328-19C. HCPCS: E0950
For use with 20"-22" back width.

2. Mounting Hardware

- Mounting Bracket, Chest Bar/Lap Tray \$75
Part: MEC2104886

26.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

- Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950
- Simple Solutions Medium Tray with Cut Out \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950
- Simple Solutions Large Tray with Cut Out \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

- Swing Away Mounting Hardware \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	S	M	M	M	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	S	M	M	M	L	L	
Drop-In, T-Arm	S	S	S	S	S	S	M	M	M	L	L	
Two-Post, Flip-Back	S	S	S	S	S	S	M	M	M	L	L	
Cantilever	S	S	S	S	S	M	M	M	M	L	L	
Q-Captain, Flip-Up						S	M			L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	M	M	M	M	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q
:Compatible :Armrest configuration may not hold tray position
 *Some armrests listed may not be available on this order form.

27. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).

27.1. Enhanced Display

- Q-Logic 3 Enhanced Display Kit \$1,980
Part: CTL167726
Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

- | | | |
|--|--|---|
| <ul style="list-style-type: none"> Black, Soft Touch Standard
Part: DGN137728
Also matches Night Train shrouds. Inferno Red \$63
Part: DGN137729-PRD5 Back in Black \$63
Part: DGN137729-PBK1 Ghost Grey \$63
Part: DGN137729-PGY6 | <ul style="list-style-type: none"> Cosmic Purple \$63
Part: DGN137729-PPR8 Ice Pink \$63
Part: DGN137729-PPK5 Pretty in Pink \$63
Part: DGN137729-PPK3 Evolution Orange \$63
Part: DGN137729-POR7 Steel Blue \$63
Part: DGN137729-PBU13 | <ul style="list-style-type: none"> Royale Blue \$63
Part: DGN137729-PBU14 Bermuda Blue \$63
Part: DGN137729-PBU15 Blue Jean \$63
Part: DGN137729-PBU11 Lemon Crush \$63
Part: DGN137729-PYL4 Forest Green \$63
Part: DGN137729-PGN13 |
| <ul style="list-style-type: none"> Military Digital Camo \$63
Part: DGN137729-WCM7 Faded Glory \$63
Part: DGN137729-WFL2 Black Camouflage \$63
Part: DGN137729-WCM9-PCL2 | | |

Alternative controls and electronics section continued on following page.

27. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

27.2. SCIM Module

Optional.

- SCIM Module for Q-Logic 3 \$525
Part: CTL168055
 Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



27.3. Q-Logic 3 SAJ

- Q-Logic 3 EX SAJ \$840
Part: CTL143980
 Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



27.4. Head Array

- i-Drive 4.0 Bluetooth Tri-array Series, \$4,774
 3 Switch with Egg Switch
Part: ST-IDHBT200-U-Q, HCPCS: E2330
- Junior, Curved, i-Drive Tri-Array with \$4,774
 Mounting Hardware
Part: ST-IDHBT450-U-Q, HCPCS: E2330
- i-Drive 4.0 Bluetooth Pro Series, \$4,825
 3 Switch with Egg Switch
Part: ST-IDHBT300-U-Q, HCPCS: E2330
- Adult, Curved, i-Drive Tri-Array with \$4,774
 Mounting Hardware
Part: ST-IDHBT400-U-Q, HCPCS: E2330



27.5. Sip and Puff Drive

1. Sip and Puff Module

- Sip and Puff Module for Q-Logic 3 \$1,990
Part: CTL168056, HCPCS: E2325
 Includes module and mounting hardware.

2. Straw and Mount

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

- Stealth Sip and Puff for Comfort Series, \$318
 Left Mount (SP625-CL)
Part: ST-SP625M-CL-Q, HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.
- Stealth Sip and Puff for Stealth i-Drive, \$318
 Head Array (SP625M-U)
Part: ST-SP625M-U-Q, HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.
- Stealth Sip and Puff for Comfort Series, \$318
 Right Mount (SP622-CR)
Part: ST-SP625M-CR-Q, HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



3. Sip & Puff Accessories

- Sip & Puff Straws, 12 pack (SP625-S12) \$35
Part: ST-SP625-S12-Q

27.6. Egg Switches

- Red Egg Switch (ESRED) \$87
Part: SWTMCRO1077
- Blue Egg Switch (ESBLU) \$87
Part: SWTMCRO1079
- Black Egg Switch (ESBLK) \$87
Part: SWTMCRO1081
- Green Egg Switch (ESGRN) \$87
Part: SWTMCRO1078
- Yellow Egg Switch (ESYEL) \$87
Part: SWTMCRO1080



27.7. Switch Adapter

- 2 to 1 Switch Adapter \$15
Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

Alternative controls and electronics section continued on following pages.

27. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

27.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520
Part: CTL168057

27.9. USB Mobile Device Charger

USB Mobile Device Charger \$184

Part: ELE168069

Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

