



Q6 Edge® 4 Order Form with TRU Balance® 4 Positioning

Account Number: _____ Date: _____
 Provider Name: _____
 Contact: _____
 Phone: _____ Fax: _____
 Email: _____
 PO Number: _____
 Quote Name/Marked For: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
 B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
 C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 7 mph programmable motors and off board charger. XLR USB charger is included. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 4 3SP-SS \$8,559	Quantum Q6 Edge 4 3MP-SS \$8,373
Part: Q6 EDGE 4 3SP-SS. HCPCS: K0856	Part: Q6 EDGE 4 3MP-SS. HCPCS: K0861
Single power option.	Multiple actuator power option.

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 7 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology \$1,360
Coming soon. Part: XRF000268. HCPCS: K0108

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

- | | | | |
|---|--|---|---|
|  Inferno Red
Part: DGN2401066-PRD5 |  Cosmic Purple
Part: DGN2401066-PPR8 |  Steel Blue
Part: DGN2401066-PBU13 |  Lemon Crush
Part: DGN2401066-PYL4 |
|  Back in Black
Part: DGN2401066-PBK1 |  Ice Pink
Part: DGN2401066-PPK5 |  Royale Blue
Part: DGN2401066-PBU14 |  Forest Green
Part: DGN2401066-PGN13 |
|  Night Train
Part: DGN2401066-PBK3 |  Pretty in Pink
Part: DGN2401066-PPK3 |  Bermuda Blue
Part: DGN2401066-PBU15 | |
|  Ghost Grey
Part: DGN2401066-PGY6 |  Evolution Orange
Part: DGN2401066-POR7 |  Blue Jean
Part: DGN2401066-PBU11 | |

3.1b. Hydrographic Shrouds

- | | | |
|--|--|---|
| <input type="checkbox"/>  Military Digital Camo \$263
Part: DGN2401066-WCM7 | <input type="checkbox"/>  Faded Glory \$263
Part: DGN2401066-WFL2 | <input type="checkbox"/>  Black Camouflage \$263
Part: DGN2401066-WCM9-PCL2
<i>Coming soon.</i> |
|--|--|---|

4. BASE MOUNTED LIGHTS

Optional. Select one option.

- | | |
|--|--|
| Accent Light Package Standard
Part: ELE2404196 | Full Light Package \$1,045
Part: ELE2404199
Requires QLNES 10 Key, Q-Logic 3 EX Joystick for Full Lights, or Q-Logic 3e Joystick. |
|--|--|

5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

- Quantum Backup Camera** \$368
Part: MEC1910365

1. Optional Mounting Bracket

- | | |
|---|---|
| Backup Camera Headrest Bracket \$21
Part: MEC1910373
<i>Works with .5" x .5" vertical supports only.</i> | Backup Camera Track Mount Bracket \$27
Part: MEC2104227 |
|---|---|

2. Optional Display Mounting Bracket

- Backup Camera Q-Logic Joystick Mount Bracket** \$105
Part: MEC2106123



6. DRIVE WHEELS

REQUIRED. Select a type of drive wheel.

- | | |
|---|--|
| 3"x14" Flat Free Tires Standard
Part: PTOINDV3525 | 4"x14" Knobby, Flat Free Black Tires \$580
Part: XRF000368 |
| 3"x14" Pneumatic Tires No Charge
Part: WHL1713925 | |

7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

- | | |
|--|--|
| Group 24, Gel Cell, Interceptor 2 Battery \$399 ea.
Part: BAT2404214. HCPCS: E2361
Battery tray configuration. | NF-22, AGM, MK Flat Top Battery \$432 ea.
<i>Coming soon.</i> Part: BAT2501912. HCPCS: E2361
Without battery tray configuration |
| NF-22, Gel Cell, Interceptor 2 Battery \$355 ea.
<i>Coming soon.</i> Part: BAT2500310. HCPCS: E2361
Without battery tray configuration. | |

8. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

8.1. Expandable Controls

1. Expandable Controls Required Components

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected elsewhere, Customer Service will quote the order with these items.

- | | |
|---|---|
| Q-Logic 3 Controller Module \$730
Part: XRF000002. HCPCS: E2377 | Q-Logic 3 Expandable Controller Harness \$458
Part: XRF000264. HCPCS: E2313 |
|---|---|

2. Q-Logic 3 Joystick

- | | |
|--|---|
| Q-Logic 3 EX Joystick \$840
Part: CTL167724 | Q-Logic 3 EX Joystick for Full Lights \$840
Part: CTL167725 |
| Q-Logic 3e Joystick \$735
Part: CTL1910600
<i>Can operate full light package. Not available with alternative or specialty controls.</i> | No Joystick, Setup for Q-Logic
Part: XRF000010 |



Electronics section continued on following page.

8. ELECTRONICS CONTINUED

8.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2.

Black, Soft TouchStandard Part: DGN124307 Also matches Night Train shrouds.	Pretty in Pink \$63 Part: DGN124336-PPK3	Lemon Crush \$63 Part: DGN124336-PYL4
Inferno Red \$63 Part: DGN124336-PRD5	Evolution Orange \$63 Part: DGN124336-POR7	Forest Green \$63 Part: DGN124336-PGN13
Back in Black \$63 Part: DGN124336-PBK1	Steel Blue \$63 Part: DGN124336-PBU13	Military Digital Camo \$63 Part: DGN124336-WCM7
Ghost Grey \$63 Part: DGN124336-PGY6	Royale Blue \$63 Part: DGN124336-PBU14	Faded Glory \$63 Part: DGN124336-WFL2
Cosmic Purple \$63 Part: DGN124336-PPR8	Bermuda Blue \$63 Part: DGN124336-PBU15	Black Camouflage \$63 Part: DGN124336-WCM9-PCL2
Ice Pink \$63 Part: DGN124336-PPK5	Blue Jean \$63 Part: DGN124336-PBU11	

9. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90 Part: ACCASMB1860. HCPCS: E2323	Extended I-Handle (G) \$90 Part: ACCASMB1857. HCPCS: E2323
Chin Cup (C) \$90 Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)	Bodypoint 3" Goal Post Handle \$104 Part: PLSKNCP1057. HCPCS: E2323
T-Bar (D) \$126 Part: ACCASMB1858. HCPCS: E2323	Bodypoint 4" Goal Post Handle \$104 Part: PLSKNCP1058. HCPCS: E2323
Medium Ball Handle (E) \$69 Part: ACCASMB1859. HCPCS: E2323	Q-Logic 3 Carrot Handle (A) No Charge Part: DGN124308
Soft Ball Handle (F) \$80 Part: ACCASMB1856. HCPCS: E2323	



10. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

10.1a. Fixed, In-line Joystick Bracket


Fixed, In-line, Bracket for Q-Logic 3Standard Part: FRMASMB10617
--



10.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319 Part: ACC1908434. HCPCS: E1032
--



10.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side \$245 Part: ACC2201041. HCPCS: E1032	Basic Swing Away, Q-Logic, Right Side \$245 Part: ACC2201042. HCPCS: E1032
---	--

10.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount \$368 Part: XRF000003. HCPCS: E1032 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.	Multi-Axis Mount \$126 Part: XRF000004 Allows height adjustment of joystick.
---	---



10.2. Joystick Mount Location

Left Side Mount Part: PTOINDV1225	Right Side Mount Part: PTOINDV1226
---	--

11. TRU BALANCE® 4 SEATING

Maximum seat size is 22"W x 22"D.

11.1. TRU Balance 4 Configuration

REQUIRED.

TRU-Balance 4 Power Tilt Only \$8,595 Part: PTOINDV4014. HCPCS: E1002 50° Center of Gravity (CG) power tilt seating system.	TRU-Balance 4 Power Tilt and Recline \$12,790 Part: PTOINDV3954. HCPCS: E1007 50° of posterior tilt and 165° of recline.
--	---

TRU Balance 4 seating section continued on following page.

11. TRU BALANCE 4 SEATING CONTINUED

11.2. 12" Power Lift

Requires combined selection of TRU-Balance 4 Power Tilt and Recline and Articulating Foot Platform. Requires selection in 12.3 Anterior Tilt. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

Power Lift with iLevel \$3,725
 Part: PTOINDV3956. HCPCS: E2298/E2300^
 12" scissor lift with positional feedback.

11.3. Anterior Tilt

Chest support is required with anterior tilt 20 and 30 degree functions. See chest support section. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.

1. Anterior Tilt Angle

Only available with Power Lift.

10° Anterior Tilt Angle No Charge
 Part: MEC2009356
 20° Anterior Tilt Angle \$1,995
 Part: MEC2009399. HCPCS: K0108
 Requires chest support and knee block.
 30° Anterior Tilt Angle \$2,595
 Part: MEC2009400. HCPCS: K0108
 Requires chest support and knee block.
 0° Anterior Tilt No Charge
 Part: PTOINDV3980

11.4. Knee Block Assembly

REQUIRED with 20° or 30° anterior tilt. Optional with 10° and 0° anterior tilt.

Anterior Knee Block Assembly (ST-KBATB4ATC-1-Q) \$314
 Part: ST-KBATB4ATC-1-Q. HCPCS: K0108
 Anterior Knee Block Assembly, Curved (ST-KBATB4ATC-U-Q) .. \$354
 Part: ST-KBATB4ATC-U-Q. HCPCS: K0108

11.5. Occupied Transit Kit

Optional. The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

Occupied Transit Option \$289
 Part: KIT152047
 ISO standard occupied transit.

11.6. Seat Width

REQUIRED with seat selection. 12"-15" widths not available with Power Lift. Maximum seat width is 22". Part number is SETWDTH plus number following # below.

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

11.7. Seat Width Package

REQUIRED with seat selection.

12"-16" Width Package Part: XRF000006 Available with Power Tilt Only.	16"-20" Width Package Part: XRF000008	20"-24" Width Package Part: XRF000158 Maximum width is 22".
---	--	---

11.8. Seat Depth

REQUIRED with seat selection. 12"-13" depth not available with Power Lift. Maximum seat depth is 22". Part number is SETDPH plus number following # below.

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

11.9. Seat Depth Package

REQUIRED with seat selection.

1a. Power Tilt Depth Package

12"-16" Depth Package Part: XRF000011	16"-20" Depth Package Part: XRF000013	20"-24" Depth Package Part: XRF000159 Maximum depth is 22".
--	--	---

1b. Power Lift Depth Package

16"-22" Depth Package
Part: XRF000309

11.10. Seat to Floor Height

REQUIRED.

Low Seat to Floor Part: PTOINDV3976	Middle Seat to Floor Part: PTOINDV3559	Highest Seat to Floor Part: PTOINDV3560
--	---	--


Low	Med	High
18.0"	18.5"	19.0"
TRU Balance 4 Seat to Floor Heights		

TRU Balance 4 Seating section continued on following pages.

11. TRU BALANCE 4 SEATING CONTINUED

11.11. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.



1. Cushion Type

- TRU-Comfort 2 SPP**
Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608
- Premiere P**
Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606
- Essence SPP**
Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608
- Glacial SP**
Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604
- Solution SPP**
Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608
- Zen SP**
Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604
- Spectrum Foam SPP**
Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608
- Simplicity G**
Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602
- Spectrum Gel SPP**
Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

- Small Wedge Positioning Kit**\$39
Part: POP1909006-S
14"W and smaller
- Medium Wedge Positioning Kit**.....\$72
Part: POP1909006-M
14"W and smaller
- Large Wedge Positioning Kit**..... \$105
Part: POP1909006-L
14"W and smaller
- Small Wedge Positioning Kit**\$39
Part: POP1712340-S
15"W and larger
- Medium Wedge Positioning Kit**.....\$72
Part: POP1712340-M
15"W and larger
- Large Wedge Positioning Kit**..... \$105
Part: POP1712340-L
15"W and larger

3. Cushion Top Foam

- Visco Top Foam**
Reference part: XRF000380
- No Top Foam**
Reference part: XRF000382
- FRG Top Foam**
Reference part: XRF000381

4. Cushion Liner

- Stealth Urethane Liner**
Reference part: XRF000383
- No Liner**
Reference part: XRF000385
- Stealth Polyurethane Liner**
Reference part: XRF000384

5. Cushion Cover

- CoolCore Cover**
Reference part: XRF000386
- Stealth Tek Cover**
Reference part: XRF000388
- Mesh Cover**
Reference part: XRF000387

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

- ABS Flat Rigidizer**
Reference part: XRF000389
MSRP: \$48-\$69

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam	Liner		Cover					
TRU-Comfort 2	-	✓	-	✓	-	✓	✓	-	-
Essence	-	✓	-	✓	-	✓	✓	✓	-
Solution	-	✓	-	-	✓	✓	✓	✓	✓
Spectrum Foam	-	✓	✓	-	✓	✓	✓	✓	✓
Spectrum Gel	-	✓	✓	-	✓	✓	✓	✓	✓
Premiere	-	✓	✓	✓	✓	✓	✓	-	-
Glacial	✓	-	-	✓	✓	✓	✓	✓	✓
Zen	-	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	-	-	✓	-	-	-	-	✓

TRU Balance 4 Seating section continued on following pages.

11. TRU BALANCE 4 SEATING CONTINUED

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
 Reference Part: XRF000390

 QTY

Additional Stealth Polyurethane Liner
 Reference Part: XRF000391

 QTY

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
 Reference Part: XRF000392

 QTY

Additional Stealth Tek Cover
 Reference Part: XRF000394

 QTY

Additional Mesh Cover
 Reference Part: XRF000393

 QTY

11.12. Back Type Selection

REQUIRED. Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back
 Part: XRF000014
 Select a back cushion configuration in the subsequent section.
 Back Cane Option
 Part: XRF000015
 Available only with TRU-Balance 3 Power Tilt (1 actuator) selections.
 Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back \$184
 Part: XRF000016
 Not available with power recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

11.13. Back Cane Selection

REQUIRED when Back Cane Option is selected above. Select type and height of back canes.

1. Back Cane Type

Angle Adjustable Back Cane No Charge
 Part: MEC140266
 Flip Forward Back Cane \$53
 Part: XRF000273

2" Endomorph (1" per side) Back Cane \$74
 Part: XRF000269
 4" Endomorph (2" per side) Back Cane \$74
 Part: XRF000271

2. Back Cane Height

16"H Cane
 Part: FRMASMB15834
 17"H Cane
 Part: FRMASMB15835

20"H Cane
 Part: FRMASMB15838
 21"H Cane
 Part: FRMASMB15839

24"H Cane
 Part: FRMASMB15842
 25"H Cane
 Part: FRMASMB15843

11.14. Back Cushion Configuration

REQUIRED. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

1a. Sport Back/Omit/After-Market Interface

Sport Back No Charge
 Part: XRF000017
 Complete following width package and size.
 Omit Back Pan/Cushion No Charge
 Part: PTOINDV1287
 Select a following width package.

After-market Back Interface \$535
 Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.

1b. TRU-Comfort Back

TRU-Comfort Back (Static/Flip Forward), 12" to 20" Width \$744
 Part: TC4V10W21W. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$855
 Part: TC4V22W32W. HCPCS: E2621
 Complete following width package and size.

1c. Q-Back 2.0, AL Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL14AS. HCPCS: E2613
 Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL14AT. HCPCS: E2613
 Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL16AS. HCPCS: E2613
 Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL16AT. HCPCS: E2613
 Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568
 Part: QB2AL16AX. HCPCS: E2613
 Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL16AXX. HCPCS: E2613
 Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL18AT. HCPCS: E2613
 Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$568
 Part: QB2AL18AX. HCPCS: E2613
 Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL18AXX. HCPCS: E2613
 Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL20AS. HCPCS: E2613
 Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL20AT. HCPCS: E2613
 Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$568
 Part: QB2AL20AX. HCPCS: E2613
 Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL20AXX. HCPCS: E2613

TRU Balance 4 Seating section continued on following pages.

11. TRU BALANCE 4 SEATING CONTINUED

1d. Q-Back 2.0, AL Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL14DS. HCPCS: E2620	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL18DT. HCPCS: E2620
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL14DT. HCPCS: E2620	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL18DX. HCPCS: E2620
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL16DS. HCPCS: E2620	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18DXX. HCPCS: E2620
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL16DT. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL16DX. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL20DT. HCPCS: E2620
Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL16DXX. HCPCS: E2620	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL20DX. HCPCS: E2620
Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL18DS. HCPCS: E2620	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20DXX. HCPCS: E2620

1e. Q-Back 2.0, AL Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL14ES. HCPCS: E2620	Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL18EX. HCPCS: E2620
Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$703 Part: QB2AL16ES. HCPCS: E2620	Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18EXX. HCPCS: E2620
Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$703 Part: QB2AL16ET. HCPCS: E2620	Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL20ES. HCPCS: E2620
Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL16EX. HCPCS: E2620	Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703 Part: QB2AL20ET. HCPCS: E2620
Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL18ES. HCPCS: E2620	Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL20EX. HCPCS: E2620
Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$703 Part: QB2AL18ET. HCPCS: E2620	Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20EXX. HCPCS: E2620

1f. Q-Back 2.0, AFT Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT16AS. HCPCS: E2613	AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT18AS. HCPCS: E2613
--	--

1g. Q-Back 2.0, AFT Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$643 Part: QB2AFT16DS. HCPCS: E2620	AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$643 Part: QB2AFT18DT. HCPCS: E2620
AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$643 Part: QB2AFT16DT. HCPCS: E2620	AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AFT20DT. HCPCS: E2620
AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AFT18DS. HCPCS: E2620	

1h. Q-Back 2.0, AFT Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AFT16EXX. HCPCS: E2620	
---	--

11.15. TRU Balance 4 Back Shroud

REQUIRED. The back shroud is a feature of the TRU Balance 4 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Inferno Red Part: DGN1714085-PRD5	Ice Pink Part: DGN1714085-PPK5	Bermuda Blue Part: DGN1714085-PBU15	Faded Glory Part: DGN1714085-WFL2
Back in Black Part: XRF000021	Pretty in Pink Part: XRF000027	Blue Jean Part: XRF000031	Black Camouflage Part: DGN1714085-WCM9-PCL2
Night Train Part: XRF000022	Evolution Orange Part: XRF000028	Lemon Crush Part: XRF000032	
Ghost Grey Part: DGN1714085-PGY6	Steel Blue Part: DGN1714085-PBU13	Forest Green Part: DGN1714085-PGN13	
Cosmic Purple Part: DGN1714085-PPR8	Royale Blue Part: DGN1714085-PBU14	Military Digital Camo Part: DGN1714085-WCM7	

TRU Balance 4 Seating section continued on following pages.

12. ARMRESTS CONTINUED

12.1b. Two-Post, Height Adjustable, Removable Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

L R Full Length, Two-Post Armrest for Power Seating No Charge
Left part: XRF000047. Right part: XRF000048

L R Desk Length, Two-Post Armrest for Power Seating.... No Charge
Left part: XRF000049. Right part: XRF000050

12.1c. Pediatric, Removable, Height Adjustable Armrests

Height adjustment from 7.75" to 9.75" in .5" increments. **Not available with full length armpads or waterfall armpads.**

L R Pediatric, Removable Armrests \$79 ea.
Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973

12.2. Armrest Height

REQUIRED. Set left and right armrest height. Note the available height range in the option description.

1a. Flip-Back and Two-Post Left Armrest Height

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

1b. Pediatric, Left Armrest Height

7.75" H PTOINDV3722	8.25" H PTOINDV3724	8.75" H PTOINDV3726	9.25" H PTOINDV3728	9.75" H PTOINDV3730
------------------------	------------------------	------------------------	------------------------	------------------------

2a. Flip-Back and Two-Post Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

2b. Pediatric, Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	8.75" H PTOINDV3725	9.25" H PTOINDV3727	9.75" H PTOINDV3729
------------------------	------------------------	------------------------	------------------------	------------------------

12.3. Armrest Size

REQUIRED. Select a left (L) and right (R) armrest size. This also affects armpad selection.

L R 10" Length
L part: XRF000057. R part: XRF000123

L R 14" Length
L part: XRF000058. R part: XRF000124

L R 17" Length
L part: XRF000059. R part: XRF000125
Only available with Flip-Back, 8"-16" H
Armrests and straight arm pads.

13. ARPADS

REQUIRED.

13.1a. Omit Armpad

Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection..... No Charge
Part: XRF000060
Not available with Flip-back, 8"-16" H armrests.

Interface for Aftermarket Pad No Charge
Part: XRF000308
Only available with Flip-back, 8"-16" armrests. Not available with chest bar.

13.1b. Foam Armpads

L R Straight, Armpad Standard
Left part: XRF000061. Right part: XRF000126

L R Waterfall, Armpad No Charge
Left part: XRF000062. Right part: XRF000127
Not available with Flip-back, 8"-16" armrests.

13.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right\$195
Part: ARM1806807
Falcon Gel, Full Length Left, Desk Length Right.....\$171
Part: ARM1806813

Falcon Gel, Desk Length Left, Full Length Right.....\$171
Part: ARM1806810
Falcon Gel, Desk Length Left and Right\$147
Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right\$311.22
Part: ARM136438
Falcon Gel, Full Length Left, Desk Length Right.....\$287.93
Part: ARM136442

Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
Part: ARM136439
Falcon Gel, Desk Length Left and Right \$265.02
Part: ARM136437

14. ARMREST BRACKETS

Standard Fixed Receiver Standard
Part: XRF000063

Stealth Height Adjustable Receiver \$263
Part: XRF000064. HCPCS: K0108

15. HEADRESTS

Optional. Select a headrest pad and bracket.

15.1. Headrest Pads, Redesign Back Assembly

Only available with 14" to 22" back widths of Static Type or Flip Forward back, without occupied transit.

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955

Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q. HCPCS: E0955

Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955

Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033

Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033

15.2. Q Cam Lock

Optional.

Q Cam Lock..... \$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

16. LAP BELTS

REQUIRED. Select a lap belt below.

16.1a. Standard Lap Belts

- 38" Lap Belt Part: ACCBELT1006 50" Lap Belt Part: ACCBELT1000 60" Lap Belt Part: ACCBELT1002 70" Lap Belt Part: ACCBELT1003 80" Lap Belt Part: ACCBELT1004 90" Lap Belt Part: ACCBELT1001

16.1b. Optional Lap Belts

Stealth Positioning, 60" Single Pull, Push Button, Medium..... \$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q

Stealth Positioning, 60" Single Pull, Push Button, Extra Large.....\$50
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q

Stealth Positioning, 60" Single Pull, Push Button, Large \$50
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q

60" Sure-Lok Lap Belt (Non-retractable)..... \$120
Part: ACCBELT1016
Required with ISO 7176-19 Compliant Transit Option.

17. LEGRESTS

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

17.1a. Basic Center Mount Foot Platform

1. Basic Center Mount Foot Platform Assembly

Basic Center Mount..... No Charge
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates \$180
Part: XRF000312. HCPCS: K0040

Long Extension, 5"x9" Footplates \$180
Part: XRF000314. HCPCS: K0040

Medium Extension, 5"x9" Footplates \$180
Part: XRF000313. HCPCS: K0040

3. Basic Center Mount Calf Pads

Small Calf Pads, Basic Center Mount..... \$158
Part: RIG2104369

Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997
Only available with Long Extension.

4. Basic Center Mount Additional Accessory

Basic Center Mount Footplate Connector \$16
Part: RIG2104240

17.1b. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Center Mount Foot Platform Type

Dual Actuator, Power Articulating Foot Platform..... \$3,595
Part: RIG2009246. HCPCS: E1012+K0040
Positional feedback with independent control of articulation.

Quick Adjustable Center Mount Foot Platform..... \$263
Part: XRF000311. HCPCS: K0040
Not available with power lift.

2. Extension and Angle Adjustable Footplates

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x8" Footplates
Part: RIG2404015

Long Extension, 5"x6" Footplates
Part: RIG151505

Medium Extension, 5"x8" Footplates
Part: RIG2404016

Long Extension, 5"x8" Footplates
Part: RIG2404017

Medium Extension, 5"x11" Footplates
Part: RIG2404018

Long Extension, 5"x11" Footplates
Part: RIG2404019

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

Legrests section continued on following page.

17. LEGRESTS CONTINUED

3a. Center Mount Calf Pad

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
 Part: RIGASMB7120018
 Small Calf Pads, Quick Adjustable Center Mount Only, \$158
 6"Wx6"H (pair)
 Part: XRF000156

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
 Part: RIGASMB7120019
 Large Calf Pads, Quick Adjustable Center Mount Only, \$158
 8" taper to 6"Wx8"H (pair)
 Part: XRF000157

3b. Center Mount Calf Panel

Available with Quick Adjustable center mount.

Small Calf Panel, 18"Wx8"H \$184
 Part: FRMASMB7110012

Large Calf Panel, 24"Wx8"H \$237
 Part: RIG123825
 Only available on widths of 20" or greater.

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
 Part: RIG152519
 Only available with the 5x6 footplates.
 Width Extension for 8" \$69
 Part: RIG2404919
 Width Extension for 11" \$69
 Part: RIG2404920

Outside Edge \$69
 Part: RIG2404921
 Available with 5"x8" and 5"x11" footplates.
 Front and Rear Cap \$16
 Part: RIG2404923
 Available with 5"x8" and 5"x11" footplates.

5. Center Mount Additional Accessories

Small Foam Calf Pad Cover (pair) \$67
 Part: ST-CSCS0775-Q
 Small Foam Panel Cover (pair) \$176
 Part: ST-CSCP0920-Q
 Footplate Connector \$16
 Part: RIGASMB7120023

Footplate Riser \$75
 Part: XRF000348
 Not available with long extension. Select legrest length between 2" to 9" at 1" increments.
 Footplate Riser Calf Panel \$50
 Part: FRMASMB7110011
 Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222
 Part: ST-FPG5008T-Q. HCPCS: K0108
 5"x11" Tapered Gel Footplate Covers \$270
 Part: ST-FPG5011T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension \$247
 Part: ST-FPG5008E-Q. HCPCS: K0108
 5"x11" Tapered Gel Footplate Cover for Width Extension \$270
 Part: ST-FPG5011E-Q. HCPCS: K0108

17.2. Legrest Accessories

Wheel Bumper \$37
 Part: XRF000100

17.3. Legrest Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)
11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	

18. POWER POSITIONING ELECTRONICS

REQUIRED.

18.1a. Q-Logic 3 Controls through Joystick

Single Actuator Through Controller and Interface \$2,163
 Part: XRF000335. HCPCS: E2310
 Single Actuator Through Controller and Interface (AAM) \$3,040
 Part: XRF000068. HCPCS: E2310

Multiple Actuator Through Controller and Interface \$3,040
 Part: XRF000336. HCPCS: E2311

18.1b. Q-Logic 3 Controls with iAccess Switch Control

Single Actuator Control Through Flat Button \$2,547
 iAccess and Controller and Interface
 Part: XRF000337. HCPCS: E2310
 Single Actuator Control Through Toggle Button \$2,547
 iAccess and Controller and Interface
 Part: XRF000338. HCPCS: E2310
 Single Actuator Control Through Flat Button \$3,460
 iAccess and Controller and Interface (AAM)
 Part: XRF000070. HCPCS: E2310

Single Actuator Control Through Toggle Button \$3,460
 iAccess and Controller and Interface (AAM)
 Part: XRF000074. HCPCS: E2310
 Multiple Actuator Control Through Flat Button \$3,460
 iAccess and Controller and Interface
 Part: XRF000339. HCPCS: E2311
 Multiple Actuator Control Through Toggle Button \$3,460
 iAccess and Controller and Interface
 Part: XRF000340. HCPCS: E2311

Power positioning electronics section continued on following page.

18. POWER POSITIONING ELECTRONICS CONTINUED

18.2. iAccess Mount

REQUIRED with iAccess options.

1. iAccess Mount Selection

Joystick Mount No Charge Part: MEC1907500	Attendant Mount \$26 Part: MEC1906564
Under Arm Mount, Flip-Back Armrest..... No Charge Part: MEC1906531	Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$26 Part: MEC2001392
Under Arm Mount, Tube Style Armrest No Charge Part: MEC1906530	Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$26 Part: MEC2001391

2. Under Arm Mount Side

Left Side Under Arm Mount Part: PTOINDV3678	Right Side Under Arm Mount Part: PTOINDV3679
--	---

19. ATTENDANT CONTROL


Optional. Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control \$525 Part: XRF000078. HCPCS: E2331 Select this option if no charge option criteria is not met.	Attendant Control with Tilt, Recline, Alternative No Charge Controls and Q-Logic Display or SCIM Module Part: XRF000079
---	---

20. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.

Accessory Bar \$132 Part: XRF000128 Horizontal rail attached to rear of seat. Standard with the depth adjustable vent tray and oxygen holder.	Oxygen Holder..... \$210 Part: ACC125006. HCPCS: E2208 Includes accessory bar, a credit will be applied PTOIN-DV3769.
Backpack Holder (Pair) (A)..... \$69 Part: ACC147774	Personal Item Hook..... \$69 Part: ACC125008
Cane/Crutch Holder \$105 Part: ACC152515. HCPCS: E2207	Personal Item Pouch \$84 Part: ACC2101017
Clothing Guard (Pair)..... \$237 Part: ACC125007	Privacy Bag..... \$42 Part: ACC2005723
Compact Vent Tray \$1,465 Part: ACC133889. HCPCS: E1029	Push Handles (pair)..... \$195 Part: ACC2005027 Not available with oxygen holder or vent tray.
Cup Holder No Charge Part: ACC125003	Quantum Backpack \$168 Part: ACC1905853
Flag Kit \$69 Part: ACC1700443	Rear View Mirror \$69 Part: ACC155833
Flip-Back Armrest, Side Track Accessory Bracket ... \$11 Part: MEC1809178	Reflector Kit \$53 Part: ACC124670
Front Mounted, Cell Phone/ Media Holder, \$69 X-Grip (E) Part: ACC155835	Stealth Large Tablet Holder \$420 Part: ACC141818
Glove Box (C)..... \$158 Part: ACC125011	Transfer Handle (D)..... \$105 ea. Part: MEC2205680 x2. HCPCS: E0705
Grip Drink Holder (20085) \$15.50 Part: ACC152993	Depth Adjustable Vent Tray \$595 Part: ACC2302500. HCPCS: E1029 Includes accessory bar, a credit will be applied PTOIN-DV3769.
Grip Lap Board (20005)..... \$38.95 Part: ACC152992	Remote Stop Switch \$1,295 Part: XRF000379
Hydration System (B)..... \$237 Part: ACC125010	
Large, Black, Weather Cover \$202 Part: ACCCOVR1021	



Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea. Part: ACC2407046	Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea. Part: ACC2405496
--	---

21. BRANCH PATCHES

Optional.

Air Force Patch\$21 Part: ACC140240	Coast Guard Patch\$21 Part: ACC140239	Navy Patch\$21 Part: ACC140244
Army Patch\$21 Part: ACC140241	Marines Patch\$21 Part: ACC140243	Purple Heart Patch\$21 Part: ACC140238



22. POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

22.1. Stealth Unilink Pelvic/Thigh Guide

1. UniLink Pelvic/Thigh Hardware

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea. Removable Mounting Hardware QTY Part: ST-UTTP5-Q, HCPCS: E1028/E1034^	UniLink Flip Down Hip/Thigh Support, Left..... \$332 ea. Part: ST-UTS25100-L-Q, HCPCS: E1028/E1034^
UniLink, 7" Height Adjustable, \$276 ea. Removable Mounting Hardware QTY Part: ST-UTTP7-Q, HCPCS: E1028/E1034^	UniLink Flip Down Hip/Thigh Support, Right \$332 ea. Part: ST-UTS25100-R-Q, HCPCS: E1028/E1034^
UniLink Multi-Angle Thigh Guide, Left \$358.87 ea. Part: ST-UTT8A11602-L-Q, HCPCS: E1028/E1034^	Flip Down with LEFT Full Surface Contact \$281 ea. Feature (TWBADD-L & LPHW407) Part: POS140428, HCPCS: E1028/E1034^
UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea. Part: ST-UTT8A11602-R-Q, HCPCS: E1028/E1034^	Flip Down with RIGHT Full Surface Contact..... \$281 ea. Feature (TWBADD-R & LPHW407) Part: POS140429, HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea. QTY Part: ST-UTP544056100-Q, HCPCS: E0956/E0953^	UniLink, 4" x 12" Pad \$106 ea. QTY Part: ST-UTP54D056100-Q, HCPCS: E0956/E0953^
UniLink, 4" x 8" Pad..... \$106 ea. QTY Part: ST-UTP548056100-Q, HCPCS: E0956/E0953^	

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link\$81 ea. QTY Part: ST-UIL001150A000-Q	UniLink, 2" Link.....\$81 ea. QTY Part: ST-UIL001200A000-Q
--	---

4. UniLink Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

22.2. Stealth Thoracic Laterals

Swing-away laterals may interfere with 2-post, flip up armrests.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$664 Part: ST-ULS4042000-Q, HCPCS: E1034 Only available with Redesigned back assemblies.	Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$332 Part: ST-ULPK404103L00-Q, HCPCS: E1034 Only available with Active Contour Q-Backs.
Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664 Part: ST-ULPK404108B00-Q, HCPCS: E1034 Only available with Sport or TRU Comfort backs.	Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$332 Part: ST-ULPK404103R00-Q, HCPCS: E1034 Only available with Active Contour Q-Backs.
Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664 Part: ST-ULPK404104B00-Q, HCPCS: E1034 Only available with back cane and omit back configurations.	

2. Pads

Curved Pads, 3"x5", pair (TWBL3X5)..... \$250 Part: INDPART2663, HCPCS: E0956	Flat Pads, 3"x5", pair (ULP846050320)..... \$250 Part: ST-ULP846050320-Q, HCPCS: E0956
Curved Pads, 4"x6", pair (TWBL4X6)..... \$250 Part: INDPART2664, HCPCS: E0956	Flat Pads, 4"x6", pair (ULP835050320) \$250 Part: ST-ULP835050320-Q, HCPCS: E0956

3. Installation

Install Laterals Part: PTOINDV3499	Do Not Install Part: PTOINDV3500
---------------------------------------	-------------------------------------

22.3. Chest Bar Assembly

Chest support is required with anterior tilt 20 and 30 degree functions. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.

Small Chest Bar, 14" Foam Pad (ST-CBATB4ATT-S-Q)..... \$368 Part: ST-CBA-TB4-S, HCPCS: K0108 Available with 14" to 20" back widths.	Medium Chest Bar, 16" Foam Pad (ST-CBATB4ATT-M-Q)..... \$368 Part: ST-CBA-TB4-M, HCPCS: K0108 Available with 17" to 22" back widths.
---	--

Positioning components section continued on following page.

22. POSITIONING COMPONENTS CONTINUED

22.4. Stealth Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$154 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$154 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$154 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960



Contoured and Classic Chest Harness

22.5. Stealth Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor, (IFCS4113)..... \$97 ea.
Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor, (IFCS4114)..... \$97 ea.
Part: ST-IFCS4114-Q. HCPCS: E0960



Chest Strap

22.6. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right..... \$318
Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left \$318
Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$434
2 Upper
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028



Left and Right Facial Lateral with Pads

2. Facial Lateral Pads and Accessories


Egg Switch Mount (ESM) (A)\$57 ea.
Part: INDPART2815

Swing-Away Joystick Mount, 18" Rod (C)..... \$276 ea.
Part: ST-SUS18-Q

Facial Pad, 4".5x2", Cool Core (900M-Q)\$61 ea.
Part: ST-900M-Q. HCPCS: K0108

Allows mounting of facial pads and attaches to the headrest

Facial Pad, 6"x2", Cool Core (920M-Q).....\$61 ea
Part: ST-920M-Q. HCPCS:K0108



Headrest Accessories

22.7. Stealth Elbow Support


1. Elbow Support Bracket

Left, Swing-Away, Elbow Support \$199 ea.
Bracket Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support..... \$199 ea.
Bracket Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket..... \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
Part: ST-UETB4021FL-Q

Right, UniLink, Fixed, Elbow Support Bracket.....\$211
for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
Part: ST-UETB4021FR-Q



Elbow Supports

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$57 ea.
Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4).....\$57 ea.
Part: POS124067. HCPCS: K0108

22.8. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

22.9. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

23. LAP TRAY

23.1. Quantum Lap Tray

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

- Black Lap Tray, 11" W \$200
Part: MEC1905328-11B. HCPCS: E0950
For use with 12"-13" back width.
- Black Lap Tray, 13" W \$200
Part: MEC1905328-13B. HCPCS: E0950
For use with 14"-15" back width.
- Black Lap Tray, 15" W \$200
Part: MEC1905328-15B. HCPCS: E0950
For use with 16"-17" back width.

- Black Lap Tray, 17" W \$200
Part: MEC1905328-17B. HCPCS: E0950
For use with 18"-19" back width.
- Black Lap Tray, 19" W \$200
Part: MEC1905328-19B. HCPCS: E0950
For use with 20"-22" back width.

1b. Clear Lap Tray

- Clear Lap Tray, 11" W \$200
Part: MEC1905328-11C. HCPCS: E0950
For use with 12"-13" back width.
- Clear Lap Tray, 13" W \$200
Part: MEC1905328-13C. HCPCS: E0950
For use with 14"-15" back width.
- Clear Lap Tray, 15" W \$200
Part: MEC1905328-15C. HCPCS: E0950
For use with 16"-17" back width.

- Clear Lap Tray, 17" W \$200
Part: MEC1905328-17C. HCPCS: E0950
For use with 18"-19" back width.
- Clear Lap Tray, 19" W \$200
Part: MEC1905328-19C. HCPCS: E0950
For use with 20"-22" back width.

2. Mounting Hardware

- Mounting Bracket, Chest Bar/Lap Tray \$75
Part: MEC2104886

23.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

- Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950
- Simple Solutions Medium Tray with Cut Out \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950
- Simple Solutions Large Tray with Cut Out \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

- Swing Away Mounting Hardware \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	S	M	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	S	M	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	S	M	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	S	M	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						S	M	ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	M	S	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q
 :Compatible :Armrest configuration may not hold tray position
 *Some armrests listed may not be available on this order form.

24. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).

24.1. Enhanced Display

Q-Logic 3 Enhanced Display Kit..... \$1,980
 Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard Part: DGN137728	Cosmic Purple \$63 Part: DGN137729-PPR8	Royale Blue \$63 Part: DGN137729-PBU14	Military Digital Camo \$63 Part: DGN137729-WCM7
Inferno Red \$63 Part: DGN137729-PRD5	Ice Pink \$63 Part: DGN137729-PPK5	Bermuda Blue \$63 Part: DGN137729-PBU15	Faded Glory \$63 Part: DGN137729-WFL2
Back in Black \$63 Part: DGN137729-PBK1	Pretty in Pink \$63 Part: DGN137729-PPK3	Blue Jean \$63 Part: DGN137729-PBU11	Black Camouflage \$63 Part: DGN137729-WCM9-PCL2
Night Train \$63 Part: DGN137729-PBK3	Evolution Orange \$63 Part: DGN137729-POR7	Lemon Crush \$63 Part: DGN137729-PYL4	
Ghost Grey \$63 Part: DGN137729-PGY6	Steel Blue \$63 Part: DGN137729-PBU13	Forest Green \$63 Part: DGN137729-PGN13	

24.2. SCIM Module

Optional.

SCIM Module for Q-Logic 3 \$525
 Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

24.3. Q-Logic 3 SAJ

Q-Logic 3 EX SAJ..... \$840
 Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

24.4. Sip and Puff Drive

1. Sip and Puff Module

Sip and Puff Module for Q-Logic 3..... \$1,990
 Part: CTL168056. HCPCS: E2325

Includes module and mounting hardware.



Sip & Puff Module

2. Sip and Puff Straw and Mount

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

Stealth Sip and Puff for Comfort Series, Left Mount (SP625M-CL) \$318
 Part: ST-SP625M-CL-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, Head Array (SP625M-U) \$318
 Part: ST-SP625M-U-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series, Right Mount (SP622M-CR) \$318
 Part: ST-SP625M-CR-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Stealth Sip & Puff

3. Sip & Puff Accessories

Sip & Puff Straws, 12 pack (SP625-S12)..... \$35
 Part: ST-SP625-S12-Q

24.5. Egg Switches

Red Egg Switch (ESRED)..... \$87
 Part: SWTMCR01077

Blue Egg Switch (ESBLU) \$87
 Part: SWTMCR01079

Black Egg Switch (ESBLK)..... \$87
 Part: SWTMCR01081

Green Egg Switch (ESGRN) \$87
 Part: SWTMCR01078

Yellow Egg Switch (ESYEL) \$87
 Part: SWTMCR01080



Egg Switch

Alternative controls and electronics section continued on following page.

24. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

24.6. Switch Adapter

2 to 1 Switch Adapter \$15
Part: DWR1361H037
Connects 2 single switches to a stereo jack to act as a dual switch.

24.7. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520
Part: CTL168057

24.8. USB Mobile Device Charger

USB Mobile Device Charger \$184
Part: ELE168069
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

