



Q6 Edge® 4 Order Form with TRU Balance® Flex Seating

Account Number: _____ Date: _____
 Provider Name: _____
 Contact: _____
 Phone: _____ Fax: _____
 Email: _____
 PO Number: _____
 Quote Name/Marked For: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
 B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
 C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 7mph programmable motors and off board charger. XLR USB charger is included. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 4 3S-SS \$8,870
 Part: Q6 EDGE 4 3S-SS. HCPCS: K0848, Pending PDAC Code Verification
 Solid seat pan standard with no power option.

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 7 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology \$1,360
 Coming soon. Part: XRF000268. HCPCS: K0108

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

- | | | | |
|---|--|---|---|
|  Inferno Red
Part: DGN2401066-PRD5 |  Cosmic Purple
Part: DGN2401066-PPR8 |  Steel Blue
Part: DGN2401066-PBU13 |  Lemon Crush
Part: DGN2401066-PYL4 |
|  Back in Black
Part: DGN2401066-PBK1 |  Ice Pink
Part: DGN2401066-PPK5 |  Royale Blue
Part: DGN2401066-PBU14 |  Forest Green
Part: DGN2401066-PGN13 |
|  Night Train
Part: DGN2401066-PBK3 |  Pretty in Pink
Part: DGN2401066-PPK3 |  Bermuda Blue
Part: DGN2401066-PBU15 | |
|  Ghost Grey
Part: DGN2401066-PGY6 |  Evolution Orange
Part: DGN2401066-POR7 |  Blue Jean
Part: DGN2401066-PBU11 | |

3.1b. Hydrographic Shrouds

- | | | |
|--|--|---|
| <input type="checkbox"/>  Military Digital Camo \$263
Part: DGN2401066-WCM7 | <input type="checkbox"/>  Faded Glory \$263
Part: DGN2401066-WFL2 | <input type="checkbox"/>  Black Camouflage \$263
Part: DGN2401066-WCM9-PCL2
<i>Coming soon.</i> |
|--|--|---|

4. BASE MOUNTED LIGHTS

Optional. Select one option.

- | | |
|--|--|
| Accent Light Package Standard
Part: ELE2404196 | Full Light Package \$1,045
Part: ELE2404199
Requires QLNES 10 Key, Q-Logic 3 EX Joystick for Full Lights, or Q-Logic 3e Joystick. |
|--|--|

5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

- Quantum Backup Camera** \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

- | | |
|--|---|
| Backup Camera Headrest Bracket \$21
Part: MEC1910373
Works with .5" x .5" vertical supports only. | Backup Camera Track Mount Bracket \$27
Part: MEC2104227 |
|--|---|

2. Optional Display Mounting Bracket

- Backup Camera Q-Logic Joystick Mount Bracket** \$105
Part: MEC2106123



6. DRIVE WHEELS

REQUIRED. Select a type of drive wheel.

- | | |
|---|--|
| 3"x14" Flat Free Tires Standard
Part: PTOINDV3525 | 4"x14" Knobby, Flat Free Black Tires \$580
Part: XRF000368 |
| 3"x14" Pneumatic Tires No Charge
Part: WHL1713925 | |

7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered. Uses 8 amp charger.

- | | |
|--|--|
| Group 24, Gel Cell, Interceptor 2 Battery \$399 ea.
Part: BAT2404214. HCPCS: E2361
Battery tray configuration. | NF-22, AGM, MK Flat Top Battery \$432 ea.
<i>Coming soon.</i> Part: BAT2501912. HCPCS: E2361
Without battery tray configuration |
| NF-22, Gel Cell, Interceptor 2 Battery \$355 ea.
<i>Coming soon.</i> Part: BAT2500310. HCPCS: E2361
Without battery tray configuration. | |

8. ELECTRONICS

REQUIRED. Select one type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control found on the [Alternative Controls and Electronics Order Form](#).

8.1a. Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

- | | |
|--|--|
| QLNES 6 Key JoystickStandard
Part: CTL2205379
Controls up to 2 actuators. | QLNES 10 Key Joystick \$420
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package. |
|--|--|



Electronics section continued on following page.

8. ELECTRONICS CONTINUED

8.1b. Expandable Controls

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

- Q-Logic 3 EX Controller Module \$730
Part: XRF000002. HCPCS: E2377
- Q-Logic 3 Expandable Controller Harness..... \$458
Part: XRF000264. HCPCS: E2313



2. Q-Logic 3 Joystick

- Q-Logic 3 EX Joystick \$840
Part: CTL167724
- Q-Logic 3 EX Joystick for Full Lights \$840
Part: CTL167725
- Q-Logic 3e Joystick \$735
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.
- No Joystick, Setup for Q-Logic 3
Part: XRF000010
Requires selection of Enhanced Display

8.2. Q-Logic Magnetic Charger

Available with QLNES. Included with Q-Logic 3e Joystick.

- Magnetic to XLR Charger Adapter with Wall Holder \$60
Part: ELE2103635

8.3. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

- | | | |
|--|---|--|
| Black, Soft Touch Standard
Part: DGN124307
Also matches Night Train shrouds. | Pretty in Pink \$63
Part: DGN124336-PPK3 | Lemon Crush \$63
Part: DGN124336-PYL4 |
| Inferno Red \$63
Part: DGN124336-PRD5 | Evolution Orange \$63
Part: DGN124336-POR7 | Forest Green \$63
Part: DGN124336-PGN13 |
| Back in Black \$63
Part: DGN124336-PBK1 | Steel Blue \$63
Part: DGN124336-PBU13 | Military Digital Camo \$63
Part: DGN124336-WCM7 |
| Ghost Grey \$63
Part: DGN124336-PGY6 | Royale Blue \$63
Part: DGN124336-PBU14 | Faded Glory \$63
Part: DGN124336-WFL2 |
| Cosmic Purple \$63
Part: DGN124336-PPR8 | Bermuda Blue \$63
Part: DGN124336-PBU15 | Black Camouflage \$63
Part: DGN124336-WCM9-PCL2 |
| Ice Pink \$63
Part: DGN124336-PPK5 | Blue Jean \$63
Part: DGN124336-PBU11 | |

9. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

- Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323
- Chin Cup (C) \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)
- T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323
- Medium Ball Handle (E) \$69
Part: ACCASMB1859. HCPCS: E2323
- Soft Ball Handle (F) \$80
Part: ACCASMB1856. HCPCS: E2323
- Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323
- Bodypoint 3" Goal Post Handle \$104
Part: PLSKNCP1057. HCPCS: E2323
- Bodypoint 4" Goal Post Handle \$104
Part: PLSKNCP1058. HCPCS: E2323
- Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308



10. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.


10.1a. Fixed, In-line Joystick Bracket

- Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617

10.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

- Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1032



10.1c. Basic Swing Away Joystick Mount

- Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1032
- Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1032

Joystick mounting brackets section continued on following pages.


10. JOYSTICK MOUNTING BRACKETS CONTINUED

10.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount.....\$368
 Part: XRF000003. HCPCS: E1032
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount.....\$126
 Part: XRF000004
 Allows height adjustment of joystick.



Retract4



Multi-Axis Joystick Mount

10.2. Joystick Mount Location

Left Side Mount
 Part: PTOINDV1225

Right Side Mount
 Part: PTOINDV1226

11. TRU BALANCE FLEX SEATING SYSTEM

Seat width and depth can be adjusted in 1" increments from 15"W-20"W and 16"D-20"D.

11.1. TRU Balance Flex Frame Selection

REQUIRED.

TRU Balance Flex Seating SystemNo Charge
 Part: PTOINDV3853

11.2. Occupied Transit Kit

Optional. The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

Occupied Transit Option\$289
 Part: KIT152047
 ISO standard occupied transit.

11.3. Seat Width

REQUIRED. Maximum seat width is 20". Part number is SETWDTH plus number following # below.

10" W (#1000)	12" W (#1002)	14" W (#1004)	16" W (#1006)	18" W (#1008)	20" W (#1010)	22" W (#1012)
11" W (#1001)	13" W (#1003)	15" W (#1005)	17" W (#1007)	19" W (#1009)	21" W (#1011)	

11.4. Seat Depth

REQUIRED. Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

10" D (#1000)	12" D (#1002)	14" D (#1004)	16" D (#1006)	18" D (#1008)	20" D (#1010)
11" D (#1001)	13" D (#1003)	15" D (#1005)	17" D (#1007)	19" D (#1009)	

High	18.5" at front with 4° posterior angle
Med	17.5" at front with 2° posterior angle
TRU Balance Flex Seat to Floor Heights	

11.5. Seat to Floor Height


Lowest Seat to Floor
 Part: PTOINDV3558

Middle Seat to Floor
 Part: PTOINDV3559

Highest Seat to Floor
 Part: PTOINDV3560

11.6. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.



TRU Balance Flex seating section continued on following pages.

11. TRU BALANCE FLEX SEATING SYSTEM CONTINUED

1. Cushion Type

TRU-Comfort 2 SPP
 Reference part: XRF000116. *See cushion order form.*
 Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPP
 Reference part: XRF000118. *See cushion order form.*
 Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPP
 Reference part: XRF000113. *See cushion order form.*
 Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPP
 Reference part: XRF000112. *See cushion order form.*
 Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPP
 Reference part: XRF000117. *See cushion order form.*
 Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere P
 Reference part: XRF000122. *See cushion order form.*
 Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SP
 Reference part: XRF000114. *See cushion order form.*
 Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SP
 Reference part: XRF000115. *See cushion order form.*
 Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity G
 Reference part: XRF000111. *See cushion order form.*
 Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39 Part: POP1909006-S 14"W and smaller	Medium Wedge Positioning Kit\$72 Part: POP1909006-M 14"W and smaller	Large Wedge Positioning Kit \$105 Part: POP1909006-L 14"W and smaller
Small Wedge Positioning Kit\$39 Part: POP1712340-S 15"W and larger	Medium Wedge Positioning Kit\$72 Part: POP1712340-M 15"W and larger	Large Wedge Positioning Kit \$105 Part: POP1712340-L 15"W and larger

3. Cushion Top Foam

Visco Top Foam Reference part: XRF000380	No Top Foam Reference part: XRF000382
FRG Top Foam Reference part: XRF000381	

4. Cushion Liner

Stealth Urethane Liner Reference part: XRF000383	No Liner Reference part: XRF000385
Stealth Polyurethane Liner Reference part: XRF000384	

5. Cushion Cover

CoolCore Cover Reference part: XRF000386	Stealth Tek Cover Reference part: XRF000388
Mesh Cover Reference part: XRF000387	

	None	Visco/V. Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam	Liner	Cover						
TRU-Comfort 2	-	✓	-	✓	-	✓	✓	-	-
Essence	-	✓	-	✓	-	✓	✓	-	-
Solution	-	✓	-	✓	✓	✓	✓	✓	✓
Spectrum Foam	-	✓	✓	-	✓	✓	✓	✓	✓
Spectrum Gel	-	✓	✓	-	✓	✓	✓	✓	✓
Premiere	-	✓	✓	✓	✓	✓	✓	-	-
Glacial	✓	-	-	✓	✓	✓	✓	✓	✓
Zen	-	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	-	-	✓	-	-	-	-	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
 Reference part: XRF000389
 MSRP: \$48-\$69

11.7. TRU Balance Flex Back Cane Selection

REQUIRED.

Angle Adjustable Back Cane No Charge Part: MEC140266	2" Endomorph (1" per side) Back Cane \$205 Part: XRF000270
Flip Forward Back Cane \$184 Part: XRF000274 Not available with oxygen holder, or the flip back, cane mount armrest.	4" Endomorph (2" per side) Back Cane \$205 Part: XRF000272
2" Mesomorph (1" per side) Back Cane \$205 Part: MEC1804736	

11.8. TRU Balance Flex Back Cane Height Selection

REQUIRED. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests.

14" High Canes Part: FRMASMB15832	17" High Canes Part: FRMASMB15835	20" High Canes Part: FRMASMB15838	23" High Canes Part: FRMASMB15841
15" High Canes Part: FRMASMB15833	18" High Canes Part: FRMASMB15836	21" High Canes Part: FRMASMB15839	24" High Canes Part: FRMASMB15842
16" High Canes Part: FRMASMB15834	19" High Canes Part: FRMASMB15837	22" High Canes Part: FRMASMB15840	25" High Canes Part: FRMASMB15843

TRU Balance Flex seating section continued on following pages.

11. TRU BALANCE FLEX SEATING SYSTEM CONTINUED

11.9. TRU Balance Flex Back Cushion Options

Select one type of back. See the seat back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

1a. Sport Back Selection

Sport Back Cushion, 21" Wide or Less..... \$394
Part: XRF000361

Sport Back Cushion, 22" Wide or More..... \$604
Part: XRF000362

1b. TRU-Comfort Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width..... \$744
Part: TC4V10W21WBC. HCPCS: E2620

TRU-Comfort Back, Back Cane Option, 22" Width \$855
Part: TC4V22W32WBC. HCPCS: E2621

		Back Width												
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"			
Overall Height	14"	S	T	-	-	-	-	-	-	-	-	-	-	
	15"	S	T	-	-	-	-	-	-	-	-	-	-	
	16"	S	T	S	T	S	T	S	T	S	T	S	T	
	17"	S	T	S	T	S	T	S	T	S	T	S	T	
	18"	-	S	T	S	T	S	T	S	T	S	T	S	T
	19"	-	S	T	S	T	S	T	S	T	S	T	S	T
	20"	-	S	T	S	T	S	T	S	T	S	T	S	T
	21"	-	S	T	S	T	S	T	S	T	S	T	S	T
	22"	-	S	T	S	T	S	T	S	T	S	T	S	T
	23"	-	S	T	S	T	S	T	S	T	S	T	S	T
	24"	-	S	T	S	T	S	T	S	T	S	T	S	T
	25"	-	S	T	S	T	S	T	S	T	S	T	S	T
26"	-	S	T	S	T	S	T	S	T	S	T	S	T	
27"	-	S	T	S	T	S	T	S	T	S	T	S	T	
28"	-	S	T	S	T	S	T	S	T	S	T	S	T	

Sport Back TRU-Comfort Back

TRU Balance Flex Seat Back Availability Chart

1f. Omit Back Selection

Omit Back..... No Charge
Part: PTOINDV1287

11.10. Back Width Package

16"-20" Back Width Package
Part: XRF000038

18"-22" Back Width Package
Part: XRF000039

11.11. Overall Back Height

REQUIRED with back selection. Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges.

14" Height Part: SETHGHT1027	17" Height Part: SETHGHT1030	20" Height Part: SETHGHT1033	23" Height Part: SETHGHT1036	26" Height Part: SETHGHT1039
15" Height Part: SETHGHT1028	18" Height Part: SETHGHT1031	21" Height Part: SETHGHT1034	24" Height Part: SETHGHT1037	27" Height Part: SETHGHT1040
16" Height Part: SETHGHT1029	19" Height Part: SETHGHT1032	22" Height Part: SETHGHT1035	25" Height Part: SETHGHT1038	28" Height Part: SETHGHT1041

11.12. Armrests

REQUIRED. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. TRU Balance Flex Seat, Drop-In Armrests

Includes full length, flat armpads. Available with 14" deep seats and greater.

Full Length, Drop-In Armrests..... No Charge
 Left part: ARM1803253. Right part: ARM156742

Desk Length, Drop-In Armrests No Charge
 Left part: ARM1803254. Right part: ARM156743

1b. Pediatric, Height Adjustable, Removable Armrests

Available with straight armpads only.

Desk Length, Pediatric Armrests \$79
 Left part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973

1c. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

Full Length, Flip Back Armrests..... \$63
 Left part: XRF000371. Right part: XRF000372

Desk Length, Flip Back Armrests \$63
 Left part: XRF000373. Right part: XRF000374

TRU Balance Flex seating section continued on following pages.

11. TRU BALANCE FLEX SEATING SYSTEM CONTINUED

1d. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

- Full Length, Two-Post Armrests.....\$158 ea. Desk Length, Two-Post Armrests\$158 ea.
- Left part: XRF000051. Right part: XRF000052. HCPCS: E0973* *Left part: XRF000053. Right part: XRF000054. HCPCS: E0973*

1e. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

- Full Length, Single Post Armrests \$79 ea. Desk Length, Single Post Armrests \$79 ea.
- Left part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973* *Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973*

11.13. Armpads

REQUIRED with armrest selection except Drop-In Armrests. Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

1a. Straight Armpads

- Straight, Full Length Left and Right No Charge Straight, Desk Length Left, Full Length Right No Charge
- Part: ARM134275* *Part: ARM134280*
- Straight, Full Length Left, Desk Length Right No Charge Straight, Desk Length Left and Right..... No Charge
- Part: ARM134281* *Part: ARM134254*

1b. Waterfall Armpads

- Waterfall, Full Length Left and Right..... No Charge Waterfall, Desk Length Left, Full Length Right No Charge
- Part: ARM134285* *Part: ARM134288*
- Waterfall, Full Length Left, Desk Length Right No Charge Waterfall, Desk Length Left and Right No Charge
- Part: ARM134289* *Part: ARM134282*

1c. Falcon Gel Armpads

- Falcon Gel, Full Length Left and Right\$311.22 Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
- Part: ARM136438* *Part: ARM136439*
- Falcon Gel, Full Length Left, Desk Length Right.....\$287.93 Falcon Gel, Desk Length Left and Right \$265.02
- Part: ARM136442* *Part: ARM136437*

12. HEADRESTS

Optional. Select a headrest pad and bracket.

12.1. Headrests, Comfort Plus with Cool Core

1. Headrest Pads

- Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179 Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
- Part: ST-CP08CCFPM-Q. HCPCS: E0955* *Part: ST-CP14CCFPM-Q. HCPCS: E0955*
- Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179 Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
- Part: ST-CP10CCFPM-Q. HCPCS: E0955* *Part: ST-CP06CCFPM-Q. HCPCS: E0955*

2. Headrest Brackets

- Comfort Plus Headrest Mounting Bracket, UniLink.....\$311 Multi-Axis Removable Headrest Bracket.....\$174
- Part: ST-UHB2D33200-Q. HCPCS: E1033* *Part: ST-HMO475-17M-Q. HCPCS: E1033*

12.2. Q Cam Lock

Optional.

- Q Cam Lock.....\$50
- Part: XRF000376*

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

13. LAP BELTS

REQUIRED. Select one lap belt below.

13.1a. Standard Lap Belts

- 38" Lap Belt *Part: ACCBELT1006*
- 50" Lap Belt *Part: ACCBELT1000*
- 60" Lap Belt *Part: ACCBELT1002*
- 70" Lap Belt *Part: ACCBELT1003*
- 80" Lap Belt *Part: ACCBELT1004*
- 90" Lap Belt *Part: ACCBELT1001*

13.1b. Optional Lap Belts

- Stealth Retractable Lap Belt (IFRB5350).....\$181 Stealth Positioning, 60" Single Pull, Push Button,\$50
- Part: ST-IFRB5350* *Extra Large, Padded (11.5"), 2" Belt (IFPB5315)*
- Stealth Positioning, 60" Single Pull, Push Button, Medium.....\$50 *Part: ST-IFPB5315-Q*
- (7.5") Padded, 1.5" Belt (IFPB5213)*
- Part: ST-IFPB5213-Q*
- Stealth Positioning, 60" Single Pull, Push Button, Large\$50 60" Sure-Lok Lap Belt (Non-retractable)..... \$120
- Padded (9.5"), 2" Belt (IFPB5314)* *Part: ACCBELT1016*
- Part: ST-IFPB5314-Q* *Required with ISO 7176-19 Compliant Transit Option.*

14. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

14.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

1. Center Mount Foot Platform Type

Basic Center Mount..... **No Charge**
 Part: XRF000323

2. Basic Center Mount Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates \$180 Long Extension, 5"x9" Footplates \$180
 Part: XRF000312. HCPCS: K0040 Part: XRF000314. HCPCS: K0040

Medium Extension, 5"x9" Footplates \$180
 Part: XRF000313. HCPCS: K0040

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount..... \$158
 Part: RIG2104369

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16
 Part: RIG2104240

14.1b. Center Mount Foot Platform

Select legrest length in following section.

1. Foot Platform Assembly

Quick Adjustable Center Mount Foot Platform \$263 Power Articulating Foot Platform \$3,145
 Part: XRF000311. HCPCS: K0040 Part: PTOINDV3754. HCPCS: E1012+K0040

2. Extension and Angle Adjustable Footplates

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates Part: RIG152577	Long Extension, 5"x6" Footplates Part: RIG151505
Pediatric Extension, 5"x8" Footplates Part: RIG2404015	Long Extension, 5"x8" Footplates Part: RIG2404017
Medium Extension, 5"x8" Footplates Part: RIG2404016	Long Extension, 5"x11" Footplates Part: RIG2404019
Medium Extension, 5"x11" Footplates Part: RIG2404018	

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard Part: RIGASMB7120018	Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard Part: RIGASMB7120019
Small Calf Pads, Quick Adjustable Center Mount Only, \$158 6"Wx6"H (pair) Part: XRF000156	Large Calf Pads, Quick Adjustable Center Mount Only, \$158 8" taper to 6"Wx8"H (pair) Part: XRF000157

3b. Center Mount Calf Panel

Small Calf Panel, 18"Wx8"H..... \$184 Part: FRMASMB7110012	Large Calf Panel, 24"Wx8"H \$237 Part: RIG123825 Only available on widths of 20" or greater.
---	--

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69 Part: RIG152519 Only available with the 5x6 footplates.	Outside Edge \$69 Part: RIG2404921 Available with 5"x8" and 5"x11" footplates.
Width Extension for 8" \$69 Part: RIG2404919	Front and Rear Cap \$16 Part: RIG2404923 Available with 5"x8" and 5"x11" footplates.
Width Extension for 11" \$69 Part: RIG2404920	

5. Center Mount Additional Accessories

Small Foam Calf Pad Cover (pair)..... \$67 Part: ST-CSCS0775-Q	Footplate Riser \$75 Part: XRF000348 Select legrest length between 4" to 8" for Power Articulating and 4"-7" with Quick Adjustable, at 1" increments.
Small Foam Panel Cover (pair) \$176 Part: ST-CSCP0920-Q	Footplate Riser Calf Panel \$50 Part: FRMASMB7110011 Only available and required with Footplate Riser selection and extension lengths of 6" to 9".
Footplate Connector \$16 Part: RIGASMB7120023	

Legrests section continued on following pages.

14. LEGRESTS CONTINUED

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers..... \$222 Part: ST-FPG5008T-Q. HCPCS: K0108	5"x8" Tapered Gel Footplate Covers for Width Extension \$247 Part: ST-AFPG5008E-Q. HCPCS: K0108
5"x11" Tapered Gel Footplate Covers..... \$270 Part: ST-AFPG5011T-Q. HCPCS: K0108	5"x11" Tapered Gel Footplate Cover for Width Extension \$270 Part: ST-AFPG5011E-Q. HCPCS: K0108

14.1c. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 60°-90° Swing-Away Legrest, Static Seat \$86 ea. Pair part: XRF000082. Left: XRF000083. Right: XRF000084 Not available with 11"-13" extensions.	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea. Pair part: XRF000085. Left: XRF000086. Right: XRF000087 For use only with 11"-13" lower extensions.
--	--

2. Lower Extension

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 11"-13" Length Pair part: XRF000088. Left: XRF000089. Right: XRF000090.	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 13"-16.5" Length Pair part: XRF000094. Left: XRF000095. Right: XRF000096.	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 15.5"-19.5" Length Pair part: XRF000097. Left: XRF000098. Right: XRF000099.
---	---	---

3a. Extruded Footplate(s)

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 7.25"Wx6"D Extruded Footplate..... Standard Part pair: FRMASMB7847. Left: FRMASMB7861 Right: FRMASMB7862	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8.5"Wx6"D Extruded Footplate Standard Part pair: FRMASMB7848. Left: FRMASMB7863 Right: FRMASMB7864 May cause interference on seat widths less than 18."
---	---

3b. Swing-Away Legrest, Angle Adjustable Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 4"Wx6"D \$90 ea. Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560 HCPCS: K0040+E0951	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx11"D ¹ \$90 ea. Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564 HCPCS: K0040+E0951
<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx6"D..... \$90 ea. Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561 HCPCS: K0040+E0951	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8"Wx11"D ² \$90 ea. Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565 HCPCS: K0040+E0951

4. Swing-Away Legrest Accessories

Gel Neoprene Footrest Sleeve (pair) \$84 Part: INDPART2243. HCPCS: K0108

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea. <small>QTY</small> Part: INDPART2822. HCPCS: E1020	Therafin 32387 Amputee Support Gel Cover \$155 ea. <small>QTY</small> Part: INDPART2823. HCPCS: K0108
---	--

14.1d. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> Elevating Legrest Assembly \$171 ea. Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990 Select lower extension and footplates

2a. Lower Extension

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 11"-14" Length Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030 For reinforced only	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 14"-17" Length Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035 For reinforced only
<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103 For angle adjustable only	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106 For angle adjustable only

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad..... \$342 Part: FRMASMB5301. HCPCS: E1020	Right Residual Limb Support with Calf Pad \$342 Part: FRMASMB5300. HCPCS: E1020
--	--

3a. Reinforced Cast Footplate(s)

<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 7.25"Wx6"D Reinforced Footplate..... Standard Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8.5"Wx6"D Reinforced Footplate Standard Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062 May cause interference on seat widths less than 18."
---	---

Legrests section continued on following pages.

14. LEGRESTS CONTINUED

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

<table border="0"> <tr> <td style="text-align: center;">P</td> <td style="text-align: center;">L</td> <td style="text-align: center;">R</td> <td>4"Wx6"D \$90 ea. Pair parts: FRMASMB8812 x2. HCPCS: K0040</td> </tr> <tr> <td style="text-align: center;">P</td> <td style="text-align: center;">L</td> <td style="text-align: center;">R</td> <td>5"Wx6"D..... \$90 ea. Pair parts: FRMASMB8811 x2. HCPCS: K0040</td> </tr> <tr> <td style="text-align: center;">P</td> <td style="text-align: center;">L</td> <td style="text-align: center;">R</td> <td>5"Wx8"D with Heel Loops..... \$90 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951</td> </tr> </table>	P	L	R	4"Wx6"D \$90 ea. Pair parts: FRMASMB8812 x2. HCPCS: K0040	P	L	R	5"Wx6"D..... \$90 ea. Pair parts: FRMASMB8811 x2. HCPCS: K0040	P	L	R	5"Wx8"D with Heel Loops..... \$90 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951	<table border="0"> <tr> <td style="text-align: center;">P</td> <td style="text-align: center;">L</td> <td style="text-align: center;">R</td> <td>6"Wx8"D¹ with Heel Loops \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951</td> </tr> <tr> <td style="text-align: center;">P</td> <td style="text-align: center;">L</td> <td style="text-align: center;">R</td> <td>6"Wx11"D¹ with Heel Loops..... \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951</td> </tr> <tr> <td style="text-align: center;">P</td> <td style="text-align: center;">L</td> <td style="text-align: center;">R</td> <td>8"Wx11"D² with Heel Loops..... \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951</td> </tr> </table>	P	L	R	6"Wx8"D ¹ with Heel Loops \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951	P	L	R	6"Wx11"D ¹ with Heel Loops..... \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951	P	L	R	8"Wx11"D ² with Heel Loops..... \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951
P	L	R	4"Wx6"D \$90 ea. Pair parts: FRMASMB8812 x2. HCPCS: K0040																						
P	L	R	5"Wx6"D..... \$90 ea. Pair parts: FRMASMB8811 x2. HCPCS: K0040																						
P	L	R	5"Wx8"D with Heel Loops..... \$90 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951																						
P	L	R	6"Wx8"D ¹ with Heel Loops \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951																						
P	L	R	6"Wx11"D ¹ with Heel Loops..... \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951																						
P	L	R	8"Wx11"D ² with Heel Loops..... \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951																						

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
Part: INDPART2242. HCPCS: E0953

14.2. Legrest Accessories

Heel Loops (pair)..... \$53 Part: FRMASMB7873. HCPCS: E0951 Standard on angle adjustable and multi axis footplates when noted.	Wheel Bumper..... \$19 ea. Part: XRF000100
--	---

14.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	4" L (#3991)	6" L (#3993)	8" L (#3680)	9" L (#3682)	10" L (#3684)	11" L (#3686)	12" L (#3688)	13" L (#3690)	14" L (#3692)	15" L (#3694)
3" L (#3990)	5" L (#3992)	7" L (#3994)	8.5" L (#3681)	9.5" L (#3683)	10.5" L (#3685)	11.5" L (#3687)	12.5" L (#3689)	13.5" L (#3691)	14.5" L (#3693)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)
7.5" L (#3471)	8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	

15. POWER OPTION ELECTRONICS

REQUIRED for articulating foot platform.

Single Actuator Through Controller Interface \$2,163 Part: XRF000335. HCPCS: E2310	Single Actuator Through Controller and Interface (AAM).....\$3,040 Part: XRF000068. HCPCS: E2310
---	---

16. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Flex Seat Cup Holder..... No Charge Part: ACCASMB2634 Not available on TRU Balance Flex seating without two-post armrests.	Grip Drink Holder (20085) \$15.50 Part: ACC152993
Oxygen Holder for Flex Seat \$210 Part: ACC125006. HCPCS: E2207	Accessory Bar\$132 Part: XRF000128
Cane/Crutch Holder \$105 Flex part: ACC152515. HCPCS: E2207	Large, Black, Weather Cover \$202 Part: ACCCOVR1021
Cell Phone Holder for 12" Lift Captain Seat..... \$70 Part: ACC1611300	TRU Balance 3 Transfer Handle.....\$105 ea. Part: MEC2205680. HCPCS: E0705 Not available with Heavy Duty armrests.
Flag Holder \$55 Part: ACCASMB2415	Privacy Bag.....\$42 Part: ACC2005723
Saddle Bag.....\$40 Part: DWR1020D010	Personal Item Pouch\$84 Part: ACC2101017
Stealth Cell Phone/Media Device Holder \$215 Part: ACC7120006	Remote Stop Switch \$1,295 Part: XRF000379
Quantum Backpack..... \$168 Part: ACC1905853	
Grip Lap Board (20005).....\$38.95 Part: ACC152992	

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea. Part: ACC2407046	Aftermarket Accessory Mounting Bracket, Non-slotted.. \$53 ea. Part: ACC2405496
--	--

17. BRANCH PATCHES

Optional.

Air Force Patch\$21 Part: ACC140240	Coast Guard Patch\$21 Part: ACC140239	Navy Patch\$21 Part: ACC140244
Army Patch\$21 Part: ACC140241	Marines Patch\$21 Part: ACC140243	Purple Heart Patch\$21 Part: ACC140238



18. POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

18.1. Stealth Thoracic Laterals

Swing-away laterals may interfere with 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, Rail Mounts, pair (TWBSTL)..... \$562 Part: INDPART2896. HCPCS: E1034 For use with Synergy Back.	Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664 Part: ST-ULPK404104B00-Q. HCPCS: E1034 Only available with back cane and omit back configurations.
--	---

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$250 Part: INDPART2663. HCPCS: E0956	Flat Pads, 3"x5", pair (ULP846050320)..... \$250 Part: ST-ULP846050320-Q. HCPCS: E0956
Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250 Part: INDPART2664. HCPCS: E0956	Flat Pads, 4"x6", pair (ULP835050320) \$250 Part: ST-ULP835050320-Q. HCPCS: E0956

3. Installation

Install Laterals Part: PTOINDV3499	Do Not Install Part: PTOINDV3500
---------------------------------------	-------------------------------------

4. Lateral Extensions

REQUIRED with Synergy Back and thoracic lateral selections.

Lateral Extension for Synergy Back (pair)..... \$28 ea. Part: XRF000150
--

18.2. Stealth UniLink Pelvic/Thigh Guide

1. UniLink Pelvic/Thigh Hardware

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea. Removable Mounting Hardware Part: ST-UTTP5-Q. HCPCS: E1028/E1034^ QTY	Flip Down with LEFT Full Surface Contact \$281 ea. Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028/E1034^
UniLink, 7" Height Adjustable, \$276 ea. Removable Mounting Hardware Part: ST-UTTP7-Q. HCPCS: E1028/E1034^ QTY	Flip Down with RIGHT Full Surface Contact..... \$281 ea. Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028/E1034^
UniLink Multi-Angle Thigh Guide, Left \$358.87 ea. Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^	UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea. Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea. Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^ QTY	UniLink, 4" x 12" Pad \$106 ea. Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^ QTY
UniLink, 4" x 8" Pad..... \$106 ea. Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^ QTY	

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea. Part: ST-UIL001150A000-Q QTY	UniLink, 2" Link..... \$81 ea. Part: ST-UIL001200A000-Q QTY
--	---

4. UniLink Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

18.3. Stealth Position Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea. Part: ST-IFPS1103-Q. HCPCS: E0960	Upper Contour, Medium (IFPS2103)..... \$154 ea. Part: ST-IFPS2103-Q. HCPCS: E0960
Upper Classic, Large (IFPS1104) \$154 ea. Part: ST-IFPS1104-Q. HCPCS: E0960	Upper Contour, Large (IFPS2104) \$154 ea. Part: ST-IFPS2104-Q. HCPCS: E0960



18.4. Stealth Position Chest Straps

Chest Strap, Medium\$97 ea. w/ Belt Anchor, (IFCS4113) Part: ST-IFCS4113-Q. HCPCS: E0960	Chest Strap, Large\$97 ea. w/ Belt Anchor, (IFCS4114) Part: ST-IFCS4114-Q. HCPCS: E0960
--	---



Positioning components section continued on following page.

18. POSITIONING COMPONENTS CONTINUED

18.5. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right..... \$318 ea.
Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left \$318 ea.
Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$434 ea.
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028



2. Facial Lateral Pads and Accessories


Egg Switch Mount (ESM) (A) \$57 ea.
Part: INDPART2815

Facial Pad, 4".5x2", Cool Core (900M-Q) \$61 ea.
Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2", Cool Core (920M-Q)..... \$61 ea.
Part: ST-920M-Q. HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C) \$276 ea.
Part: ST-SUS18-Q

Allows mounting of facial pads and attaches to the headrest



18.6. Stealth Elbow Support

1. Elbow Support Bracket

Left, Swing-Away, Elbow Support Bracket, Tube \$199 ea.
Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket, Tube \$199 ea.
Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$57 ea.
Part: POS124067. HCPCS: K0108

18.7. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

18.8. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

19. LAP TRAY

19.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
 Part: ST-SSLXS010001171-Q, HCPCS: E0950

Simple Solutions Large Tray with Cut Out..... \$309
 Part: ST-SSLXL010001171-Q, HCPCS: E0950

Simple Solutions Medium Tray with Cut Out..... \$309
 Part: ST-SSLXM010001171-Q, HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware..... \$302
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	S	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q
 :Compatible  :Armrest configuration may not hold tray position
 *Some armrests listed may not be available on this order form.