



Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® HD Order Form with Contoured and Synergy Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering. Options noted with "XRef" have multiple possible part numbers based on system configuration. Please see [cross-reference spreadsheet](#) for the XRef part numbers.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 450 lb. weight capacity with standard 4.7 mph programmable motors, 8 amp off board charger, and unoccupied transit loops. Dimensions of base without legrests: 26"Wx36.75"L.

Quantum Q6 Edge 3HD-SS.....\$8,595

Part: Q6 EDGEHD 3HD-SS. HCPCS: K0850

Solid seat pan standard with no power option.

Quantum Q6 Edge 3HD-C.....\$8,595

Part: Q6 EDGEHD 3HD-C. HCPCS: K0851

Contoured Seat (Captain) Standard.

Account Number: _____ Date: _____ Marked for: _____

2. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

2.1. Match Colors

Match All Color Options to Selection Below

Part: N/A. Any subsequent color selections, if available, will be matched to color chosen below when the order is placed. Please note that some color selections, such as joystick armor, may have a MSRP.

2.2. Shroud Color

 Red Corvette Part: DGN169178-PRD1	 Night Train Part: DGN169178-PBK3	 Deep Purple Part: DGN169178-PPR5	 Green Machine Part: DGN169178-PGN11
 Electric Blue Part: DGN169178-PBU9	 Pretty in Pink Part: DGN169178-PPK3	 Grey Street Part: DGN169178-PSL3	To Be Determined Part: N/A
 Back in Black Part: DGN169178-PBK1	 Evolution Orange Part: DGN169178-POR7	 Midnight Blue Part: DGN169178-PBU1	

3. FENDER LIGHTS

Optional. Fender lights operate by pressing included single switch button. Fender lighting is not available with NE 4-Key Joystick or Single Actuator Control through Toggle Box.

Switch Fender Lights \$250
Part: DGN154117

4. FULL LIGHT PACKAGE

Optional. Not available with contoured seating. Requires NE+ 8 Key or Q-Logic 3 Joystick for Full Light Package.

Seat Mounted, LED, Full Light Package \$995
Part: ELE131477

5. BASE MOUNTED FOOT PLATFORM

REQUIRED. Select a foot platform option. May interfere with seat mounted legrests.

Base Mounted Foot Platform Standard Part: FRMASMB16077	Omit base mounted foot platform No Charge Part: PTOINDV3674
Base Mounted, Bell Shaped Foot Platform No Charge Part: RIGASMB7120041	

6. BATTERIES

REQUIRED. The powerchair requires 2 batteries and will be installed if ordered.

Group 24, 70 Amp Gel Cell, Interceptor Battery \$365 ea. Part: BAT1704614. HCPCS: E2363	Group 24, 70 Amp Gel Cell, MK Battery \$495 ea. Part: BATGEL1009. HCPCS: E2363
--	---

No batteries. Provider will assume responsibility of batteries.
Part: Configuration selection

7. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

7.1. Drive Wheels

3" x 14" Run Flat Tires Standard Part: PTOINDV3510	4" x 14" Flat Free, Knobby, Extended \$275.50 ea. Traction Tires Part: WHL121698 x2
3" x 14" Flat Free Tires No Charge Part: WHL124116 x2	4" x 14" Pneumatic, Knobby, Extended \$275.50 ea. Traction Tires Part: WHL121703 x2
3" x 14" Pneumatic Tires No Charge Part: WHL124117 x2	

7.2. Drive Wheel Accent Color

Optional. Not available with knobby tires. Sold in Pairs (parts below will be quantity of 2).

Match to Base Color Part: Configuration selection, see specific part below	Non-matching Color Part: Configuration selection, see specific part below
Back in Black \$24.98 ea. Part: DGN164245-PBK1 x2	Pretty in Pink \$24.98 ea. Part: DGN164245-PPK3 x2
Red Corvette \$24.98 ea. Part: DGN164245-PRD1 x2	Evolution Orange \$24.98 ea. Part: DGN164245-POR7 x2
Electric Blue \$24.98 ea. Part: DGN164245-PBU9 x2	Deep Purple \$24.98 ea. Part: DGN164245-PPR5 x2
Night Train \$24.98 ea. Part: DGN164245-PBK3 x2	Grey Street \$24.98 ea. Part: DGN164245-PSL3 x2
	Midnight Blue \$24.98 ea. Part: DGN164245-PBU10 x2
	Green Machine \$24.98 ea. Part: DGN164245-PGN11 x2
	To Be Determined Follow-up required.

8. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Forms](#).

8.1a. Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick.....Standard Part: CTLDC1554 No actuator control.	NE+ 8-Key Joystick (2 Actuators).....\$400 Part: CTLDC1561 Control up to 2 actuators and LED, Full Light Package through joystick.
NE+ 6-Key Joystick (2 Actuators)\$200 Part: CTLDC1560 Control up to 2 actuators through joystick.	No Joystick, Setup for NE+ Joystick No Charge Part: Configuration selection

**8.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module\$695 Part: See XRef. HCPCS: E2377	Q-Logic 3 Expandable Controller\$425 Harness Part: HAR144005. HCPCS: E2313
---	--

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick.....\$700 Part: CTL167724	No Joystick, Setup for Q-Logic 3 Part: Configuration selection.
Q-Logic 3 EX Joystick for LED, Full Light Package...\$700 Part: CTL167725	

**8.2. Q-Logic 3 Joystick Armor**

REQUIRED with Q-Logic 3 EX Joystick selections above. Reference colors shown on page 2.

Black, Soft TouchStandard Part: DGN124307	Non-matching Shroud\$50 Part: PTOINDV3422 Select color below	Match Shroud\$50 Part: PTOINDV3421
--	--	---

Red Corvette Part: DGN124336-PRD1	Back in Black Part: DGN124336-PBK1	Pretty in Pink Part: DGN124336-PPK3	Deep Purple Part: DGN124336-PPR5	Midnight Blue Part: DGN124336-PBU10
Electric Blue Part: DGN124336-PBU9	Night Train Part: DGN124336-PBK3	Evolution Orange Part: DGN124336-POR7	Grey Street Part: DGN124336-PSL3	Green Machine Part: DGN124336-PGN11

9. JOYSTICK HANDLES

Optional. The NE/NE Plus joysticks are standard with a carrot handle (See A in image). The standard Q-Logic 3 joystick knob can be seen in the image in section 8.1b. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B)\$85 Part: ACCASMB1860. HCPCS: E2323	Extended I-Handle (G)\$85 Part: ACCASMB1857. HCPCS: E2323
Chin Cup (C)\$85 Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)	Bodypoint 3" Goal Post Handle\$99 Part: PLSKNCP1057. HCPCS: E2323
T-Bar (D)\$120 Part: ACCASMB1858. HCPCS: E2323	Bodypoint 4" Goal Post Handle\$99 Part: PLSKNCP1058. HCPCS: E2323
Medium Ball Handle (E)\$65 Part: ACCASMB1859. HCPCS: E2323	Q-Logic 3 Carrot Handle (A) No Charge Part: DGN124308
Soft Ball Handle (F)\$65 Part: ACCASMB1856. HCPCS: E2323	


**10. JOYSTICK MOUNTING BRACKETS**

Optional. Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket


10.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, BracketStandard Q-Logic Part: FRMASMB10617. NE/NE+ Part: FRMASMB12744	
Left Side Mount Part: PTOINDV1225	Right Side Mount Part: PTOINDV1226

**10.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" is used with track/TRU-Balance 3 style arms when the joystick is inline with armrest. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armrest on tube style armrests.

Swing-away, Inline, Left Side Mount.....\$245 Part: See XRef. HCPCS: E1028	Swing-away, 1" Offset, Left Side Mount\$245 Part: See XRef. HCPCS: E1028
Swing-away, Inline, Right Side Mount\$245 Part: See XRef. HCPCS: E1028	Swing-away, 1" Offset, Right Side Mount.....\$245 Part: See XRef. HCPCS: E1028



Joystick mounting brackets section continued on following page.

10. JOYSTICK MOUNTING BRACKETS CONTINUED**10.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$350

NE/NE+ Part: MEC138876. HCPCS: E1028

Q-Logic Part: MEC138860. HCPCS: E1028

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-axis Mount \$100

NE/NE+ Part: ACC142528. HCPCS: E1028

Q-Logic Part: ACC131216. HCPCS: E1028


Allows height adjustment of joystick.

Stealth Swing-Away/Flip Down Bracket \$294.25

SAJ Part: ST-TWBM480CJ-Q. HCPCS: E1028

Q-Logic Part: ST-TWBM480QL-Q. HCPCS: E1028

Not available with contoured seat. Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.



Retract4



Multi-Axis Joystick Mount

Specialty Mount Side

Left Side Mount

Part: PTOINDV1225

Right Side Mount

Part: PTOINDV1226

11. CONTOURED SEATING**Optional.** Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.**11.1a. Captain Style Seating****1a. Depth Adjustable High Back with Limited Recline**

Full length armrests must be selected with a depth adjustable high back.

20"Wx18"-20"D with 60" Lap Belt No Charge

Part: SETASMB6896

22"Wx20"-22"D with 70" Lap Belt No Charge

Part: SETASMB6897

1b. High Back with Limited Recline

20"Wx16"D No Charge

Part: SETASMB2780

20"Wx22"D No Charge

Part: SETASMB4889

1c. Pinchless Hinge High Back

20"Wx24"D No Charge

Part: SETASMB2380

22"Wx22"D No Charge

Part: SETASMB2383

24"Wx20"D No Charge

Part: SETASMB2387

22"Wx16"D No Charge

Part: SETASMB2381

22"Wx24"D No Charge

Part: SETASMB2384

24"Wx22"D No Charge

Part: SETASMB2388

22"Wx18"D No Charge

Part: SETASMB2171

24"Wx16"D No Charge

Part: SETASMB2385

24"Wx24"D No Charge

Part: SETASMB2389

22"Wx20"D No Charge

Part: SETASMB2382

24"Wx18"D No Charge

Part: SETASMB2386

1d. Heavy Duty High Back

Seat includes a lumber cutout, static hinges, and a Synergy Seat frame.

26"Wx18"D \$450

Part: SETASMB4565

26"Wx26"D \$450

Part: SETASMB4569

28"Wx22"D \$450

Part: SETASMB4573

26"Wx20"D \$450

Part: SETASMB4566

26"Wx28"D \$450

Part: SETASMB4570

28"Wx24"D \$450

Part: SETASMB4574

26"Wx22"D \$450

Part: SETASMB4567

28"Wx18"D \$450

Part: SETASMB4571

28"Wx26"D \$450

Part: SETASMB4575

26"Wx24"D \$450

Part: SETASMB4568

28"Wx20"D \$450

Part: SETASMB4572

28"Wx28"D \$450

Part: SETASMB4576

11.1b. Contoured Seat with Solid Seat Pan**1a. High Back with Limited Recline and Solid Seat Pan**

20"Wx16"D No Charge

Part: SETASMB3183

20"Wx20"D No Charge

Part: SETASMB3185

20"Wx18"D No Charge

Part: SETASMB3184

20"Wx22"D No Charge

Part: SETASMB5008

1b. Pinchless Hinge High Back with Solid Seat Pan

20"Wx24"D No Charge

Part: SETASMB2982

22"Wx22"D No Charge

Part: SETASMB2986

24"Wx20"D No Charge

Part: SETASMB2990

22"Wx16"D No Charge

Part: SETASMB2983

22"Wx24"D No Charge

Part: SETASMB2987

24"Wx22"D No Charge

Part: SETASMB2991

22"Wx18"D No Charge

Part: SETASMB2984

24"Wx16"D No Charge

Part: SETASMB2988

24"Wx24"D No Charge

Part: SETASMB2992

22"Wx20"D No Charge

Part: SETASMB2985

24"Wx18"D No Charge

Part: SETASMB2989

Contoured seating section continued on following page.

11. CONTOURED SEATING CONTINUED**1c. Heavy Duty High Back with Solid Seat Pan**

Seat includes a lumber cutout, static hinges, and a Synergy Seat frame.

26"Wx18"D..... \$1000 Part: SETASMB5478	26"Wx26"D..... \$1000 Part: SETASMB5482	28"Wx22"D \$1000 Part: SETASMB8486
26"Wx20"D..... \$1000 Part: SETASMB5479	26"Wx28"D..... \$1000 Part: SETASMB5483	28"Wx24"D \$1000 Part: SETASMB5487
26"Wx22"D..... \$1000 Part: SETASMB5480	28"Wx18"D..... \$1000 Part: SETASMB8484	28"Wx26"D..... \$1000 Part: SETASMB5488
26"Wx24"D..... \$1000 Part: SETASMB5481	28"Wx20"D \$1000 Part: SETASMB8485	28"Wx28"D \$1000 Part: SETASMB5489

11.1c. Special Order Contoured Seat

Complete information below. Contact inside sales specialist for pricing.

Type of Seat: _____ Width: _____ Depth: _____

11.2. Seat to Floor Height

Lowest Seat to Floor Part: PTOINDV3558	Middle Seat to Floor Part: PTOINDV3559	High Seat to Floor Part: PTOINDV3560
---	---	---

11.3. Contoured Seat Flip-Up, Height Adjustable, Removable Armrests**REQUIRED and only available with non-heavy duty, High-Back seat. Depth Adjustable High-Back must select full length armrests.** Select a Left and Right Armrest or to Omit them.

Full Length, Left Armrest \$85 Part: FRMASMB8614. HCPCS: E0973	Full Length, Right Armrest..... \$85 Part: FRMASMB8615. HCPCS: E0973
Desk Length, Left Armrest..... \$85 Part: FRMASMB8612. HCPCS: E0973	Desk Length, Right Armrest \$85 Part: FRMASMB8613. HCPCS: E0973
No Armrest Standard Part: N/A	

11.4. Heavy Duty, Contoured Seat, Height Adjustable, Removable Armrests**REQUIRED and only available with heavy duty High-Back seats.** Select armrests and armpads. Heavy duty armrests are sold in pairs.**1. Bariatric/Heavy Duty Armrests**

Bariatric, HD, Height Adjustable, Full Length, Standard Removable 8"H-12"H Left part: FRMASMB6728. Right part: FRMASMB6729	Heavy Duty, Height Adjustable, Full Length, \$85 ea. Removable 8"H-12"H Left part: FRMASMB9982. Right part: FRMASMB9983. HCPCS: E0973
Bariatric, HD, Height Adjustable, Full Length, Standard Removable 10"H-14"H Left part: FRMASMB6730. Right part: FRMASMB6731	Heavy Duty, Height Adjustable, Full Length, \$85 ea. Removable 10"H-14"H Left part: FRMASMB9986. Right part: FRMASMB9987. HCPCS: E0973

2. Armpads for Bariatric/Heavy Duty Armrests

Required with bariatric or heavy duty armrests

Bariatric Armpads No Charge Part: UPSASMB1027 Only for use with Bariatric Armrests.	Straight Armpads..... No Charge Part: UPSASMB1021 Only for use with Heavy Duty Armrests
No Armpad No Charge Part: N/A	Waterfall Armpads..... No Charge Part: UPSASMB1024 Only for use with Heavy Duty Armrests

12. SYNERGY SEATING SYSTEM**Optional.** Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 28"W x 32"D.**12.1. Synergy Frame Selection**

Synergy Seating Frame \$434 Part: PTOINDV2995
--

12.2. Occupied Transit, WC-19**Optional.** Synergy Seating only with widths 20"-24". To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

WC-19 Occupied Transit Option..... \$275 Part: KITASMB1830

12.3. Seat Dimension - Width**REQUIRED with Synergy seating.** Part number is SETWDTH plus number following # below.

20" Wide (#1010)	21" Wide (#1011)	22" Wide (#1012)	23" Wide (#1013)	24" Wide (#1014)	25" Wide (#1056)	26" Wide (#1054)	27" Wide (#1057)	28" Wide (#1055)
---------------------	---------------------	---------------------	---------------------	---------------------	---------------------	---------------------	---------------------	---------------------

12.4. Seat Dimension - Width Package**REQUIRED with Synergy seating.**

20"-24" Width Package Part: MEC120262	24"-26" Width Package Part: MEC120263	26"-28" Width Package Part: MEC120264
--	--	--

Synergy seating section continued on following page.

12. SYNERGY SEATING SYSTEM CONTINUED**12.5. Seat Dimension - Depth****REQUIRED with Synergy seating.** Part number is SETDPH plus number following # below.

15" Deep (#1005)	17" Deep (#1007)	19" Deep (#1009)	21" Deep (#1011)	23" Deep (#1013)	25" Deep (#1060)	27" Deep (#1062)
16" Deep (#1006)	18" Deep (#1008)	20" Deep (#1010)	22" Deep (#1012)	24" Deep (#1014)	26" Deep (#1061)	28" Deep (#1063)

12.6. Seat Dimension - Depth Package**REQUIRED with Synergy seating.**

15"-20" Depth Package Part: FRMASMB7732	20"-24" Depth Package Part: FRMASMB7733	24"-28" Depth Package Part: FRMASMB6605
--	--	--

12.7. Synergy Seat to Floor**REQUIRED with Synergy seating.**

16.5" Part: PTOINDV1888	17.5" Part: PTOINDV1889	18.5" Part: PTOINDV1890
----------------------------	----------------------------	----------------------------

12.8. Synergy Back Options**REQUIRED with Synergy seating.** See the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Back widths larger or smaller than seat widths on the availability chart are allowed with endomorph or mesomorph back canes.**1. Back Cane Selection**

Synergy Back Cane No Charge Part: FRMASMB15844. Sling/Transit Part: FRMASMB15845 Select a back cane height. Back width will be the same as seat width.	2" Endomorph (1" per side) Back Cane \$195 Part: FRMASMB16393. Transit part: FRMASMB12115 Select a cane height. Back width will be 2" smaller than seat width.
2" Mesomorph (1" per side) Back Cane \$195 Part: FRMASMB16387. Transit part: FRMASMB12067 Select a cane height. Back width will be 2" larger than seat width.	4" Endomorph (2" per side) Back Cane \$195 Part: FRMASMB16396. Transit part: FRMASMB12139 Select a cane height. Back width will be 4" smaller than seat width.
4" Mesomorph (2" per side) Back Cane \$195 Part: FRMASMB16390. Transit part: FRMASMB12091 Select a cane height. Back width will be 4" larger than seat width.	

2. Synergy Back Cane Height

Omit this selection if a TRU-Comfort back is selected. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles. Parts: See XRef, configuration selection with push handles selection.

14" High Canes Part: FRMASMB15832	17" High Canes Part: FRMASMB15835	20" High Canes Part: FRMASMB15838	23" High Canes Part: FRMASMB15841
15" High Canes Part: FRMASMB15833	18" High Canes Part: FRMASMB15836	21" High Canes Part: FRMASMB15839	24" High Canes Part: FRMASMB15842
16" High Canes Part: FRMASMB15834	19" High Canes Part: FRMASMB15837	22" High Canes Part: FRMASMB15840	25" High Canes Part: FRMASMB15843

3. Back Selection

Synergy Back \$275 Part: See XRef Select a back cane type and height and back height.	HD Sling Back \$98 Part: See XRef Select a back height. Available with 24"-28" Width and 21"-25" overall heights. Back cane height must match back height. Not available with flip up cantilever arms.
TRU-Comfort 4-way Stretch, 20" Width or less \$675 Part: See XRef. HCPCS: E2620 Select back height.	Omit Back No Charge Part: PTOINDV1287 Select a back cane type and height in following section.
TRU-Comfort 4-way Stretch, 22" Width or greater \$775 Part: See XRef. HCPCS: E2621 Select back height.	

		Back Width																
		16"	17"	18"	19"	20"	21"	22"	23"	24"	25"	26"	27"	28"	29"	30"	31"	32"
Overall Height	17"	S	S	S	S	S	-	S	-	S	-	-	-	-	-	-	-	-
	18"	S	S	S	S	S	-	S	-	S	-	-	-	-	-	-	-	-
	19"	S	-	S	S	S	S	S	S	S	-	S	S	S	S	-	S	-
	20"	STC	-	STC	STC	STC	S	S	S	S	-	S	S	S	S	-	S	-
	21"	STC	S	STC	STC	STC	S	S	S	S	S	-	S	S	-	S	-	-
	22"	STC	STC	STC	STC	STC	S	S	STC	S	S	-	S	S	-	S	-	-
	23"	STC	STC	STC	STC	STC	S	S	STC	S	S	S	S	S	-	S	S	S
	24"	STC	STC	STC	STC	STC	S	STC	STC	STC	S	S	S	S	-	S	S	S
	25"	STC	STC	STC	STC	STC	S	STC	STC	STC	S	S	S	S	-	S	S	S
	27"	TC	-	TC	TC	TC	-	TC	-	TC	-	-	-	-	-	-	-	-
	28"	TC	-	TC	TC	TC	-	TC	-	TC	-	-	-	-	-	-	-	-
S Synergy Back TC TRU-Comfort Back STC Synergy/TRU-Comfort Back																		

Synergy and TRU-Comfort Back Availability Chart

4. Synergy Overall Back Height

Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Parts: See XRef, configuration selection with back selection.

17" H	18" H	19" H	20" H	21" H	22" H	23" H	24" H	25" H	27" H	28" H
-------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------

Synergy seating section continued on following page.

12. SYNERGY SEATING SYSTEM CONTINUED**12.9. Lateral Extensions****REQUIRED with Synergy Back and thoracic lateral selections.**

Lateral Extension for Synergy Back (pair)..... \$26.25 ea.
 16"OH back part: FRMASMB16274. 17"OH back part: FRMASMB16280. 20"&21"OH back part: FRMASMB16298. 24"&25"OH back part: FRMASMB16281

12.10. Synergy Armrests**REQUIRED with Synergy seating.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. Bariatric, Single Post, Heavy Duty, Height Adjustable, Removable Armrests**

Left and right armrests must be the same armrest height. Only available with bariatric armpads.

L R Full Length, Bariatric Armrests, 8"-12" H No Charge Left part: FRMASMB6728. Right part: FRMASMB6729	L R Full Length, Bariatric Armrests, 10"-14"H No Charge Left part: FRMASMB6730. Right part: FRMASMB6731.
---	--

1b. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

L R Full Length, Heavy Duty Armrests, 8"-12" H \$85 ea. Left part: FRMASMB9982. Right part: FRMASMB9983. HCPCS: E0973	L R Desk Length, Heavy Duty Armrests, 8"-12"H \$85 ea. Left part: FRMASMB9984. Right part: FRMASMB9985. HCPCS: E0973
--	---

1c. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

L R Full Length, Heavy Duty Armrests, 10"-14"H \$85 ea. Left part: FRMASMB9986. Right part: FRMASMB9987. HCPCS: E0973	L R Desk Length, Heavy Duty Armrests, 10"-14"H \$85 ea. Left part: FRMASMB9988. Right part: FRMASMB9989. HCPCS: E0973
--	--

1d. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

L R Full Length, Cantilever Armrests \$97.50 ea. Left part: FRMASMB16188. Right part: FRMASMB16189	L R Desk Length, Cantilever Armrests \$97.50 ea. Left part: FRMASMB16190. Right part: FRMASMB16191
--	--

12.11. Synergy Armpads**REQUIRED with armrest selection above.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

Bariatric Armpads No Charge
 Part: UPSASMB1025
 Only available with bariatric armrests.

Straight Armpads No Charge
 Desk/Desk part: UPSASMB1022. Full/Full part: UPSASMB1021.
 Desk/Full part: UPSASMB1023
 Not available with bariatric armrests.

Waterfall Armpads No Charge
 Desk/Desk part: UPSASMB1025. Full/Full part: UPSASMB1024.
 Desk/Full part: UPSASMB1026
 May cause interference with seat back. Not available with bariatric armrests.

Omit Synergy Armpads No Charge
 Part: PTOINDV1373

12.12. Synergy Seat Joystick Receiver Bracket**REQUIRED.**

Standard Fixed Receiver Standard
 Part: FRMASMB16062

Stealth Height Adjustable Receiver \$214
 Tube Arm Part: ST-ARM250C-Q. HCPCS: K0108

13. HEADRESTS**Optional.** Select a headrest pad and bracket. **Not available with contoured seating or sling back.****13.1. Headrest Pads**

Stealth 8" Comfort Plus (CP450)..... \$151.94 Part: POS148484. HCPCS: E0955 Stealth 10" Comfort Plus (CP250) (A) \$151.94 Part: SETHEAD1057. HCPCS: E0955 Stealth 14" Comfort Plus (CP150). \$151.94 Part: INDPART2327. HCPCS: E0955	Stealth Small Contoured (9.75"Wx3"H) \$64.20 (STL-OTS250) Part: SETHEAD1052. HCPCS: E0955 Stealth Large Contoured (11"Wx6"H) (B) \$80.25 (STL-OTLG150) Part: SETHEAD1053. HCPCS: E0955
---	---


13.2. Headrest Brackets

Stealth's Worlds Best Removable.....\$220.50 Headrest Bracket (TWB485-1) (C) Part: FRMASMB15999. HCPCS: E1028	Stealth Multi-Axis Removable Headrest \$160.50 Bracket (HMO475P) (D) Part: HDWASMB1084. HCPCS: E1028
---	--

13.3. Flip-down Headrest Bracket

Stealth Flip-down Headrest Bracket (FDM350)..... \$224.70
 Part: INDPART2219

The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.



14. LAP BELTS**REQUIRED except with select contoured seats with included lap belt.** Select one lap belt below.**14.1a. Standard Lap Belts**

38" Lap Belt Part: ACCBELT1006	50" Lap Belt Part: ACCBELT1000	60" Lap Belt Part: ACCBELT1002	70" Lap Belt Part: ACCBELT1003	80" Lap Belt Part: ACCBELT1004	90" Lap Belt Part: ACCBELT1001
-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------

14.1b. Optional Lap Belts

Not available with contoured seating.

i-Fit 60" Single Pull, Push Button, Medium\$86 (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q	i-Fit 60" Single Pull, Push Button, Extra Large\$86 Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q
i-Fit 60" Single Pull, Push Button, Large\$86 Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q	60" Sure-Lok Lap Belt (Non-retractable)..... \$120 Required with WC-19 Transit Option. Part: ACCBELT1016

15. SYNERGY SEAT WITH POWER AFP OPTION ELECTRONICS**REQUIRED for Synergy Seat with articulating foot platform.**

One Actuator Control (Articulating Foot Platform)..... \$1,895 through NE+ Joystick Part: ELEASMB6430. HCPCS: E2310	One Actuator Control (Articulating Foot Platform)..... \$1,895 Through Q-Logic 3 Electronics No fender lights part: ELE1705592. With fender lights part: ELE1705562. HCPCS: E2310
---	--

16. LEGRESTS**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**16.1. Heavy Duty, Elevating Legrests**

Select a legrest with extension and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Heavy Duty, Elevating Legrest Selection

☐ P ☐ L ☐ R	Heavy Duty, Elevating Legrest(s) with Calf Pad(s)..... \$162.50 ea. Pair part: FRMASMB6702. Left: FRMASMB8444. Right: FRMASMB8445. HCPCS: E0990
--	--

2. Lower Extension

☐ P ☐ L ☐ R	11"-14" Length Part: FRMASMB4755. Left: RIG138026. Right: RIG138030 For reinforced footplates only.	☐ P ☐ L ☐ R	14"-17" Length Part: FRMASMB4756. Left: RIG138033. Right: RIG138035 For reinforced footplates only.
☐ P ☐ L ☐ R	12"-14.5" Length Part: FRMASMB4752. Left: FRMASMB5391. Right: FRMASMB4723 For angle adjustable footplates only.	☐ P ☐ L ☐ R	14.5"-18" Length Part: FRMASMB4753. Left: FRMASMB4724. Right: FRMASMB5392 For angle adjustable footplates only.

3. Reinforced Footplates(s)

☐ P ☐ L ☐ R	7.25"Wx6"D Extruded Footplate..... Standard Part pair: FRMASMB4748. Left: RIG138055. Right: RIG138056	☐ P ☐ L ☐ R	8.5"Wx6"D Extruded Footplate Standard Part pair: FRMASMB4747. Left: RIG138060. Right: RIG138062
--	---	--	---

4. Angle Adjustable Footplates(s) with Heel Loops

☐ P ☐ L ☐ R	6"Wx8"D \$85 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795. HCPCS: K0040+E0951	☐ P ☐ L ☐ R	8"Wx11"D \$85 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799. HCPCS: K0040+E0951 May cause interference with seat widths less than 20"
☐ P ☐ L ☐ R	6"Wx11"D \$85 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797. HCPCS: K0040+E0951		

16.2. Center Mount Foot Platform**Not available with contoured seating.** If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.**1. Foot Platform Assembly**

Power Articulating Foot Platform \$2,995 Part: PTOINDV3754. HCPCS: E1012+K0040	Manually Adjustable Center Mount Foot Platform..... \$250 Part: RIG163640. HCPCS: K0040
---	--

2. Lower Extension and Angle Adjustable Footplates

10"-16" Extension, 5"x6" Footplates Part: RIG150206	10"-16" Extension, 5"x11" Footplates Part: RIG150212	13"-19" Extension, 5"x8" Footplates Part: RIG151506
10"-16" Extension, 5"x8" Footplates Part: RIG150203	13"-19" Extension, 5"x6" Footplates Part: RIG151505	13"-19" Extension, 5"x11" Footplates Part: RIG151508

Center Mount Foot Platform Compatibility Matrix	Synergy Seat			
	Low	Medium	High	
10"-16" Length Extension	!	✓	✓	☒ No interference issues at this seat-to-floor height. ☐ May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
13"-19" Length Extension	!	!	!	



Headrest section continued on following page.

16. LEGRESTS CONTINUED**3. Center Mount Calf Support**

*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads (pair) 6"Wx6"H..... \$150*	Small Calf Panel 18"Wx8"H..... \$150
Part: RIGASMB7120018	Part: FRMASMB7110012
Large Calf Pads (pair) 8" taper to 6"Wx8"H \$150*	Large Calf Panel 24"Wx8"H \$200
Part: RIGASMB7120019	Part: RIG123824
	Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$50	4" Outside Edge (pair) \$50
Part: RIG152517	Part: RIG152510
Width Extension for 11"D Footplates \$50	Only available with 6"D footplates and 8"D footplate extension.
Part: RIG152518	6" Outside Edge (pair) \$50
2" Bell-Shaped Length Extension \$50	Part: RIG152516
Part: RIG152519	Available with 8"D and 11"D footplates and 11"D footplate extension.
Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.	

5. Additional Accessories

Reversed Dartex Calf Pad Cover (pair) \$150	Wheel Bumpers (pair)..... \$12.50 ea.
Small part: RIG133717. Large part: RIG133719	Part: WHLASMB1578
Reversed Dartex Panel Cover (pair) \$150	Footplate Connector \$15
Small part: RIG134673. Large part: RIG134744	Part: RIGASMB7120023

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers..... \$205	6.5"x8" Tapered Gel Footplate Covers for Width Extension..... \$225
Part: ST-FPG5008T-Q	Part: ST-FPG6508T-Q
5"x11" Tapered Gel Footplate Covers..... \$225	7"x11" Tapered Gel Footplate Cover for Width Extension \$250
Part: ST-FPG5011T-Q	Part: ST-FPG7011T-Q

16.3. Heavy Duty, Drop-In, Swing-Away Legrests

Select a legrest with extension and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Drop-In Legrest Selection

☐ P ☐ L ☐ R 70° HD Drop-in with 12"-14.5" No Charge	☐ P ☐ L ☐ R 70° HD Drop-in with 14.5-18" No Charge
Length Lower Extension	Length Lower Extension
Pair part: FRMASMB4762. Left: RIG137706. Right: RIG137707	Pair part: FRMASMB4853. Left: RIG137708. Right: RIG137709

2. Angle Adjustable Footplate(s) with Heel Loops

☐ P ☐ L ☐ R 6"Wx8"D ¹ \$85 ea.	☐ P ☐ L ☐ R 8"Wx11"D ² \$85 ea.
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.	Pair part: FRMASMB4744. Left: RIG137708. Right: RIG137709.
HCPCS: K0040+E0951	HCPCS: K0040+E0951
☐ P ☐ L ☐ R 6"Wx11"D ¹ \$85 ea.	May cause interference on seat widths less than 20".
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.	
HCPCS: K0040+E0951	
May cause interference on seat widths less than 18".	

3. Accessories

☐ P ☐ L ☐ R Wheel Bumper(s) \$12.50 ea.
Pair part: WHLASMB1578. Left part: WHL169510. Right part: WHL169511

17. SET LEGREST LENGTH

REQUIRED with matching left and right legrest configuration. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Non-matching or individual legrest configurations lengths are not set by the manufacturer. The part number will be PTOINDV plus the number after # below.

8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)
8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	

18. SEAT CUSHION


Optional. Specific cushion/component sizes available, part numbers, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

18.1. Seat Cushion Model Selection**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover. Weight capacity for 16" to 20" W is 350 lbs, over 20" W is 400 lbs.

Part series: ST-SIMGXXYY

Simplicity G, 21"W, E2601 \$110
Simplicity G, 22"W to 24"W, E2602 \$165

**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner. Weight capacity for 20" W is 350 lbs, over 20" W is 650 lbs.

Part series: ST-SPCFSPXXYY


Spectrum Foam SPP, 20"W, E2607 \$460
Spectrum Foam SPP, 22"W to 32"W, E2608 \$530

**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner. Weight capacity for 20" W is 350 lbs, over 20" W is 550 lbs.

Part series: ST-SOLSPPXXYY


Solution SPP, 20"W, E2607 \$375
Solution SPP, 22"W to 24"W, E2608 \$430

**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover. Weight capacity is for 20"W is 350 lbs, over 20"W is 450 lbs.

Part series: ST-GLSPXXYY


Glacial SP, 21"W, E2603 \$330
Glacial SP, 22"W to 24"W, E2604 \$440

**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover. Weight capacity for 20"W-21"W is 350 lbs, 22"W and over is 400 lbs.

Part series: ST-ZENSPXXYY


Zen SP, 20"W to 21" W, E2603 \$315
Zen SP, 22"W to 24" W, E2604 \$425

**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Weight capacity for 16" to 20" W is 350 lbs, over 20" W is 850 lbs.

Part series: ST-TC2SPPXXYY


TRU-Comfort 2 SPP, 20"W, E2607 \$400
TRU-Comfort 2 SPP, 22"W, E2608 \$500

**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner. Weight capacity for 20" W is 350 lbs, over 20" W is 550 lbs.

Part series: ST-SPCGSPPXXYY

Spectrum Gel SPP, 20"W, E2607 \$460
Spectrum Gel SPP, 22"W to 24"W, E2608 \$530

**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Weight capacity for 20" W is 350 lbs, over 20" W is 550 lbs.

Part series: ST-ESSSPPXXYY

Essence SPP, No Wedge, 20" Wide, E2607 \$400
Essence SPP, No Wedge, 22" Wide, E2608 \$500
Essence SPP, Small Wedge, 20" Wide, E2607 \$435
Essence SPP, Small Wedge, 22" Wide, E2608 \$535
Essence SPP, Medium Wedge, 20" Wide, E2607 \$465
Essence SPP, Medium Wedge, 22" Wide, E2608 \$565
Essence SPP, Large Wedge, 20" Wide, E2607 \$495
Essence SPP, Large Wedge, 22" Wide, E2608 \$595

**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover. Weight capacity for 20"W is 350 lbs, 22"W is 500 lbs..

Part series: ST-PREMPXXYY

Premiere P, 20"W, E2605 \$325
Premiere P, 22"W, E2606 \$475

**18.2. Seat Cushion Dimensions**

REQUIRED with cushion model selection. Select to match or pick non-matching dimensions. See respective cushion order form for available sizes.

Match Seat Width and Depth

Cushion Width

20" W 22" W 26" W 30" W
21" W 24" W 28" W 32" W

Cushion Depth

20" D 21" D 22" D 23" D 24" D 26" D 28" D

18.3. Seat Cushion Substitute/Add Fluid Proof Liners

Fluid Proof Urethane Liner \$50-\$70

Part: URL1709070-XXYY

Not available with Simplicity, Essence, or TRU-Comfort 2.

Fluid Proof Polyurethane Liner \$55-\$75

Essence/TC2 part: PUL1709080-XXYY. Others part: PUL1709030-XXYY

Not available with Simplicity.

18.4. Seat Cushion Substitute Cover

An additional cover may be ordered or the standard cover may be substituted with a cover by checking the option in subsection 2. The ""*prefix" of the part number will be the part prefix of the cushion model selected above. For example, SPCGSPP-NBLK1816 will be an additional coolcore cover for a Spectrum Gel cushion.

Stealth Tek 2 Cover No Charge

Part: *prefix-TBLKXXYY

Not available with Simplicity, Essence, TRU-Comfort 2, and Premiere.

Mesh Cover No Charge

Part: *prefix-MBLKXXYY

Not available with Simplicity or TRU-Comfort 2.

19. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be mounted by the provider.

Contoured Seat Oxygen Holder \$118 Part: FRMASMB2292. HCPCS: E2208	Saddle Bag \$35 Part: DWR1020D010
Contoured Seat Cup Holder \$20 Non-bariatric part: ACCASMB2287. Bariatric arm part: ACCASMB2632	Stealth Cell Phone/Media Device Holder \$195 Part: ACC7120006
Synergy Seat Cup Holder No Charge Part: ACCASMB2634	Stealth Large Tablet Holder \$375 Part: ACC7120005
Synergy Seat Oxygen Holder \$150 Not available with Synergy Manual Recline. Part: FRMASMB7030. HCPCS: E2208	Quantum Backpack \$150 Part: ACC1905853
Walker Holder \$103 Synergy part: FRMASMB7031	Grip Lap Board (20005) \$38.95 Part: ACC152992
Cane/Crutch Holder \$110 Contoured Seat part: ACCASMB2294. HCPCS: E2207 Synergy part: FRMASMB4822. HCPCS: E2207	Grip Drink Holder (20085) \$15.50 Part: ACC152993
Rear Basket \$118 Contoured seat part: ACCASMB1914 Synergy part: FRMASMB4821	Accessory Bar \$100 20"W-24"W part: MEC146018. 24"W-28"W part: MEC146050. 28"W-32"W part: MEC146051
Flag Holder \$25 Contoured seat part: ACCASMB1930 Synergy part: ACCASMB2415	Large, Black, Weather Cover \$202 Part: ACCCOVR1021

20. BRANCH PATCHES**Optional.**

Air Force Patch \$20 Part: ACC140240	Coast Guard Patch \$20 Part: ACC140239	Navy Patch \$20 Part: ACC140244
Army Patch \$20 Part: ACC140241	Marines Patch \$20 Part: ACC140243	Purple Heart Patch \$20 Part: ACC140238

**21. ELECTRONICS ACCESSORIES****Optional.**

Q-Logic EX/NE Programmer \$545 Part: ELEASMB6976	Q-Logic PC Programming Disk with Cable \$750 Part: ELEASMB5216
---	---

22. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional. Only available with the Synergy Seating system.** Select a type of specialty controls below or to refer to the specialty controls order form.**22.1. Enhanced Display**

Q-Logic 3 Enhanced Display Kit \$1,720
Part: CTL167726
Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud**Required with Q-Logic 3 Enhanced Display** selection above. Reference colors shown on page 2.

Black, Soft Touch Standard Part: DGN137728	Match Joystick Shroud to Base Color \$50 Part: See color below
Non-matching Joystick Shroud \$50 Select color below	

Red Corvette Part: DGN137729-PRD1	Back in Black Part: DGN137729-PBK1	Pretty in Pink Part: DGN137729-PPK3	Deep Purple Part: DGN137729-PPR5	Midnight Blue Part: DGN137729-PBU10
Electric Blue Part: DGN137729-PBU9	Night Train Part: DGN137729-PBK3	Evolution Orange Part: DGN137729-POR7	Grey Street Part: DGN137729-PSL3	Green Machine Part: DGN137729-PGN11

**22.2. SCIM Module****Optional.**

SCIM Module Kit for Q-Logic 3 \$525
Part: CTL168055
Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.

**22.3. Q-Logic 3 SAJ**

Q-Logic 3 EX SAJ \$700
Part: CTL143980
Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



Alternate controls and electronics section continued on following page.

22. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**22.4. Head Array**

Adult, i-Drive Tri-Array Series, 3-Switch \$4,414
 Head Array (IDHBT200-1)
 Part: ST-IDHBT200-1-Q. HCPCS: E2330
 Features bluetooth, head support, 3 non-proportional
 proximity sensors for drive control, an egg switch for re-
 set mode change, and i-Drive smart box.

Adult, i-Drive Pro Series Head Array \$4,461
 (IDHBT300-1)
 Part: ST-IDHBT300-1-Q. HCPCS: E2330
 Features bluetooth, head support with occipital pad and
 2 swing-away spot pads for left and right side, 3 non-pro-
 portional proximity sensors, an egg switch for reset
 mode change, and i-Drive smart box.

**22.5. Sip and Puff Drive****1. Sip and Puff Module**

Sip and Puff Module for Q-Logic 3 \$1,895
 Part: CTL168056. HCPCS: E2325
 Includes module and mounting hardware.

2. Sip and Puff Straw and Mount

Stealth Sip and Puff for Comfort Series, \$294.25
 Left Mount (SP625-CL)
 Part: ST-SP625-CL-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth
 Comfort Plus headrests.

Stealth Sip and Puff for Stealth Ultra, \$294.25
 3-Switch Head Array (SP625-U)
 Part: ST-SP625-U-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth
 i-Drive Head Arrays.



Stealth Sip and Puff for Comfort Series, \$294.25
 Right Mount (SP622-CR)
 Part: ST-SP625-CR-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth
 Comfort Plus headrests.

**22.6. Egg Switches**

Red Egg Switch (ESRED) \$80.85 ea.
 Part: SWTMCRO1077

Green Egg Switch (ESGRN) \$80.85 ea.
 Part: SWTMCRO1078

Blue Egg Switch (ESBLU) \$80.85 ea.
 Part: SWTMCRO1079

Yellow Egg Switch (ESYEL) \$80.85 ea.
 Part: SWTMCRO1080

Black Egg Switch (ESBLK) \$80.85 ea.
 Part: SWTMCRO1081

**22.7. Switch Adapter**

2 to 1 Switch Adapter \$13.39
 Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

22.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,400
 Part: CTL168057

23. POSITIONING COMPONENTS

Optional. Not available with contoured seating. More options are available on the [Positioning Components Order Form](#).

23.1. Stealth Thoracic Laterals

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, Rail Mounts, pair (TWBTLTC) \$472.50
 Part: INDPART2896. HCPCS: E1028

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5) \$210
 Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6) \$210
 Part: INDPART2664. HCPCS: E0956

3. Installation

Install Laterals
 Part: PTOINDV3499

Do Not Install
 Part: PTOINDV3500

Positioning components section continued on following pages.

23. POSITIONING COMPONENTS CONTINUED**23.2. Stealth UniLink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable & Removable Hardware... \$245 ea.
 QTY Part: ST-UTTP5-Q

UniLink, 7" Height Adjustable & Removable Hardware... \$255 ea.
 QTY Part: ST-UTTP7-Q

2. UniLink Pelvic/Thigh Pads

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad\$95 ea.
 QTY Part: ST-UTP544056100-Q

UniLink, 4" x 12" Pad \$125 ea.
 QTY Part: ST-UTP54D056100-Q

UniLink, 4" x 8" Pad..... \$105 ea.
 QTY Part: ST-UTP548056100-Q

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$75 ea.
 QTY Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$75 ea.
 QTY Part: ST-UIL001200A00-Q



UniLink Pad & Hardware

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

23.3. Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW) \$236.25 ea.
 QTY Part: INDPART2804. HCPCS: E1028

Removable Bracket (LP430)..... \$262.50 ea.
 QTY Part: INDPART2805. HCPCS: E1028

Removable Bracket with 1" Offset \$400 ea.
 QTY (LP430 & LHW-121)
 Part: INDPART2854. HCPCS: E1028

Flip Down with LEFT Full Surface Contact \$236.25 ea.
 Feature (TWBADD-L & LPHW407)
 Part: POS140428. HCPCS: E1028

Flip Down with RIGHT Full Surface Contact..... \$236.25 ea.
 Feature (TWBADD-R & LPHW407)
 Part: POS140429. HCPCS: E1028

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122).....\$26.25
 Part: INDPART2821. HCPCS: K0108
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129)\$53.50
 Part: POS140430. HCPCS: K0108
 Only available with the flip down full surface bracket.

3. Pelvic/Thigh Guide Pad

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Short Pad, 4"x4" (TWBADD-4X4)..... \$75 ea.
 QTY Part: INDPART2801. HCPCS: E0953^/E0956

Medium Pad, 4"x8" (TWBADD-4X8)..... \$75 ea.
 QTY Part: INDPART2802. HCPCS: E0953^/E0956

Long Pad, 4"x12" (TWBADD-4X12)..... \$75 ea.
 QTY Part: INDPART2803. HCPCS: E0953^/E0956

4. Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

23.4. Stealth i-Fit Chest Harness

Upper Classic, Medium (IFPS1103) \$142 ea.
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$142 ea.
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$142 ea.
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$142 ea.
 Part: ST-IFPS2104-Q. HCPCS: E0960



Contoured and Classic Chest Harness

23.5. Stealth i-Fit Chest Straps

Chest Strap, Medium w/ Belt Anchor \$90 ea.
 (IFCS4113)
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor \$90 ea.
 (IFCS4114)
 Part: ST-IFCS4114-Q. HCPCS: E0960



Chest Strap

23.6. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. The 1" ball kit must be selected to mount the facial lateral to the headrest bracket. Also select a facial pad(s) with the facial lateral(s).

1. Facial Laterals

Facial Lateral for Comfort Series.....\$294.25 ea.
 Right Side (CPS9-R-HW)
 Part: INDPART2647. HCPCS: E1028

Facial Laterals for Comfort Series \$347.75 ea.
 Left and Right Side (CPS9-B-HW)
 Part: INDPART2649. HCPCS: E1028

Facial Lateral for Comfort Series.....\$294.25 ea.
 Left Side (CPS9-L-HW)
 Part: INDPART2648. HCPCS: E1028



Left and Right Facial Lateral with Pads

Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**2. Facial Lateral Pads and Accessories**

1" Ball Kit for Facial Laterals (TWB-1BK) (D) \$80.25
 Part: POS150532
 Required with facial laterals to mount to comfort headrest brackets.

Egg Switch Mount (ESM) (A) \$52.50 ea.
 Part: INDPART2815


Figure 8 Bracket with 13" Rod \$267.50 ea.
 and Egg Switch Mount (B)
 Part: ST-F8675-13-Q

Facial Pad, 4"x2" (900-Q) \$53.50 ea.
 Part: ST-900-Q. HCPCS: K0108

Facial Pad, 6"x2" (920-Q) \$53.50 ea.
 Part: ST-920-Q. HCPCS: K0108

Swing-Away Joystick Mount, \$250 ea.
 18" Rod (C)
 Part: ST-SUS18-Q

Allows mounting of facial pads and attaches to the headrest

**23.7. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket \$183.75 ea.
 Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket \$183.75 ea.
 Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$52.50 ea.
 Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$52.50 ea.
 Part: POS124067. HCPCS: K0108

23.8. Stealth Foot Boxes

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Stealth Small Foot Boxes, 5"Wx8"Dx7"H \$498.75
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954^

Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$498.75
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954^

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$498.75
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954^

23.9. Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D \$266
 (SPLE111) (pair)
 Part: INDPART2783. HCPCS: E0951+E0952

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D \$266
 (SPLE114) (pair)
 Part: INDPART2785. HCPCS: E0951+E0952

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D \$266
 (SPLE112) (pair)
 Part: INDPART2784. HCPCS: E0951+E0952

24. TRU-BALANCE 3 MULTIPLIER HARNESS

REQUIRED. Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

No Additional Multiplier (NA) No Charge
 Part: Configuration selection

Secondary Multiplier (SM) \$221.21
 Part: HAR144004

Special Order Configuration (SOF)
 Part: Configuration selection

Arm Mounts
 Joystick (Hand Control)
 Enhanced Display
 iAccess
 USB Charger
 SAJ

Back Mounts
 Attendant Control
 SCIM
 Sip n Puff
 ECU

	1 Arm	2 Arm	3+ Arm
0 Back	NA	SM	SOF
1 Back	NA	SM	SOF
2 Back	SM	SM	SOF
3+ Back	SOF	SOF	SOF

Q-Logic 3 Multiplier with Synergy Seat Chart

25. REFERENCE ORDER FORMS

Optional. Select below to note that additional order forms will accompany this order form. This will also be a link in interactive forms.

[Alternative Controls & Electronics Order Form](#)

[Positioning Components Order Form](#)

Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2019 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.

