



Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® HD Order Form with TRU-Balance® Bariatric Tilt

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

Q6 edge[®]
HD



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____



A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

1. BASE MODEL

REQUIRED. Select a group 3 power base. 450 lb. weight capacity with standard 4.7 mph programmable motors, 8 amp off board charger, and unoccupied transit loops. Dimensions of base without legrests: 26"Wx36.75"L.

Quantum Q6 Edge HD 3SPHD-SS \$9,025

Part: Q6 EDGEHD 3SPHD-SS. HCPCS: K0858

Single power option.

2. BASE COLOR**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

 Red Corvette Part: DGN169178-PRD1	 Night Train Part: DGN169178-PBK3	 Deep Purple Part: DGN169178-PPR5	 Green Machine Part: DGN169178-PGN11
 Electric Blue Part: DGN169178-PBU9	 Pretty in Pink Part: DGN169178-PPK3	 Grey Street Part: DGN169178-PSL3	
 Back in Black Part: DGN169178-PBK1	 Evolution Orange Part: DGN169178-POR7	 Midnight Blue Part: DGN169178-PBU10	

3. FENDER LIGHTS**Optional.** Fender lights operate by pressing included single switch button. Fender lighting is not available with NE 4-Key Joystick or Single Actuator Control through Toggle Box.Switch Fender Lights \$263
Part: XRF000177**4. QUANTUM BACKUP CAMERA****Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.Quantum Backup Camera \$368
Part: MEC1910365
Available for Q-Logic 3 electronics only. NE/NE+ coming soon.**1. Optional Camera Mounting Bracket**Backup Camera Headrest Bracket \$21
Part: MEC1910373
Works with .5" x .5" vertical supports only.Backup Camera Track Mount Bracket \$27
Part: MEC2104227**2. Optional Display Mounting Bracket**Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123Backup Camera NE/NE+ Joystick Mount Bracket \$105
Part: MEC2106132**5. DRIVE WHEELS****REQUIRED.** Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".**5.1. Drive Wheels**

3" x 14" Run Flat Tires Standard Part: PTOINDV3510	4" x 14" Flat Free, Knobby, Extended \$290 ea. Traction Tires Part: WHL121698 x2
3" x 14" Flat Free Tires No Charge Part: WHL124116 x2	4" x 14" Pneumatic, Knobby, Extended \$290 ea. Traction Tires Part: WHL121703 x2
3" x 14" Pneumatic Tires No Charge Part: WHL124117 x2	

5.2. Drive Wheel Accent Color**Optional. Not available with knobby tires.**

Back in Black Accent Color \$63 Part: XRF000179	Pretty in Pink Accent Color \$63 Part: XRF000185	Midnight Blue Accent Color \$63 Part: XRF000187
Red Corvette Accent Color \$63 Part: XRF000180	Evolution Orange Accent Color \$63 Part: XRF000186	Green Machine Accent Color \$63 Part: XRF000190
Electric Blue Accent Color \$63 Part: XRF000188	Deep Purple Accent Color \$63 Part: XRF000184	
Night Train Accent Color \$63 Part: XRF000181	Grey Street Accent Color \$63 Part: XRF000182	

6. BATTERIES**Optional.** The powerchair requires 2 batteries and will be installed if ordered.Group 24, 70 Amp Gel Cell, Interceptor 2 Battery \$384 ea.
Part: BAT2100315. HCPCS: E2363Group 24, 70 Amp Gel Cell, MK Battery \$520 ea.
Part: BATGEL1009. HCPCS: E2363**7. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Forms](#).**7.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick Standard Part: CTLDC1554 Available with Static Seat or Tilt Only through Toggle.	NE+ 6-Key Joystick \$210 Part: XRF000001 Control up to 2 actuators through joystick.
NE+ 6-Key Joystick No Charge Part: CTLDC1560 Control of one actuator through Joystick on single power base model (3SPHD-SS) only.	NE+ 8-Key Joystick \$420 Part: CTLDC1561 Control up to 2 actuators and LED, Full Light Package through joystick.
No Joystick, Setup for NE+ Joystick No Charge Part: XRF000009	



Electronics section continued on following page.

7. ELECTRONICS CONTINUED**7.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers. This sections allows a selection of common configurations.

1. Q-Logic 2 Required Components

Q-Logic 2 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

Q-Logic 2 EX Controller Module\$730
Part: XRF000197. HCPCS: E2377

Q-Logic 2 Expandable Controller Harness.....\$447
Part: DWR1111H065. HCPCS: E2313

2. Q-Logic 2 Joystick

Q-Logic 2 EX Joystick\$735
Part: CTL123217

Q-Logic 2 EX Joystick for Full Lights\$840
Part: CTL124653

No Joystick, Setup for Q-Logic 2

Part: XRF000010

Requires selection of Enhanced Display.

Stand Alone Joystick with\$840
External Charger Port

Part: ELEASMB5205

Up to 3 drive profiles and 1 seating profile when used alone. Use in conjunction with Enhanced Display for up to 5 profiles and auxiliary functions such as Bluetooth Mouse and IR. Includes charger port mounting bracket and multiplier box harness.

3. Q-Logic 2 Joystick Armor

REQUIRED with Q-Logic 2 EX Joystick selection above. Reference colors shown in section 2.

Black, Soft Touch Standard
Part: DGN124307

Night Train.....\$63
Part: DGN124336-PBK3

Red Corvette.....\$63
Part: DGN124336-PRD1

Pretty in Pink\$63
Part: DGN124336-PPK3

Electric Blue\$63
Part: DGN124336-PBU9

Evolution Orange.....\$63
Part: DGN124336-POR7

Back in Black.....\$63
Part: DGN124336-PBK1

Deep Purple.....\$63
Part: DGN124336-PPR5

Grey Street\$63
Part: DGN124336-PSL3

Midnight Blue\$63
Part: DGN124336-PBU10

Green Machine\$63
Part: DGN124336-PGN11



Stand Alone Joystick



Q-Logic 2 EX Joystick

8. JOYSTICK HANDLES

Optional. The Q-Logic 2 EX joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323

Soft Ball Handle (F) \$80
Part: ACCASMB1856. HCPCS: E2323

Chin Cup (C) \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323

T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323

Bodypoint 3" Goal Post Handle\$104
Part: PLSKNCP1057. HCPCS: E2323

Medium Ball Handle (E) \$69
Part: ACCASMB1859. HCPCS: E2323

Bodypoint 4" Goal Post Handle\$104
Part: PLSKNCP1058. HCPCS: E2323



Joystick Handles

9. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

9.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket for Q-Logic.....Standard
Part: FRMASMB10617

Fixed, In-line, Bracket for NE/NE+Standard
Part: FRMASMB12744

9.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic \$294
Part: ACC1908434. HCPCS: E1028

Swing-away, NE/NE+.....\$294
Part: ACC1908437. HCPCS: E1028


9.1c. Basic Swing Away Joystick Mount

Basic Swing Away, NE/NE+, Left Side.....\$245
Part: ACC2201047. HCPCS: E1028

Basic Swing Away, NE/NE+, Right Side.....\$245
Part: ACC2201048. HCPCS: E1028

Basic Swing Away, Q-Logic, Left Side\$245
Part: ACC2201041. HCPCS: E1028

Basic Swing Away, Q-Logic, Right Side\$245
Part: ACC2201042. HCPCS: E1028



Swing Away Joystick

Joystick mounting brackets section continued on following page.

9. JOYSTICK MOUNTING BRACKETS CONTINUED**9.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount\$368

Part: XRF000003. HCPCS: E1028

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount \$126


Part: XRF000004

Allows height adjustment of joystick.

Stealth Swing-Away/Flip Down Bracket\$250

Part: XRF000005. HCPCS: E1028

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.



Retract4



Multi-Axis Joystick Mount

9.2. Joystick Mount Location**Left Side Mount**

Part: PTOINDV1225

Right Side Mount

Part: PTOINDV1226

10. TRU-BALANCE BARIATRIC POWER TILT**REQUIRED.****10.1. Bariatric Tilt Seating****TRU-Balance Bariatric Power Tilt**\$7,613

Part: PTOINDV3263. HCPCS: E1002

40° center of gravity power tilt.

10.2. Seat Dimension - Width**REQUIRED.** Maximum seat width is 28". 15" depth is not available on widths over 23"W. Part number is SETWDTH plus number following # below.

20" W (#1010)	21" W (#1011)	22" W (#1012)	23" W (#1013)	24" W (#1014)	25" W (#1056)	26" W (#1054)	27" W (#1057)	28" W (#1055)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

10.3. Seat Dimension - Width Package**REQUIRED.****20"-24" Width Package**

Part: MEC120322

24"-28" Width Package

Part: MEC120323

10.4. Seat Dimension - Depth**REQUIRED.** Maximum seat depth is 28". Part number is SETDPTH plus number following # below.

15" D (#1005)	17" D (#1007)	19" D (#1009)	21" D (#1011)	23" D (#1013)	25" D (#1060)	27" D (#1062)
16" D (#1006)	18" D (#1008)	20" D (#1010)	22" D (#1012)	24" D (#1014)	26" D (#1061)	28" D (#1063)

10.5. Seat Dimension - Depth Package**REQUIRED.****15"-20" Depth Package**

Part: FRMASMB10521

Not available with power articulating foot platform.

20"-24" Depth Package

Part: FRMASMB10524

24"-28" Depth Package

Part: MEC147637

10.6. Seat to Floor Height**REQUIRED.****Middle Seat to Floor**

Part: PTOINDV3559

Highest Seat to Floor

Part: PTOINDV3560

Med	High
Bariatric Tilt	
19.1875"	20.1875"
Seat to Floor Height	

10.7. Matching Seat CushionOptional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion that matches the seating width and depth. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Weight capacity for 16" to 20" W is 350 lbs, over 20"W is 450 lbs.

Part series: ST-TC2SPPXXYY. HCPCS: E2607/E2608

TRU-Comfort 2 SPP, Matching Seat Size

Reference part: XRF000116. MSRP: \$420/\$525

**Stealth Solution SPP Cushion**


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Includes urethane liner. Weight capacity for 20" W is 350 lbs, over 20"W is 550 lbs.

Part series: ST-SOLSPPXXYY. HCPCS: E2607/E2608

Solution SPP, Matching Seat Size

Reference part: XRF000113. MSRP: \$394/\$452

**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Includes urethane liner. Weight capacity for 20" W is 350 lbs, over 20"W is 550 lbs.

Part series: ST-SPCGSPPXXYY. HCPCS: E2607/E2608

Spectrum Gel SPP, Matching Seat Size

Reference part: XRF000117. MSRP: \$483/\$557

**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Includes urethane liner. Weight capacity for 20" W is 350 lbs, over 20"W is 550 lbs.

Part series: ST-SPCFSPXXYY. HCPCS: E2607/E2608

Spectrum Foam SPP, Matching Seat Size

Reference part: XRF000112. MSRP: \$483/\$557



TRU-Balance Bariatric Power Tilt section continued on following page.

10. TRU-BALANCE BARIATRIC POWER TILT CONTINUED**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover. Weight capacity for 20"W is 350 lbs, 22"W is 500 lbs.

Part series: ST-PREMPXXYY. HCPCS: E2605/E2606

Premiere P, Matching Seat Size

Reference part: XRF000122. MSRP: \$341/\$499

**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover. Weight capacity for 20"W-21"W is 350 lbs, 22"W and over is 400 lbs.

Part series: ST-ZENSPXXYY. HCPCS: E2603/E2604

Zen SP, Matching Seat Size

Reference part: XRF000115. MSRP: \$331/\$446


**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover. Weight capacity for 16" to 20" W is 350 lbs, over 20"W is 450 lbs.

Part series: ST-GLSPXXYY. HCPCS: E2603/E2604

Glacial SP, Matching Seat Size

Reference part: XRF000114. MSRP: \$347/\$462

**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Weight capacity for 20" W is 350 lbs, over 20"W is 550 lbs.

Part series: ST-ESSSPXXYY. HCPCS: E2607/E2608

Essence SPP, No Wedge, Matching Seat Size

Reference part: XRF000118. MSRP: \$420/\$525

Essence SPP, Small Wedge, Matching Seat Size

Reference part: XRF000119. MSRP: \$457/\$562

Essence SPP, Medium Wedge, Matching Seat Size

Reference part: XRF000120. MSRP: \$488/\$593

Essence SPP, Large Wedge, Matching Seat Size

Reference part: XRF000121. MSRP: \$520/\$625


**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover. Weight capacity for 16" to 20" W is 350 lbs, over 20"W is 400 lbs.

Part series: ST-SIMGXXYY. HCPCS: E2601/E2602

Simplicity G, Matching Seat Size

Reference part: XRF000111. MSRP: \$116/\$173

**10.8. Back Cane Type****REQUIRED.**

Angle Adjustable Back Canes..... No Charge

Part: FRMASMB15845

2" Mesomorph (1" per side) Back Canes..... \$205

Part: FRMASMB16387

2" Endomorph (1" per side) Back Canes..... \$205

Part: FRMASMB16393

4" Endomorph (2" per side) Back Canes \$205

Part: FRMASMB16396

10.9. Back Cane Height**REQUIRED.**

18" High Canes

Part: FRMASMB15836

19" High Canes

Part: FRMASMB15837

20" High Canes

Part: FRMASMB15838

21" High Canes

Part: FRMASMB15839

22" High Canes

Part: FRMASMB15840

23" High Canes

Part: FRMASMB15841

24" High Canes

Part: FRMASMB15842

25" High Canes

Part: FRMASMB15843

10.10. Back Cushion Selection**REQUIRED.** See the TRU-Balance Bariatric Power Tilt back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1. Synergy Back Selection**

Synergy Back..... No Charge

Part: XRF000199

Select overall back height in later section. Back cushion width is determined by seat width and back cane type.

2. TRU-Comfort Back Selection

TRU-Comfort Back Cushion, 20" Width or Less..... \$709

HCPCS: E2620. Select size below. Height is overall. Part number is SETASMB plus number following #.

16"Wx21"H

Part: #7151

16"Wx25"H

Part: #7159

17"Wx23"H

Part: #7155

17"Wx28"H

Part: #7167

18"Wx23"H

Part: #7156

18"Wx28"H

Part: #7168

19"Wx25"H

Part: #7162

20"Wx21"H

Part: #7153

20"Wx25"H

Part: #7163

16"Wx23"H

Part: #7154

16"Wx28"H

Part: #7166

17"Wx25"H

Part: #7160

18"Wx21"H

Part: #7152

18"Wx25"H

Part: #7161

19"Wx23"H

Part: #7157

19"Wx28"H

Part: #7169

20"Wx23"H

Part: #7158

20"Wx28"H

Part: #7170

TRU-Comfort Back Cushion, 22" Width..... \$814

HCPCS: E2621. Select size below. Height is overall. Part number is SETASMB plus number following #.

22"Wx25"H

Part: #7164

22"Wx28"H

Part: #7171

24"Wx25"H

Part: #7165

24"Wx28"H

Part: #7172

3. Omit Back Selection

Omit Back..... No Charge

Part: PTOINDV1287

4. Overall Back Height

Refer to the TRU-Balance Bariatric Power Tilt back availability chart for stock order back size ranges.

17" Overall Back Height

Part: XRF000139

18" Overall Back Height

Part: XRF000140

19" Overall Back Height

Part: XRF000141

20" Overall Back Height

Part: XRF000142

21" Overall Back Height

Part: XRF000143

22" Overall Back Height

Part: XRF000144

23" Overall Back Height

Part: XRF000145

24" Overall Back Height

Part: XRF000146

25" Overall Back Height

Part: XRF000147

26" Overall Back Height

Part: XRF000200

27" Overall Back Height

Part: XRF000201

28" Overall Back Height

Part: XRF000202

		Back Width										
		20"	21"	22"	23"	24"	25"	26"	27"	28"		
Overall Height	17"	S	S	S	S	S	-	-	-	-	Legend	
	18"	S	S	S	S	S	-	-	-	-		
	19"	S	S	S	S	S	S	S	S	S		
	20"	S	S	S	S	S	S	S	S	S		
	21"	S	S	S	S	S	S	S	S	S		
	22"	S	S	S	S	S	S	S	S	S		
	23"	S	S	S	S	S	S	S	S	S		
	24"	S	T	S	S	S	S	T	S	S		
	25"	S	T	S	S	S	S	T	S	S		
	26"	S	S	S	S	S	S	-	-	-		
	27"	-	-	-	-	-	-	-	-	-		
	28"	T	-	T	-	T	-	-	-	-		

Synergy and TRU-Comfort Back Availability Chart

TRU-Balance Bariatric Power Tilt section continued on following page.

10. TRU-BALANCE BARIATRIC POWER TILT CONTINUED**10.11. Lateral Extensions****REQUIRED with Synergy Back and thoracic lateral selections.**Lateral Extension for Synergy Back (pair)..... \$56
Part: XRF000150**10.12. Bariatric Tilt Armrests****REQUIRED.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. Bariatric, Single Post, Heavy Duty, Height Adjustable, Removable Armrests**

Left and right armrest heights must be the same armrest height. May interfere with swing-away laterals.

☐ ☐	Full Length, Bariatric Armrests, 8"-12" H No Charge	☐ ☐	Full Length, Bariatric Armrests, 10"-14"H No Charge
☐ ☐	Left part: FRMASMB6728. Right part: FRMASMB6729	☐ ☐	Left part: FRMASMB6730. Right part: FRMASMB6731.

1b. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

☐ ☐	Full Length, Heavy Duty Armrests, 8"-12" H No Charge	☐ ☐	Desk Length, Heavy Duty Armrests, 8"-12"H No Charge
☐ ☐	Left part: FRMASMB9982. Right part: FRMASMB9983.	☐ ☐	Left part: FRMASMB9984. Right part: FRMASMB9985.

1c. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

☐ ☐	Full Length, Heavy Duty Armrests, 10"-14"H No Charge	☐ ☐	Desk Length, Heavy Duty Armrests, 10"-14"H No Charge
☐ ☐	Left part: FRMASMB9986. Right part: FRMASMB9987.	☐ ☐	Left part: FRMASMB9988. Right part: FRMASMB9989.

1d. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Angle adjustable from -30° to +30°.

☐ ☐	Full Length, Cantilever Armrests \$103 ea.	☐ ☐	Desk Length, Cantilever Armrests \$103 ea.
☐ ☐	Left part: FRMASMB16188. Right part: FRMASMB16189	☐ ☐	Left part: FRMASMB16190. Right part: FRMASMB16191

10.13. Armrest Height**REQUIRED.** Set left and right armrest height. Note the available height range in the option description.**1. Left Armrest Height**

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

2. Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

10.14. Bariatric Tilt Armpads**REQUIRED with armrest selection above.** Armpads are configured as pairs.

Bariatric Armpads No Charge	No Armpad Selection No Charge
Part: UPSASMB1025	Part: PTOINDV1373
Only available with bariatric armrests.	
Straight Armpads Standard	Waterfall Armpads No Charge
Part: XRF000148	Part: XRF000149

11. JOYSTICK BRACKET RECEIVER**REQUIRED.**

Standard Fixed Receiver Standard	Stealth Height Adjustable Receiver \$225
Part: XRF000063	Part: XRF000064. HCPCS: K0108

12. HEADRESTS**Optional.** Select a headrest pad and bracket.**12.1a. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$170.25	Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) \$170.25
Part: ST-CP08CCFPM-Q. HCPCS: E0955	Part: ST-CP14CCFPM-Q. HCPCS: E0955
Stealth 10" Comfort Plus, Cool Core (CP10CCFPM) \$170.25	Stealth 6", Comfort Plus, Cool Core (CP06CCFPM) \$170.25
Part: ST-CP10CCFPM-Q. HCPCS: E0955	Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink..... \$296.50
Part: ST-UHB2D33200-Q. HCPCS: E1028

Headrests section continued on following page.

12. HEADRESTS CONTINUED**12.1b. Headrests, Comfort Plus and Stealth Contoured****1. Headrest Pads**Stealth 8" Comfort Plus (CP450)..... \$170.25
Part: POS148484. HCPCS: E0955Stealth 14" Comfort Plus (CP150). \$170.25
Part: INDPART2327. HCPCS: E0955Stealth 10" Comfort Plus (CP250) (A) \$170.25
Part: SETHEAD1057. HCPCS: E0955**2. Headrest Brackets**Stealth Multi-Axis Removable Headrest Bracket (HMO475P) (D) \$165.32
Part: HDWASMB1084. HCPCS: E1028**13. LAP BELTS****REQUIRED.** Select a lap belt below.**13.1a. Standard Lap Belts**38" Lap Belt
Part: ACCBELT100650" Lap Belt
Part: ACCBELT100060" Lap Belt
Part: ACCBELT100270" Lap Belt
Part: ACCBELT100380" Lap Belt
Part: ACCBELT100490" Lap Belt
Part: ACCBELT1001**13.1b. Optional Lap Belts**Stealth Retractable Lap Belt (IFRB5350)..... \$172.53
Part: ST-IFRB5350Stealth Positioning 60" Single Pull, Push Button, Large \$47.38
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-QStealth Positioning 60" Single Pull, Push Button, Medium..... \$47.38
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-QStealth Positioning 60" Single Pull, Push Button, \$47.38
Extra Large, Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q**14. BASE MOUNTED FOOT PLATFORM****REQUIRED.** A base mounted foot platform is not available with center mount foot platforms. May interfere with seat mounted legrests.Base Mounted Foot Platform Standard
Part: FRMASMB16077Omit Base Mounted Foot Platform No Charge
Part: PTOINDV3674Base Mounted, Bell Shaped Foot Platform No Charge
Part: RIGASMB7120041**15. LEGRESTS****Optional. Please note the set lower leg length may vary based on your seat type and seat to floor height selections.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.**15.1a. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform AssemblyPower Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040Manually Adjustable Center Mount Foot Platform..... \$263
Part: XRF000155. HCPCS: K0040**2. Lower Extension and Angle Adjustable Footplates**Medium Extension, 5"x6" Footplates
Part: RIG150206Long Extension, 5"x6" Footplates
Part: RIG151505Medium Extension, 5"x8" Footplates
Part: RIG150203Long Extension, 5"x8" Footplates
Part: RIG151506Medium Extension, 5"x11" Footplates
Part: RIG150212Long Extension, 5"x11" Footplates
Part: RIG151508

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Medium	Long
10"-16"	13.5"-19.5"
Manually Adjustable Extension Range	

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019Small Calf Pads, Adjustable Center Mount Only, \$158
6"Wx6"H (pair)
Part: XRF000156Large Calf Pads, Adjustable Center Mount Only, \$158
8" taper to 6"Wx8"H (pair)
Part: XRF000157**3b. Center Mount Calf Panel**

Calf pads are standard with a power articulating foot platform.

Small Calf Panel for AFP, 18"Wx8"H..... \$184
Part: FRMASMB7110007Small Calf Panel for Adjustable Center Mount, 18"Wx8"H \$184
Part: FRMASMB7110012Large Calf Panel for AFP, 24"Wx8"H \$237
Part: RIG123824Large Calf Panel for Adjustable Center Mount, 24"Wx8"H..... \$237
Part: RIG123825

Only available on widths of 20" or greater.

Only available on widths of 20" or greater.

Legrests section continued on following page.

15. LEGRESTS CONTINUED**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$69 Part: RIG152517	4" Outside Edge (pair) \$69 Part: RIG152510 Only available with 6"D footplates and 8"D footplate extension.
Width Extension for 11"D Footplates \$69 Part: RIG152518	6" Outside Edge (pair) \$69 Part: RIG152516 Available with 8"D and 11"D footplates and 11"D footplate extension.
2" Bell-Shaped Length Extension \$69 Part: RIG152519 Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.	

5. Center Mount Additional Accessories

Small Reversed Dartex Calf Pad Cover (pair) \$158 Part: RIG133717	Large Foam Calf Pad Cover (pair) \$78.72 Part: ST-CSCS0985-Q
Large Reversed Dartex Calf Pad Cover (pair) \$210 Part: RIG133719	Small Foam Panel Cover (pair) \$167.57 Part: ST-CSCP0920-Q
Small Reversed Dartex Panel Cover (pair) \$210 Part: RIG134673	Large Foam Panel Cover (pair) \$268.77 Part: ST-CSCP0927-Q
Large Reversed Dartex Panel Cover (pair) \$210 Part: RIG134744	Footplate Connector \$16 Part: RIGASMB7120023
Small Foam Calf Pad Cover (pair) \$63.74 Part: ST-CSCS0775-Q	

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$211.15 Part: ST-FPG5008T-Q. HCPCS: K0108	5"x8" Tapered Gel Footplate Covers for Width Extension ... \$231.75 Part: ST-FPG6508T-Q. HCPCS: K0108
5"x11" Tapered Gel Footplate Covers \$231.75 Part: ST-FPG5011T-Q. HCPCS: K0108	5"x11" Tapered Gel Footplate Cover for Width Extension ... \$257.50 Part: ST-FPG7011T-Q. HCPCS: K0108

15.1b. Heavy Duty, Drop-In, Swing-Away Legrests

Select a configuration, a legrest with extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Drop-In Legrest Selection

☐ P ☐ L ☐ R 70° HD Drop-in for 8.5"-12" Lower Extension No Charge Pair part: FRMASMB4762. Left: RIG137706. Right: RIG137707	☐ P ☐ L ☐ R 70° HD Drop-in for 13"-18" or 16"-21" Lower No Charge Extension Pair part: FRMASMB4853. Left: RIG137708. Right: RIG137709
---	--

2. Lower Extension

☐ P ☐ L ☐ R 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103.	☐ P ☐ L ☐ R 16"-21" HD Drop-Ins/18.5"-22.5" ELR, Ext Length Pair part: FRMASMB4754. L: FRMASMB4725. R: FRMASMB5393
☐ P ☐ L ☐ R 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106	

3. Angle Adjustable Footplate(s) with Heel Loops

☐ P ☐ L ☐ R 6"Wx8"D \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795. HCPCS: K0040+E0951	☐ P ☐ L ☐ R 8"Wx11"D \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799. HCPCS: K0040+E0951
☐ P ☐ L ☐ R 6"Wx11"D \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797. HCPCS: K0040+E0951	

15.1c. Heavy Duty, Elevating Legrests

Select a configuration, a legrest, extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Heavy Duty, Elevating Legrest Selection

☐ P ☐ L ☐ R Heavy Duty, Elevating Legrest(s) with Calf Pad(s) \$171 ea. Pair part: FRMASMB6702. Left: FRMASMB8444. Right: FRMASMB8445. HCPCS: E0990

2. Lower Extension

☐ P ☐ L ☐ R 11"-14" Length Part: FRMASMB4755. Left: RIG138026. Right: RIG138030 For reinforced footplates only.	☐ P ☐ L ☐ R 14"-17" Length Part: FRMASMB4756. Left: RIG138033. Right: RIG138035 For reinforced footplates only.
☐ P ☐ L ☐ R 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103. For angle adjustable footplates only.	☐ P ☐ L ☐ R 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106 For angle adjustable footplates only.

3. Reinforced Footplates(s)

☐ P ☐ L ☐ R 7.25"Wx6"D Extruded Footplate Standard Part pair: FRMASMB4748. Left: RIG138055. Right: RIG138056	☐ P ☐ L ☐ R 8.5"Wx6"D Extruded Footplate Standard Part pair: FRMASMB4747. Left: RIG138060. Right: RIG138062
---	--

Legrests section continued on following page.

15. LEGRESTS CONTINUED**4. Angle Adjustable Footplates(s) with Heel Loops**

P L R 6"Wx8"D \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.
 HCPCS: K0040+E0951

P L R 8"Wx11"D \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799.
 HCPCS: K0040+E0951
 May cause interference with seat widths less than 20"

P L R 6"Wx11"D \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.
 HCPCS: K0040+E0951

15.2. Legrest Accessories

Wheel Bumper \$37
 Part: XRF000100

15.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

8.5" L (#3681)	9.5" L (#3683)	10.5" L (#3685)	11.5" L (#3687)	12.5" L (#3689)	13.5" L (#3691)	14.5" L (#3693)	15.5" L (#3695)	16.5" L (#3697)	17.5" L (#3699)	18.5" L (#3701)
9" L (#3682)	10" L (#3684)	11" L (#3686)	12" L (#3688)	13" L (#3690)	14" L (#3692)	15" L (#3694)	16" L (#3696)	17" L (#3698)	18" L (#3700)	19" L (#3702)

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	17.5" L (#3491)
9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)	18" L (#3492)

16. POWER POSITIONING ELECTRONICS

REQUIRED.

16.1a. NE, Tilt Only through Toggle

Single Actuator Control, Power Tilt Through Toggle No Charge

Part: ELEASMB6395

Select a side to mount toggle. Only available with NE 4-Key joystick.

Left Side Toggle

Part: ELEASMB5156

Right Side Toggle

Part: ELEASMB5157

16.1b. Single Function NE+/Q-Logic 2 through Joystick

Single Actuator Control Through NE+ Joystick or Q-Logic 2 Controller and Interface (AM1) \$2,163

Part: ELEASMB7218. HCPCS: E2310

16.1c. Two Function NE+ through Joystick

Two Actuator Control Through NE+ Joystick \$3,040

Part: ELE122740. HCPCS: E2311

16.1d. Q-Logic 2 Controls through Joystick

Single Actuator Through Q-Logic 2 \$3,040

Controller and Interface (AAM)

Part: ELEASMB7217. HCPCS: E2310

Two Actuator Through Q-Logic 2 \$3,040

Controller and Interface (AAM)

Part: ELE122739. HCPCS: E2311

Q-Logic 2 Optional Switch Control

Not available with AM1 electronics. Switch will be mounted to same side of joystick.

4 Way Toggle

Part: ELEASMB5925

4 Way Button Toggle

Part: ELEASMB5926

Quad Push Button

Part: ELEASMB5924

**17. PROVIDER INSTALLED ACCESSORIES**

Optional. Accessories may need to be installed by the provider.

Oxygen Holder \$210

Part: ACC125006 HCPCS: E2208

Hydration System \$237

Part: ACC125010

Deluxe Cell Phone/Media Holder \$205

Not available with desk length armrests.

Part: ACC7120006

Cup Holder \$21

Part: ACCASMB2634

Quantum Backpack \$168

Part: ACC1905853

Grip Lap Board (20005) \$38.95

Part: ACC152992

Grip Drink Holder (20085) \$15.50

Part: ACC152993

Accessory Bar \$132

Part: XRF000128

Large, Black, Weather Cover \$202

Part: ACCCOVR1021

Q-Logic EX/NE Programmer \$573

Part: ELEASMB6976

Q-Logic PC Programming Disk with Cable \$445

Part: ELEASMB5216

Personal Item Pouch \$84

Part: ACC2101017

18. BRANCH PATCHES

Optional.

Air Force Patch\$21
Part: ACC140240

Army Patch\$21
Part: ACC140241

Coast Guard Patch\$21
Part: ACC140239

Marines Patch\$21
Part: ACC140243

Navy Patch\$21
Part: ACC140244

Purple Heart Patch\$21
Part: ACC140238

**19. POSITIONING COMPONENTS**Optional. More options are available on the [Positioning Components Order Form](#).**19.1. Stealth Thoracic Laterals**

Laterals are not available with the omit back option.

1. Thoracic Laterals Bracket

Swing-away, Rail Mounts, pair (TWBLTLC)\$535
Part: INDPART2896. HCPCS: E1028

2. Thoracic Laterals Pads

Curved Pads, 3"x5", pair (TWBL-3X5) \$238
Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6) \$238
Part: INDPART2664. HCPCS: E0956

Flat Pads, 3"x5", pair (TWBLP3X5) \$238
Part: INDPART2676. HCPCS: E0956

Flat Pads, 4"x6", pair (TWBLP4X6) \$238
Part: INDPART2675. HCPCS: E0956

3. Thoracic Laterals Installation

Install Laterals
Part: PTOINDV3499

Do Not Install
Part: PTOINDV3500

19.2. Stealth Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$146.25 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$146.25 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103) \$146.25 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$146.25 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960

19.3. Stealth Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor (IFCS4113) .. \$92.70 ea.
Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor (IFCS4114) .. \$902.70ea.
Part: ST-IFCS4114-Q. HCPCS: E0960

19.4. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Stealth, Swing-Away, Headrest Mount Hardware

Facial Lateral for Comfort Series\$303.08 ea.
Right Side (CPS9-R-HW)
Part: INDPART2647. HCPCS: E1028

Facial Lateral for Comfort Series\$303.08 ea.
Left Side (CPS9-L-HW)
Part: INDPART2648. HCPCS: E1028

Facial Laterals for Comfort Series \$358.18 ea.
Left and Right Side (CPS9-B-HW)
Part: INDPART2649. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Right\$303.08
Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left\$303.08
Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$468.39
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$468.39
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$413.29
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028

2. Stealth, Swing-Away, Headrest Mount, Facial Lateral Pads and Accessories


Egg Switch Mount (ESM) (A) \$54.08 ea.
Part: INDPART2815

Figure 8 Bracket with 13" Rod \$275.53 ea.
and Egg Switch Mount (B)
Part: ST-F8675-13-Q

Facial Pad, 4.5"x2" (900M-Q) \$55.11 ea.
Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2" (920M-Q) \$55.11 ea.
Part: ST-920M-Q. HCPCS: K0108

Swing-Away Joystick Mount, \$262.50 ea.
18" Rod (C)
Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest

**19.5. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket\$189.26 ea.
Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket\$189.26 ea.
Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$54.08 ea.
Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)\$54.08 ea.
Part: POS124067. HCPCS: K0108

Positioning components section continued on following page.

19. POSITIONING COMPONENTS CONTINUED**19.6. Stealth Foot Boxes**

Stealth Small Foot Boxes, 5"Wx8"Dx7"H.....\$513.71
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H\$513.71
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H.....\$513.71
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

19.7. Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D \$279
(SPLE111) (pair)
Part: INDPART2783. HCPCS: E0951+E0952

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D..... \$279
(SPLE112) (pair)
Part: INDPART2784. HCPCS: E0951+E0952

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D \$279
(SPLE114) (pair)
Part: INDPART2785. HCPCS: E0951+E0952

20. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Select a type of specialty controls below or to refer to the specialty controls order form.

20.1. Enhanced Display

Q-Logic 2 Enhanced Display Kit.....\$1,980
Part: KIT140039

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 2 Enhanced Display Shroud

Required with Q-Logic 2 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard
Part: DGN137728

Red Corvette..... \$63
Part: DGN137729-PRD1

Back in Black..... \$63
Part: DGN137729-PBK1

Night Train..... \$63
Part: DGN137729-PBK3

Grey Street \$63
Part: DGN137729-PSL3

Deep Purple..... \$63
Part: DGN137729-PPR5

Pretty in Pink \$63
Part: DGN137729-PPK3

Midnight Blue \$63
Part: DGN137729-PBU10

Electric Blue \$63
Part: DGN137729-PBU9

Green Machine \$63
Part: DGN137729-PGN11



Enhanced Display

20.2. SCIM Module

SCIM Module Kit for Q-Logic 2.....\$683
Part: KIT138989

Available with Q-Logic 2 EX Joystick selection above. Allows use of joystick display for specialty devices and offers up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

20.3. Q-Logic 2 Stand Alone Joystick

Q-Logic 2 EX Stand Alone Joystick.....\$840
Part: ELEASMB5205

Up to 3 drive profiles and 1 seating profile when used alone. Use in conjunction with Enhanced Display for up to 5 profiles and auxiliary functions such as Bluetooth Mouse and IR. Includes charger port mounting bracket and multiplier box harness.



EX Stand Alone Joystick

20.4. Head Array

i-Drive 4.0 Bluetooth Tri-array Series,..... \$4,546.42
3 Switch with Egg Switch
Part: ST-IDHBT200-U-Q. HCPCS: E2330

Adult, Curved, i-Drive Tri-Array with \$4,546.42
Mounting Hardware
Part: ST-IDHBT400-U-Q. HCPCS: E2330

Junior, Curved, i-Drive Tri-Array with \$4,546.42
Mounting Hardware
Part: ST-IDHBT450-U-Q. HCPCS: E2330

Adult, i-Drive Pro Series Head Array \$4,546.42
(IDHBT300-1)

Part: ST-IDHBT300-1-Q. HCPCS: E2330
Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

20.5. Sip and Puff Drive**1. Sip and Puff Module**

Sip and Puff Module for Q-Logic 2.....\$2,185
Part: KITASMB1697. HCPCS: E2325
Includes module and bus cable.



Sip & Puff Module

Alternative controls and electronics section continued on following page.

20. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**2. Sip and Puff Straw and Mount**

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

Stealth Sip and Puff for Comfort Series,\$303.08
Left Mount (SP625M-CL)

Part: ST-SP625M-CL-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Comfort Series,\$303.08
Right Mount (SP622M-CR)

Part: ST-SP625M-CR-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive,\$303.08
Head Array (SP625-U)

Part: ST-SP625-U-Q. E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.



Stealth Sip & Puff

3. Sip & Puff Accessories

Sip & Puff Straws, 12 pack (SP625-S12)..... \$33.06
Part: ST-SP625-S12-Q

Sip & Puff Bus Cable\$447
Part: DWR1111H046



Egg Switch

20.6. Egg Switches

Red Egg Switch (ESRED)..... \$83.28
Part: SWTMCRO1077

Blue Egg Switch (ESBLU) \$83.28
Part: SWTMCRO1079

Black Egg Switch (ESBLK) \$83.28
Part: SWTMCRO1081

Green Egg Switch (ESGRN) \$83.28
Part: SWTMCRO1078

Yellow Egg Switch (ESYEL) \$83.28
Part: SWTMCRO1080

20.7. Switch Adapter

2 to 1 Switch Adapter \$15

Part: DWR1361H037

Connects 2 single switches to a stereo jack to act as a dual switch.