



Quantum Rehab®


401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® HD Order Form with TRU Balance® 3 Positioning

Account Number: _____ Date: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
PO Number: _____
Quote Name/Marked for: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____

Q6 edge®
HD



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 450 lb. weight capacity with standard 4.7 mph programmable motors, 8 amp off board charger, and unoccupied transit loops. Dimensions of base without legrests: 26"Wx36.75"L.

Quantum Q6 Edge HD 3HD-SS \$9,657

Part: Q6 EDGEHD 3HD-SS. HCPCS: K0850

Solid seat pan standard with no power option.

Quantum Q6 Edge HD 3SPHD-SS \$9,318

Part: Q6 EDGEHD 3SPHD-SS. HCPCS: K0858

Single power option.

Quantum Q6 Edge HD 3MPHD-SS \$9,115

Part: Q6 EDGEHD 3MPHD-SS. HCPCS: K0862

Multiple power option.

Account Number: _____ Date: _____ Quote Name/Marked for: _____

2. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.


 Red Corvette Part: DGN169178-PRD1	 Night Train Part: DGN169178-PBK3	 Deep Purple Part: DGN169178-PPR5	 Green Machine Part: DGN169178-PGN11
 Electric Blue Part: DGN169178-PBU9	 Pretty in Pink Part: DGN169178-PPK3	 Grey Street Part: DGN169178-PSL3	
 Back in Black Part: DGN169178-PBK1	 Evolution Orange Part: DGN169178-POR7	 Midnight Blue Part: DGN169178-PBU10	

3. FENDER LIGHTS

Optional. Fender lights operate by pressing included single switch button. Fender lighting is not available with NE 4-Key Joystick or Single Actuator Control through Toggle Box.

Switch Fender Lights without Lift System..... \$263
Part: XRF000177

Switch Fender Lights with Lift System No Charge
Part: DGN154117



4. QUANTUM BACKUP CAMERA

Optional. Requires Q-Logic electronics. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

Available for Q-Logic 3 electronics only. NE/NE+ coming soon.

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket..... \$21
Part: MEC1910373

Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$27
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123

Backup Camera NE/NE+ Joystick Mount Bracket \$105
Part: MEC2106132



5. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

5.1. Drive Wheels

3" x 14" Run Flat Tires..... Standard
Part: PTOINDV3510

4" x 14" Flat Free, Knobby, Extended Traction Tires..... \$290 ea.
Part: WHL121698 x2

3" x 14" Flat Free Tires..... No Charge
Part: WHL124116 x2

4" x 14" Pneumatic, Knobby, Extended Traction Tires \$290 ea.
Part: WHL121703 x2

3" x 14" Pneumatic Tires..... No Charge
Part: WHL124117 x2

5.2. Drive Wheel Accent Color

Optional. Not available with knobby tires.

Back in Black Accent Color..... \$63
Part: XRF000179

Pretty in Pink Accent Color \$63
Part: XRF000185

Midnight Blue Accent Color \$63
Part: XRF000187

Red Corvette Accent Color \$63
Part: XRF000180

Evolution Orange Accent Color \$63
Part: XRF000186

Green Machine Accent Color \$63
Part: XRF000190

Electric Blue Accent Color \$63
Part: XRF000188

Deep Purple Accent Color..... \$63
Part: XRF000184

Night Train Accent Color \$63
Part: XRF000181

Grey Street Accent Color..... \$63
Part: XRF000182

6. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

Group 24, 70 Amp Gel Cell, Interceptor 2 Battery \$399 ea.
Part: BAT2100315. HCPCS: E2363

Group 24, 70 Amp Gel Cell, MK Battery \$541 ea.
Part: BATGEL1009. HCPCS: E2363

7. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Forms](#).

7.1a. Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick.....Standard
Part: CTLDC1554

Available with Static Seat or Tilt Only through Toggle.

NE+ 6-Key Joystick No Charge
Part: CTLDC1560

Control of one actuator through Joystick on single power base model (3SPHD-SS) only.

NE+ 6-Key Joystick \$210
Part: XRF000001

Control up to 2 actuators through joystick.



NE+ 6 Key Joystick

7.1b. Expandable Controls

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module\$730
Part: XRF000002. HCPCS: E2377

Q-Logic 3 Expandable Controller Harness.....\$458
Part: XRF000264. HCPCS: E2313

2. Q-Logic 3 Joystick

Q-Logic 3e Joystick\$735
Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.

Q-Logic 3 EX Joystick.....\$840
Part: CTL167724

Q-Logic 3 EX Joystick for Full Lights\$840
Part: CTL167725

No Joystick, Setup for Q-Logic

Part: XRF000010

Requires selection of Enhanced Display.



Q-Logic 3 Joysticks

7.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Reference colors shown on page 2.

Black, Soft Touch Standard
Part: DGN124307

Red Corvette..... \$63
Part: DGN124336-PRD1

Electric Blue \$63
Part: DGN124336-PBU9

Back in Black..... \$63
Part: DGN124336-PBK1

Night Train..... \$63
Part: DGN124336-PBK3

Pretty in Pink \$63
Part: DGN124336-PPK3

Evolution Orange..... \$63
Part: DGN124336-POR7

Deep Purple..... \$63
Part: DGN124336-PPR5

Grey Street \$63
Part: DGN124336-PSL3

Midnight Blue \$63
Part: DGN124336-PBU10

Green Machine \$63
Part: DGN124336-PGN11

8. JOYSTICK HANDLES

Optional. The NE/NE Plus joysticks are standard with a carrot handle (See A in image). The standard Q-Logic 3 joystick knob can be seen in the image in section 7.1b. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323

Chin Cup (C) \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323

Medium Ball Handle (E)..... \$69
Part: ACCASMB1859. HCPCS: E2323

Soft Ball Handle (F)..... \$80
Part: ACCASMB1856. HCPCS: E2323

Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323

Bodypoint 3" Goal Post Handle\$104
Part: PLSKNCP1057. HCPCS: E2323

Bodypoint 4" Goal Post Handle\$104
Part: PLSKNCP1058. HCPCS: E2323

Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308



Joystick Handles

9. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount and side to be used with a joystick selected in previous section.

9.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket for Q-Logic 3Standard
Part: FRMASMB10617

Fixed, In-line, Bracket for NE/NE+Standard
Part: FRMASMB12744

9.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1028

Swing-away, NE/NE+..... \$319
Part: ACC1908437. HCPCS: E1028


9.1c. Basic Swing Away Joystick Mount

Basic Swing Away, NE/NE+, Left Side\$245
Part: ACC2201047. HCPCS: E1028


Basic Swing Away, NE/NE+, Right Side.....\$245
Part: ACC2201048. HCPCS: E1028

Basic Swing Away, Q-Logic, Left Side\$245
Part: ACC2201041. HCPCS: E1028

Basic Swing Away, Q-Logic, Right Side\$245
Part: ACC2201042. HCPCS: E1028



Swing Away Joystick



Retract4

Joystick mounting brackets section continued on following page.

9. JOYSTICK MOUNTING BRACKETS CONTINUED**9.1d. Specialty Joystick Mounts**

Retract4 Mount..... \$368
 Part: XRF000003. HCPCS: E1028
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount..... \$126
 Part: XRF000004
 Allows height adjustment of joystick.

**9.2. Joystick Mount Location**

Left Side Mount
 Part: PTOINDV1225

Right Side Mount
 Part: PTOINDV1226

10. TRU BALANCE® 3 SEATING - QUANTUM INTELLIGENT SYSTEMS™

REQUIRED. Maximum seat size is 32"W x 25"D.

10.1. TRU Balance 3 Configuration

TRU Balance 3 Static Seat..... \$2,410
 Part: PTOINDV3442
 20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU Balance 3 Heavy Duty, Power Recline \$9,318
 Part: PTOINDV3645. HCPCS: E1004
 175° of power recline.

TRU Balance 3, Heavy Duty, Power Tilt..... \$8,668
 Part: PTOINDV3643. HCPCS: E1002
 50° Center of Gravity (CG) power tilt seating system.

TRU Balance 3 Heavy Duty, Power Tilt & Recline \$13,675
 Part: PTOINDV3644. HCPCS: E1007
 50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

10.2. 10" Power Lift

Optional. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

TRU Balance 3, Heavy Duty, 10" Power Lift with iLevel® and Fender Lights \$3,725
 Part: PTOINDV3797. HCPCS: E2298/E2300
 400 lb weight capacity. 24"W x 24"D maximum seat size. Includes fender lights found in section 5.

**10.3. Seat Width**

REQUIRED. Maximum non-lift seat width is 32"W. *Widths over 24" are not available with 10" power lift. Part number is SETWDTH plus number following # below.

20" W (#1010)	22" W (#1012)	24" W (#1014)	26" W* (#1054)	28" W* (#1055)	30" W* (#1059)	32" W* (#1061)
21" W (#1011)	23" W (#1013)	25" W* (#1056)	27" W* (#1057)	29" W* (#1058)	31" W* (#1060)	

10.4. Seat Width Package

REQUIRED.

20"-24" Width Package
 Part: XRF000158

24"-28" Width Package
 Part: XRF000192

28"-32" Width Package
 Part: XRF000193

10.5. Seat Depth

REQUIRED. Maximum seat depth is 25". Part number is SETDPTH plus number following # below.

20" D (#1010)	21" D (#1011)	22" D (#1012)	23" D (#1013)	24" D (#1014)	25" D (#1060)
------------------	------------------	------------------	------------------	------------------	------------------

10.6. Seat Depth Package

REQUIRED.

16"-20" Depth Package
 Part: XRF000013
 Minimum Seat Depth is 20"

20"-24" Depth Package
 Part: XRF000159
 Use with 25" Depth

Low	Med	High
Static Seat, Recline		
16.5"	17.5"	18.5"
Tilt, Tilt & Recline, Lift		
18.5"	19.5"	20.5"
TRU Balance 3 Seat to Floor Heights		

10.7. Seat to Floor Height

REQUIRED.


Lowest Seat to Floor
 Part: PTOINDV3558

Middle Seat to Floor
 Part: PTOINDV3559

Highest Seat to Floor
 Part: PTOINDV3560

10.8. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the **respective complete Stealth Cushion Order Forms**. Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.



TRU Balance 3 Seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED**1. Cushion Type**

TRU-Comfort 2 SPP
Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPP
Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPP
Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPP
Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPP
Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere P
Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SP
Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SP
Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity G
Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller

Medium Wedge Positioning Kit.....\$72
Part: POP1909006-M
14"W and smaller

Large Wedge Positioning Kit..... \$105
Part: POP1909006-L
14"W and smaller

Small Wedge Positioning Kit\$39
Part: POP1712340-S
15"W and larger

Medium Wedge Positioning Kit.....\$72
Part: POP1712340-M
15"W and larger

Large Wedge Positioning Kit..... \$105
Part: POP1712340-L
15"W and larger

3. Cushion Top Foam

Visco Top Foam
Reference part: XRF000380

No Top Foam
Reference part: XRF000382

FRG Top Foam
Reference part: XRF000381

4. Cushion Liner

Stealth Urethane Liner
Reference part: XRF000383

No Liner
Reference part: XRF000385

Stealth Polyurethane Liner
Reference part: XRF000384

5. Cushion Cover

CoolCore Cover
Reference part: XRF000386

Stealth Tek Cover
Reference part: XRF000388

Mesh Cover
Reference part: XRF000387

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389
MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390

Additional Stealth Polyurethane Liner
Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392

Additional Stealth Tek Cover
Reference Part: XRF000394

Additional Mesh Cover
Reference Part: XRF000393

10.9. Back Type Selection**REQUIRED.** Match choice in this section with the options in the following back cushion configuration section.

TRU Balance 3 Back Assembly (Static Back)
Part: XRF000194
Select a back cushion configuration in the subsequent section.

Back Cane Option
Part: XRF000015
Available only with TRU Balance 3 Static Seat and Power Tilt selections. Select back cane type, back cane height, and back cushion configuration.

Flip Forward Back..... \$184
Part: XRF000220

Available with Static Seat and Power Tilt (1 actuator). Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

	None	Visco/V Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

TRU Balance 3 Seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED**10.10. Back Cane Selection (Back Cane Option Only)****REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type**

Angle Adjustable No Charge Part: MEC140266	1" Endo..... \$74 Part: XRF000269	2" Endo \$74 Part: XRF000271
---	--------------------------------------	---------------------------------------

2. Back Cane Height

16" High Cane Part: FRMASMB15834	20" High Cane Part: FRMASMB15838	24" High Cane Part: FRMASMB15842
17" High Cane Part: FRMASMB15835	21" High Cane Part: FRMASMB15839	25" High Cane Part: FRMASMB15843

10.11. Back Cushion Configuration**REQUIRED.** Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1a. Sport Back/Omit/After-Market Interface**

Sport Back No Charge Part: XRF000017 Complete following width package and size.	After-market Back Interface \$210 Part: XRF000018 Select a back width within a range of 16"-24"W. Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.
Omit Back Pan/Cushion No Charge Part: PTOINDV1287 Select a following width package.	

1b. TRU-Comfort Back

TRU-Comfort Back (Static/Flip Forward), 14" to 20" Width..... \$744 Part: TC4V10W21W. HCPCS: E2620 Complete following width package and size.	TRU-Comfort Back (Static/Flip Forward), 22" Width or Greater.... \$855 Part: TC4V22W32W. HCPCS: E2621 Complete following width package and size.
TRU-Comfort Back, Back Cane Option, 16" to 20" Width..... \$744 Part: TC4V10W21WBC. HCPCS: E2620 Complete following width package and size.	TRU-Comfort Back, Back Cane Option, 22" Width or Greater.. \$855 Part: TC4V22W32WBC. HCPCS: E2621 Complete following width package and size.

1c. Q-Back 2.0, AL Series, Active Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline..... \$568 Part: QB2AL18AS. HCPCS: E2613	Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$568 Part: QB2AL20AS. HCPCS: E2613
Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$568 Part: QB2AL18AT. HCPCS: E2613	Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$568 Part: QB2AL20AT. HCPCS: E2613
Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline..... \$568 Part: QB2AL18AX. HCPCS: E2613	Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline..... \$568 Part: QB2AL20AX. HCPCS: E2613
Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline..... \$703 Part: QB2AL18AXX. HCPCS: E2613	Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline..... \$703 Part: QB2AL20AXX. HCPCS: E2613

1d. Q-Back 2.0, AL Series, Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL18DS. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL18DT. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL20DT. HCPCS: E2620
Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL18DX. HCPCS: E2620	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL20DX. HCPCS: E2620
Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18DXX. HCPCS: E2620	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20DXX. HCPCS: E2620

1e. Q-Back 2.0, AL Series, Extra Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL18ES. HCPCS: E2620	Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL20ES. HCPCS: E2620
Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$703 Part: QB2AL18ET. HCPCS: E2620	Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703 Part: QB2AL20ET. HCPCS: E2620
Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL18EX. HCPCS: E2620	Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL20EX. HCPCS: E2620
Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18EXX. HCPCS: E2620	Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20EXX. HCPCS: E2620

1f. Q-Back 2.0, AFT Series, Active Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT16AS. HCPCS: E2613	AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT18AS. HCPCS: E2613
--	--

TRU Balance 3 Seating section continued on following page.

10. TRU BALANCE 3 SEATING CONTINUED**1g. Q-Back 2.0, AFT Series, Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$643
Part: QB2AFT16DS. HCPCS: E2620AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$643
Part: QB2AFT16DT. HCPCS: E2620AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AFT18DS. HCPCS: E2620AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$643
Part: QB2AFT18DT. HCPCS: E2620AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AFT20DT. HCPCS: E2620**1h. Q-Back 2.0, AFT Series, Extra Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AFT16EXX. HCPCS: E2620**10.12. TRU Balance 3 Back Shroud****REQUIRED.** The back shroud is a feature of the TRU Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown in section 2.Red Corvette
Part: XRF000019Night Train
Part: XRF000022Deep Purple
Part: XRF000025Evolution Orange
Part: XRF000028Electric Blue
Part: XRF000030Back in Black
Part: XRF000021Grey Street
Part: XRF000023Pretty in Pink
Part: XRF000027Midnight Blue
Part: XRF000029Green Machine
Part: XRF000033**10.13. Lateral Positioning Wedge****Optional.**Lateral Positioning Wedge, pair \$16 ea.
Part: SET123465

Available with Sport and TRU-Comfort 2 backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

10.14. Back Width Package**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.Match Seat Width Selections
Part: PTOINDV346820"W-24"W Package
Part: XRF00004028"W-32"W Package
Part: XRF00019616"W-20"W Package
Part: XRF00003824"W-28"W Package
Part: XRF000195**10.15. Back Width****REQUIRED** except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.

18" W (#3566)	19" W (#3567)	20" W (#3568)	21" W (#3569)	22" W (#3570)	24" W (#3571)	26" W (#3572)	28" W (#3573)	30" W (#3574)	32" W (#3575)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

10.16. Overall Back Height**Static Seat and Power Tilt available height range:** 16"-25"H. **Odd heights include 4" gap and even heights include 3" gap** between seat pan and bottom of back cushion. **Power Recline available height range:** 19"-28"H. **Odd heights include 6" gap and even heights include 7" gap.** See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

16" Overall Height Part: SETHGHT1029	23" Overall Height Part: SETHGHT1036
17" Overall Height Part: SETHGHT1030	24" Overall Height Part: SETHGHT1037
18" Overall Height Part: SETHGHT1031	25" Overall Height Part: SETHGHT1038
19" Overall Height Part: SETHGHT1032	26" Overall Height Part: SETHGHT1039
20" Overall Height Part: SETHGHT1033	27" Overall Height Part: SETHGHT1040
21" Overall Height Part: SETHGHT1034	28" Overall Height Part: SETHGHT1041
22" Overall Height Part: SETHGHT1035	

	Tilt or Static Back Width										Recline Back Width									
	18"	19"	20"	21"	22"	24"	26"	28"	30"	32"	18"	19"	20"	21"	22"	24"	26"	28"	30"	32"
Overall Back Height	16"	ST	ST	ST	ST	S	ST	ST	S	S	S	S	-	-	-	-	-	-	-	-
	17"	ST	ST	ST	ST	S	ST	ST	S	S	S	S	-	-	-	-	-	-	-	-
	18"	ST	ST	ST	ST	S	ST	ST	ST	ST	S	S	-	-	-	-	-	-	-	-
	19"	ST	ST	ST	ST	S	ST	ST	ST	ST	S	S	ST	ST	ST	S	ST	S	S	S
	20"	ST	ST	ST	ST	S	ST	ST	ST	ST	S	S	ST	ST	ST	S	ST	S	S	S
	21"	ST	ST	ST	ST	S	ST	ST	ST	ST	S	S	ST	ST	ST	S	ST	ST	ST	ST
	22"	ST	ST	ST	ST	S	ST	ST	S	S	S	S	ST	ST	ST	S	ST	ST	ST	S
	23"	ST	ST	ST	ST	S	ST	ST	S	S	S	S	ST	ST	ST	S	ST	ST	ST	S
	24"	ST	S	ST	S	S	S	S	S	S	S	S	ST	ST	ST	S	ST	ST	ST	S
	25"	ST	S	ST	S	S	S	S	S	S	S	S	ST	ST	ST	S	ST	ST	S	S
	26"	-	-	-	-	-	-	-	-	-	ST	ST	ST	S	ST	ST	S	S	S	S
	27"	-	-	-	-	-	-	-	-	-	ST	S	ST	S	S	S	S	S	S	S
	28"	-	-	-	-	-	-	-	-	-	ST	S	ST	S	S	S	S	S	S	S

S Sport Back ST Sport or TRU-Comfort 2 Back

Sport and TRU-Comfort 2 Back Availability Chart

10.17. Power Tilt Privacy Flap**Optional.** A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.Privacy Flap \$48
Part: XRF000351

11. TRU BALANCE 3 ARMRESTS**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.**11.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Flip-back armrests include an armrest downpost on the Q6 Edge HD. Not available with waterfall armpads. Not available with Back Cane Option.

L R Flip-Back, 22"-26"Widths, 8"-16"H Armrests Standard
 Left part: ARM2004694. Right part: ARM2004693

L R Flip-Back, 26"-32"Widths, 8"-16"H Armrests Standard
 Left part: ARM2004696. Right part: ARM2004695

11.1b. Height Adjustable, Removable, Heavy Duty Armrests

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable, \$158
 Heavy Duty (Sold in pairs)
 Part: XRF000055. HCPCS w/ Static Seat: E0973

10"-14"H, Single Post, Height Adjustable, Removable, \$158
 Heavy Duty (Sold in pairs)
 Part: XRF000056. HCPCS w/ Static Seat: E0973

11.2. Armrest Height

Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

2. Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

11.3. Armrest Size**REQUIRED with armrest selection above.** Select a left (L) and right (R) armrest size. This also affects armpad selection.

L R 10" Length
 L part: XRF000057. R part: XRF000123

L R 14" Length
 L part: XRF000058. R part: XRF000124

L R 17" Length
 L part: XRF000059. R part: XRF000125
 Only available with Flip-Back Armrests
 and straight arm pads.

12. TRU BALANCE 3 ARMPADS**REQUIRED with armrest selection above except for Drop-In, T-Arm.** Please note that waterfall and gel armpads may cause interference with the seat back.**12.1a. Omit Armpad**Select if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection.....No Charge
 Part: XRF000060
 Not available with Flip-back, 8"-16" H armrests.

12.1b. Foam Armpads

L R Straight, Armpad Standard
 Left part: XRF000061. Right part: XRF000126

L R Waterfall, Armpad No Charge
 Left part: XRF000062. Right part: XRF000127
 Not available with Flip-back, 8"-16" armrests.

12.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right.....\$195 Part: ARM1806807	Falcon Gel, Desk Length Left, Full Length Right.....\$171 Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right.....\$171 Part: ARM1806813	Falcon Gel, Desk Length Left and Right\$147 Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right.....\$311.22 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right.....\$287.93 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right.....\$287.93 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$265.02 Part: ARM136437

13. TRU BALANCE 3 JOYSTICK BRACKET RECEIVER**REQUIRED.**

Standard Fixed Receiver Standard
 Part: XRF000063

Stealth Height Adjustable Receiver..... \$236
 Part: XRF000064. HCPCS: K0108

14. TRU BALANCE 3 HEADRESTS**Optional.** Select a headrest pad and bracket.**14.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q. HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q. HCPCS: E0955**2. Headrest Brackets**Comfort Plus Headrest Mounting Bracket, UniLink\$311
Part: ST-UHB2D33200-Q. HCPCS: E1034Multi-Axis Removable Headrest Bracket\$174
Part: ST-HMO475-17M-Q. HCPCS: E1034**14.2. Q Cam Lock****Optional.**Q Cam Lock\$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

15. LAP BELTS**REQUIRED.** Select a lap belt below.**15.1a. Standard Lap Belts**38" Lap Belt
Part: ACCBELT100650" Lap Belt
Part: ACCBELT100060" Lap Belt
Part: ACCBELT100270" Lap Belt
Part: ACCBELT100380" Lap Belt
Part: ACCBELT100490" Lap Belt
Part: ACCBELT1001**15.1b. Optional Lap Belts**Stealth Retractable Lap Belt (IFRB5350)\$181
Part: ST-IFRB5350Stealth Positioning, 60" Single Pull, Push Button, Medium\$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-QStealth Positioning, 60" Single Pull, Push Button, Large\$50
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-QStealth Positioning, 60" Single Pull, Push Button,\$50
Extra Large, Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q**16. BASE MOUNTED FOOT PLATFORM****REQUIRED.** Select a foot platform option. A base mounted foot platform is not available with center mount foot platforms. May interfere with seat mounted legrests.Base Mounted Foot Platform Standard
Part: FRMASMB16077Omit Base Mounted Foot Platform No Charge
Part: PTOINDV3674Base Mounted, Bell Shaped Foot Platform No Charge
Part: RIGASMB7120041**17. LEGRESTS****Optional.** Please note the set lower leg length may vary based on your seat type and seat to floor height selections. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**17.1a. Basic Center Mount Foot Platform****1. Basic Center Mount Foot Platform Assembly**Basic Center Mount No Charge
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Extension and Angle Adjustable FootplatesShort Extension, 5"x9" Footplates \$180
Part: XRF000312. HCPCS: K0040Medium Extension, 5"x9" Footplates \$180
Part: XRF000313. HCPCS: K0040Long Extension, 5"x9" Footplates \$180
Part: XRF000314. HCPCS: K0040Medium Lower Extension, 6.25"W Tapered 5"W x 12"D \$180
Part: XRF000349. HCPCS: K0040Long Lower Extension, 6.25"W Tapered 5"W x 12"D \$180
Part: XRF000350. HCPCS: K0040**3. Basic Center Mount Calf Pads**Small Calf Pads, Basic Center Mount\$158
Part: RIG2104369**4. Basic Center Mount Additional Accessory**Basic Center Mount Footplate Connector \$16
Part: RIG2104240

Legrests section continued on following page.

17. LEGRESTS CONTINUED**17.1b. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform AssemblyPower Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040Manually Adjustable Center Mount Foot Platform \$263
Part: XRF000155. HCPCS: K0040**2. Lower Extension and Angle Adjustable Footplates**Medium Extension, 5"x6" Footplates
Part: RIG150206Long Extension, 5"x6" Footplates
Part: RIG151505Medium Extension, 5"x8" Footplates
Part: RIG2404016Long Extension, 5"x8" Footplates
Part: RIG2404017Medium Extension, 5"x11" Footplates
Part: RIG2404018Long Extension, 5"x11" Footplates
Part: RIG2404019

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Medium	Long
10"-16"	13.5"-19.5"
Manually Adjustable Extension Range	

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019Small Calf Pads, Adjustable Center Mount Only, \$158
6"Wx6"H (pair)
Part: XRF000156Large Calf Pads, Adjustable Center Mount Only, \$158
8" taper to 6"Wx8"H (pair)
Part: XRF000157**3b. Center Mount Calf Panel**

Calf pads are standard with a power articulating foot platform.

Small Calf Panel for AFP, 18"Wx8"H \$184
Part: FRMASMB7110007Small Calf Panel for Adjustable Center Mount, 18"Wx8"H \$184
Part: FRMASMB7110012Large Calf Panel for AFP, 24"Wx8"H \$237
Part: RIG123824
Only available on widths of 20" or greater.Large Calf Panel for Adjustable Center Mount, 24"Wx8"H \$237
Part: RIG123825
Only available on widths of 20" or greater.**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
Part: RIG152519
Only available with the 5x6 footplates.Outside Edge \$69
Part: RIG2404921
Available with 5"x8" footplates.Width Extension for 8" \$69
Part: RIG2404919Front and Rear Cap \$16
Part: RIG2404923Width Extension for 11" \$69
Part: RIG2404920

Available with 5"x8" footplates.

5. Center Mount Additional AccessoriesSmall Foam Calf Pad Cover (pair) \$67
Part: ST-CSCS0775-QLarge Foam Panel Cover (pair) \$282
Part: ST-CSCP0927-QLarge Foam Calf Pad Cover (pair) \$83
Part: ST-CSCS0985-QFootplate Connector \$16
Part: RIGASMB7120023Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-Q**6. Stealth Gel Footplate Covers**5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K01085"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K01085"x11" Tapered Gel Footplate Covers \$270
Part: ST-AFPG5011T-Q. HCPCS: K01085"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-AFPG5011E-Q. HCPCS: K0108**17.1c. Heavy Duty, Drop-In, Swing-Away Legrests**

Select a configuration, a legrest, extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Drop-In Legrest Selection

☐ P	☐ L	☐ R	70° HD Drop-in for 8.5"-12" Lower Extension No Charge Pair part: FRMASMB4762. Left: RIG137706. Right: RIG137707
----------------------------	----------------------------	----------------------------	--

☐ P	☐ L	☐ R	70° HD Drop-in 13"-18" or 16"-21" Lower No Charge Extension Pair part: FRMASMB4853. Left: RIG137708. Right: RIG137709
----------------------------	----------------------------	----------------------------	---

2. Lower Extension

☐ P	☐ L	☐ R	8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103.
----------------------------	----------------------------	----------------------------	--

☐ P	☐ L	☐ R	16"-21" HD Drop-Ins/18.5"-22.5" ELR, Ext Length Pair part: FRMASMB4754. L: FRMASMB4725. R: FRMASMB5393
----------------------------	----------------------------	----------------------------	---

☐ P	☐ L	☐ R	13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106
----------------------------	----------------------------	----------------------------	--

Legrests section continued on following page.

17. LEGRESTS CONTINUED**3. Angle Adjustable Footplate(s) with Heel Loops**

☐ ☐ ☐ 6"Wx8"D¹ \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx11"D¹ \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.
 HCPCS: K0040+E0951
 May cause interference on seat widths less than 18".

☐ ☐ ☐ 8"Wx11"D² \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799.
 HCPCS: K0040+E0951
 May cause interference on seat widths less than 20".

17.1d. Heavy Duty, Elevating Legrests

Select a configuration, a legrest, extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Heavy Duty, Elevating Legrest Selection

☐ ☐ ☐ Heavy Duty, Elevating Legrest(s) with Calf Pad(s) \$171 ea.
 Pair part: FRMASMB6702. Left: FRMASMB8444. Right: FRMASMB8445. HCPCS: E0990

2. Lower Extension

☐ ☐ ☐ 11"-14" Length
 Part: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced footplates only.

☐ ☐ ☐ 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103.
 For angle adjustable footplates only.

☐ ☐ ☐ 14"-17" Length
 Part: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced footplates only.

☐ ☐ ☐ 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable footplates only.

3. Reinforced Footplates(s)

☐ ☐ ☐ 7.25"Wx6"D Extruded Footplate Standard
 Part pair: FRMASMB4748. Left: RIG138055.
 Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Extruded Footplate Standard
 Part pair: FRMASMB4747. Left: RIG138060.
 Right: RIG138062

4. Angle Adjustable Footplates(s) with Heel Loops

☐ ☐ ☐ 6"Wx8"D \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx11"D \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.
 HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799.
 HCPCS: K0040+E0951
 May cause interference with seat widths less than 20"

17.2. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

8.5" L (#3681)	9.5" L (#3683)	10.5" L (#3685)	11.5" L (#3687)	12.5" L (#3689)	13.5" L (#3691)	14.5" L (#3693)	15.5" L (#3695)	16.5" L (#3697)	17.5" L (#3699)	18.5" L (#3701)
9" L (#3682)	10" L (#3684)	11" L (#3686)	12" L (#3688)	13" L (#3690)	14" L (#3692)	15" L (#3694)	16" L (#3696)	17" L (#3698)	18" L (#3700)	19" L (#3702)

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	17.5" L (#3491)
9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)	18" L (#3492)

18. TRU BALANCE 3 POWER POSITIONING ELECTRONICS

REQUIRED with power seating or legrest options.

18.1a. Tilt Only through Toggle

Single Actuator Control, Power Tilt Only, Through Toggle (NE 4-Key Joystick Only) No Charge
 Part: ELE1905907
 Select toggle side. May require extended lead time when selected with Fender Lights.

1. Toggle Mount Side Selection

Left Side Toggle
 Part: PTOINDV1839

Right Side Toggle
 Part: PTOINDV1838

18.1b. Actuator Control through Joystick**1. Actuator Control Selection**

Single Actuator Through Controller \$2,163
 Interface
 Part: XRF000335. HCPCS: E2310

Multiple Actuator Through Controller and \$3,040
 Interface
 Part: XRF000336. HCPCS: E2311

Single Actuator Through Controller and \$3,040
 Interface (AAM)
 Part: XRF000068. HCPCS: E2310



TRU Balance 3 power positioning electronics section continued on following page.

18. TRU BALANCE 3 POWER POSITIONING ELECTRONICS CONTINUED**2. Optional Switch Control**

Not available with AM1 electronics. 4-Way switches will be mounted to same side of joystick.

4-Function Push Button

Part: CTL123218

See following section for switch mounting.

4-Way Toggle

Part: ELEASMB5925

4-Way Button Toggle

Part: ELEASMB5926

3. User Controls Side Selection**Required with 4-Way Toggle and 4-Way Button Toggle.**

Mount User Controls on Left

Part: PTOINDV3996

Mount User Controls on Right

Part: PTOINDV3997

18.1c. Actuator Control with iAccess Switch

Single Actuator Control Through Flat Button\$2,547

iAccess, Controller, and Interface

Part: XRF000337. HCPCS: E2310

Single Actuator Control Through Toggle Button\$2,547

iAccess, Controller, and Interface

Part: XRF000338. HCPCS: E2310

Single Actuator Control Through Flat Button\$3,460

iAccess and Controller and Interface (AAM)

Part: XRF000070. HCPCS: E2310

Single Actuator Control Through Toggle Button\$3,460

iAccess and Controller and Interface (AAM)

Part: XRF000074. HCPCS: E2310

Multiple Actuator Control Through Flat Button\$3,460

iAccess, Controller, and Interface

Part: XRF000339. HCPCS: E2311

Multiple Actuator Control Through Toggle Button.....\$3,460

iAccess, Controller, and Interface

Part: XRF000340. HCPCS: E2311

18.2. 4-Function Push Button/iAccess Mount**REQUIRED** with 4-function push button and iAccess options. Under arm mount is not available with TRU Balance 3 HD armrests.

Joystick Mount No Charge

Part: MEC1907500

Attendant Mount \$26

Part: MEC1906564

Under Arm Mount, Flip-Back Armrest..... No Charge

Part: MEC1906531

Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$26

Part: MEC2001392

Under Arm Mount, Tube Style Armrest No Charge

Part: MEC1906530

Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$26

Part: MEC2001391

1. Under Arm Mount Side Selection

Left Side Under Arm Mount

Part: PTOINDV3678

Right Side Under Arm Mount

Part: PTOINDV3679

19. TRU BALANCE 3 ATTENDANT CONTROL**Optional.** *Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control \$525

Part: XRF000078. HCPCS: E2331

Select this option if no charge option criteria is not met.

Attendant Control with Tilt, Recline, Alternative No Charge


Controls and Q-Logic Display or SCIM Module

Part: XRF000079

20 FULL LIGHT PACKAGE**Optional.** Requires Q-Logic 3 Joystick for Full Light Package.

Seat Mounted, LED, Full Light Package\$1,045

Part: XRF000130

Seat Mounted LED
Lights (One Shown)

21. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be mounted by the provider.

TRU Balance 3 Cane/Crutch Holder \$105
 Part: ACC152515. HCPCS: E2207

Quantum Backpack \$168
 Part: ACC1905853

Grip Lap Board (20005)..... \$38.95
 Part: ACC152992

Grip Drink Holder (20085) \$15.50
 Part: ACC152993

TRU Balance 3 Hydration System (B) \$237
 Part: ACC125010

TRU Balance 3 Transfer Handle (D) \$105 ea.
 Part: MEC2205680. HCPCS: E0705
 Not available with Heavy Duty armrests.

TRU Balance 3 Cup Holder No Charge
 Part: ACC125003
 For use with TRU Balance 3 (track style) armrests.

Cup Holder (Tube-style Armrest)..... No Charge
 Part: ACCASMB2634
 For use with tube-style armrests.

TRU Balance 3 Personal Item Hook \$69
 Part: ACC125008
 Not available with back cane option.

TRU Balance 3 Reflector Kit..... \$53
 Part: ACC124670

TRU Balance 3 Backpack Holder (Pair) (A) \$69
 Part: ACC147774
 Not available with back cane option.

TRU Balance 3 Front Mounted \$69
 Cell Phone/ Media Holder, X-Grip (E)
 Part: ACC155835
 Not available with Heavy Duty armrests.

TRU Balance 3 Rear View Mirror..... \$69
 Part: ACC155833
 Not available with Heavy Duty armrests.

TRU Balance 3 Oxygen Holder \$210
 Part: ACC125006. HCPCS: E2208
 Includes accessory bar, a credit will be applied PTOIN-DV3769.

TRU Balance 3 Clothing Guard (Pair) \$237
 Part: ACC125007
 Not available with Heavy Duty armrests.

TRU Balance 3 Glove Box (C) \$158
 Part: ACC125011

Stealth Large Tablet Holder \$625.28
 Part: ACC141818
 Not available with Heavy Duty armrests.

Compact Vent Tray \$1,465
 Part: ACC133889. HCPCS: E1029

Push Handles (pair)..... \$195
 Part: ACC2005027
 Not available with oxygen holder, back cane option, or vent tray. Required with attendant manual recline.

Accessory Bar \$132
 Part: XRF000128
 Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder. Not available with back cane option.

Large, Black, Weather Cover \$202
 Part: ACCCOVR1021

Q-Logic EX/NE Programmer \$573
 Part: ELEASMB6976


Privacy Bag \$42
 Part: ACC2005723

Flip-Back Armrest, Side Track Accessory Bracket ... \$11
 Part: MEC1809178
 Only available with flip-back, 8"-16" H armrests. This adds a 47mm channel to the side of the armrest to allow mounting of accessories with t-nuts.

Personal Item Pouch \$84
 Part: ACC2101017

Depth Adjustable Vent Tray \$595
 Part: ACC2302500. HCPCS: E1029
 Includes accessory bar, a credit will be applied PTOIN-DV3769.

Remote Stop Switch \$1,295
 Part: XRF000379

**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
 QTY Part: ACC2407046

Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
 QTY Part: ACC2405496

22. BRANCH PATCHES**Optional.**

Air Force Patch \$21
 Part: ACC140240

Coast Guard Patch \$21
 Part: ACC140239

Navy Patch \$21
 Part: ACC140244

Army Patch \$21
 Part: ACC140241

Marines Patch \$21
 Part: ACC140243

Purple Heart Patch \$21
 Part: ACC140238

**23. POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#).**23.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Thoracic Lateral Bracket

Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664
 Part: ST-ULPK404108B00-Q. HCPCS: E1034
 Only available with Sport or TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
 Part: ST-ULPK404104B00-Q. HCPCS: E1034
 Only available with back cane and omit back configurations.

Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$332
 Part: ST-ULPK404103L00-Q. HCPCS: E1034
 Only available with Active Contour Q-Backs.

Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$332
 Part: ST-ULPK404103R00-Q. HCPCS: E1034
 Only available with Active Contour Q-Backs.

Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**2. Thoracic Lateral Pads**

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$250
 Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250
 Part: INDPART2664. HCPCS: E0956

3. Thoracic Lateral Installation

Install Laterals
 Part: PTOINDV3499

Do Not Install
 Part: PTOINDV3500

23.2. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea.
 Removable Mounting Hardware
 Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Left..... \$332 ea.
 Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^

UniLink, 7" Height Adjustable, \$276 ea.
 Removable Mounting Hardware
 Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Right..... \$332 ea.
 Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
 Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$281 ea.
 Feature (TWBADD-L & LPHW407)
 Part: POS140428. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea.
 Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact..... \$281 ea.
 Feature (TWBADD-R & LPHW407)
 Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

UniLink, 4" x 4" Pad \$106 ea.
 Part: ST-UTP544056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 12" Pad \$106 ea.
 Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 8" Pad..... \$106 ea.
 Part: ST-UTP548056100-Q. HCPCS: E0953/E0956

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea.
 Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$81 ea.
 Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

23.4. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea.
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$154 ea.
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$154 ea.
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$154 ea.
 Part: ST-IFPS2104-Q. HCPCS: E0960

23.5. Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor, (IFCS4113)..... \$97 ea.
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor, (IFCS4114)..... \$97 ea.
 Part: ST-IFCS4114-Q. HCPCS: E0960

23.6. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right..... \$318
 Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492
 2 Upper & 1 Lower Left
 Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left \$318
 Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$434
 2 Upper
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492
 2 Upper & 1 Lower Right
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

2. Facial Lateral Pads and Accessories


Egg Switch Mount (ESM) (A) \$57 ea.
 Part: INDPART2815

Swing-Away Joystick Mount, 18" Rod (C)..... \$276 ea.
 Part: ST-SUS18-Q

Facial Pad, 4.5"x2" (900M-Q) \$61 ea.
 Part: ST-900M-Q. HCPCS: K0108

Allows mounting of facial pads and attaches to the headrest

Facial Pad, 6"x2" (920M-Q)..... \$61 ea.
 Part: ST-920M-Q. HCPCS: K0108



Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**23.7. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support \$199 ea.
Bracket Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support \$199 ea.
Bracket Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
Part: ST-UETB4021FL-Q

Right, UniLink, Fixed, Elbow Support Bracket \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
Part: ST-UETB4021FR-Q

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$57 ea.
Part: POS124067. HCPCS: K0108

**23.8. Stealth Foot Boxes**

Stealth Small Foot Boxes, 5"Wx8"Dx7"H \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

23.9. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211) \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

24. LAP TRAY**24.1. Quantum Lap Tray**

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

Black Lap Tray, 11" W \$200
Part: MEC1905328-11B. HCPCS: E0950
For use with 12"-13" back width.

Black Lap Tray, 17" W \$200
Part: MEC1905328-17B. HCPCS: E0950
For use with 18"-19" back width.

Black Lap Tray, 13" W \$200
Part: MEC1905328-13B. HCPCS: E0950
For use with 14"-15" back width.

Black Lap Tray, 19" W \$200
Part: MEC1905328-19B. HCPCS: E0950
For use with 20"-22" back width.

Black Lap Tray, 15" W \$200
Part: MEC1905328-15B. HCPCS: E0950
For use with 16"-17" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$200
Part: MEC1905328-11C. HCPCS: E0950
For use with 12"-13" back width.

Clear Lap Tray, 17" W \$200
Part: MEC1905328-17C. HCPCS: E0950
For use with 18"-19" back width.

Clear Lap Tray, 13" W \$200
Part: MEC1905328-13C. HCPCS: E0950
For use with 14"-15" back width.

Clear Lap Tray, 19" W \$200
Part: MEC1905328-19C. HCPCS: E0950
For use with 20"-22" back width.

Clear Lap Tray, 15" W \$200
Part: MEC1905328-15C. HCPCS: E0950
For use with 16"-17" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray \$75
Part: MEC2104886

24.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
Part: ST-SSLXHD1000060-Q

Lap tray section continued on following page.

23. LAP TRAY CONTINUED

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	SM	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

25. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).

25.1. Enhanced Display

Q-Logic 3 Enhanced Display Kit.....\$1,980
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard
Part: DGN137728

Night Train.....\$63
Part: DGN137729-PBK3

Pretty in Pink\$63
Part: DGN137729-PPK3

Green Machine\$63
Part: DGN137729-PGN11

Red Corvette.....\$63
Part: DGN137729-PRD1

Grey Street.....\$63
Part: DGN137729-PSL3

Midnight Blue\$63
Part: DGN137729-PBU10

Back in Black.....\$63
Part: DGN137729-PBK1

Deep Purple.....\$63
Part: DGN137729-PPR5

Electric Blue\$63
Part: DGN137729-PBU9



Enhanced Display

25.2. SCIM Module

Optional.

SCIM Module Kit for Q-Logic 3.....\$525
Part: CTL168055

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

25.3. Q-Logic 3 SAJ

Q-Logic 3 EX SAJ.....\$840
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

25.4. Head Array

i-Drive 4.0 Bluetooth Tri-array Series,.....\$4,774
3 Switch with Egg Switch
Part: ST-IDHBT200-U-Q. HCPCS: E2330

Junior, Curved, i-Drive Tri-Array with\$4,774
Mounting Hardware
Part: ST-IDHBT450-U-Q. HCPCS: E2330

i-Drive 4.0 Bluetooth Pro Series,.....\$4,825
3 Switch with Egg Switch
Part: ST-IDHBT300-U-Q. HCPCS: E2330

Adult, Curved, i-Drive Tri-Array with\$4,774
Mounting Hardware
Part: ST-IDHBT400-U-Q. HCPCS: E2330



Adult Pro-Series

Alternative controls and electronics section continued on following page.

25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**25.5. Sip and Puff Drive****1. Sip and Puff Module**

Sip and Puff Module for Q-Logic 3.....\$1,990
 Part: CTL168056. HCPCS: E2325
 Includes module and mounting hardware.

**2. Sip and Puff Straw and Mount**

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

Stealth Sip and Puff for Comfort Series, \$318
 Left Mount (SP625M-CL)
 Part: ST-SP625M-CL-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, \$318
 Head Array (SP625-U)
 Part: ST-SP625-U-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series, \$318
 Right Mount (SP622M-CR)
 Part: ST-SP625M-CR-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

**3. Sip & Puff Accessories**

Sip & Puff Straws, 12 pack (SP625-S12)..... \$35
 Part: ST-SP625-S12-Q

25.6. Egg Switches

Red Egg Switch (ESRED)..... \$87
 Part: SWTMCRO1077
 QTY
 Blue Egg Switch (ESBLU)..... \$87
 Part: SWTMCRO1079
 QTY
 Black Egg Switch (ESBLK)..... \$87
 Part: SWTMCRO1081
 QTY

Green Egg Switch (ESGRN) \$87
 Part: SWTMCRO1078
 QTY
 Yellow Egg Switch (ESYEL) \$87
 Part: SWTMCRO1080
 QTY

**25.7. Switch Adapter**

2 to 1 Switch Adapter..... \$15
 Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

25.8. Environmental Controls

Q-Logic 3 ECU Module Kit\$2,520
 Part: CTL168057

25.9. USB Mobile Device Charger

USB Mobile Device Charger \$184
 Part: ELE168069
 Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

