

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: [quantumorders@pridemobility.com](mailto:quantumorders@pridemobility.com)**Q6 Edge® Z Order Form with Contoured and Synergy Seating**

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_  
Provider Name: \_\_\_\_\_  
Contact: \_\_\_\_\_  
Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
Email: \_\_\_\_\_  
PO Number: \_\_\_\_\_  
Marked for: \_\_\_\_\_  
Ship to Address: \_\_\_\_\_  
City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

**Q6 edge<sup>®</sup> Z****INTRODUCTION**


This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email ([quantumorders@pridemobility.com](mailto:quantumorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering. Options noted with "XRef" have multiple possible part numbers based on system configuration. Please see [cross-reference spreadsheet](#) for the XRef part numbers.

**Thank you for choosing Quantum!****PATIENT INFORMATION**

Weight: \_\_\_\_\_ Height: \_\_\_\_\_



A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_  
B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_  
C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_

**1. BASE MODEL**

**REQUIRED.** Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors, 8 amp off board charger, and switch operated fender lights. Dimensions of base without legrests: 26.25"W x 35.5"L.

Quantum Q6 Edge Z 3S-SS..... \$7,895

Part: Q6 EDGE Z 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

Quantum Q6 Edge Z 3S-C..... \$7,895

Part: Q6 EDGE Z 3S-C. HCPCS: K0849

Captain seat standard.

**2. BASE COLOR****REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.**2.1. Match Colors**

Match All Color Options to Selection Below

Part: N/A. Any subsequent color selections, if available, will be matched to color chosen below when the order is placed. Please note that some color selections, such as joystick armor, may have a MSRP.

**2.2. Power Base Shroud**

|                                                                                                                                 |                                                                                                                                  |                                                                                                                                    |                                                                                                                                    |
|---------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|
|  <b>Red Corvette</b><br>Part: DGN1610108-PRD1  |  <b>Night Train</b><br>Part: DGN1610108-PBK3    |  <b>Evolution Orange</b><br>Part: DGN1610108-POR7 |  <b>Green Machine</b><br>Part: DGN1610108-PGN11 |
|  <b>Electric Blue</b><br>Part: DGN1610108-PBU9 |  <b>Pretty in Pink</b><br>Part: DGN1610108-PPK3 |  <b>Grey Street</b><br>Part: DGN1610108-PSL3      | <b>To Be Determined</b><br>Part: N/A                                                                                               |
|  <b>Back in Black</b><br>Part: DGN1610108-PBK1 |  <b>Deep Purple</b><br>Part: DGN1610108-PPR5    |  <b>Midnight Blue</b><br>Part: DGN1610108-PBU10   |                                                                                                                                    |

**3. DRIVE MOTORS****REQUIRED.**6 mph Drive Motors ..... Standard  
Part: MOT150186Community Use Package ..... \$2,125  
Part: PTOINDV3262

300 lbs. weight capacity, 4-pole high speed motor package, 4" wide flat-free or pneumatic extended traction tires. Increases overall base width by 1.5".

High Speed Motor Package ..... \$1,574  
Part: MOT150141


Programmable speed up to 8 mph. High speed motors decrease per-charge range. Not available with iLevel seating.

Accu-Trac Motor Technology ..... \$1,295  
Part: MOT150181. HCPCS: K0108

Requires Q-Logic 2 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. Programmable speed up to 6 mph.

**4. BASE MOUNTED FOOT PLATFORM****REQUIRED.** Select a foot platform option. May interfere with seat mounted legrests.Base Mounted Foot Platform ..... Standard  
Part: FRMASMB16077

Not available with center mount foot platforms.

Omit Base Mounted Foot Platform ..... No Charge  
Part: Configuration selectionBell Shaped, Base Mounted Foot Platform ..... \$235  
Part: RIGASMB7120041**5. BATTERIES****REQUIRED.** The powerchair requires 2 batteries and will be installed if ordered.Group 24, 70 Amp Gel Cell, Interceptor Battery ..... \$365 ea.  
Part: BAT1704614. HCPCS: E2363No batteries. Provider will assume responsibility of batteries.  
Part: Configuration selectionGroup 24, 70 Amp Gel Cell, MK Battery ..... \$495 ea.  
Part: BATGEL1009. HCPCS: E2363**6. DRIVE WHEELS****REQUIRED.** Select a type of drive wheels. Knobby tires increase overall base width by 1.5".**6.1. Drive Wheels**14" Run Flat Tires ..... Standard  
Part: PTOINDV351014" Run Flat Tires for Drive Wheel Accents ..... Standard  
Part: WHL1703421 x24"W Flat Free, Knobby, Extended Traction Tires ..... \$275.50 ea.  
Part: WHL121698 x24"W Pneumatic, Knobby, Extended Traction Tires ..... \$275.50 ea.  
Part: WHL153737 x2**6.2. Drive Wheel Accent Color****Optional.** Not available with knobby tires. Select with WHL1703421 above. Sold in Pairs (parts below will be quantity of 2).Match to Base Color ..... \$24.98 ea.  
Part: Configuration selection, see specific part belowNon-matching Color ..... \$24.98 ea.  
Part: Configuration selection, see specific part belowRed Corvette  
Part: DGN164245-PRD1Night Train  
Part: DGN164245-PBK3Deep Purple  
Part: DGN164245-PPR5Grey Street  
Part: DGN164245-PSL3Green Machine  
Part: DGN164245-PGN11Electric Blue  
Part: DGN164245-PBU9Pretty in Pink  
Part: DGN164245-PPK3Evolution Orange  
Part: DGN164245-POR7Midnight Blue  
Part: DGN164245-PBU10To Be Determined  
Follow-up required.Back in Black  
Part: DGN164245-PBK1**7. FULL LIGHT PACKAGE****Optional.** Requires NE+ 8 Key or Q-Logic 3 Joystick for Full Light Package.Seat Mounted, LED, Full Light Package ..... \$995  
Part: ELE131477

**8. ELECTRONICS**

**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Forms](#).

**8.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE+ 6-Key Joystick ..... Standard

Part: CTLDC1560

Control up to 2 actuators through joystick.

No Joystick, Setup for NE+ Joystick ..... No Charge

Part: Configuration selection

NE+ 8-Key Joystick (2 Actuators).....\$400

Part: CTLDC1561

Control up to 2 actuators and Seat Mounted LED Light Package through joystick.



NE+ 6 Key Joystick

**8.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

**1. Q-Logic 3 Required Components**

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module ..... \$695

Part: See XRef. HCPCS: E2377

Q-Logic 3 Expandable Controller Harness.....\$425

Part: HAR144005. HCPCS: E2313

**2. Q-Logic 3 Joystick**

Q-Logic 3 EX Joystick .....\$700

Part: CTL167724

Q-Logic 3 EX Joystick for Full Lights .....\$700

Part: CTL167725

No Joystick, Setup for Q-Logic

Part: Configuration selection.



Q-Logic 3 EX Joystick

**8.2. Q-Logic 3 Joystick Armor**

**REQUIRED with Q-Logic 3 Joystick selections above.** Reference colors shown on page 2.

Black, Soft Touch ..... Standard

Part: DGN124307

Non-matching Armor (Select below) ..... \$50

Part: PTOINDV3422

Match Armor to Base Color ..... \$50

Part: PTOINDV3421

Red Corvette

Part: DGN124336-PRD1

Night Train

See Black, Soft Touch

Deep Purple

Part: DGN124336-PPR5

Midnight Blue

Part: DGN124336-PBU10

Electric Blue

Part: DGN124336-PBU9

Pretty in Pink

Part: DGN124336-PPK3

Grey Street

Part: DGN124336-PSL3

Green Machine

Part: DGN124336-PGN11

Back in Black

Part: DGN124336-PBK1

Evolution Orange

Part: DGN124336-POR7

**9. JOYSTICK HANDLES**

**REQUIRED.** The NE Plus are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) ..... \$85

Part: ACCASMB1860. HCPCS: E2323

Extended I-Handle (G) ..... \$85

Part: ACCASMB1857. HCPCS: E2323

Chin Cup (C) ..... \$85

Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

Bodypoint 3" Goal Post Handle ..... \$99

Part: PLSKNCP1057. HCPCS: E2323

T-Bar (D) ..... \$120

Part: ACCASMB1858. HCPCS: E2323

Bodypoint 4" Goal Post Handle ..... \$99

Part: PLSKNCP1058. HCPCS: E2323

Medium Ball Handle (E) ..... \$65

Part: ACCASMB1859. HCPCS: E2323

Q-Logic 3 Carrot Handle (A) ..... No Charge

Part: DGN124308

Soft Ball Handle (F) ..... \$65

Part: ACCASMB1856. HCPCS: E2323



Joystick Handles

**10. JOYSTICK MOUNTING BRACKETS**

**Optional.** Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket

**10.1a. Fixed, In-line Joystick Bracket**

Fixed, In-line, Bracket ..... Standard


Q-Logic Part: FRMASMB10617. NE/NE+ Part: FRMASMB12744

Left Side Mount

Part: PTOINDV1225

Right Side Mount

Part: PTOINDV1226



Swing-Away (Retractable) (right side)

**10.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

Swing-away, Inline, Left Side Mount ..... \$245

NE/NE+ Part: ACC142526. Q-Logic: ACC142524. HCPCS: E1028

Swing-away, 1" Offset, Left Side Mount ..... \$245


NE/NE+ Part: ACC142747. Q-Logic: ACC142742. HCPCS: E1028

Swing-away, Inline, Right Side Mount ..... \$245

NE/NE+ Part: ACC142527. Q-Logic: ACC142525. HCPCS: E1028

Swing-away, 1" Offset, Right Side Mount ..... \$245

NE/NE+ Part: ACC142748. Q-Logic: ACC142743. HCPCS: E1028



Retract4

Joystick mounting brackets continued on following page.

**10. JOYSTICK MOUNTING BRACKETS CONTINUED****10.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

**Retract4 Mount** ..... \$350

NE/NE+ Part: MEC138876. HCPCS: E1028

Q-Logic Part: MEC138860. HCPCS: E1028

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

**Multi-Axis Mount** ..... \$100

NE/NE+ Part: ACC142528. HCPCS: E1028

Q-Logic Part: ACC131216. HCPCS: E1028

SAJ Part: ACC142740. HCPCS: E1028

Allows height adjustment of joystick.


**Stealth Swing-Away/Flip Down Bracket** ..... \$294.25

NE/NE+ Part: ST-TWBM480QL-Q. HCPCS: E1028

Q-Logic Part: ST-TWBM480QL-Q. HCPCS: E1028

SAJ Part: ST-TWBM480CJ-Q. HCPCS: E1028

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.

**Specialty Mount Side****Left Side Mount**

Part: PTOINDV1225

**Right Side Mount**

Part: PTOINDV1226

**11. CONTOURED SEATING****Optional.** Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.**11.1a. Captain Seating****1a. Depth Adjustable High Back with Limited Recline**

Full length armrests must be selected with a depth adjustable high back.

18"Wx18"-20"D with 60" Lap Belt ..... No Charge

Part: SETASMB6895

22"Wx20"-22"D with 70" Lap Belt ..... No Charge

Part: SETASMB6897

20"Wx18"-20"D with 60" Lap Belt ..... No Charge

Part: SETASMB6896

**1b. High Back with Limited Recline**

16"Wx16"D ..... No Charge

Part: SETASMB2770

16"Wx22"D ..... No Charge

Part: SETASMB4887

20"Wx16"D ..... No Charge

Part: SETASMB2780

16"Wx18"D ..... No Charge

Part: SETASMB2771

18"Wx16"D ..... No Charge

Part: SETASMB2775

20"Wx22"D ..... No Charge

Part: SETASMB4889

16"Wx20"D ..... No Charge

Part: SETASMB2772

18"Wx22"D ..... No Charge

Part: SETASMB4888

**1c. Medium Back**

16"Wx16"D ..... No Charge

Part: SETASMB2488

16"Wx18"D ..... No Charge

Part: SETASMB2560

18"Wx16"D ..... No Charge

Part: SETASMB2159

18"Wx18"D ..... No Charge

Part: SETASMB2551

**11.1b. Contoured Seat with Solid Seat Pan****1a. High Back with Limited Recline and Solid Seat Pan**

16"Wx16"D ..... No Charge

Part: SETASMB3178

18"Wx16"D ..... No Charge

Part: SETASMB3181

20"Wx16"D ..... No Charge

Part: SETASMB3183

16"Wx18"D ..... No Charge

Part: SETASMB3179

18"Wx18"D ..... No Charge

Part: SETASMB2848

20"Wx18"D ..... No Charge

Part: SETASMB3184

16"Wx20"D ..... No Charge

Part: SETASMB3180

18"Wx20"D ..... No Charge

Part: SETASMB3182

20"Wx20"D ..... No Charge

Part: SETASMB3185

16"Wx22"D ..... No Charge

Part: SETASMB5006

18"Wx22"D ..... No Charge

Part: SETASMB5007

20"Wx22"D ..... No Charge

Part: SETASMB5008

**1b. Medium Back with Solid Seat Pan**

16"Wx16"D ..... No Charge

Part: SETASMB3112

16"Wx18"D ..... No Charge

Part: SETASMB3113

18"Wx16"D ..... No Charge

Part: SETASMB2492

18"Wx18"D ..... No Charge

Part: SETASMB3114

**11.1c. Special Order Contoured Seat**

Complete information below. Contact inside sales specialist for pricing.

Type of Seat: \_\_\_\_\_

Width: \_\_\_\_\_

Depth: \_\_\_\_\_

Contoured seating section continued on following page.

**11. CONTOURED SEATING CONTINUED****11.2. Contoured Seat Flip-Up, Height Adjustable, Removable Armrests****REQUIRED with High Back Seat except with Depth Adjustable High Back.** Select a Left and Right Armrest or to Omit them.

|                                                                       |                                                                        |
|-----------------------------------------------------------------------|------------------------------------------------------------------------|
| Full Length, Left Armrest.....\$85<br>Part: FRMASMB8614. HCPCS: E0973 | Full Length, Right Armrest.....\$85<br>Part: FRMASMB8615. HCPCS: E0973 |
| Desk Length, Left Armrest.....\$85<br>Part: FRMASMB8612. HCPCS: E0973 | Desk Length, Right Armrest.....\$85<br>Part: FRMASMB8613. HCPCS: E0973 |
| No Armrest ..... Standard<br>Part: N/A                                |                                                                        |

**12. SYNERGY SEATING SYSTEM****Optional.** Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 24"W x 24"D.**12.1. Synergy Frame Selection**

|                                                            |
|------------------------------------------------------------|
| Synergy Seating System .....No Charge<br>Part: PTOINDV2995 |
|------------------------------------------------------------|

**12.2. Synergy Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

|                                                              |
|--------------------------------------------------------------|
| WC-19 Occupied Transit Option.....\$275<br>Part: KITASMB1830 |
|--------------------------------------------------------------|

**12.3. Seat Dimension - Width****REQUIRED with Synergy seating.** Maximum seat width is 24". Part number is SETWDTH plus number following # below.

|                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|
| 12" W<br>(#1002) | 13" W<br>(#1003) | 14" W<br>(#1004) | 15" W<br>(#1005) | 16" W<br>(#1006) | 17" W<br>(#1007) | 18" W<br>(#1008) | 19" W<br>(#1009) | 20" W<br>(#1010) | 21" W<br>(#1011) | 22" W<br>(#1012) | 23" W<br>(#1013) | 24" W<br>(#1014) |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|

**12.4. Seat Dimension - Width Package****REQUIRED with Synergy seating.**

|                                     |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| 10"-14" Width<br>Part: FRMASMB18011 | 14"-18" Width<br>Part: FRMASMB18012 | 16"-20" Width<br>Part: FRMASMB18013 | 20"-24" Width<br>Part: FRMASMB18014 |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|

**12.5. Seat Dimension - Depth****REQUIRED with Synergy seating.** Maximum seat depth is 24". Part number is SETDPTH plus number following # below.

|                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|
| 12" D<br>(#1002) | 13" D<br>(#1003) | 14" D<br>(#1004) | 15" D<br>(#1005) | 16" D<br>(#1006) | 17" D<br>(#1007) | 18" D<br>(#1008) | 19" D<br>(#1009) | 20" D<br>(#1010) | 21" D<br>(#1011) | 22" D<br>(#1012) | 23" D<br>(#1013) | 24" D<br>(#1014) |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|

**12.6. Seat Dimension - Depth Package****REQUIRED with Synergy seating.** Width and depth package part: See XRef: TB3SPACK

|                                            |                                             |                                            |                                            |
|--------------------------------------------|---------------------------------------------|--------------------------------------------|--------------------------------------------|
| 10"-14" Depth Package<br>Part: FRMASMB6981 | 14"-18" Depth Package<br>Part: FRMASMB17823 | 15"-20" Depth Package<br>Part: FRMASMB7732 | 20"-24" Depth Package<br>Part: FRMASMB7733 |
|--------------------------------------------|---------------------------------------------|--------------------------------------------|--------------------------------------------|

**12.7. Synergy Seat to Floor****REQUIRED with Synergy seating.**

|                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|
| 16.50" STF<br>Part: PTOINDV1888 | 17.50" STF<br>Part: PTOINDV1889 | 18.50" STF<br>Part: PTOINDV1890 |
|---------------------------------|---------------------------------|---------------------------------|

**12.8. Synergy Back Options****REQUIRED with Synergy seating.** See the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1. Back Cane Selection**

|                                                                                                                                                                   |                                                                                                                                                                                                                                                                |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Synergy Back Cane ..... No Charge<br>Part: FRMASMB15844<br>Select a back cane height.                                                                             | 2" Endomorph (1" per side) Back Cane.....\$195<br>Part: FRMASMB12115<br>Select a cane height.                                                                                                                                                                  |
| Flip Forward Back Cane ..... \$150<br>Part: FRMASMB15846<br>Select a back cane height. Not available with oxygen holder, or cane mounted and cantilever armrests. | 4" Endomorph (2" per side) Back Cane.....\$195<br>Part: FRMASMB12139<br>Select a cane height.                                                                                                                                                                  |
| 2" Mesomorph (1" per side) Back Cane.....\$195<br>Part: FRMASMB12067<br>Select a cane height.                                                                     | Manual Recline Back ..... \$750<br>Part: FRMASMB12910. HCPCS: E1226<br>Select an overall back height. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, rear accessories, or articulating foot platform. |
| 4" Mesomorph (2" per side) Back Cane.....\$195<br>Part: FRMASMB12091<br>Select a cane height.                                                                     |                                                                                                                                                                                                                                                                |

**2. Synergy Back Cane Height**

Omit this selection if a TRU-Comfort back is desired. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles.

|                                      |                                      |                                      |                                      |
|--------------------------------------|--------------------------------------|--------------------------------------|--------------------------------------|
| 14" High Canes<br>Part: FRMASMB15832 | 17" High Canes<br>Part: FRMASMB15835 | 20" High Canes<br>Part: FRMASMB15838 | 23" High Canes<br>Part: FRMASMB15841 |
| 15" High Canes<br>Part: FRMASMB15833 | 18" High Canes<br>Part: FRMASMB15836 | 21" High Canes<br>Part: FRMASMB15839 | 24" High Canes<br>Part: FRMASMB15842 |
| 16" High Canes<br>Part: FRMASMB15834 | 19" High Canes<br>Part: FRMASMB15837 | 22" High Canes<br>Part: FRMASMB15840 | 25" High Canes<br>Part: FRMASMB15843 |

Synergy seating continued on following page.



**12. SYNERGY SEATING SYSTEM CONTINUED****3. Back Selection**

Synergy Back Cushion..... \$275  
Part: See XRef

TRU-Comfort, 20" Width or less..... \$675  
Part: See XRef. HCPCS: E2620

TRU-Comfort, 22" Width or greater..... \$775  
Part: See XRef. HCPCS: E2621

Sling Back..... \$98  
Part: See XRef

Back cane height must match back height. Not available with flip up cantilever arms.

Omit Back..... No Charge  
Part: PTOINDV1287

**4. Synergy Overall Back Height**

Parts: See XRef, configuration selection with back selection.

|          |          |          |          |          |
|----------|----------|----------|----------|----------|
| 12" High | 15" High | 18" High | 21" High | 24" High |
| 13" High | 16" High | 19" High | 22" High | 25" High |
| 14" High | 17" High | 20" High | 23" High |          |

|                |     | Width |     |     |     |     |     |     |     |     |     |     |     |
|----------------|-----|-------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
|                |     | 12"   | 14" | 15" | 16" | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" |
| Overall Height | 15" | S     | S   | S   | -   | -   | -   | -   | -   | -   | -   | -   | -   |
|                | 16" | S     | S   | S   | -   | -   | -   | -   | -   | -   | -   | -   | -   |
|                | 17" | -     | S   | S   | S   | S   | S   | S   | S   | S   | S   | S   | S   |
|                | 18" | -     | S   | S   | S   | S   | S   | S   | S   | S   | S   | S   | S   |
|                | 19" | -     | S   | S   | S   | S   | S   | S   | S   | S   | S   | S   | S   |
|                | 20" | -     | S   | S   | STC | S   | STC | STC | STC | S   | S   | S   | STC |
|                | 21" | -     | S   | S   | STC | S   | STC | STC | STC | S   | S   | S   | STC |
|                | 22" | -     | S   | S   | STC | STC | STC | STC | STC | S   | S   | STC | STC |
|                | 23" | -     | S   | S   | STC | STC | STC | STC | STC | S   | S   | STC | STC |
|                | 24" | -     | S   | S   | STC | STC | STC | STC | STC | S   | STC | STC | STC |
|                | 25" | -     | S   | S   | STC | STC | STC | STC | STC | S   | STC | STC | STC |

**S** Synergy Back    **STC** Synergy or TRU-Comfort Back

**Synergy and TRU-Comfort Back Availability Chart**

**12.9. Synergy Armrests**

**REQUIRED with Synergy seating.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

**1a. Flip Back, Cane Mounted Armrests**

May interfere with swing-away laterals.

☐ ☐ Full Length, Flip Back Armrests..... Standard  
L R Left part: FRMASMB7038. Right part: FRMASMB7039

☐ ☐ Desk Length, Flip Back Armrests..... Standard  
L R Left part: FRMASMB7040. Right part: FRMASMB7041

**1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests**

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐ Full Length, Two-Post Armrests.....\$150 ea.  
L R Left part: FRMASMB7062. Right part: FRMASMB7063. HCPCS: E0973

☐ ☐ Desk Length, Two-Post Armrests.....\$150 ea.  
L R Left part: FRMASMB7060. Right part: FRMASMB7061. HCPCS: E0973

**1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests**

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H.....\$75 ea.  
L R Left part: FRMASMB9982. Right part: FRMASMB9983. HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H.....\$75 ea.  
L R Left part: FRMASMB9984. Right part: FRMASMB9985. HCPCS: E0973

**1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests**

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H.....\$75 ea.  
L R Left part: FRMASMB9986. Right part: FRMASMB9987. HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14"H.....\$75 ea.  
L R Left part: FRMASMB9988. Right part: FRMASMB9989. HCPCS: E0973

**1e. Pediatric, Removable, Height Adjustable Armrests**

Only available with desk length armpads.

☐ ☐ Pediatric, Removable Armrests.....\$75 ea.  
L R Left part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973

**1f. Height Adjustable, Cantilever Armrests**

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

☐ ☐ Full Length, Cantilever Armrests.....\$97.50 ea.  
L R Left part: FRMASMB16188. Right part: FRMASMB16189

☐ ☐ Desk Length, Cantilever Armrests.....\$97.50 ea.  
L R Left part: FRMASMB16190. Right part: FRMASMB16191

**1g. Single Post, Height Adjustable, Removable Armrests**

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests.....\$75 ea.  
L R Left part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973

☐ ☐ Desk Length, Single Post Armrests.....\$75 ea.  
L R Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973

**12.10. Synergy Armpads**

**REQUIRED with armrest selection above.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

Straight Armpads..... Standard  
Desk/Desk part: UPSASMB1022. Full/Full part: UPSASMB1021.  
Desk/Full part: UPSASMB1023

Waterfall Armpads..... No Charge  
Desk/Desk part: UPSASMB1025. Full/Full part: UPSASMB1024.  
Desk/Full part: UPSASMB1026

No Armpad Selection..... No Charge  
Part: PTOINDV1373

**12.11. Joystick Bracket Receiver**

**REQUIRED.**

Standard Fixed Receiver..... Standard  
Part: FRMASMB16062


Stealth Height Adjustable Receiver..... \$214  
Tube Arm Part: ST-ARM250C-Q. HCPCS: K0108

**13. HEADRESTS****Optional.** Select a headrest pad and bracket. **Not available with contoured seating.****13.1. Headrest Pads**Stealth 8" Comfort Plus (CP450)..... \$151.94  
Part: POS148484. HCPCS: E0955Stealth 10" Comfort Plus (CP250)(A) ..... \$151.94  
Part: SETHEAD1057. HCPCS: E0955Stealth 14" Comfort Plus (CP150) ..... \$151.94  
Part: INDPART2327. HCPCS: E0955Pediatric Comfort Plus (6"Wx4.5"H) ..... \$151.94  
(CPN350)  
Part: SETHEAD1069. HCPCS: E0955Stealth Small Contoured (9.75"Wx3"H)..... \$64.20  
(STL-OTSM250)

Part: SETHEAD1052. HCPCS: E0955

Stealth Large Contoured (11"Wx6"H) (B)..... \$80.25  
(STL-OTLG150)

Part: SETHEAD1053. HCPCS: E0955

**13.2. Headrest Brackets**Stealth's Worlds Best Removable..... \$220.50  
Headrest Bracket (TWB485-1) (C)  
Part: FRMASMB15999. HCPCS: E1028Stealth Multi-Axis Removable Headrest..... \$160.50  
Bracket (HMO475P) (D)  
Part: HDWASMB1084. HCPCS: E1028**13.3. Flip-down Headrest Bracket**Stealth Flip-down Headrest Bracket (FDM350)..... \$224.70  
Part: INDPART2219

The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

**14. LAP BELTS****REQUIRED.** Select a lap belt below.**14.1a. Standard Lap Belts**38" Lap Belt  
Part: ACCBELT100650" Lap Belt  
Part: ACCBELT100060" Lap Belt  
Part: ACCBELT100270" Lap Belt  
Part: ACCBELT100380" Lap Belt  
Part: ACCBELT100490" Lap Belt  
Part: ACCBELT1001**14.1b. Optional Lap Belts**i-Fit 60" Single Pull, Push Button, Medium ..... \$86  
(7.5") Padded, 1.5" Belt (IFPB5213)  
Part: ST-IFPB5213-Qi-Fit 60" Single Pull, Push Button, Extra Large ..... \$86  
Padded (11.5"), 2" Belt (IFPB5315)  
Part: ST-IFPB5315-Qi-Fit 60" Single Pull, Push Button, Large ..... \$86  
Padded (9.5"), 2" Belt (IFPB5314)  
Part: ST-IFPB5314-Q**15. POWER LEGREST CONTROL****Optional.** Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.Independent Legrests (Requires 2 actuator controls)  
Part: Configuration selectionArticulating Foot Platform (Requires 1 Actuator Control)  
Part: Configuration selectionCombined Legrests (Requires 1 actuator control)  
Part: Configuration selection**16. SYNERGY POWER OPTION ELECTRONICS****REQUIRED for Synergy Seat with Power Legrests.** Select the desired control type for the seating system.One Actuator Control through NE+ Joystick ..... \$1,895  
Part: See XRef. HCPCS: E2310One Actuator Control through Q-Logic 3 Electronics ..... \$2,895  
Part: See XRef. HCPCS: E2310Two Actuator Control through NE+ Joystick ..... \$2,895  
Part: See XRef. HCPCS: E2311Multiple Actuator Control through Q-Logic 3 Electronics..... \$2,895  
Part: See XRef. HCPCS: E2311**17. LEGRESTS****Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**17.1. Center Mount Foot Platform with Angle Adjustable Footplates**If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section. **Not available with contoured seating.****1. Foot Platform Assembly**Power Articulating Foot Platform ..... \$2,995  
Part: PTOINDV3754. HCPCS: E1012+K0040Manually Adjustable Center Mount Foot Platform... \$250  
12"-16" width part: FRMASMB17944. 14"-18" width part:  
FRMASMB17945. 16"-20" width part: FRMASMB17946.  
HCPCS: K0040.**2. Lower Extension and Angle Adjustable Footplates**8"-12" Extension, 5"x6" A.A. Footplates  
Part: RIG152577  
Not available with power articulating platform.10"-16" Extension, 5"x6" A.A. Footplates  
Part: RIG15020613"-19" Extension, 5"x6" A.A. Footplates  
Part: RIG1515058"-12" Extension, 5"x8" A.A. Footplates  
Part: RIG152578  
Not available with power articulating platform.10"-16" Extension, 5"x8" A.A. Footplates  
Part: RIG15020313"-19" Extension, 5"x8" A.A. Footplates  
Part: RIG15150610"-16" Extension, 5"x11" A.A. Footplates  
Part: RIG15021213"-19" Extension, 5"x11" A.A. Footplates  
Part: RIG151508

Legrests section continued on following pages.

**17. LEGRESTS CONTINUED**

| Center Mount Foot Platform<br>Compatibility Matrix | Synergy Seat |        |      |                                                                                                                                                                                                              |
|----------------------------------------------------|--------------|--------|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                    | Low          | Medium | High |                                                                                                                                                                                                              |
| 10"-16" Length Extension                           | !            | !      | ✓    | <div>✓ No interference issues at this seat-to-floor height.</div> <div>! May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.</div> |
| 13"-19" Length Extension                           | !            | !      | !    |                                                                                                                                                                                                              |

**3. Center Mount Calf Support**

\*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

|                                                         |                                             |
|---------------------------------------------------------|---------------------------------------------|
| Small Calf Pads (pair) 6"Wx6"H..... \$150*              | Small Calf Panel 18"Wx8"H..... \$150        |
| Part: RIGASMB7120018                                    | Part: FRMASMB7110012                        |
| Large Calf Pads (pair) 8" taper to 6"Wx8"H ..... \$150* | Large Calf Panel 24"Wx8"H ..... \$200       |
| Part: RIGASMB7120019                                    | Part: RIG123824                             |
|                                                         | Only available on widths of 20" or greater. |

**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

|                                                                                             |                                                                      |
|---------------------------------------------------------------------------------------------|----------------------------------------------------------------------|
| Width Extension for 8"D Footplates ..... \$50                                               | 4" Outside Edge (pair) ..... \$50                                    |
| Part: RIG152517                                                                             | Part: RIG152510                                                      |
| Width Extension for 11"D Footplates ..... \$50                                              | Only available with 6"D footplates and 8"D footplate extension.      |
| Part: RIG152518                                                                             | 6" Outside Edge (pair) ..... \$50                                    |
| 2" Bell-Shaped Length Extension ..... \$50                                                  | Part: RIG152516                                                      |
| Part: RIG152519                                                                             | Available with 8"D and 11"D footplates and 11"D footplate extension. |
| Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates. |                                                                      |

**5. Additional Accessories**

|                                                   |                                        |
|---------------------------------------------------|----------------------------------------|
| Reversed Dartex Calf Pad Cover (pair) ..... \$150 | Wheel Bumpers (pair) ..... \$12.50 ea. |
| Small part: RIG133717. Large part: RIG133719      | Part: WHLASMB1578 x2                   |
| Reversed Dartex Panel Cover (pair) ..... \$150    | Footplate Connector ..... \$15         |
| Small part: RIG134673. Large part: RIG134744      | Part: RIGASMB7120023                   |

**6. Stealth Gel Footplate Covers**

|                                                |                                                                     |
|------------------------------------------------|---------------------------------------------------------------------|
| 5"x8" Tapered Gel Footplate Covers ..... \$205 | 6.5"x8" Tapered Gel Footplate Covers for Width Extension..... \$225 |
| Part: ST-FPG5008T-Q                            | Part: ST-FPG6508T-Q                                                 |
| 5"x11" Tapered Gel Footplate Covers..... \$225 | 7"x11" Tapered Gel Footplate Cover for Width Extension ..... \$250  |
| Part: ST-FPG5011T-Q                            | Part: ST-FPG7011T-Q                                                 |

**17.2. 70° Power Articulating, Swing Away Legrests with Angle Adjustable Footplates and Heel Loops**

4" of articulating range. Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Articulating Legrest Upper Assembly**

|                                                       |                                                       |
|-------------------------------------------------------|-------------------------------------------------------|
| P L R 70° Articulating, Swing-Away.....\$1,197.50 ea. | P L R 70° Articulating, Swing-Away.....\$1,197.50 ea. |
| Legrest with 14.5"-18" Extension                      | Legrest with 11"-13" Extension                        |
| Pair Part: FRMASMB9073. Left Part: FRMASMB9074.       | Pair part: FRMASMB11444. Left Part: FRMASMB11366.     |
| Right Part: FRMASMB9075. HCPCS: E1010+K0040+E0951     | Right Part: FRMASMB11367. HCPCS: E1010+K0040+E0951    |

| 70° Power Articulating, Swing-Away<br>Legrest Compatibility Matrix | Synergy Seat |        |      |                                                                                                                                                                                                                                                                           |
|--------------------------------------------------------------------|--------------|--------|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                    | Low          | Medium | High |                                                                                                                                                                                                                                                                           |
| Angle Adjustable, 11" - 13"                                        | ✓            | ✓      | ✓    | <div>✓ This rigging configuration causes no interference issues at this seat-to-floor height.</div> <div>! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.</div> |
| Angle Adjustable, 14.5" - 18"                                      | !            | !      | !    |                                                                                                                                                                                                                                                                           |

**2. Angle Adjustable Footplate with Heel Loop and Wheel Bumper**

\*The 6"x11" footplates may cause interference with seat widths less than 18". \*The 8"x11" footplates may cause interference with seat widths less than 20".

|                                                                      |                                                                      |                                                                      |
|----------------------------------------------------------------------|----------------------------------------------------------------------|----------------------------------------------------------------------|
| P L R 4"Wx6"D                                                        | P L R 5"Wx8"D                                                        | P L R 6"Wx11"D <sup>1</sup>                                          |
| Pair part: FRMASMB7448. Left part: RIG133543. Right part: RIG133566. | Pair part: FRMASMB7450. Left part: RIG133545. Right part: RIG133568. | Pair part: FRMASMB7452. Left part: RIG133548. Right part: RIG133570. |
| P L R 5"Wx6"D                                                        | P L R 6"Wx8"D <sup>1</sup>                                           | P L R 8"Wx11"D <sup>2</sup>                                          |
| Pair part: FRMASMB7449. Left part: RIG133544. Right part: RIG133567. | Pair part: FRMASMB7451. Left part: RIG133547. Right part: RIG133569. | Pair part: FRMASMB7453. Left part: RIG133549. Right part: RIG133571. |

**3. Power Articulating, Swing-Away Legrest Accessories**

|                                        |                                                                                                  |
|----------------------------------------|--------------------------------------------------------------------------------------------------|
| Calf Pad (pair) 6"Wx6"H..... No Charge | Falcon Gel Padded Knee Adductor Buttons (KBE)(pair) .....\$125                                   |
| Part: MEC153735                        | Part: INDPART2242. HCPCS: E0953                                                                  |
|                                        | Please verify the appropriate billing code to be used for this item with all third-party payors. |

Legrests section continued on following pages.



**17. LEGRESTS CONTINUED****17.3. 60°-90° Swing-Away Legrests**

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Swing-Away Legrest Upper**

**60°-90° Swing-Away Legrest**..... No Charge  
 Pair part: RIG120062. Left: FRMASMB17558.  
 Right: FRMASMB17559

**60°-90° Swing-Away Legrest, Short**..... No Charge  
 Pair part: RIG136481. Left: RIG132607.  
 Right: RIG132609.

**2. Lower Extension**

**11"-13" Length**  
 Part: See XRef  
 N/A with Multi-axis.

**13"-16.5" Length**  
 Part: See XRef  
 N/A with Multi-Axis.

**11"-16.5" Length**  
 Part: See XRef  
 N/A with Angle Adjustable.

**15.5"-19.5" Length**  
 Part: See XRef

| 60°-90° Swing-Away Legrest Compatibility Matrix | Contoured Seat |        |      | Synergy Seat |        |      |
|-------------------------------------------------|----------------|--------|------|--------------|--------|------|
|                                                 | Low            | Medium | High | Low          | Medium | High |
| Extruded, 13" - 16.5"                           | ✓              | ✓      | ✓    | !            | !      | ✓    |
| Extruded, 15.5" - 19.5"                         | !              | ✓      | ✓    | X            | !      | !    |
| Angle Adjustable, 13" - 16.5"                   | ✓              | ✓      | ✓    | !            | !      | ✓    |
| Angle Adjustable, 15.5" - 19.5"                 | !              | ✓      | ✓    | X            | !      | !    |
| Multi-Axis, 11" - 16.5"                         | ✓              | ✓      | ✓    | !            | !      | ✓    |
| Multi-Axis, 15.5" - 19.5"                       | !              | ✓      | ✓    | X            | !      | !    |

11"-13" Extensions cause no interference at all seat-to-floor heights.

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

X This rigging configuration does not work with this seat-to-floor height.

**3a. Extruded Footplate(s)**

**7.25"Wx6"D Extruded Footplate**..... Standard  
 Part pair: FRMASMB7847. Left: FRMASMB7861  
 Right: FRMASMB7862

**8.5"Wx6"D Extruded Footplate**..... Standard  
 Part pair: FRMASMB7848. Left: FRMASMB7863  
 Right: FRMASMB7864  
 May cause interference on seat widths less than 18."

**3b. Angle Adjustable Footplate(s) with Heel Loop(s)**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

**4"Wx6"D** ..... \$85 ea.  
 Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560  
 HCPCS: K0040+E0951

**6"Wx8"D<sup>1</sup>** ..... \$85 ea.  
 Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563  
 HCPCS: K0040+E0951

**5"Wx6"D** ..... \$85 ea.  
 Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561  
 HCPCS: K0040+E0951

**6"Wx11"D<sup>1</sup>** ..... \$85 ea.  
 Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564  
 HCPCS: K0040+E0951

**5"Wx8"D** ..... \$85 ea.  
 Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562  
 HCPCS: K0040+E0951

**8"Wx11"D<sup>2</sup>** ..... \$85 ea.  
 Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565  
 HCPCS: K0040+E0951

**3c. Multi-Axis Footplate(s) with Heel Loops(s)**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

**4"Wx6"D** ..... \$247.50 ea.  
 Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552  
 HCPCS: K0040+E0951

**6"Wx8"D<sup>1</sup>** ..... \$247.50 ea.  
 Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555  
 HCPCS: K0040+E0951

**5"Wx6"D** ..... \$247.50 ea.  
 Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553  
 HCPCS: K0040+E0951

**6"Wx11"D<sup>1</sup>** ..... \$247.50 ea.  
 Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558  
 HCPCS: K0040+E0951

**5"Wx8"D** ..... \$247.50 ea.  
 Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554  
 HCPCS: K0040+E0951

**8"Wx11"D<sup>2</sup>** ..... \$247.50 ea.  
 Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559  
 HCPCS: K0040+E0951

**4. 60°-90° Swing-Away Legrest Accessories**

**Heel Loops (pair)**..... \$50  
 Part: FRMASMB7873. HCPCS: E0951  
 Standard on angle adjustable and multi-axis footplates.

**Wheel Bumper**..... \$12.50 ea.  
 Part: See XRef

**Gel Neoprene Footrest Sleeve (pair)** ..... \$80  
 Part: INDPART2243. HCPCS: K0108

**5. Amputee Supports with Swing-Away Upper**

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

**Therafin 31515 Multi-Axis Amputee Support** ..... \$230 ea.  
 Part: INDPART2822. HCPCS: E1020

**Therafin 32387 Amputee Support Gel Cover** ..... \$155 ea.  
 Part: INDPART2823. HCPCS: K0108

Legrests section continued on following page.

**17. LEGRESTS CONTINUED****17.4. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Elevating Legrests Upper**

☐ ☐ ☐ Elevating Legrest Assembly ..... \$162.50 ea.  
 Pair part: FRMASMB483. Left: RIG138012  
 Right: RIG138013. HCPCS: E0990  
 Select lower extension and footplates

**2. Lower Extension**

☐ ☐ ☐ 11"-14" Length .....  
 Pair: FRMASMB4755  
 Left: RIG138026  
 Right: RIG138030  
 For reinforced only

☐ ☐ ☐ 14"-17" Length .....  
 Pair: FRMASMB4756  
 Left: RIG138033  
 Right: RIG138035  
 For reinforced only

☐ ☐ ☐ 12"-14.5" Length .....  
 Part: See XRef  
 For angle adjustable only

☐ ☐ ☐ 14.5"-18" Length .....  
 Part: See XRef  
 For angle adjustable only

| Elevating Legrest<br>Compatibility Matrix | Contoured Seat |        |      | Synergy Seat |        |      |
|-------------------------------------------|----------------|--------|------|--------------|--------|------|
|                                           | Low            | Medium | High | Low          | Medium | High |
| Standard, 11"-14"                         | ✓              | ✓      | ✓    | ✓            | ✓      | ✓    |
| Standard, 14"-17"                         | ✓              | ✓      | ✓    | !            | !      | !    |
| Angle Adjustable, 12"-14.5"               | ✓              | ✓      | ✓    | ✓            | ✓      | ✓    |
| Angle Adjustable, 14.5"-18"               | ✓              | ✓      | ✓    | !            | !      | !    |

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

**3a. Reinforced Cast Footplate(s)**

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate..... Standard  
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate ..... Standard  
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062  
 May cause interference on seat widths less than 18."

**3b. Angle Adjustable Footplate(s)**

\*The 6"x11" footplates may cause interference with seat widths less than 18". \*The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

☐ ☐ ☐ 4"Wx6"D ..... \$85 ea.  
 Pair parts: FRMASMB8812 x2.  
 HCPCS: K0040

☐ ☐ ☐ 5"Wx6"D..... \$85 ea.  
 Pair parts: FRMASMB8811 x2.  
 HCPCS: K0040

☐ ☐ ☐ 5"Wx8"D with Heel Loops..... \$85 ea.  
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792  
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx8"D<sup>1</sup> with Heel Loops..... \$85 ea.  
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795  
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx11"D<sup>1</sup> with Heel Loops ..... \$85 ea.  
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797  
 HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D<sup>2</sup> with Heel Loops ..... \$85 ea.  
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799  
 HCPCS: K0040+E0951

**4. Elevating Legrest Accessories**

Heel Loops (pair)..... \$50  
 Part: FRMASMB7873. HCPCS: E0951  
 Standard on angle adjustable footplates.

Wheel Bumper..... \$12.50 ea.  
 Part: See XRef  
 Available on angle adjustable footplates.

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) .....\$125  
 Part: INDPART2242. HCPCS: E0953  
 Please verify the appropriate billing code to be used for this item with all third-party payors.

**18. SET LEGREST LENGTH**

**REQUIRED with matching left and right legrest configuration.** Non-matching or individual legrest configurations lengths are not set by the manufacturer.

|                   |                   |                    |                    |                    |                    |                    |                    |                    |                    |                    |
|-------------------|-------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|
| 7" L<br>(#3470)   | 8.5" L<br>(#3473) | 10" L<br>(#3476)   | 11.5" L<br>(#3479) | 13" L<br>(#3482)   | 14.5" L<br>(#3485) | 16" L<br>(#3488)   | 17.5" L<br>(#3491) | 19" L<br>(#3494)   | 20.5" L<br>(#3497) | 22" L<br>(#3669)   |
| 7.5" L<br>(#3471) | 9" L<br>(#3474)   | 10.5" L<br>(#3477) | 12" L<br>(#3480)   | 13.5" L<br>(#3483) | 15" L<br>(#3486)   | 16.5" L<br>(#3489) | 18" L<br>(#3492)   | 19.5" L<br>(#3495) | 21" L<br>(#3498)   | 22.5" L<br>(#3670) |
| 8" L<br>(#3472)   | 9.5" L<br>(#3475) | 11" L<br>(#3478)   | 12.5" L<br>(#3481) | 14" L<br>(#3484)   | 15.5" L<br>(#3487) | 17" L<br>(#3490)   | 18.5" L<br>(#3493) | 20" L<br>(#3496)   | 21.5" L<br>(#3600) |                    |

**19. SEAT CUSHION**


**Optional. Specific cushion/component sizes available, part numbers, weight capacities, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#).** This section allows a quick selection of a single cushion. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

**19.1. Seat Cushion Model Selection****Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.  
 Part series: ST-SIMGXXYY

Simplicity G, 12"W to 21" W, E2601..... \$110


Simplicity G, 22"W to 24"W, E2602..... \$165

**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
 Part series: ST-SPCFSPXXYY

Spectrum Foam SPP, 16"W to 20" W, E2607 .....\$460


Spectrum Foam SPP, 22"W to 24"W, E2608..... \$530

**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
 Part series: ST-TC2SPPXXYY

TRU-Comfort 2 SPP, 16"W to 20" W, E2607 .....\$400

TRU-Comfort 2 SPP, 22"W, E2608..... \$500

**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
 Part series: ST-SPCGSPXXYY

Spectrum Gel SPP, 14"W to 20" W, E2607 .....\$460

Spectrum Gel SPP, 22"W to 24"W, E2608..... \$530



Seat cushion section continued on following page.


**19. SEAT CUSHION CONTINUED****Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SOLSPXXYY

Solution SPP, 12"W to 20" W, E2607 ..... \$375

Solution SPP, 22"W to 24"W, E2608 ..... \$430


**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.

Part series: ST-GLSPXXYY

Glacial SP, 12"W to 20"W, E2603 ..... \$330

Glacial SP, 22"W to 24"W, E2604 ..... \$440


**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.

Part series: ST-ZENSPXXYY

Zen SP, 12"W to 21" W, E2603 ..... \$315

Zen SP, 22"W to 24"W, E2604 ..... \$425

**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Part series: ST-ESSSPXXYY

Essence SPP, No Wedge, 16"-20" Wide, E2607 ..... \$400

Essence SPP, No Wedge, 22" Wide, E2608 ..... \$500

Essence SPP, Small Wedge, 16"-20" Wide, E2607 .... \$435


Essence SPP, Small Wedge, 22" Wide, E2608 ..... \$535

Essence SPP, Medium Wedge, 16"-20" Wide, E2607 . \$465

Essence SPP, Medium Wedge, 22" Wide, E2608 ..... \$565

Essence SPP, Large Wedge, 16"-20" Wide, E2607 .... \$495

Essence SPP, Large Wedge, 22" Wide, E2608 ..... \$595

**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.

Part series: ST-PREMPXXYY

Premiere P, 16"W to 20" W, E2605 ..... \$325

Premiere P, E2606, 22" W ..... \$475

**19.2. Seat Cushion Dimensions****REQUIRED with cushion model selection.** Select to match seating or pick non-matching dimensions. See respective cushion order form for pricing and available sizes.

Match Seat Width and Depth

**Cushion Width**

|       |       |       |       |       |       |
|-------|-------|-------|-------|-------|-------|
| 12" W | 15" W | 17" W | 19" W | 21" W | 24" W |
| 14" W | 16" W | 18" W | 20" W | 22" W |       |

**Cushion Depth**

|       |       |       |       |       |       |
|-------|-------|-------|-------|-------|-------|
| 12" D | 15" D | 17" D | 19" D | 21" D | 24" D |
| 14" D | 16" D | 18" D | 20" D | 22" D |       |

**19.3. Seat Cushion Substitute Fluid Proof Liners**

Fluid Proof Urethane Liner ..... \$50-\$70

Part: URL1709070-XXYY

Not available with Simplicity, Essence, or TRU-Comfort 2.

Fluid Proof Polyurethane Liner ..... \$55-\$75

Essence/TC2 part: PUL1709080-XXYY. Others part: PUL1709030-XXYY

Not available with Simplicity.

**19.4. Seat Cushion Substitute Cover**

An additional cover may be ordered or the standard cover may be substituted with a cover by checking the option in subsection 2. The "\*"prefix" of the part number will be the part prefix of the cushion model selected above. For example, SPCGSP-NBLK1816 will be an additional coolcore cover for a Spectrum Gel cushion.

Stealth Tek 2 Cover ..... No Charge

Part: \*prefix-TBLKXXYY

Not available with Simplicity, Essence, Premiere, or TRU-Comfort 2.

Mesh Cover ..... No Charge

Part: \*prefix-MBLKXXYY

Not available with Simplicity or TRU-Comfort 2.

**20. ACCESSORIES****Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Contoured Seat Cup Holder ..... \$20

Part: ACCASMB2287

Synergy Seat Cup Holder ..... No Charge

Part: ACCASMB2634

Oxygen Holder for Contoured Seat ..... \$118

Part: ACCASMB2292. HCPCS: E2208

Oxygen Holder for Synergy Seat ..... \$150

Not available with Synergy Manual Recline.

Part: FRMASMB7030. HCPCS: E2208

Walker Holder ..... \$103

Contoured seat part: ACCASMB1912. Synergy part: FRMASMB7031

Cane/Crutcher Holder ..... \$110

Contoured seat part: ACCASMB2294. HCPCS: E2207

Synergy part: FRMASMB4822. HCPCS: E2207

Cell Phone Holder for 12" Lift Captain Seat ..... \$70

Part: ACC1611300

Rear Basket ..... \$118

Contoured seat part: ACCASMB1914. Synergy part: FRMASMB4821

Flag Holder ..... \$25

Contoured seat part: ACCASMB1930. Synergy part: ACCASMB2415.

Saddle Bag ..... \$35

Part: DWR1020D010

Stealth Cell Phone/Media Device Holder ..... \$195

Part: ACC7120006

Stealth Large Tablet Holder ..... \$375

Part: ACC7120005

Quantum Backpack ..... \$150

Part: ACC1905853

Grip Lap Board (20005) ..... \$38.95

Part: ACC152992

Grip Drink Holder (20085) ..... \$15.50

Part: ACC152993

Accessory Bar ..... \$100

12"W-16"W part: MEC145955. 16"W-20"W part: MEC143975.

20"W-22"W part: MEC146018

Large, Black, Weather Cover ..... \$202

Part: ACCCOVR1021

**21. BRANCH PATCHES****Optional.**

Air Force Patch ..... \$20

Part: ACC140240

Army Patch ..... \$20

Part: ACC140241

Coast Guard Patch ..... \$20

Part: ACC140239

Marines Patch ..... \$20

Part: ACC140243

Navy Patch ..... \$20

Part: ACC140244

Purple Heart Patch ..... \$20

Part: ACC140238



Military Patches

**22. ELECTRONICS ACCESSORIES****Optional.**Q-Logic EX/NE Programmer ..... \$545  
Part: ELEASMB6976Q-Logic PC Programming Disk with Cable ..... \$750  
Part: ELEASMB5216**23. ALTERNATIVE CONTROLS AND ELECTRONICS****Optional. Only available with the Synergy Seating system.** Select a type of specialty controls below or to refer to the specialty controls order form.**23.1. Enhanced Display**Q-Logic 3 Enhanced Display Kit.....\$1,720  
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

**Q-Logic 3 Enhanced Display Shroud**

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch ..... Standard Non-matching Shroud ..... \$50 Match Shroud to Base Color .... \$50  
Part: DGN137728 Select color below Part: See color belowRed Corvette  
Part: DGN137729-PRD1Night Train  
Part: DGN137729-PBK3Deep Purple  
Part: DGN137729-PPR5Midnight Blue  
Part: DGN137729-PBU10Electric Blue  
Part: DGN137729-PBU9Pretty in Pink  
Part: DGN137729-PPK3Grey Street  
Part: DGN137729-PSL3Green Machine  
Part: DGN137729-PGN11Back in Black  
Part: DGN137729-PBK1Evolution Orange  
Part: DGN137729-POR7

Enhanced Display

**23.2. SCIM Module****Optional.**☐ SCIM Module for Q-Logic 3 ..... \$525  
Part: CTL168055

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

**23.3. Q-Logic 3 SAJ**☐ Q-Logic 3 EX SAJ ..... \$700  
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

**23.4. Head Array**Adult, i-Drive Tri-Array Series, 3-Switch ..... \$4,414  
Head Array (IDHBT200-1)

Part: ST-IDHBT200-1-Q. HCPCS: E2330

Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.

Adult, i-Drive Pro Series Head Array ..... \$4,461  
Head Array (IDHBT300-1)

Part: ST-IDHBT300-1-Q. HCPCS: E2330

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

**23.5. Sip and Puff Drive****1. Sip and Puff Module**Sip and Puff Module for Q-Logic 3 ..... \$1,895  
Part: CTL168056. HCPCS: E2325

Includes module and mounting hardware.

**2. Sip and Puff Straw and Mount**Stealth Sip and Puff for Comfort Series, ..... \$294.25  
Left Mount (SP625-CL)

Part: ST-SP625-CL-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, ..... \$294.25  
Head Array (SP625-U)

Part: ST-SP625-U-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.



Sip &amp; Puff Module

Stealth Sip and Puff for Comfort Series, ..... \$294.25  
Right Mount (SP622-CR)

Part: ST-SP625-CR-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Stealth Sip &amp; Puff

**23.6. Egg Switches**Red Egg Switch (ESRED) ..... \$80.85 ea.  
Part: SWTMCRO1077Blue Egg Switch (ESBLU) ..... \$80.85 ea.  
Part: SWTMCRO1079Black Egg Switch (ESBLK) ..... \$80.85 ea.  
Part: SWTMCRO1081Green Egg Switch (ESGRN) ..... \$80.85 ea.  
Part: SWTMCRO1078Yellow Egg Switch (ESYEL) ..... \$80.85 ea.  
Part: SWTMCRO1080

Egg Switch

Alternative controls and electronics section continued on following page.

**23. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED****23.7. Switch Adapter**

2 to 1 Switch Adapter ..... \$13.39  
 Part: DWR1361H037  
 Connects 2 single switches to a stereo jack to act as a dual switch.

**23.8. Environmental Controls**

Q-Logic 3 ECU Module Kit ..... \$2,400  
 Part: CTL168057

**24. POSITIONING COMPONENTS**

**Optional. Only available with the Synergy Seating system.** More options are available on the [Positioning Components Order Form](#).

**24.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

**1. Bracket**

Swing-away, Rail Mounts, pair (TWBTLTC) ..... \$472.50  
 Part: INDPART2677. HCPCS: E1028

**2. Pads**

Curved Pads, 3"x5", pair (TWBL-3X5) ..... \$210  
 Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6) ..... \$210  
 Part: INDPART2664. HCPCS: E0956

**3. Installation**

Install Laterals .....  
 Part: PTOINDV3499

Do Not Install .....  
 Part: PTOINDV3500

**24.2. Stealth Unilink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable, ..... \$245 ea.  
 Removable Mounting Hardware  
 Part: ST-UTTP5-Q

UniLink, 7" Height Adjustable, ..... \$255 ea.  
 Removable Mounting Hardware  
 Part: ST-UTTP7-Q

**2. UniLink Pelvic/Thigh Pads**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad ..... \$95 ea.  
 Part: ST-UTP544056100-Q

UniLink, 4" x 12" Pad ..... \$125 ea.  
 Part: ST-UTP54D056100-Q

UniLink, 4" x 8" Pad ..... \$105 ea.  
 Part: ST-UTP548056100-Q



UniLink Pad &amp; Hardware

**3. UniLink Pelvic/Thigh Additional Links**

UniLink, 1.5" Link ..... \$75 ea.  
 Part: ST-UIL001150A000-Q

UniLink, 2" Link ..... \$75 ea.  
 Part: ST-UIL001200A00-Q

**4. UniLink Pelvic/Thigh Guide Installation**

Install on Left .....  
 Part: PTOINDV3585

Install on Right .....  
 Part: PTOINDV3586

Do Not Install .....  
 Part: PTOINDV3719

**24.3. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

**1. Pelvic/Thigh Guide Bracket**

Flip Down Bracket (TWBTLADD-HW) ..... \$236.25 ea.  
 Part: INDPART2804. HCPCS: E1028

Removable Bracket (LP430) ..... \$262.50 ea.  
 Part: INDPART2805. HCPCS: E1028

Removable Bracket with 1" Offset ..... \$400 ea.  
 (LP430 & LHW-121)  
 Part: INDPART2854. HCPCS: E1028

Flip Down with LEFT Full Surface Contact ..... \$236.25 ea.  
 Feature (TWBADD-L & LPHW407)  
 Part: POS140428. HCPCS: E1028

Flip Down with RIGHT Full Surface Contact ..... \$236.25 ea.  
 Feature (TWBADD-R & LPHW407)  
 Part: POS140429. HCPCS: E1028

**2. Pelvic/Thigh Guide Adapter**

2" Height Adapter Plate (LHW-122) ..... \$26.25  
 Part: INDPART2821. HCPCS: K0108  
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129) ..... \$53.50  
 Part: POS140430. HCPCS: K0108  
 Only available with the flip down full surface bracket.

**3. Pelvic/Thigh Guide Pad**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Short Pad, 4"x4" (TWBADD-4X4) ..... \$75 ea.  
 Part: INDPART2801. HCPCS: E0953^/E0956

Long Pad, 4"x12" (TWBADD-4X12) ..... \$75 ea.  
 Part: INDPART2803. HCPCS: E0953^/E0956

Medium Pad, 4"x8" (TWBADD-4X8) ..... \$75 ea.  
 Part: INDPART2802. HCPCS: E0953^/E0956

**4. Pelvic/Thigh Guide Installation**

Install on Left .....  
 Part: PTOINDV3585

Install on Right .....  
 Part: PTOINDV3586

Do Not Install .....  
 Part: PTOINDV3719

Manufacturer installed positioning components section continued on following page.



**24. POSITIONING COMPONENTS CONTINUED****24.4. Stealth i-Fit Chest Harness**

Upper Classic, Medium (IFPS1103) ..... \$142 ea.  
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) ..... \$142 ea.  
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$142 ea.  
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) ..... \$142 ea.  
 Part: ST-IFPS2104-Q. HCPCS: E0960

**24.5. Stealth i-Fit Chest Straps**

Chest Strap, Medium w/ Belt Anchor ..... \$90 ea.  
 (IFCS4113)  
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor ..... \$90 ea.  
 (IFCS4114)  
 Part: ST-IFCS4114-Q. HCPCS: E0960

**24.6. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. The 1" ball kit must be selected to mount the facial lateral to the headrest bracket. Also select a facial pad(s) with the facial lateral(s).

Facial Lateral for Comfort Series.....\$294.25 ea.  
 Right Side (CPS9-R-HW)  
 Part: INDPART2647. HCPCS: E1028

Facial Laterals for Comfort Series ..... \$347.75 ea.  
 Left and Right Side (CPS9-B-HW)  
 Part: INDPART2649. HCPCS: E1028

Facial Lateral for Comfort Series.....\$294.25 ea.  
 Left Side (CPS9-L-HW)  
 Part: INDPART2648. HCPCS: E1028

**Facial Lateral Pads and Accessories**

1" Ball Kit for Facial Laterals (TWB-1BK) (D) ..... \$80.25  
 Part: POS150532  
 Required with facial laterals to mount to comfort headrest brackets.


Facial Pad, 4"x2" (900-Q)..... \$53.50 ea.  
 Part: ST-900-Q. HCPCS: K0108

Egg Switch Mount (ESM) (A) ..... \$52.50 ea.  
 Part: INDPART2815

Facial Pad, 6"x2" (920-Q)..... \$53.50 ea.  
 Part: ST-920-Q. HCPCS: K0108

Figure 8 Bracket with 13" Rod..... \$267.50 ea.  
 and Egg Switch Mount (B)  
 Part: ST-F8675-13-Q

Swing-Away Joystick Mount,..... \$250 ea.  
 18" Rod (C)  
 Part: ST-SUS18-Q  
 Allows mounting of facial pads and attaches to the headrest

**24.7. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Tube Style Arms (ARMS-ES-L) .....\$183.75 ea.  
 Part: POS124064. HCPCS: E1028

Right, Swing-Away, Tube Style Arms (ARMS-ES-R).....\$183.75 ea.  
 Part: POS124066. HCPCS: E1028

**2. Elbow Support Pad**

3"Wx3"H Elbow Pad (ARMS-ES-3X3) .....\$52.50 ea.  
 Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4).....\$52.50 ea.  
 Part: POS124067. HCPCS: K0108

**24.8. Stealth Foot Boxes**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Stealth Small Foot Boxes, 5"Wx8"Dx7"H.....\$498.75  
 (SPLE107) (pair)  
 Part: INDPART2786. HCPCS: E0954^

Stealth Large Foot Boxes, 7"Wx12"Dx11"H.....\$498.75  
 (SPLE109) (pair)  
 Part: INDPART2788. HCPCS: E0954^

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H ....\$498.75  
 (SPLE108) (pair)  
 Part: INDPART2787. HCPCS: E0954^

**24.9. Stealth Padded Shoe Holder**

Stealth Small Padded Shoe Holder, .....\$266  
 3.75"Wx6.25"D (SPLE111) (pair)  
 Part: INDPART2783. HCPCS: E0951+E0952

Stealth X-Large Padded Shoe Holder, .....\$266  
 5"Wx10.75"D (SPLE114) (pair)  
 Part: INDPART2785. HCPCS: E0951+E0952

Stealth Medium Padded Shoe Holder, .....\$266  
 4"Wx8.25"D (SPLE112) (pair)  
 Part: INDPART2784. HCPCS: E0951+E0952



## 25. Q-LOGIC 3 MULTIPLIER HARNESS

**REQUIRED.** Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

No Additional Multiplier (NA)..... No Charge  
Part: Configuration selection

Secondary Multiplier (SM)..... \$221.21  
Part: HAR144004

Special Order Configuration (SOF)  
Part: Configuration selection

|             |                         |
|-------------|-------------------------|
| Arm Mounts  | Joystick (Hand Control) |
|             | Enhanced Display        |
|             | iAccess                 |
|             | USB Charger             |
|             | SAJ                     |
| Back Mounts | Attendant Control       |
|             | SCIM                    |
|             | Sip n Puff              |
|             | ECU                     |

|                             |                         |                            |         |       |       |        |
|-----------------------------|-------------------------|----------------------------|---------|-------|-------|--------|
| NA No additional Multiplier | SM Secondary Multiplier | SOF Special Order Required |         | 1 Arm | 2 Arm | 3+ Arm |
|                             |                         |                            | 0 Back  | NA    | SM    | SOF    |
|                             |                         |                            | 1 Back  | NA    | SM    | SOF    |
|                             |                         |                            | 2 Back  | SM    | SM    | SOF    |
|                             |                         |                            | 3+ Back | SOF   | SOF   | SOF    |

Q-Logic 3 Multiplier with Synergy Seat Chart

## 26. REFERENCE ORDER FORMS

**Optional.** Select below to note that additional order forms will accompany this order form. This will also be a link in interactive forms.

[Alternative Controls & Electronics Order Form](#)

[Positioning Components Order Form](#)

### Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2019 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.



youtube.com/  
QuantumRehab



facebook.com/  
QuantumRehab



twitter.com/  
quantumrehab



instagram.com/  
quantum\_rehab