

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com**Q6 Edge® Z Order Form with TRU-Balance® 3 Positioning**

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

Q6 edge[®] Z**INTRODUCTION**


This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!**PATIENT INFORMATION**

Weight: _____ Height: _____



A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors, 8 amp off board charger, and switch operated fender lights. Dimensions of base without legrests: 26.25"W x 35.5"L.

Quantum Q6 Edge Z 3S-SS.....\$8,290

Part: Q6 EDGE Z 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

Quantum Q6 Edge Z 3MP-SS.....\$8,290

Part: Q6 EDGE Z 3MP-SS. HCPCS: K0861

Multiple power option.

Quantum Q6 Edge Z 3SP-SS.....\$8,290

Part: Q6 EDGE Z 3SP-SS. HCPCS: K0856

Single power option.

Account Number: _____ Date: _____ Quote Name/Marked for: _____

2. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

 Red Corvette Part: DGN1610108-PRD1	 Night Train Part: DGN1610108-PBK3	 Evolution Orange Part: DGN1610108-POR7	 Green Machine Part: DGN1610108-PGN11
 Electric Blue Part: DGN1610108-PBU9	 Pretty in Pink Part: DGN1610108-PPK3	 Grey Street Part: DGN1610108-PSL3	
 Back in Black Part: DGN1610108-PBK1	 Deep Purple Part: DGN1610108-PPR5	 Midnight Blue Part: DGN1610108-PBU10	

3. QUANTUM BACKUP CAMERA

Optional. Requires Q-Logic electronics. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368

Part: MEC1910365

Available for Q-Logic 3 electronics only. NE/NE+ coming soon.

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket..... \$21

Part: MEC1910373

Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$27

Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket\$105

Part: MEC2106123

Backup Camera NE/NE+ Joystick Mount Bracket\$105

Part: MEC2106132



4. DRIVE MOTORS

REQUIRED.

6 mph Drive Motors Standard

Part: MOT150186

5. DRIVE WHEELS

REQUIRED. Select a type of drive wheels. Knobby tires increase overall base width by 1.5".

5.1. Drive Wheels

14" Run Flat Tires Standard

Part: PTOINDV3510

4"W Flat Free, Knobby, Extended Traction Tires \$290 ea.

Part: WHL121698 x2

14" Run Flat Tires for Drive Wheel Accents Standard

Part: WHL1703421 x2

4"W Pneumatic, Knobby, Extended Traction Tires \$290 ea.

Part: WHL153737 x2

5.2. Drive Wheel Accent Color

Optional. Only available WHL1703421 above. Sold in Pairs (parts below will be quantity of 2).

Back in Black Accent Color..... \$63

Part: XRF000179

Deep Purple Accent Color..... \$63

Part: XRF000184

Electric Blue Accent Color \$63

Part: XRF000188

Red Corvette Accent Color \$63

Part: XRF000180

Pretty in Pink Accent Color \$63

Part: XRF000185

Green Machine Accent Color \$63

Part: XRF000190

Night Train Accent Color \$63

Part: XRF000181

Evolution Orange Accent Color \$63

Part: XRF000186

Grey Street Accent Color..... \$63

Part: XRF000182

Midnight Blue Accent Color \$63

Part: XRF000187

6. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor 2 Battery \$341 ea.

Part: BAT2004376. HCPCS: E2361

Group 24, 70 Amp Gel Cell, MK Battery \$520 ea.

Part: BATGEL1009. HCPCS: E2363

7. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Forms](#).

7.1a. Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE+ 6-Key JoystickStandard

Part: CTLDC1560

Control up to 2 actuators through joystick

No Joystick, Setup for NE+ Joystick No Charge

Part: XRF000009

NE+ 8-Key Joystick\$420

Part: CTLDC1561

Control up to 2 actuators and Seat Mounted LED Light Package through joystick.



NE+ 6 Key Joystick

7.1b. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module\$730

Part: XRF000002. HCPCS: E2377

Q-Logic 3 Expandable Controller Harness.....\$447

Part: XRF000264. HCPCS: E2313

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick\$840

Part: CTL167724

Q-Logic 3e Joystick\$735

Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.

Q-Logic 3 EX Joystick for Full Lights\$840

Part: CTL167725

No Joystick, Setup for Q-Logic

Part: XRF000010

Requires selection of Enhanced Display.



Q-Logic 3 EX Joystick

7.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Reference colors shown on page 2.

Black, Soft TouchStandard

Part: DGN124307

Also matches Night Train shrouds.

Red Corvette\$63

Part: DGN124336-PRD1

Electric Blue\$63

Part: DGN124336-PBU9

Back in Black\$63

Part: DGN124336-PBK1

Pretty in Pink\$63

Part: DGN124336-PPK3

Evolution Orange.....\$63

Part: DGN124336-POR7

Deep Purple.....\$63

Part: DGN124336-PPR5

Grey Street.....\$63

Part: DGN124336-PSL3

Midnight Blue\$63

Part: DGN124336-PBU10

Green Machine\$63

Part: DGN124336-PGN11

8. JOYSTICK HANDLES

REQUIRED. The NE Plus joystick is standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B)\$90

Part: ACCASMB1860. HCPCS: E2323

Chin Cup (C).....\$90

Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

T-Bar (D)\$126

Part: ACCASMB1858. HCPCS: E2323

Medium Ball Handle (E).....\$69

Part: ACCASMB1859. HCPCS: E2323

Soft Ball Handle (F).....\$80

Part: ACCASMB1856. HCPCS: E2323

Extended I-Handle (G)\$90

Part: ACCASMB1857. HCPCS: E2323

Bodypoint 3" Goal Post Handle\$104


Part: PLSKNCP1057. HCPCS: E2323

Bodypoint 4" Goal Post Handle\$104

Part: PLSKNCP1058. HCPCS: E2323

Q-Logic 3 Carrot Handle (A) No Charge

Part: DGN124308



Joystick Handles

9. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

9.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket for Q-Logic 3Standard

Part: FRMASMB10617

Fixed, In-line, Bracket for NE/NE+Standard

Part: FRMASMB12744

Left Side Mount

Part: PTOINDV1225

Right Side Mount

Part: PTOINDV1226

9.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3\$294

Part: ACC1908434. HCPCS: E1028

Swing-away, NE/NE+.....\$294

Part: ACC1908437. HCPCS: E1028

9.1c. Basic Swing Away Joystick Mount

Basic Swing Away, NE/NE+, Left Side\$245

Part: ACC2201047. HCPCS: E1028

Basic Swing Away, NE/NE+, Right Side.....\$245


Part: ACC2201048. HCPCS: E1028

Basic Swing Away, Q-Logic, Left Side\$245

Part: ACC2201041. HCPCS: E1028

Basic Swing Away, Q-Logic, Right Side\$245

Part: ACC2201042. HCPCS: E1028



Swing Away Joystick

Joystick mounting brackets continued on following page.

9. JOYSTICK MOUNTING BRACKETS CONTINUED**9.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount.....\$368

Part: XRF000003. HCPCS: E1028

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount\$126


Part: XRF000004

Allows height adjustment of joystick.

Stealth Swing-Away/Flip Down Bracket\$250

Part: XRF000005. HCPCS: E1028

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.



Retract4



Multi-Axis Joystick Mount

Specialty Mount Side
☐ Left Side Mount
Part: PTOINDV1225

☐ Right Side Mount
Part: PTOINDV1226
10. TRU-BALANCE® 3 SEATING - QUANTUM INTELLIGENT SYSTEMS™**REQUIRED.** Maximum seat size is 24"W x 24"D.**10.1. TRU-Balance 3 Configuration**

TRU-Balance 3 Static Seat\$2,410

Part: PTOINDV3442

20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU-Balance 3 Power Tilt Only\$7,944

Part: PTOINDV3443. HCPCS: E1002

50° Center of Gravity (CG) power tilt seating system.

TRU-Balance 3 Power Recline Only\$7,975

Part: PTOINDV3444. HCPCS: E1004

175° of power recline.

TRU-Balance 3 Power Tilt & Recline.....\$12,490

Part: PTOINDV3446. HCPCS: E1007

50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

10.2. 12" Power Lift**Optional.** Select one option below.

TRU-Balance 3, 12" Power Lift with iLevel®\$2,730

Part: PTOINDV3773. HCPCS: E2300

12" of scissor lift with iLevel® System. Not available with high speed motors or manual recline.

**10.3. Occupied Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. Only available with 16"-20" width and depth packages. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

WC-19 Occupied Transit Option.....\$289

Part: KIT152047

10.4. Seat Width**REQUIRED.** Maximum seat width is 24". Part number is SETWDTH plus number following # below.

12" W	13" W	14" W	15" W	16" W	17" W	18" W	19" W	20" W	21" W	22" W	23" W	24" W
(#1002)	(#1003)	(#1004)	(#1005)	(#1006)	(#1007)	(#1008)	(#1009)	(#1010)	(#1011)	(#1012)	(#1013)	(#1014)

10.5. Seat Width Package**REQUIRED.**

12"-16" Width Package

Part: XRF000006

14"-16" Width Package

Part: XRF000007

16"-20" Width Package

Part: XRF000008

20"-24" Width Package

Part: XRF000158

10.6. Seat Depth**REQUIRED.** Maximum seat depth is 24". Part number is SETDPTH plus number following # below.

12" D	13" D	14" D	15" D	16" D	17" D	18" D	19" D	20" D	21" D	22" D	23" D	24" D
(#1002)	(#1003)	(#1004)	(#1005)	(#1006)	(#1007)	(#1008)	(#1009)	(#1010)	(#1011)	(#1012)	(#1013)	(#1014)

10.7. Seat Depth Package**REQUIRED.**

12"-16" Depth Package

Part: XRF000011

14"-16" Depth Package

Part: XRF000012

16"-20" Depth Package

Part: XRF000013

20"-24" Depth Package

Part: XRF000159

10.8. Seat to Floor Height**REQUIRED.**

Lowest Seat to Floor

Part: PTOINDV3558

Middle Seat to Floor

Part: PTOINDV3559

Highest Seat to Floor

Part: PTOINDV3560

Low	Med	High
Static Seat, Recline		
16.5"	17.5"	18.5"
Tilt, Tilt & Recline, Lift		
18.5"	19.5"	20.5"
TRU-Balance 3 Seat to Floor Heights		

TRU-Balance 3 Seating continued on following pages.


10. TRU-BALANCE 3 SEATING CONTINUED**10.9. Matching Seat Cushion**

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, and special order cushions can be found on the respective complete Stealth Cushion Order Forms. This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

Stealth TRU-Comfort 2 SPP Cushion


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-TC2SPPXXYY. HCPCS: E2607/E2608

TRU-Comfort 2 SPP, Matching Seat Size
Reference part: XRF000116. MSRP: \$420/\$525

**Stealth Solution SPP Cushion**


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Includes urethane liner.
Part series: ST-SOLSPPXXYY. HCPCS: E2607/E2608

Solution SPP, Matching Seat Size
Reference part: XRF000113. MSRP: \$394/\$452

**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Includes urethane liner.
Part series: ST-SPCFSPXXYY. HCPCS: E2607/E2608

Spectrum Foam SPP, Matching Seat Size
Reference part: XRF000112. MSRP: \$483/\$557

**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-PREMPXXYY. HCPCS: E2605/E2606

Premiere P, Matching Seat Size
Reference part: XRF000122. MSRP: \$341/\$499

**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ZENSPXXYY. HCPCS: E2603/E2604

Zen SP, Matching Seat Size
Reference part: XRF000115. MSRP: \$331/\$446

**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Includes urethane liner.
Part series: ST-SPCGSPPXXYY. HCPCS: E2607/E2608

Spectrum Gel SPP, Matching Seat Size
Reference part: XRF000117. MSRP: \$483/\$557

**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ESSSPXXYY. HCPCS: E2607/E2608

Essence SPP, No Wedge, Matching Seat Size
Reference part: XRF000118. MSRP: \$420/\$525

Essence SPP, Small Wedge, Matching Seat Size
Reference part: XRF000119. MSRP: \$457/\$562


Essence SPP, Medium Wedge, Matching Seat Size
Reference part: XRF000120. MSRP: \$488/\$593

Essence SPP, Large Wedge, Matching Seat Size
Reference part: XRF000121. MSRP: \$520/\$625

**Stealth Glacial SP Cushion**


Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-GLSPXXYY. HCPCS: E2603/E2604

Glacial SP, Matching Seat Size
Reference part: XRF000114. MSRP: \$347/\$462

**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.
Part series: ST-SIMGXXYY. HCPCS: E2601/E2602

Simplicity G, Matching Seat Size
Reference part: XRF000111. MSRP: \$116/\$173

**10.10. Back Type Selection**

REQUIRED. Match choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back

Part: XRF000014
12" and 24" back width will use the TRU-Balance 3 Back assembly.
Select a back cushion configuration in the subsequent section.

Back Cane Option

Part: XRF000015
Available only with TRU-Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back..... \$184

Part: XRF000016
12" and 24" back width will use the TRU-Balance 3 Back assembly. Only available with 2-post or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

10.11. Back Cane Selection (Back Cane Option Only)

REQUIRED when Back Cane Option is selected above. Select type and height of back canes.

1. Back Cane Type

Angle Adjustable Back Cane..... No Charge
Part: MEC140266

Flip Forward Back Cane \$53
Part: XRF000273

2" Endomorph (1" per side) Back Cane \$74
Part: XRF000269

4" Endomorph (2" per side) Back Cane..... \$74
Part: XRF000271

2. Back Cane Height

16"H Cane
Part: FRMASMB15834

17"H Cane
Part: FRMASMB15835

20"H Cane
Part: FRMASMB15838

21"H Cane
Part: FRMASMB15839

24"H Cane
Part: FRMASMB15842

25"H Cane
Part: FRMASMB15843

10.12. Back Cushion Configuration

REQUIRED. Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

1a. Sport Back/Omit/After-Market Interface

Sport Back..... No Charge
Part: XRF000017
Complete following width package and size.

Omit Back Pan/Cushion..... No Charge
Part: PTOINDV1287
Select a following width package.

After-market Back Interface \$509.85
Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm trough. Does not include brackets to mount after-market backs.

TRU-Balance 3 Seating continued on following pages.

10. TRU-BALANCE 3 SEATING CONTINUED**1b. TRU-Comfort Back**

TRU-Comfort Back (Static/Flip Forward), 12" Width \$709
 Part: SET123622. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 14" to 20" Width \$709
 Part: TC4V10W21W. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width \$709
 Part: TC4V10W21WBC. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$814
 Part: TC4V22W32W. HCPCS: E2621
 Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 22" Width \$814
 Part: TC4V22W32WBC. HCPCS: E2621
 Complete following width package and size.

1c. ADI Back, AL Series, 4-Point, Active Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 14" Wide x 13" or 14" Overall Height \$592.25
 Part: ALBR14LA4. HCPCS: E2613

Active Contour, 14" Wide x 16" or 17" Overall Height \$592.25
 Part: ALBR14MA4. HCPCS: E2613

Active Contour, 16" Wide x 19" or 20" Overall Height \$592.25
 Part: ALBR16SA4. HCPCS: E2613

Active Contour, 16" Wide x 21" or 22" Overall Height \$592.25
 Part: ALBR16TA4. HCPCS: E2613

Active Contour, 16" Wide x 23" or 24" Overall Height \$592.25
 Part: ALBR16XA4. HCPCS: E2613

Active Contour, 18" Wide x 19" or 20" Overall Height \$592.25
 Part: ALBR18SA4. HCPCS: E2613

Active Contour, 18" Wide x 21" or 22" Overall Height \$592.25
 Part: ALBR18TA4. HCPCS: E2613

Active Contour, 18" Wide x 23" or 24" Overall Height \$592.25
 Part: ALBR18XA4. HCPCS: E2613

Active Contour, 20" Wide x 19" or 20" Overall Height \$592.25
 Part: ALBR20SA4. HCPCS: E2613

Active Contour, 20" Wide x 21" or 22" Overall Height \$592.25
 Part: ALBR20TA4. HCPCS: E2613

Active Contour, 20" Wide x 23" or 24" Overall Height \$592.25
 Part: ALBR20XA4. HCPCS: E2613

1d. ADI Back, AL Series, 4-Point, Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 14" Wide x 13" or 14" Overall Height \$700.40
 Part: ALBR14LD4. HCPCS: E2615

Deep Contour, 14" Wide x 16" or 17" Overall Height \$700.40
 Part: ALBR14MD4. HCPCS: E2615

Deep Contour, 16" Wide x 19" or 20" Overall Height \$700.40
 Part: ALBR16SD4. HCPCS: E2615

Deep Contour, 16" Wide x 21" or 22" Overall Height \$700.40
 Part: ALBR16TD4. HCPCS: E2615

Deep Contour, 16" Wide x 23" or 24" Overall Height \$700.40
 Part: ALBR16XD4. HCPCS: E2615

Deep Contour, 18" Wide x 19" or 20" Overall Height \$700.40
 Part: ALBR18SD4. HCPCS: E2615

Deep Contour, 18" Wide x 21" or 22" Overall Height \$700.40
 Part: ALBR18TD4. HCPCS: E2615

Deep Contour, 20" Wide x 19" or 20" Overall Height \$700.40
 Part: ALBR20SD4. HCPCS: E2615

Deep Contour, 20" Wide x 21" or 22" Overall Height \$700.40
 Part: ALBR20TD4. HCPCS: E2615

Deep Contour, 20" Wide x 23" or 24" Overall Height \$700.40
 Part: ALBR20XD4. HCPCS: E2615

1e. ADI Back, AL Series, 4-Point, Extra Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 16" Wide x 21" or 22" Overall Height .. \$782.80
 Part: ALBR16T64. HCPCS: E2615

Extra Deep Contour, 16" Wide x 23" or 24" Overall Height .. \$782.80
 Part: ALBR16X64. HCPCS: E2615

Extra Deep Contour, 18" Wide x 19" or 20" Overall Height .. \$782.80
 Part: ALBR18S64. HCPCS: E2615

Extra Deep Contour, 18" Wide x 21" or 22" Overall Height .. \$782.80
 Part: ALBR18T64. HCPCS: E2615

Extra Deep Contour, 18" Wide x 23" or 24" Overall Height .. \$782.80
 Part: ALBR18X64. HCPCS: E2615

Extra Deep Contour, 20" Wide x 19" or 20" Overall Height .. \$782.80
 Part: ALBR20S64. HCPCS: E2615

Extra Deep Contour, 20" Wide x 21" or 22" Overall Height .. \$782.80
 Part: ALBR20T64. HCPCS: E2615

Extra Deep Contour, 20" Wide x 23" or 24" Overall Height .. \$782.80
 Part: ALBR20X64. HCPCS: E2615

1f. Q-Back 2.0, AL Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$540.75
 Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$540.75
 Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$540.75
 Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$540.75
 Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$540.75
 Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$669.50
 Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$540.75
 Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$540.75
 Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$540.75
 Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$669.50
 Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$540.75
 Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$540.75
 Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$540.75
 Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$669.50
 Part: QB2AL20AXX. HCPCS: E2613

TRU-Balance 3 Seating section continued on following pages.

10. TRU-BALANCE 3 SEATING CONTINUED**1g. Q-Back 2.0, AL Series, Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$612.85
Part: QB2AL14DS. HCPCS: E2620Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$612.85
Part: QB2AL14DT. HCPCS: E2620Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$612.85
Part: QB2AL16DS. HCPCS: E2620Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$612.85
Part: QB2AL16DT. HCPCS: E2620Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$612.85
Part: QB2AL16DX. HCPCS: E2620Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$669.50
Part: QB2AL16DXX. HCPCS: E2620Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$612.85
Part: QB2AL18DS. HCPCS: E2620Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$612.85
Part: QB2AL18DT. HCPCS: E2620Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$612.85
Part: QB2AL18DX. HCPCS: E2620Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$669.50
Part: QB2AL18DXX. HCPCS: E2620Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$612.85
Part: QB2AL20DS. HCPCS: E2620Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$612.85
Part: QB2AL20DT. HCPCS: E2620Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$612.85
Part: QB2AL20DX. HCPCS: E2620Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$669.50
Part: QB2AL20DXX. HCPCS: E2620**1h. Q-Back 2.0, AL Series, Extra Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline . \$669.50
Part: QB2AL14ES. HCPCS: E2620Extra Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline.. \$669.50
Part: QB2AL14ET. HCPCS: E2620Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline.. \$669.50
Part: QB2AL16ES. HCPCS: E2620Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline.. \$669.50
Part: QB2AL16ET. HCPCS: E2620Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline.. \$669.50
Part: QB2AL16EX. HCPCS: E2620Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline.. \$669.50
Part: QB2AL16EXX. HCPCS: E2620Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline . \$669.50
Part: QB2AL18ES. HCPCS: E2620Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline.. \$669.50
Part: QB2AL18ET. HCPCS: E2620Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline . \$669.50
Part: QB2AL18EX. HCPCS: E2620Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline . \$669.50
Part: QB2AL18EXX. HCPCS: E2620Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline . \$669.50
Part: QB2AL20ES. HCPCS: E2620Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline . \$669.50
Part: QB2AL20ET. HCPCS: E2620Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline . \$669.50
Part: QB2AL20EX. HCPCS: E2620Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline . \$669.50
Part: QB2AL20EXX. HCPCS: E2620**1i. Q-Back 2.0, AFT Series, Active Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$540.75
Part: QB2AFT16AS. HCPCS: E2613AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$540.75
Part: QB2AFT18AS. HCPCS: E2613**1j. Q-Back 2.0, AFT Series, Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline... \$612.85
Part: QB2AFT16DS. HCPCS: E2620AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline... \$612.85
Part: QB2AFT16DT. HCPCS: E2620AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline .. \$612.85
Part: QB2AFT18DS. HCPCS: E2620AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline... \$612.85
Part: QB2AFT18DT. HCPCS: E2620AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline .. \$612.85
Part: QB2AFT20DT. HCPCS: E2620**1k. Q-Back 2.0, AFT Series, Extra Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$669.50
Part: QB2AFT16EXX. HCPCS: E2620**10.13. TRU-Balance 3 Back Shroud****REQUIRED.** The back shroud is a feature of the TRU-Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.Red Corvette
Part: XRF000019Grey Street
Part: XRF000023Evolution Orange
Part: XRF000028Green Machine
Part: XRF000033Back in Black
Part: XRF000021Deep Purple
Part: XRF000025Midnight Blue
Part: XRF000029Night Train
Part: XRF000022Pretty in Pink
Part: XRF000027Electric Blue
Part: XRF000030

TRU-Balance 3 Seating continued on following page.

10. TRU-BALANCE 3 SEATING CONTINUED**10.14. Lateral Positioning Wedge****Optional.**

Lateral Positioning Wedge, pair \$16 ea.

Part: SET123465

Available with Sport and TRU-Comfort 2 backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

10.15. Back Cushion Configuration - Back Width Package**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections

Part: PTOINDV3468

14"W-16"W Package

Part: XRF000037

18"W-22"W Package

Part: XRF000039

10"W-14"W Package

Part: XRF000036

16"W-20"W Package

Part: XRF000038

20"W-24"W Package

Part: XRF000040

10.16. Back Cushion Configuration - Back Width**REQUIRED** except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.

12" W (#3562)	14" W (#3563)	15" W (#3856)	16" W (#3564)	17" W (#3565)	18" W (#3566)	19" W (#3567)	20" W (#3568)	21" W (#3569)	22" W (#3570)	24" W (#3571)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

10.17. Back Cushion Configuration - Overall Back Height

Static Seat and Power Tilt available height range: 14"-25"H. **Odd heights include 4" gap and even heights include 3" gap** between seat pan and bottom of back cushion.

Power Recline available height range: 19"-28"H. **Odd heights include 6" gap and even heights include 7" gap.** See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

14" Height

Part: SETHGHT1027

18" Height

Part: SETHGHT1031

22" Height

Part: SETHGHT1035

26" Height

Part: SETHGHT1039

30" Height

Part: SETHGHT1043

15" Height

Part: SETHGHT1028

19" Height

Part: SETHGHT1032

23" Height

Part: SETHGHT1036

27" Height

Part: SETHGHT1040

31" Height

Part: SETHGHT1102

16" Height

Part: SETHGHT1029

20" Height

Part: SETHGHT1033

24" Height

Part: SETHGHT1037

28" Height

Part: SETHGHT1041

32" Height

Part: SETHGHT1103

17" Height

Part: SETHGHT1030

21" Height

Part: SETHGHT1034

25" Height

Part: SETHGHT1038

29" Height

Part: SETHGHT1042

		Tilt or Static Width												Recline Width											
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"		
Overall Height	14"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"		
	15"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"		
	16"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	16"		
	17"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	17"		
	18"	-	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	18"		
	19"	-	S	Q	T	S	Q	T	S	Q	T	S	Q	T	S	S	S	T	S	T	S	T	S	19"	
	20"	-	S	Q	T	S	Q	T	S	Q	T	S	Q	T	S	S	S	T	S	T	S	T	S	20"	
	21"	-	S	Q	T	S	Q	T	S	Q	T	S	Q	T	S	-	S	T	S	T	S	T	S	21"	
	22"	-	S	Q	T	S	Q	T	S	Q	T	S	Q	T	S	-	S	Q	T	S	Q	T	S	22"	
	23"	-	S	T	S	T	S	Q	T	S	T	S	Q	T	S	-	S	Q	T	S	Q	T	S	23"	
	24"	-	S	T	S	T	S	Q	T	S	T	S	Q	T	S	-	S	Q	T	S	Q	T	S	24"	
	25"	-	S	T	S	T	S	Q	T	S	T	S	Q	T	S	-	S	Q	T	S	Q	T	S	25"	
	26"	-	S	T	S	T	S	Q	T	S	T	S	Q	T	S	-	S	T	S	T	S	Q	T	S	26"
	27"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	Q	T	S	27"
	28"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	Q	T	S	28"
	29"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	Q	T	S	29"
	30"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	30"	
31"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	31"		
32"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	32"		

S

 Sport Back

Q

 Q-Back

T

 TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

S Sport Back
 Q Q-Back
 T TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

10.18. Power Tilt Privacy Flap**Optional.** A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap

Part: XRF000351

\$48

10.19. Back Height Extension**Optional.** Not available with back cane option, 12" back width, or 24" back width. Adds 2" of seat back pan height. 25"H selection includes a back height extension.

Back Height Extension

Part: XRF000041

\$48

11. ARMRESTS

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

11.1a. Flip-Back, 8"-16" H Armrests

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

☐ ☐ Flip-Back, 8"-16"H Armrests Standard	☐ ☐ Flip-Back, 8"-16"H Armrests for Static Seat \$289
☐ ☐ Left part: ARM1806662. Right part: ARM1806663	☐ ☐ Left part: XRF000320. Right part: XRF000321
Not available with TRU-Balance 3 Static Seat.	

Armrest Support

Optional. Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

☐ ☐ Flip-Back Armrest Downpost	\$92 ea.
☐ ☐ Left part: XRF000265. Right part: XRF000266	

11.1b. Drop-In, Height Adjustable, Removable, T-Arm Armrests

Only available with TRU-Balance 3 Static Seat. Includes flat armpads.

☐ ☐ Full Length, Drop-In Armrests..... No Charge	☐ ☐ Desk Length, Drop-In Armrests No Charge
☐ ☐ Left part: ARM1803253. Right part: ARM156742	☐ ☐ Left part: ARM1803254. Right part: ARM156743

11.1c. Two-Post, Flip-Back, Height Adjustable, Removable Armrests

Select a left and right armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes. Not available with a depth less than 16". Required with flip forward back and lap trays.

1a. Two-Post Armrest for Power Seating

☐ ☐ Full Length, Two-Post Armrest for Power Seating No Charge	☐ ☐ Desk Length, Two-Post Armrest for Power Seating..... No Charge
☐ ☐ Left part: XRF000047. Right part: XRF000048	☐ ☐ Left part: XRF000049. Right part: XRF000050

1b. Two-Post Armrest for Static Seating

☐ ☐ Full Length, Two-Post Armrest for Static Seating.....\$158 ea.	☐ ☐ Desk Length, Two-Post Armrest for Static Seating\$158 ea.
☐ ☐ Left part: XRF000051. Right part: XRF000052. HCPCS: E0973	☐ ☐ Left part: XRF000053. Right part: XRF000054. HCPCS: E0973

11.1d. Heavy Duty, Single Post, Height Adjustable, Removable Armrests

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration.

8"-12"H, Single Post, Height Adjustable, Removable, \$158	10"-14"H, Single Post, Height Adjustable, Removable,..... \$158
Heavy Duty (Sold in pairs)	Heavy Duty (Sold in pairs)
Part: XRF000055. HCPCS w/ Static Seat: E0973	Part: XRF000056. HCPCS w/ Static Seat: E0973

11.2. Armrest Height

REQUIRED. Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

2. Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

11.3. Armrest Size

REQUIRED. Select a Left (L) and Right (R) armrest size. This also affects the armpad selection below.

☐ ☐ 10" Desk Length L part: XRF000057. R part: XRF000123	☐ ☐ 14" Full Length L part: XRF000058. R part: XRF000124	☐ ☐ 17" Length L part: XRF000059. R part: XRF000125 Only available with Flip-Back, 8"-16" H Armrests and straight arm pads.
---	---	---

12. ARMPADS

REQUIRED. Select to omit armpads, select a left and right foam armpads, or a gel armpad configuration. Please note that waterfall and gel armpads may cause interference with the seat back.

12.1a. Omit Armpad

Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection..... No Charge	Interface for Aftermarket Pad No Charge
Part: XRF000060	Part: XRF000308
Not available with Flip-back, 8"-16" H armrests.	Only available with Flip-back, 8"-16" armrests.

12.1b. Foam Armpads

☐ ☐ Straight, Armpad Standard	☐ ☐ Waterfall, Armpad No Charge
☐ ☐ Left part: XRF000061. Right part: XRF000126	☐ ☐ Left part: XRF000062. Right part: XRF000127
	Not available with Flip-back, 8"-16" armrests.

Armpads section continued on following page.

12. ARMPADS CONTINUED**12.1c. Gel Armpads**

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip Back ArmrestsFalcon Gel, Full Length Left and Right\$195
Part: ARM1806807Falcon Gel, Desk Length Left, Full Length Right.....\$171
Part: ARM1806810Falcon Gel, Full Length Left, Desk Length Right.....\$171
Part: ARM1806813Falcon Gel, Desk Length Left and Right\$147
Part: ARM1806806**b. Falcon Gel Armpads for Tube Style Armrests**Falcon Gel, Full Length Left and Right\$311.22
Part: ARM136438Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
Part: ARM136439Falcon Gel, Full Length Left, Desk Length Right.....\$287.93
Part: ARM136442Falcon Gel, Desk Length Left and Right \$265.02
Part: ARM136437**13. TRU-BALANCE 3 ARMREST BRACKETS****REQUIRED.**Standard Fixed ReceiverStandard
Part: XRF000063Stealth Height Adjustable Receiver.....\$225
Part: XRF000064. HCPCS: K0108**14. TRU-BALANCE 3 HEADRESTS****Optional.** Select a headrest pad and bracket.**14.1a. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$170.25
Part: ST-CP08CCFPM-Q. HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) \$170.25
Part: ST-CP14CCFPM-Q. HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)..... \$170.25
Part: ST-CP10CCFPM-Q. HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM) \$170.25
Part: ST-CP06CCFPM-Q. HCPCS: E0955**2. Headrest Brackets**Comfort Plus Headrest Mounting Bracket, UniLink..... \$296.50
Part: ST-UHB2D33200-Q. HCPCS: E1028**14.1b. Headrests, Comfort Plus****1. Headrest Pads**Stealth 8" Comfort Plus (CP450)..... \$170.25
Part: POS148484. HCPCS: E0955Stealth 14" Comfort Plus (CP150) \$170.25
Part: INDPART2327. HCPCS: E0955Stealth 10" Comfort Plus (CP250)..... \$170.25
Part: SETHEAD1057. HCPCS: E0955Pediatric Comfort Plus (6"Wx4.5"H) (CPN350) \$170.25
Part: SETHEAD1069. HCPCS: E0955**2. Headrest Brackets**Multi-Axis Removable Headrest Bracket (HMO475P)\$165.32
Part: HDWASMB1084. HCPCS: E1028**15. LAP BELTS****REQUIRED.** Select a lap belt below.**15.1a. Standard Lap Belts**38" Lap Belt
Part: ACCBELT100650" Lap Belt
Part: ACCBELT100060" Lap Belt
Part: ACCBELT100270" Lap Belt
Part: ACCBELT100380" Lap Belt
Part: ACCBELT100490" Lap Belt
Part: ACCBELT1001**15.1b. Optional Lap Belts**Stealth Retractable Lap Belt (IFRB5350)..... \$172.53
Part: ST-IFRB5350Stealth Positioning 60" Single Pull, Push Button, \$47.38
Extra Large, Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-QStealth Positioning 60" Single Pull, Push Button, Medium..... \$47.38
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q60" Sure-Lok Lap Belt (Non-retractable)..... \$120
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016Stealth Positioning 60" Single Pull, Push Button, Large \$47.38
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q

16. POWER LEGREST CONTROL

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (adds 2 actuator control requirement)
Part: XRF000108

Articulating Foot Platform (adds 1 actuator control requirement)
Part: XRF000110

Combined Legrests (adds 2 actuator control requirement)
Part: XRF000109

17. BASE MOUNTED FOOT PLATFORM

REQUIRED. Base mounted foot platforms are not available with center mount foot platforms. May interfere with seat mounted legrests.

Base Mounted Foot Platform Standard
Part: FRMASMB16077

Bell Shaped, Base Mounted Foot Platform \$247
Part: RIGASMB7120041

Omit Base Mounted Foot Platform No Charge
Part: PTOINDV3674

18. LEGRESTS

Optional. Please note the set lower leg length may vary based on your seat type and seat to floor height selections. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

18.1a. Basic Center Mount Foot Platform**1. Basic Center Mount Foot Platform Assembly**

Basic Center Mount No Charge
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates \$180
Part: XRF000312. HCPCS: K0040

Long Extension, 5"x9" Footplates \$180
Part: XRF000314. HCPCS: K0040

Medium Extension, 5"x9" Footplates \$180
Part: XRF000313. HCPCS: K0040

3. Basic Center Mount Calf Pads

Small Calf Pads, Basic Center Mount \$158
Part: RIG2104369

4. Basic Center Mount Additional Accessory

Basic Center Mount Footplate Connector \$16
Part: RIG2104240

18.1b. Center Mount Foot Platform with Angle Adjustable Footplates

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform Assembly

Requires a selection of calf pad or calf panel.

Power Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040

Quick Adjustable Center Mount Foot Platform \$263
Part: XRF000311. HCPCS: K0040

2. AFP Lower Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility. Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates
Part: RIG152577

Medium Extension, 5"x11" Footplates
Part: RIG150212

Pediatric Extension, 5"x8" Footplates
Part: RIG152578

Long Extension, 5"x6" Footplates
Part: RIG151505

Medium Extension, 5"x6" Footplates
Part: RIG150206

Long Extension, 5"x8" Footplates
Part: RIG151506

Medium Extension, 5"x8" Footplates
Part: RIG150203

Long Extension, 5"x11" Footplates
Part: RIG151508

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019

Small Calf Pads, Manually Adjustable Center Mount Only, \$158
6"Wx6"H (pair)
Part: XRF000156

Large Calf Pads, Manually Adjustable Center Mount Only, \$158
8" taper to 6"Wx8"H (pair)
Part: XRF000157

3b. Center Mount Calf Panel

Calf pads are standard with a power articulating foot platform.

Small Calf Panel for AFP, 18"Wx8"H \$184
Part: FRMASMB7110007

Small Calf Panel for Manually Adjustable, 18"Wx8"H \$184
Part: FRMASMB7110012

Large Calf Panel for AFP, 24"Wx8"H \$237
Part: RIG123824
Only available on widths of 20" or greater.

Large Calf Panel for Manually Adjustable, 24"Wx8"H \$237
Part: RIG123825
Only available on widths of 20" or greater.

Legrests section continued on following pages.

18. LEGRESTS CONTINUED**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$69
Part: RIG152517Width Extension for 11"D Footplates \$69
Part: RIG1525182" Bell-Shaped Length Extension \$69
Part: RIG152519
Not compatible with width extensions.4" Outside Edge (pair) \$69
Part: RIG152510

Only available with 6"D footplates and 8"D footplate extension.

6" Outside Edge (pair) \$69
Part: RIG152516
Available with 8"D and 11"D footplates and 11"D footplate extension.**5. Additional Accessories**Small Reversed Dartex Calf Pad Cover (pair) \$158
Part: RIG133717Large Reversed Dartex Calf Pad Cover (pair) \$158
Part: RIG133719Small Reversed Dartex Panel Cover (pair) \$210
Part: RIG134673Large Reversed Dartex Panel Cover (pair) \$210
Part: RIG134744Small Foam Calf Pad Cover (pair) \$63.74
Part: ST-CSCS0775-QLarge Foam Calf Pad Cover (pair) \$78.72
Part: ST-CSCS0985-QSmall Foam Panel Cover (pair) \$167.57
Part: ST-CSCP0920-QLarge Foam Panel Cover (pair) \$268.77
Part: ST-CSCP0927-QFootplate Connector \$16
Part: RIGASMB7120023Footplate Riser \$75
Part: XRF000348
Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.Footplate Riser Calf Panel \$50
Part: FRMASMB7110011
Only available and required with Footplate Riser selection and extension lengths of 6" to 9".**6. Stealth Gel Footplate Covers**5"x8" Tapered Gel Footplate Covers \$211.15
Part: ST-FPG5008T-Q. HCPCS: K01085"x11" Tapered Gel Footplate Covers \$231.75
Part: ST-FPG5011T-Q. HCPCS: K01085"x8" Tapered Gel Footplate Covers for Width Extension ... \$231.75
Part: ST-FPG6508T-Q. HCPCS: K01085"x11" Tapered Gel Footplate Cover for Width Extension ... \$257.50
Part: ST-FPG7011T-Q. HCPCS: K0108**18.1c. Power Elevating Legrests**

Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Select legrest length in following section.

1. Power Elevating Legrest Upper Assembly

Includes calf pads.

Power Elevating Legrests, Pair \$2,516
Part: XRF000341. HCPCS: E1010+K0040+E0951

Legrest Length Matrix Center Mount	Seat to Floor Option			Legrest Length Matrix Center Mount	Seat to Floor Option		
	Low	Middle	High		Low	Middle	High
Power ELR: Static Seat, Recline Only				Power ELR: Tilt Only, Tilt and Recline, Any Power Lift			
8.5"-14.5" Lower Length	8.5"-14"	8.5"-14.5"	8.5"-14.5"	8.5"-14.5" Lower Length	8.5"-14.5"	8.5"-14.5"	8.5"-14.5"
14.5"-21" Lower Length	12.5"-14"	12.5"-15"	12.5"-16"	14.5"-21" Lower Length	12.5"-15"	12.5"-16"	12.5"-17"

2. Lower Extension8.5" to 14.5" Lower Length
Part: XRF00034412.5" to 21" Lower Length
Part: XRF000347**3. Angle Adjustable Footplate with Heel Loop and Wheel Bumper**

*The 6"W footplates may cause interference with seat widths less than 16". *The 8"W footplates may cause interference with seat widths less than 19".

4"Wx6"D, Pair
Part: FRMASMB74485"Wx8"D, Pair
Part: FRMASMB74506"Wx11"D¹, Pair
Part: FRMASMB74525"Wx6"D, Pair
Part: FRMASMB74496"Wx8"D¹, Pair
Part: FRMASMB74518"Wx11"D², Pair
Part: FRMASMB7453**4. Power Elevating Legrest Accessories**AEL Knee Adductor, 4"x3" Oval Pad (pair) \$512
Part: INDPART2433. HCPCS: E0953

Legrests section continued on following pages.

18. LEGRESTS CONTINUED**18.1d. Swing-Away Legrests**

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1a. 60°-90° Swing-Away Legrest Upper

60°-90° Swing-Away Legrest, Power Seating..... No Charge
 Pair part: RIG120062. Left: FRMASMB17558.
 Right: FRMASMB17559
 Available with Power Positioning only. Not available with 11"-13" extensions.

60°-90° Swing-Away Legrest, Short, Power Seating.... No Charge
 Pair part: RIG136481. Left: RIG132607.
 Right: RIG132609.
 Available with Power Positioning only. For use only with 11"-13" lower extensions.

60°-90° Swing-Away Legrest, Static Seat \$86 ea.
 Pair part: XRF000082. Left: XRF000083.
 Right: XRF000084
 Available with Static Seat only. Not available with 11"-13" extensions.

60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea.
 Pair part: XRF000085. Left: XRF000086.
 Right: XRF000087.
 Available with Static Seat only. For use only with 11"-13" lower extensions.

1b. 70° Swing-Away Legrest Upper

Only available with occupied transit WC-19 selections. Sold only in pairs. Select lower extension and footplates as pairs.

70° Swing-Away Legrest (pair)..... \$172
 Pair part: FRMASMB7841
 Not available with 11"-13" extensions.

2. Lower Extension

11"-13" Length
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.
 N/A with Multi-axis.

13"-16.5" Length
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.
 N/A with Multi-Axis.

11"-16.5" Length
 Pair part: XRF000091.
 Left: XRF000092.
 Right: XRF000093.
 Multi-Axis Only.

15.5"-19.5" Length
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3a. Extruded Footplate(s)

7.25"Wx6"D Extruded Footplate..... Standard
 Part pair: FRMASMB7847. Left: FRMASMB7861
 Right: FRMASMB7862

8.5"Wx6"D Extruded Footplate Standard
 Part pair: FRMASMB7848. Left: FRMASMB7863
 Right: FRMASMB7864
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s) with Heel Loop(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

4"Wx6"D \$90 ea.
 Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560
 HCPCS: K0040+E0951

5"Wx6"D..... \$90 ea.
 Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561
 HCPCS: K0040+E0951

5"Wx8"D \$90 ea.
 Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562
 HCPCS: K0040+E0951

6"Wx8"D¹ \$90 ea.
 Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563
 HCPCS: K0040+E0951

6"Wx11"D¹ \$90 ea.
 Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564
 HCPCS: K0040+E0951

8"Wx11"D² \$90 ea.
 Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565
 HCPCS: K0040+E0951

3c. Multi-Axis Footplate(s) with Heel Loop(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

4"Wx6"D \$260 ea.
 Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552
 HCPCS: K0040+E0951

5"Wx6"D..... \$260 ea.
 Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553
 HCPCS: K0040+E0951

5"Wx8"D \$260 ea.
 Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554
 HCPCS: K0040+E0951

6"Wx8"D¹ \$260 ea.
 Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555
 HCPCS: K0040+E0951

6"Wx11"D¹ \$260 ea.
 Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558
 HCPCS: K0040+E0951

8"Wx11"D² \$260 ea.
 Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559
 HCPCS: K0040+E0951

4. Swing-Away Legrest Accessories

Gel Neoprene Footrest Sleeve (pair) \$84
 Part: INDPART2243. HCPCS: K0108

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-13"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
11"-13" extension causes no interference at all seat-to-floor heights.									

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

Legrests section continued on following pages.

18. LEGRESTS CONTINUED**5. Amputee Supports with Swing-Away Upper**

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 QTY Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
 QTY Part: INDPART2823. HCPCS: K0108

18.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ Elevating Legrest Assembly \$171 ea.
 P L R Pair part: FRMASMB6483. Left: RIG138012.
 Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

☐ ☐ ☐ 11"-14" Length
 P L R Pair: FRMASMB4755.
 Left: RIG138026. Right: RIG138030
 For reinforced only

☐ ☐ ☐ 14"-17" Length
 P L R Pair: FRMASMB4756.
 Left: RIG138033. Right: RIG138035
 For reinforced only

☐ ☐ ☐ 8.5"-12" HD Drop-
 P L R Ins/12.5"-16.5" ELR, Ext
 Length
 Pair part: XRF000101.
 Left: XRF000102. Right: XRF000103
 For angle adjustable only

☐ ☐ ☐ 13"-18" HD Drop-
 P L R Ins/15.5"-19.5" ELR, Ext
 Length
 Pair part: XRF000104.
 Left: XRF000105. Right: XRF000106
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate Standard
 P L R Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate Standard
 P L R Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

'The 6"x11" footplates may cause interference with seat widths less than 18". 'The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

☐ ☐ ☐ 4"Wx6"D \$90 ea.
 P L R Pair parts: FRMASMB8812 x2.
 HCPCS: K0040

☐ ☐ ☐ 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
 P L R Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

☐ ☐ ☐ 5"Wx6"D \$90 ea.
 P L R Pair parts: FRMASMB8811 x2.
 HCPCS: K0040

☐ ☐ ☐ 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
 P L R Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

☐ ☐ ☐ 5"Wx8"D Footplates with Heel Loops \$90 ea.
 P L R Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D² Footplates with Heel Loops \$90 ea.
 P L R Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953
 Please verify the appropriate billing code to be used for this item with all third-party payors.

18.2. Legrest Accessories

Heel Loops (pair) \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper \$19 ea.
 Part: XRF000100

18.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

Legrests section continued on following page.

18. LEGRESTS CONTINUED**1b. Front Rigging Extension Length**

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

19. TRU-BALANCE 3 POWER POSITIONING ELECTRONICS**REQUIRED** with power seating or legrest options.**19.1a. Actuator Control through Joystick****1. Actuator Control Selection**Single Actuator Through Controller and\$2,163
Interface

Part: XRF000335. HCPCS: E2310

Multiple Actuator Through Controller and \$3,040
Interface

Part: XRF000336. HCPCS: E2311

2. Q-Logic 3 Optional Switch Control

Only available with Q-Logic 3, multiple actuator configurations. 4-Way switches will be mounted to same side of joystick.

4-Function Push Button

Part: CTL123218

See section 21.2 for switch mounting.

4-Way Toggle

Part: ELEASMB5925

4-Way Button Toggle

Part: ELEASMB5926

3. User Controls Side Selection**Required with 4-Way Toggle and 4-Way Button Toggle.**

Mount User Controls on Left

Part: PTOINDV3996

Mount User Controls on Right

Part: PTOINDV3997

**19.1b. Actuator Control with iAccess Switch**

Single Actuator Control Through Flat Button\$2,547

iAccess and Controller and Interface

Part: XRF000337. HCPCS: E2310

Multiple Actuator Control Through Flat Button\$3,460

iAccess and Controller and Interface

Part: XRF000339. HCPCS: E2311

Single Actuator Control Through Toggle Button\$2,547

iAccess and Controller and Interface

Part: XRF000338. HCPCS: E2310

Multiple Actuator Control Through Toggle Button\$3,460

iAccess and Controller and Interface

Part: XRF000340. HCPCS: E2311

19.2. 4-Function Push Button/iAccess Mount**REQUIRED** with 4-Function Push Button or iAccess options. Under arm mount is not available with TRU-Balance 3 HD armrests.

Joystick Mount No Charge

Part: MEC1907500

Attendant Mount \$26

Part: MEC1906564

Under Arm Mount, Flip-Back Armrest No Charge

Part: MEC1906531

Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$26

Part: MEC2001400

Under Arm Mount, Tube Style Armrest No Charge

Part: MEC1906530

Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$26

Part: MEC2001401

1. Under Arm Mount Side Selection

Left Side Under Arm Mount

Part: PTOINDV3678

Right Side Under Arm Mount

Part: PTOINDV3679

20. TRU-BALANCE 3 ATTENDANT CONTROL**Optional.** *Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control \$525

Part: XRF000078. HCPCS: E2331

Select this option if no charge option criteria is not met.

Attendant Control with Tilt & Recline, Alternative No Charge

Controls and Q-Logic Display or SCIM Module

Part: XRF000079


21. FULL LIGHT PACKAGE**Optional.** Requires NE+ 8 Key or Q-Logic 3 Joystick for Full Light Package.

Seat Mounted, LED, Full Light Package \$1,045

Part: XRF000130

22. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.Accessory Bar \$132
Part: XRF000128

Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.

Backpack Holder (Pair) (A) \$69
Part: ACC147774Cane/Crutch Holder \$105
Part: ACC152515. HCPCS: E2207Cell Phone/ Media Holder, X-Grip (E) \$69
Part: ACC155835
Not available with Heavy Duty armrests.Clothing Guard (Pair) \$237
Part: ACC125007Compact Vent Tray \$1,465
Part: ACC133889. HCPCS: E1029Cup Holder No Charge
Part: ACC125003
For use with track-style armrests.Cup Holder (Tube-style Armrest) No Charge
Part: ACCASMB2634
For use with tube-style armrests.Flag Kit \$69
Part: ACC1700443Flip-Back Armrest, Side Track Accessory Bracket ... \$11
Part: MEC1809178Front Mounted, Cell Phone/ Media Holder, \$69
X-Grip (E)
Part: ACC155835Glove Box (C) \$158
Part: ACC125011Grip Drink Holder (20085) \$15.50
Part: ACC152993Grip Lap Board (20005) \$38.95
Part: ACC152992Hydration System (B) \$237
Part: ACC125010Large, Black, Weather Cover \$202
Part: ACCCOVR1021Oxygen Holder \$210
Part: ACC125006. HCPCS: E2208
Includes accessory bar, a credit will be applied PTOIN-DV3769.Personal Item Hook \$69
Part: ACC125008Personal Item Pouch \$84
Part: ACC2101017Privacy Bag \$42
Part: ACC1907570Push Handles (pair) \$195
Part: ACC2005027
Not available with oxygen holder or vent tray.Q-Logic EX/NE Programmer \$573
Part: ELEASMB6976Quantum Backpack \$168
Part: ACC1905853Rear View Mirror \$69
Part: ACC155833Reflector Kit \$53
Part: ACC124670Stealth Large Tablet Holder \$400
Part: ACC141818Transfer Handle (D) \$105 ea.
Part: MEC2205680 x2. HCPCS: E0705Trilogy Vent Tray \$420
Part: ACC136791. HCPCS: E1029
Includes accessory bar, a credit will be applied PTOIN-DV3769. Contact custom for models other than 100 and 200.VOCSN Vent Tray for Accessory Bar Mount \$940
Part: ACC2000166. HCPCS: E1029
Includes accessory bar, a credit will be applied PTOIN-DV3769.VOCSN Vent Tray for Power Tilt, Seat Mount \$940
Part: ACC2000162. HCPCS: E1029VOCSN Vent Tray for Power Lift, Seat Mount \$940
Part: ACC2000159. HCPCS: E1029**23. BRANCH PATCHES****Optional.**Air Force Patch \$21
Part: ACC140240Coast Guard Patch \$21
Part: ACC140239Navy Patch \$21
Part: ACC140244Army Patch \$21
Part: ACC140241Marines Patch \$21
Part: ACC140243Purple Heart Patch \$21
Part: ACC140238**24. POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#).**24.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. BracketSwing-away, Rail Mounts, pair (TWBTLTC) \$535
Part: INDPART2677. HCPCS: E1028Lateral Hardware, Q-Back Right \$267.50
Part: ST-TWBTL-QL-R. HCPCS: E1028Lateral Hardware, Q-Back Left \$267.50
Part: ST-TWBTL-QL-L. HCPCS: E1028Lateral Hardware, Unilink, pair \$632
Part: ST-ULS4042000-Q. HCPCS: E1028**2. Pads**Curved Pads, 3"x5", pair (TWBL-3X5) \$238
Part: INDPART2663. HCPCS: E0956Curved Pads, 4"x6", pair (TWBL-4X6) \$238
Part: INDPART2664. HCPCS: E0956Flat Pads, 3"x5", pair (TWBLP3X5) \$238
Part: INDPART2676. HCPCS: E0956Flat Pads, 4"x6", pair (TWBLP4X6) \$238
Part: INDPART2675. HCPCS: E0956**3. Installation**Install Laterals
Part: PTOINDV3499Do Not Install
Part: PTOINDV3500

Positioning components section continued on following page.

24. POSITIONING COMPONENTS CONTINUED**24.2. Stealth UniLink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable & Removable Hardware .. \$252.21 ea.
 Part: ST-UTTP5-Q. HCPCS: E1028

UniLink, 7" Height Adjustable & Removable Hardware .. \$262.50 ea.
 Part: ST-UTTP7-Q. HCPCS: E1028

2. UniLink Pelvic/Thigh Pads

UniLink, 4" x 4" Pad \$101.25 ea.
 Part: ST-UTP544056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 12" Pad \$101.25 ea.
 Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 8" Pad..... \$101.25 ea.
 Part: ST-UTP548056100-Q. HCPCS: E0953/E0956

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$77.20 ea.
 Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$77.20 ea.
 Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

**24.3. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW) \$267.75 ea.
 Part: INDPART2804. HCPCS: E1028

Flip Down with LEFT Full Surface Contact \$267.75 ea.
 Feature (TWBADD-L & LPHW407)
 Part: POS140428. HCPCS: E1028

Removable Bracket (LP430)..... \$270.38 ea.
 Part: INDPART2805. HCPCS: E1028

Flip Down with RIGHT Full Surface Contact..... \$267.75 ea.
 Feature (TWBADD-R & LPHW407)
 Part: POS140429. HCPCS: E1028

Removable Bracket with 1" Offset \$363.61 ea.
 (LP430 & LHW-121)
 Part: INDPART2854. HCPCS: E1028

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122)..... \$27.04
 Part: INDPART2821. HCPCS: K0108
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129) \$54.08
 Part: POS140430. HCPCS: K0108
 Only available with the flip down full surface bracket.

3. Pelvic/Thigh Guide Pad

Short Pad, 4"x4" (TWBADD-4X4)..... \$77.25 ea.
 Part: INDPART2801. HCPCS: E0953/E0956

Long Pad, 4"x12" (TWBADD-4X12)..... \$77.25 ea.
 Part: INDPART2803. HCPCS: E0953/E0956

Medium Pad, 4"x8" (TWBADD-4X8)..... \$77.25 ea.
 Part: INDPART2802. HCPCS: E0953/E0956

4. Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

24.4. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$146.26 ea.
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$146.26 ea.
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$146.26 ea.
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$146.26 ea.
 Part: ST-IFPS2104-Q. HCPCS: E0960

24.5. Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor, (IFCS4113).. \$92.70 ea.
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor, (IFCS4114).. \$92.70 ea.
 Part: ST-IFCS4114-Q. HCPCS: E0960

24.6. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. The 1" ball kit must be selected to mount the facial lateral to the headrest bracket. Also select a facial pad(s) with the facial lateral(s).

1. Facial Lateral Hardware

Facial Lateral for Comfort Series..... \$303.08 ea.
 Right Side (CPS9-R-HW)
 Part: INDPART2647. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left \$303.08
 Part: ST-SAR2A-L-Q. HCPCS: E1028

Facial Lateral for Comfort Series..... \$303.08 ea.
 Left Side (CPS9-L-HW)
 Part: INDPART2648. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$468.39
 2 Upper & 1 Lower Right
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

Facial Laterals for Comfort Series \$358.18 ea.
 Left and Right Side (CPS9-B-HW)
 Part: INDPART2649. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$468.39
 2 Upper & 1 Lower Left
 Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Right..... \$303.08
 Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$413.29
 2 Upper
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028



Positioning components section continued on following page.

24. POSITIONING COMPONENTS CONTINUED**2. Facial Lateral Pads and Accessories**


Egg Switch Mount (ESM) (A) \$54.08 ea.
 Part: INDPART2815

Figure 8 Bracket with 13" Rod \$275.53 ea.
 and Egg Switch Mount (B)
 Part: ST-F8675-13-Q

Facial Pad, 4.5"x2" (900M-Q) \$57.68 ea.
 Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2" (920M-Q) \$57.68 ea.
 Part: ST-920M-Q. HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C) .. \$262.50 ea.
 Part: ST-SUS18-Q
 Allows mounting of facial pads and attaches to the headrest


**24.7. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support \$189.26 ea.
 Bracket Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket... \$200.85
 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
 Part: ST-UETB4021FL-Q

Right, Swing-Away, Elbow Support \$189.26 ea.
 Bracket Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

Right, UniLink, Fixed, Elbow Support Bracket.. \$200.85
 for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
 Part: ST-UETB4021FR-Q

**2. Elbow Support Pad**

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$54.08 ea.
 Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$54.08 ea.
 Part: POS124067. HCPCS: K0108

24.8. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H \$513.71
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$513.71
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$513.71
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954

26.9. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211) \$82.40 ea.
 Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$82.40 ea.
 Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$82.40 ea.
 Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

25. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).

25.1. Enhanced Display

Q-Logic 3 Enhanced Display Kit \$1,980
 Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard
 Part: DGN137728

Red Corvette \$63
 Part: DGN137729-PRD1

Electric Blue \$63
 Part: DGN137729-PBU9

Back in Black \$63
 Part: DGN137729-PBK1

Night Train \$63
 Part: DGN137729-PBK3

Pretty in Pink \$63
 Part: DGN137729-PPK3

Evolution Orange \$63
 Part: DGN137729-POR7

Deep Purple \$63
 Part: DGN137729-PPR5

Grey Street \$63
 Part: DGN137729-PSL3

Midnight Blue \$63
 Part: DGN137729-PBU10

Green Machine \$63
 Part: DGN137729-PGN11

**25.2. SCIM Module**

Optional.

☐ SCIM Module for Q-Logic 3 \$525
 Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.

**25.3. Q-Logic 3 SAJ**

Q-Logic 3 EX SAJ \$840
 Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



Alternative controls and electronics section continued on following page.

25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**25.4. Head Array**

i-Drive 4.0 Bluetooth Tri-array Series, \$4,546.42
3 Switch with Egg Switch
Part: ST-IDHBT200-U-Q. HCPCS: E2330

i-Drive 4.0 Bluetooth Pro Series, \$4,594.83
3 Switch with Egg Switch
Part: ST-IDHBT300-U-Q. HCPCS: E2330

Adult, Curved, i-Drive Tri-Array with \$4,546.42
Mounting Hardware
Part: ST-IDHBT400-U-Q. HCPCS: E2330

Junior, Curved, i-Drive Tri-Array with \$4,546.42
Mounting Hardware
Part: ST-IDHBT450-U-Q. HCPCS: E2330

Adult, i-Drive Pro Series Head Array \$4,594.83
(IDHBT300-1)
Part: ST-IDHBT300-1-Q. HCPCS: E2330

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

25.5. Sip and Puff Drive**1. Sip and Puff Module**

Sip and Puff Module for Q-Logic 3 \$1,990
Part: CTL168056. HCPCS: E2325
Includes module and mounting hardware.

2. Sip and Puff Straw and Mount

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

Stealth Sip and Puff for Comfort Series, \$303.08
Left Mount (SP625M-CL)
Part: ST-SP625M-CL-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, \$303.08
Head Array (SP625-U)
Part: ST-SP625-U-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series, \$303.08
Right Mount (SP622M-CR)
Part: ST-SP625M-CR-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Sip & Puff Module



Stealth Sip & Puff

3. Sip & Puff Accessories

Sip & Puff Straws, 12 pack (SP625-S12) \$33.06
Part: ST-SP625-S12-Q

Sip & Puff Bus Cable \$447
Part: DWR1111H046

25.6. Egg Switches

Red Egg Switch (ESRED) \$83.28
Part: SWTMCRO1077

Blue Egg Switch (ESBLU) \$83.28
Part: SWTMCRO1079

Black Egg Switch (ESBLK) \$83.28
Part: SWTMCRO1081

Green Egg Switch (ESGRN) \$83.28
Part: SWTMCRO1078

Yellow Egg Switch (ESYEL) \$83.28
Part: SWTMCRO1080



Egg Switch

25.7. Switch Adapter

2 to 1 Switch Adapter \$15
Part: DWR1361H037
Connects 2 single switches to a stereo jack to act as a dual switch.

25.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520
Part: CTL168057

25.9. USB Mobile Device Charger

USB Mobile Device Charger \$184
Part: ELE168069
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.



USB Mobile Device Charger