



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

J4 Order Form with TRU Balance® Flex Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____



A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

1. BASE MODEL

REQUIRED. 300 lb. weight capacity with standard 5 mph programmable motors, 8 amp off board charger, and flat free drive wheels. Dimensions of base without legrests: 24.5"Wx32.3"L.

J4 2S-SS \$7634

Part: J4 2S-SS. HCPCS: K0822

Solid seat pan standard.

2. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.



Red Corvette

Part: DGN1900533-PRD1



Electric Blue

Part: DGN1900533-PBU9



Back in Black

Part: DGN1900533-PBK1

3. BATTERIES**Optional.** The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Nano-Silica, Interceptor 2 Battery \$355 ea.
 Part: BAT2004376 x2. HCPCS: E2361

4. QUANTUM BACKUP CAMERA**Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
 Part: MEC1910365

1. Optional Mounting Bracket

Backup Camera Headrest Bracket \$21 Backup Camera Track Mount Bracket \$27
 Part: MEC1910373 Part: MEC2104227
 Works with .5" x .5" vertical supports only.

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
 Part: MEC2106123

**5. TRANSIT OPTIONS****Optional.** The power base is standard without transit securement.

Unoccupied Transit No Charge
 Part: KIT130229
 Allows the powerchair to be secured for unoccupied transit.

WC-19 Occupied Transit Option \$289
 Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

6. ELECTRONICS**REQUIRED.** Select one option.**6.1a. Non-Expandable Controls**

QLNES 6 Key Joystick Standard
 Part: CTL2205379

6.1b. Expandable Controls

Expandable controls are special order, please contact your inside sales representative.

6.2. Q-Logic Magnetic Charger

Available with QLNES Joystick only.

Magnetic to XLR Charger Adapter with Wall Holder \$60
 Part: ELE2103635

7. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**7.1a. Fixed, In-line Joystick Bracket**

Fixed, In-line, Bracket for Q-Logic Standard
 Part: FRMASMB10617

7.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.


Swing-away, Q-Logic \$319
 Part: ACC1908434. HCPCS: E1032

7.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side \$245 Basic Swing Away, Q-Logic, Right Side \$245
 Part: ACC2201041. HCPCS: E1032 Part: ACC2201042. HCPCS: E1032

7.2. Joystick Mount Location

Left Side Mount Right Side Mount
 Part: PTOINDV1225 Part: PTOINDV1226

**8. JOYSTICK HANDLES****Optional.** The joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90
 Part: ACCASMB1860. HCPCS: E2323

Chin Cup (C) \$90
 Part: ACCASMB1861.
 HCPCS: E2323. (E2324 when chin control)

T-Bar (D) \$126
 Part: ACCASMB1858. HCPCS: E2323

Medium Ball Handle (E) \$69
 Part: ACCASMB1859. HCPCS: E2323

Soft Ball Handle (F) \$80
 Part: ACCASMB1856. HCPCS: E2323

Extended I-Handle (G) \$90
 Part: ACCASMB1857. HCPCS: E2323

Bodypoint 3" Goal Post Handle \$104
 Part: PLSKNCP1057. HCPCS: E2323

Bodypoint 4" Goal Post Handle \$104
 Part: PLSKNCP1058. HCPCS: E2323



9. TRU BALANCE FLEX SEATING SYSTEM

Seat width and depth can be adjusted in 1" increments from 15"W-22"W and 16"D-20"D. Maximum seat size 22"W x 20"D.

9.1. TRU Balance Flex Frame Selection**REQUIRED.**

TRU Balance Flex Seating System\$351
 Part: XRF000262

9.2. TRU Balance Flex Seat to Floor**REQUIRED.** Select one option.

16.5" at front and rear
 Part: PTOINDV3459

17.5" at front with 2° posterior angle
 Part: PTOINDV3460

18.25" at front with 4° posterior angle
 Part: PTOINDV3805

9.3. TRU Balance Flex Seat Dimension - Width**REQUIRED.** Select one option. Maximum seat width is 22". Part number is SETWDTH plus number following # below.


15" W (#1005) 16" W (#1006) 17" W (#1007) 18" W (#1008) 19" W (#1009) 20" W (#1010) 21" W (#1011) 22" W (#1012)

9.4. TRU Balance Flex Seat Dimension - Depth**REQUIRED.** Select one option. Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

15" D (#1005) 16" D (#1006) 17" D (#1007) 18" D (#1008) 19" D (#1009) 20" D (#1010)

9.5. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type**

TRU-Comfort 2 SPP
 Reference part: XRF000116. [See cushion order form](#).
 Part series: TC2-###XYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPP
 Reference part: XRF000118. [See cushion order form](#).
 Part series: ESS-###XYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPP
 Reference part: XRF000113. [See cushion order form](#).
 Part series: SOL-###XYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPP
 Reference part: XRF000112. [See cushion order form](#).
 Part series: SPCF-###XYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPP
 Reference part: XRF000117. [See cushion order form](#).
 Part series: SPCG-###XYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere P
 Reference part: XRF000122. [See cushion order form](#).
 Part series: PREM-###XYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SP
 Reference part: XRF000114. [See cushion order form](#).
 Part series: GLA-###XYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SP
 Reference part: XRF000115. [See cushion order form](#).
 Part series: ZEN-###XYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity G
 Reference part: XRF000111. [See cushion order form](#).
 Part series: SIM-###XYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39 Part: POP1909006-S 14"W and smaller	Medium Wedge Positioning Kit.....\$72 Part: POP1909006-M 14"W and smaller	Large Wedge Positioning Kit..... \$105 Part: POP1909006-L 14"W and smaller
Small Wedge Positioning Kit\$39 Part: POP1712340-S 15"W and larger	Medium Wedge Positioning Kit.....\$72 Part: POP1712340-M 15"W and larger	Large Wedge Positioning Kit..... \$105 Part: POP1712340-L 15"W and larger

TRU Balance Flex seating system section continued on following pages.

9. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**3. Cushion Top Foam**

Visco Top Foam
Reference part: XRF000380

FRG Top Foam
Reference part: XRF000381

No Top Foam
Reference part: XRF000382

4. Cushion Liner

Stealth Urethane Liner
Reference part: XRF000383

Stealth Polyurethane Liner
Reference part: XRF000384

No Liner
Reference part: XRF000385

5. Cushion Cover

CoolCore Cover
Reference part: XRF000386

Mesh Cover
Reference part: XRF000387

Stealth Tek Cover
Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner					
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389
MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390
QTY _____

Additional Stealth Polyurethane Liner
Reference Part: XRF000391
QTY _____

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392
QTY _____

Additional Stealth Tek Cover
Reference Part: XRF000394
QTY _____

Additional Mesh Cover
Reference Part: XRF000393
QTY _____

9.6. TRU Balance Flex Back Cane Selection

REQUIRED. Select one option.

Angle Adjustable Back Cane No Charge
Part: MEC140266

4" Endomorph (2" per side) Back Cane \$205
Part: XRF000272

Flip Forward Back Cane \$184
Part: XRF000274
Not available with oxygen holder, or cane mounted and cantilever armrests.

Manual Recline Back \$788
Part: FRMASMB12910. HCPCS: E1226

2" Mesomorph (1" per side) Back Cane \$205
Part: MEC1804736

Select Synergy Back and an overall back height, or to omit back. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, Sport Back, rear accessories, or articulating foot platform. Pivot point is 2.5" above the seat pan.

2" Endomorph (1" per side) Back Cane \$205
Part: XRF000270

9.7. TRU Balance Flex Back Cane Height Selection

REQUIRED. Select one option.

14" High Canes
Part: FRMASMB15832

17" High Canes
Part: FRMASMB15835

20" High Canes
Part: FRMASMB15838

23" High Canes
Part: FRMASMB15841

15" High Canes
Part: FRMASMB15833

18" High Canes
Part: FRMASMB15836

21" High Canes
Part: FRMASMB15839

24" High Canes
Part: FRMASMB15842

16" High Canes
Part: FRMASMB15834

19" High Canes
Part: FRMASMB15837

22" High Canes
Part: FRMASMB15840

25" High Canes
Part: FRMASMB15843

9.8. TRU Balance Flex Back Cushion Options

REQUIRED. Select one option. See the seat back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

1a. Synergy Back Selection

Synergy Back \$289
Part: XRF000363
Select overall back height in later section. Back cushion width is determined by seat width and back cane type.

1b. Sport Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

Sport Back Cushion, 21" Wide or less \$394
Part: XRF000361

Sport Back Cushion, 22" Wide or more \$604
Part: XRF000362

TRU Balance Flex seating system section continued on following pages.

9. TRU BALANCE FLEX SEATING SYSTEM CONTINUED

		Back Width											
		14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"	25"
Overall Height	16"	P	T	P	T	P	T	P	T	P	T	P	T
	17"	P	T	S	P	T	S	P	T	S	P	T	P
	18"	P	T	S	P	T	S	P	T	S	P	T	P
	19"	P	T	S	P	T	S	P	T	S	P	T	P
	20"	P	T	S	P	T	S	P	T	S	P	T	P
	21"	P	T	S	P	T	S	P	T	S	P	T	P
	22"	P	T	S	P	T	S	P	T	S	P	T	P
	23"	P	T	S	P	T	S	P	T	S	P	T	P
	24"	P	T	S	P	T	S	P	T	S	P	T	P
	25"	P	T	S	P	T	S	P	T	S	P	T	P
		S Synergy Back P Sport Back T TRU-Comfort Back											
		TRU Balance Flex Seat Back Availability Chart											

1c. TRU-Comfort Back

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width..... \$744
Part: TC4V10W21WBC. HCPCS: E2620TRU-Comfort Back, Back Cane Option, 22" Width \$855
Part: TC4V22W32WBC. HCPCS: E2621**1d. Omit Back Selection**Omit Back..... No Charge
Part: PTOINDV1287**9.9. Back Width Package****REQUIRED** with Sport or TRU-Comfort Back selection. Available with +/-4" of seat width.16"W-20"W Package
Part: XRF00003818"W-22"W Package
Part: XRF000039**9.10. Overall Back Height**16" Overall Back Height
Part: SETHGHT102919" Overall Back Height
Part: SETHGHT103222" Overall Back Height
Part: SETHGHT103525" Overall Back Height
Part: SETHGHT103828" Height
Part: SETHGHT104117" Overall Back Height
Part: SETHGHT103020" Overall Back Height
Part: SETHGHT103323" Overall Back Height
Part: SETHGHT103626" Height
Part: SETHGHT103929" Height
Part: SETHGHT104218" Overall Back Height
Part: SETHGHT103121" Overall Back Height
Part: SETHGHT103424" Overall Back Height
Part: SETHGHT103727" Height
Part: SETHGHT1040**9.11. TRU Balance Flex Armrests****REQUIRED.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. Drop-In, Height Adjustable, Removable, T-Arm Armrests**☐ ☐ Full Length, Drop-In Armrests..... No Charge
☐ ☐ Left part: ARM1803253. Right part: ARM156742☐ ☐ Desk Length, Drop-In Armrests No Charge
☐ ☐ Left part: ARM1803254. Right part: ARM156743**1b. Flip Back, Cane Mounted Armrests**

May interfere with swing-away laterals.

☐ ☐ Full Length, Flip Back Armrests..... \$63
☐ ☐ Left part: FRMASMB7038. Right part: FRMASMB7039☐ ☐ Desk Length, Flip Back Armrests \$63
☐ ☐ Left part: FRMASMB7040. Right part: FRMASMB7041**1c. Two-Post, Flip Back, Height Adjustable, Removable Armrests**

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐ Full Length, Two-Post Armrests.....\$158 ea.
☐ ☐ Left part: XRF000051. Right part: XRF000052. HCPCS: E0973☐ ☐ Desk Length, Two-Post Armrests\$158 ea.
☐ ☐ Left part: XRF000053. Right part: XRF000054. HCPCS: E0973**1d. Single Post, Height Adjustable, Removable Armrests**

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests \$79 ea.
☐ ☐ Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973☐ ☐ Desk Length, Single Post Armrests \$79 ea.
☐ ☐ Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973**9.12. TRU Balance Flex Armpads****REQUIRED** with armrest selection above, except with drop-in armrests. Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.**1a. Straight Armpads**Straight, Full Length Left and Right No Charge
Part: ARM134275Straight, Desk Length Left, Full Length Right No Charge
Part: ARM134280Straight, Full Length Left, Desk Length Right No Charge
Part: ARM134281Straight, Desk Length Left and Right..... No Charge
Part: ARM134254

TRU Balance Flex seating system section continued on following pages.

9. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**1b. Waterfall Armpads**Waterfall, Full Length Left and Right..... No Charge
Part: ARM134285Waterfall, Full Length Left, Desk Length Right No Charge
Part: ARM134289Waterfall, Desk Length Left, Full Length Right No Charge
Part: ARM134288Waterfall, Desk Length Left and Right No Charge
Part: ARM134282**1c. Falcon Gel Armpads**Falcon Gel, Full Length Left and Right\$311.22
Part: ARM136438Falcon Gel, Full Length Left, Desk Length Right.....\$287.93
Part: ARM136442Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
Part: ARM136439Falcon Gel, Desk Length Left and Right \$265.02
Part: ARM136437**10. HEADRESTS****Optional.** Select a headrest pad and bracket.**10.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q. HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q. HCPCS: E0955**2. Headrest Brackets**Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033**10.2. Flip-down Headrest Bracket**Stealth Flip-down Headrest Bracket (FDM350).....\$170
Part: INDPART2219

The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

11. LAP BELT**REQUIRED.**60" Lap Belt
Part: ACCBELT100260" Sure-Lok Lap Belt (Non-retractable)..... \$120
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016**12. LEGRESTS****Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**12.1a. Base Mounted Foot Platform**

May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard
Part: FRMASMB16077Omit base mounted foot platform No Charge
Part: PTOINDV3674**12.1b. Basic Center Mount Foot Platform**

Select legrest length in following section.

1. Center Mount Foot Platform TypeBasic Center Mount..... No Charge
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility.

Short Extension, 5"x9" Footplates \$180
Part: XRF000312Medium Extension, 5"x9" Footplates \$180
Part: XRF000313Long Extension, 5"x9" Footplates \$180
Part: XRF000314**3. Basic Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount..... \$158
Part: RIG2104369Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997
Only available with Long Length Extension.**4. Basic Center Mount Additional Accessories**Basic Center Mount Footplate Connector \$16
Part: RIG2104240

Legrests section continued on following pages.

12. LEGRESTS CONTINUED**12.1c. Center Mount Foot Platform**

Select legrest length in following section.

1. Center Mount Foot Platform TypeQuick Adjustable Center Mount Foot Platform\$263
Part: XRF000311**2. Extension and Angle Adjustable Footplates**

See legrest length matrix for compatibility. Do not select with Basic Center Mount.

Medium Extension, 5"x6" Footplates Long Extension, 5"x6" Footplates
Part: RIG150206 Part: RIG151505Medium Extension, 5"x8" Footplates Long Extension, 5"x8" Footplates
Part: RIG2404016 Part: RIG2404017Medium Extension, 5"x11" Footplates Long Extension, 5"x11" Footplates
Part: RIG2404018 Part: RIG2404019

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf SupportSmall Calf Pads (pair) 6"Wx6"H \$158
Part: RIGASMB7120018Large Calf Pads (pair) 8" taper to 6"Wx8"H \$158
Part: RIGASMB7120019**3b. Center Mount Calf Panel**Small Calf Panel 18"Wx8"H \$184
Part: FRMASMB7110012Large Calf Panel 24"Wx8"H \$237
Part: RIG123825
Only available on widths of 20" or greater**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
Part: RIG152519
Only available with the 5x6 footplates.Outside Edge \$69
Part: RIG2404921
Available with 5"x8" footplates.Width Extension for 8" \$69
Part: RIG2404919Front and Rear Cap \$16
Part: RIG2404923Width Extension for 11" \$69
Part: RIG2404920

Available with 5"x8" footplates.

5. Center Mount Additional AccessoriesFootplate Connector \$16
Part: RIGASMB7121023Footplate Riser Calf Panel \$50
Part: FRMASMB7110011Footplate Riser \$75
Part: XRF000348
Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K01085"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K01085"x11" Tapered Gel Footplate Covers \$270
Part: ST-AFPG5011T-Q. HCPCS: K01085"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-AFPG5011E-Q. HCPCS: K0108**12.1d. Swing-Away Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

☐ P	☐ L	☐ R	60°-90° Swing-Away Legrest, Static Seat \$86 ea. Pair part: XRF000082. Left: XRF000083. Right: XRF000084. Not available with 11"-13" extensions.
----------------------------	----------------------------	----------------------------	---

☐ P	☐ L	☐ R	60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea. Pair part: XRF000085. Left: XRF000086. Right: XRF000087. For use only with 11"-13" lower extensions.
----------------------------	----------------------------	----------------------------	---

2. Lower Extension

☐ P	☐ L	☐ R	11"-13" Length Pair part: XRF000088. Left: XRF000089. Right: XRF000090.
----------------------------	----------------------------	----------------------------	--

☐ P	☐ L	☐ R	13"-16.5" Length Pair part: XRF000094. Left: XRF000095. Right: XRF000096.
----------------------------	----------------------------	----------------------------	--

☐ P	☐ L	☐ R	15.5"-19.5" Length Pair part: XRF000097. Left: XRF000098. Right: XRF000099.
----------------------------	----------------------------	----------------------------	--

3. Angle Adjustable Footplate(s)¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ P	☐ L	☐ R	5"Wx8"D with Heel Loops \$90 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: E0951
----------------------------	----------------------------	----------------------------	---

☐ P	☐ L	☐ R	6"Wx11"D ¹ with Heel Loops \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: E0951
----------------------------	----------------------------	----------------------------	---

☐ P	☐ L	☐ R	6"Wx8"D ¹ with Heel Loops \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: E0951
----------------------------	----------------------------	----------------------------	--

☐ P	☐ L	☐ R	8"Wx11"D ² with Heel Loops \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: E0951
----------------------------	----------------------------	----------------------------	---

Legrests section continued on following pages.

12. LEGRESTS CONTINUED**4. Amputee Supports with Swing-Away Upper**

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 QTY Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
 QTY Part: INDPART2823. HCPCS: K0108

12.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ Elevating Legrest Assembly \$171 ea.
 P L R Pair part: FRMASMB483. Left: RIG138012. Right: RIG138013. HCPCS: K0195
 Select lower extension and footplates

2a. Lower Extension

☐ ☐ ☐ 11"-14" Length
 P L R Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

☐ ☐ ☐ 14"-17" Length
 P L R Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

☐ ☐ ☐ 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 P L R Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

☐ ☐ ☐ 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 P L R Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	!	!	!	!	!	!	!
Angle Adjustable, 12.5"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	!

☒ This rigging configuration causes no interference issues at this seat-to-floor height.

☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate Standard
 P L R Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate Standard
 P L R Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ 5"Wx8"D with Heel Loops \$90 ea.
 P L R Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: E0951

☐ ☐ ☐ 6"Wx11"D¹ with Heel Loops \$90 ea.
 P L R Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: E0951

☐ ☐ ☐ 6"Wx8"D¹ with Heel Loops \$90 ea.
 P L R Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: E0951

☐ ☐ ☐ 8"Wx11"D² with Heel Loops \$90 ea.
 P L R Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

12.2. Legrest Accessories

Heel Loops (pair) \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper \$19 ea.
 Part: XRF000100

12.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

Legrests section continued on following pages.

12. LEGRESTS CONTINUED**1b. Front Rigging Extension Length**

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)	18" L (#3492)	19" L (#3494)
11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	17.5" L (#3491)	18.5" L (#3493)	19.5" L (#3495)

13. ACCESSORIES**Optional.**Large, Black, Weather Cover \$202
Part: ACCCOVR1021Oxygen Holder..... \$210
Flex part: ACC125006. HCPCS: E2208Flex Seat Cup Holder..... No Charge
Part: ACCASMB2634
Requires two-post armrests.Saddle Bag..... \$40
Part: DWR1020D010Flag Holder \$55
Flex part: ACCASMB2415.Cane/Crutch Holder \$105
Flex part: ACC152515. HCPCS: E2207Transfer Handle \$105 ea.
Part: MEC2205680. HCPCS: E0705
Not available with Heavy Duty armrests.Accessory Bar \$132
Part: XRF000128Privacy Bag \$42
Part: ACC20005723Personal Item Pouch \$84
Part: ACC2101017Remote Stop Switch \$1,295
Part: XRF000379**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
QTY Part: ACC2407046Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
QTY Part: ACC2405496**14. POSITIONING COMPONENTS****Optional.****14.1. Stealth Thoracic Laterals**

Laterals are not available with the omit back option.

1. BracketSwing-away, Rail Mounts, pair (TWBSTL) \$562
Part: INDPART2896. HCPCS: E1034Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034**2. Pads**Curved Pads, 3"x5", pair (TWBL-3X5) \$250
Part: INDPART2663. HCPCS: E0956Curved Pads, 4"x6", pair (TWBL-4X6) \$250
Part: INDPART2664. HCPCS: E0956**3. Installation**Install Laterals
Part: PTOINDV3499Do Not Install
Part: PTOINDV3500**4. Lateral Extensions****REQUIRED with Synergy Back and thoracic lateral selections.**Lateral Extension for Synergy Back (pair) \$56
Part: XRF000150**14.2. Stealth UniLink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea.
Removable Mounting Hardware
QTY Part: ST-UTTP5-Q. HCPCS: E1028/E1034^Flip Down with LEFT Full Surface Contact \$281 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^UniLink, 7" Height Adjustable, \$276 ea.
Removable Mounting Hardware
QTY Part: ST-UTTP7-Q. HCPCS: E1028/E1034^Flip Down with RIGHT Full Surface Contact..... \$281 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^**2. UniLink Pelvic/Thigh Pads**

^ Please verify the appropriate billing code to be used with all third-party payors

UniLink, 4" x 4" Pad \$106 ea.
QTY Part: ST-UTP544056100-Q. HCPCS: E0953/E0956UniLink, 4" x 12" Pad \$106 ea.
QTY Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956UniLink, 4" x 8" Pad..... \$106 ea.
QTY Part: ST-UTP548056100-Q. HCPCS: E0953/E0956

Positioning components section continued on following pages.

Account Number: _____ Date: _____ Quote Name: _____

14. POSITIONING COMPONENTS CONTINUED

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea.
 QTY Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$81 ea.
 QTY Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

17. LAP TRAY

17.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
 Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out..... \$309
 Part: ST-SSLXL010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out..... \$309
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware..... \$302
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	SM	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.