



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

J4E Order Form with Power Seating

Account Number: _____ Date: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
PO Number: _____
Quote Name/Marked for: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____



A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

1. BASE MODEL

REQUIRED. 300 lb. weight capacity with standard 5 mph programmable motors, 8 amp off board charger, and flat free drive wheels. Dimensions of base without legrests: 24.5"Wx32.3"L.

J4E 2SP-SS \$7,583
Part: J4E 2SP-SS. HCPCS: K0835

2. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.



Red Corvette

Part: DGN2400773-PRD1



Electric Blue

Part: DGN2400773-PBU9



Back in Black

Part: DGN2400773-PBK1

Account Number: _____ Date: _____ Quote Name: _____

3. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Nano-Silica, Interceptor 2 Battery \$355 ea.
Part: BAT2004376 x2. HCPCS: E2361

4. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Mounting Bracket

Backup Camera Headrest Bracket \$21 Backup Camera Track Mount Bracket \$27
Part: MEC1910373 Part: MEC2104227
Works with .5" x .5" vertical supports only.

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



5. TRANSIT OPTIONS

Optional. The power base is standard without transit securement.

Unoccupied Transit No Charge
Part: KIT130229
Allows the powerchair to be secured for unoccupied transit.

WC-19 Occupied Transit Option \$289
Part: KIT152047
To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

6. ELECTRONICS

REQUIRED.

6.1a. Non-Expandable Controls

QLNES 6 Key Joystick Standard
Part: CTL2205379
Controls up to 2 actuators.

6.1b. Expandable Controls

1. Expandable Controls Required Components

Power Tilt, Power Lift, and Power Articulating Foot Platform combined selections require the controller module and expandable harness below.

Q-Logic 3 Controller Module \$730 Q-Logic 3 Expandable Controller Harness \$458
Part: XRF000002. HCPCS: E2377 Part: XRF000264. HCPCS: E2313

2. Q-Logic 3 Joystick

Q-Logic 3e Joystick \$735
Part: CTL2301467

3. Q-Logic Magnetic Charger

Available with QLNES joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$60
Part: ELE2103635

6.1c. Specialty Controls

Specialty controls are special order, please contact your inside sales representative.

7. JOYSTICK MOUNTING BRACKETS

REQUIRED. Select a mount to be used with a joystick selected in previous section.

7.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617

7.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.


Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1032

7.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side \$245 Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201041. HCPCS: E1032 Part: ACC2201042. HCPCS: E1032

7.2. Joystick Mount Location

Left Side Mount Right Side Mount
Part: PTOINDV1225 Part: PTOINDV1226



8. JOYSTICK HANDLES

Optional. The joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323

Chin Cup (C) \$90
Part: ACCASMB1861.
HCPCS: E2323. (E2324 when chin control)

T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323


Medium Ball Handle (E) \$69
Part: ACCASMB1859. HCPCS: E2323

Soft Ball Handle (F) \$80
Part: ACCASMB1856. HCPCS: E2323

Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323

Bodypoint 3" Goal Post Handle \$104
Part: PLSKNCP1057. HCPCS: E2323

Bodypoint 4" Goal Post Handle \$104
Part: PLSKNCP1058. HCPCS: E2323

**9. Q4 SEAT****9.1. Q4 Seat**

REQUIRED. Select one option.

Q4 Static Seat \$2,115
Part: PTOINDV4008
Requires adjustable seat lift.

Q4 Seat Power Tilt \$5,632
Part: PTOINDV3882. HCPCS: E1002
55° center of gravity (CG) power tilt system.

9.2. 10" Power Lift

Optional. ^ Please verify the appropriate billing code to be used for this item with all third-party payors.

Q4 Seat 10" Adjustable Seat Lift \$3,725
Part: PTOINDV4009. HCPCS: E2298/E2300^
10" of scissor elevate. When elevated, maximum speed will be up to 1.25 mph.

9.3. Seat to Floor Height

REQUIRED.

Lowest Seat to Floor
Part: PTOINDV3558

Middle Seat to Floor
Part: PTOINDV3559

Highest Seat to Floor
Part: PTOINDV3560

Low	Med	High
17.25"	18.25"	19.25"

Seat to Floor Heights

9.4. Seat Width

REQUIRED. Maximum seat width is 22". Width package range will be 16"-20" wide.

16" Wide
Part: SETWDTH1006

18" Wide
Part: SETWDTH1008

20" Wide
Part: SETWDTH1010

22" Wide
Part: SETWDTH1012

9.5. Seat Width Package

REQUIRED.

16"-20" Width Package
Part: XRF000008

20"-24" Width Package
Part: XRF000158

9.6. Seat Depth

REQUIRED. Maximum seat depth is 22".

16" Deep
Part: SETDPTH1006

18" Deep
Part: SETDPTH1008

20" Deep
Part: SETDPTH1010

22" Deep
Part: SETDPTH1012

9.7. Seat Depth Package


REQUIRED.

16"-20" Depth Package
Part: XRF000013

20"-24" Depth Package
Part: XRF000159

9.8. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.



Q4 seat section continued on following pages.

9. Q4 SEAT CONTINUED**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)

Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPPReference part: XRF000118. [See cushion order form.](#)

Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPPReference part: XRF000113. [See cushion order form.](#)

Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPPReference part: XRF000112. [See cushion order form.](#)

Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPPReference part: XRF000117. [See cushion order form.](#)

Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere PReference part: XRF000122. [See cushion order form.](#)

Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SPReference part: XRF000114. [See cushion order form.](#)

Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SPReference part: XRF000115. [See cushion order form.](#)

Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity GReference part: XRF000111. [See cushion order form.](#)

Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39

Part: POP1909006-S

14"W and smaller

Medium Wedge Positioning Kit.....\$72

Part: POP1909006-M

14"W and smaller

Large Wedge Positioning Kit..... \$105

Part: POP1909006-L

14"W and smaller

Small Wedge Positioning Kit\$39

Part: POP1712340-S

15"W and larger

Medium Wedge Positioning Kit.....\$72

Part: POP1712340-M

15"W and larger

Large Wedge Positioning Kit..... \$105

Part: POP1712340-L

15"W and larger

3. Cushion Top Foam**Visco Top Foam**

Reference part: XRF000380

No Top Foam

Reference part: XRF000382

FRG Top Foam

Reference part: XRF000381

4. Cushion Liner**Stealth Urethane Liner**

Reference part: XRF000383

No Liner

Reference part: XRF000385

Stealth Polyurethane Liner

Reference part: XRF000384

5. Cushion Cover**CoolCore Cover**

Reference part: XRF000386

Stealth Tek Cover

Reference part: XRF000388

Mesh Cover

Reference part: XRF000387

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer

Reference part: XRF000389

MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner

Reference Part: XRF000390

QTY

Additional Stealth Polyurethane Liner

Reference Part: XRF000391

QTY

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover

Reference Part: XRF000392

QTY

Additional Stealth Tek Cover

Reference Part: XRF000394

QTY

Additional Mesh Cover

Reference Part: XRF000393

QTY

9.9. Back Cane Type**REQUIRED.**

Angle Adjustable No Charge

Part: MEC140266

2" Endomorph (1" per side) Back Cane \$74

Part: XRF000269

Flip Forward Back Cane \$184

Part: XRF000274

4" Endomorph (2" per side) Back Cane..... \$74

Part: XRF000271

Q4 seat section continued on following pages.

9. Q4 SEAT CONTINUED**9.10. Back Cane Height****REQUIRED.**16"H Cane
Part: FRMASMB1583420"H Cane
Part: FRMASMB1583824"H Cane
Part: FRMASMB1584217"H Cane
Part: FRMASMB1583521"H Cane
Part: FRMASMB1583925"H Cane
Part: FRMASMB15843**9.11. Back Cushion Type****REQUIRED.****1a. Synergy Back Selection**

Synergy Back.....No Charge

Part: XRF000199

Select overall back height in later section. Back cushion width is determined by seat width and back cane type.

1b. TRU-Comfort 4 Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

TRU-Comfort 4 Back, Back Cane Option, 16" to 20" Width \$744
Part: TC4V10W21WBC. HCPCS: E2620TRU-Comfort 4 Back, Back Cane Option, 22" Width..... \$855
Part: TC4V22W32WBC. HCPCS: E2621**1c. ADI Back, AL Series, 4-Point, Active Contour****Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.Active Contour, 16" Wide x 19" or 20" Overall Height \$622
Part: ALBR16SA4. HCPCS: E2613Active Contour, 18" Wide x 23" or 24" Overall Height \$622
Part: ALBR18XA4. HCPCS: E2613Active Contour, 16" Wide x 21" or 22" Overall Height \$622
Part: ALBR16TA4. HCPCS: E2613Active Contour, 20" Wide x 19" or 20" Overall Height \$622
Part: ALBR20SA4. HCPCS: E2613Active Contour, 16" Wide x 23" or 24" Overall Height \$622
Part: ALBR16XA4. HCPCS: E2613Active Contour, 20" Wide x 21" or 22" Overall Height \$622
Part: ALBR20TA4. HCPCS: E2613Active Contour, 18" Wide x 19" or 20" Overall Height \$622
Part: ALBR18SA4. HCPCS: E2613Active Contour, 20" Wide x 23" or 24" Overall Height \$622
Part: ALBR20XA4. HCPCS: E2613Active Contour, 18" Wide x 21" or 22" Overall Height \$622
Part: ALBR18TA4. HCPCS: E2613**1d. ADI Back, AL Series, 4-Point, Deep Contour****Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.Deep Contour, 16" Wide x 19" or 20" Overall Height..... \$735
Part: ALBR16SD4. HCPCS: E2615Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$735
Part: ALBR18TD4. HCPCS: E2615Deep Contour, 16" Wide x 21" or 22" Overall Height..... \$735
Part: ALBR16TD4. HCPCS: E2615Deep Contour, 20" Wide x 19" or 20" Overall Height \$735
Part: ALBR20SD4. HCPCS: E2615Deep Contour, 16" Wide x 23" or 24" Overall Height..... \$735
Part: ALBR16XD4. HCPCS: E2615Deep Contour, 20" Wide x 23" or 24" Overall Height \$735
Part: ALBR20XD4. HCPCS: E2615Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$735
Part: ALBR18SD4. HCPCS: E2615**1e. ADI Back, AL Series, 4-Point, Extra Deep Contour****Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.Extra Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$822
Part: ALBR18SE4. HCPCS: E2620Extra Deep Contour, 18" Wide x 23" or 24" Overall Height..... \$822
Part: ALBR18XE4. HCPCS: E2620Extra Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$822
Part: ALBR18TE4. HCPCS: E2620**1f. Omit Back Selection**Omit Back Pan/Cushion..... No Charge
Part: PTOINDV1287**9.12. Overall Back Cushion Height****REQUIRED with back selection.** Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion.16" Height
Part: SETHGHT102921" Height
Part: SETHGHT103417" Height
Part: SETHGHT103022" Height
Part: SETHGHT103518" Height
Part: SETHGHT103123" Height
Part: SETHGHT103619" Height
Part: SETHGHT103224" Height
Part: SETHGHT103720" Height
Part: SETHGHT103325" Height
Part: SETHGHT1038

		Back Width			
		14"	16"	18"	20"
Overall Height	16"	T	T	T	T
	17"	S	T	S	T
	18"	T	T	T	T
	19"	T	T	T	T
	20"	T	T	T	T
	21"	S	T	S	T
	22"	T	T	T	T
	23"	T	T	T	T
	24"	S	T	S	T
	25"	S	T	S	T

S Synergy Back
 T TRU-Comfort Back

Q4 Seat Back Availability Chart

10. Q4 SEAT ARMRESTS**REQUIRED.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**10.1a. Drop-In, Height Adjustable, Removable Armrests**

Not available with toggle switch. Includes straight armpads. Do not select options from armrest size nor armpads sections.

☐ ☐ Full Length, Drop-In Armrests..... No Charge ☐ ☐ Desk Length, Drop-In Armrests No Charge
 L R Left part: ARM1803253. Right part: ARM156742 L R Left part: ARM1803254. Right part: ARM156743

10.1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

☐ ☐ Full Length, Two-Post Armrests.....\$158 ea. ☐ ☐ Desk Length, Two-Post Armrests\$158 ea.
 L R Left part: XRF000051. Right part: XRF000052 L R Left part: XRF000053. Right part: XRF000054

10.1c. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals

☐ ☐ Full Length, Flip Back Armrests.....\$63 ☐ ☐ Desk Length, Flip Back Armrests\$63
 L R Left part: FRMASMB7038. Right part: FRMASMB7039 L R Left part: FRMASMB7040. Right part: FRMASMB7041

10.2. Armrest Height**REQUIRED with two-post armrest selection.** Set left and right armrest height.**1. Left Armrest Height**

9" H 10" H 11" H 12" H 13" H
 PTOINDV3536 PTOINDV3537 PTOINDV3538 PTOINDV3539 PTOINDV3540

2. Right Armrest Height

9" H 10" H 11" H 12" H 13" H
 PTOINDV3451 PTOINDV3452 PTOINDV3453 PTOINDV3454 PTOINDV3455

10.3. Armpads**REQUIRED with two-post armrest selection.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

No Armpad Selection..... No Charge ☐ ☐ Waterfall, Armpad No Charge
 Part: XRF000060 L R Left part: XRF000062. Right part: XRF000127
☐ ☐ Straight, Armpad Standard
 L R Left part: XRF000061. Right part: XRF000126

11. HEADRESTS**Optional.** Select a headrest pad and bracket.**11.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179 Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
 Part: ST-CP08CCFPM-Q. HCPCS: E0955 Part: ST-CP14CCFPM-Q. HCPCS: E0955
 Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)\$179 Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
 Part: ST-CP10CCFPM-Q. HCPCS: E0955 Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink.....\$311 Multi-Axis Removable Headrest Bracket.....\$174
 Part: ST-UHB2D33200-Q. HCPCS: E1033 Part: ST-HMO475-17M-Q. HCPCS: E1033

11.2. Q Cam Lock**Optional.**

Q Cam Lock.....\$50
 Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

12. LAP BELT**REQUIRED.**

60" Lap Belt 60" Sure-Lok Lap Belt (Non-retractable)..... \$120
 Part: ACCBELT1002 Required with ISO 7176-19 Compliant Transit Option.
 70" Lap Belt
 Part: ACCBELT1003

13. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

13.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Type

Basic Center Mount..... No Charge
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility.

Short Extension, 5"x9" Footplates No Charge
Part: MEC1906675

Long Extension, 5"x9" Footplates No Charge
Part: MEC1906677

Medium Extension, 5"x9" Footplates No Charge
Part: MEC1906680

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount..... \$158
Part: RIG2104369

Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997
Only available with Long Length Extension.

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16
Part: RIG2104240

13.1b. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected. Select legrest length in following section.

1. Center Mount Foot Platform Type

Power Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012

Quick Adjustable Center Mount Foot Platform \$263
Part: XRF000311

2. Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility. Do not select with Basic Center Mount.

Medium Extension, 5"x6" Footplates
Part: RIG150206

Long Extension, 5"x6" Footplates
Part: RIG151505

Medium Extension, 5"x8" Footplates
Part: RIG2404016

Long Extension, 5"x8" Footplates
Part: RIG2404017

Medium Extension, 5"x11" Footplates
Part: RIG2404018

Long Extension, 5"x11" Footplates
Part: RIG2404019

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Support

Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads (pair) 6"Wx6"H..... No Charge
Part: RIGASMB7120018

Large Calf Pads (pair) 8" taper to 6"Wx8"H No Charge
Part: RIGASMB7120019

3b. Center Mount Calf Panel

Small Calf Panel, 18"Wx8"H..... \$184
Part: FRMASMB7110012

Large Calf Panel, 24"Wx8"H \$237
Part: RIG123825
Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
Part: RIG152519
Only available with the 5x6 footplates.

Outside Edge \$69
Part: RIG2404921
Available with 5"x8" footplates.

Width Extension for 8" \$69
Part: RIG2404919

Front and Rear Cap \$16
Part: RIG2404923
Available with 5"x8" footplates.

Width Extension for 11" \$69
Part: RIG2404920

5. Additional Accessories

Small Foam Calf Pad Cover (pair) \$67
Part: ST-CSCS0775-Q

Footplate Riser \$75
Part: XRF000348

Large Foam Calf Pad Cover (pair) \$83
Part: ST-CSCS0985-Q

Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-Q

Footplate Riser Calf Panel \$50
Part: FRMASMB7110011

Large Foam Panel Cover (pair) \$282
Part: ST-CSCP0927-Q

Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

Footplate Connector \$16
Part: RIGASMB7120023

Legrests section continued on following pages.

13. LEGRESTS CONTINUED**6. Stealth Gel Footplate Covers**

5"x8" Tapered Gel Footplate Covers \$222
 Part: ST-FPG5008T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers \$270
 Part: ST-AFPG5011T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension \$247
 Part: ST-AFPG5008E-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Cover for Width Extension \$270
 Part: ST-AFPG5011E-Q. HCPCS: K0108

13.1c. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

☐ ☐ ☐ 60°-90° Swing-Away Legrest No Charge
 Pair part: RIG120062. Left: FRMASMB17558.
 Right: FRMASMB17559
 Not available with 11"-13" extensions.

☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short No Charge
 Pair part: RIG136481. Left: RIG132607.
 Right: RIG132609.
 For use only with 11"-13" lower extensions.

2. Lower Extension

☐ ☐ ☐ 11"-13" Length
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.

☐ ☐ ☐ 13"-16.5" Length
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.

☐ ☐ ☐ 15.5"-19.5" Length
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ 5"Wx8"D Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: E0951

☐ ☐ ☐ 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: E0951

☐ ☐ ☐ 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: E0951

☐ ☐ ☐ 8"Wx11"D² Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
 Part: INDPART2823. HCPCS: K0108

13.2. Legrest Accessories

Heel Loops (pair) \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper \$37
 Part: XRF000100

13.3 Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)	18" L (#3492)	19" L (#3494)
11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	17.5" L (#3491)	18.5" L (#3493)	19.5" L (#3495)

14. POWER ELECTRONICS

REQUIRED.

14.1a. Tilt Only through Toggle

Single Actuator Control, Power Tilt Only, Through Toggle No Charge
 Part: ELE1905907
 Available with QLNES Joystick. Select toggle side.

1. Toggle Mount Side Selection

Left Side Toggle
 Part: ELEASMB5156

Right Side Toggle
 Part: ELEASMB5157

14.1b. Actuator Control through Joystick

Single Actuator Through Controller and Interface \$2,163
 Part: XRF000335. HCPCS: E2310

Multiple Actuator Through Controller and Interface \$3,040
 Part: XRF000336. HCPCS: E2311

15. ACCESSORIES**Optional.**

Accessory Bar\$132
 Part: XRF000128
 Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.

Cup Holder \$21
 Part: XRF000354

Large, Black, Weather Cover \$202
 Part: ACCCOVR1021

Oxygen Holder..... \$210
 Part: ACC125006. HCPCS: E2208
 Includes accessory bar, a credit will be applied PTOINDV3769.

Personal Item Pouch\$84
 Part: ACC2101017

Privacy Bag\$84
 Part: ACC2005723

Q-Logic EX/NE Programmer \$573
 Part: ELEASMB6976

Transfer Handle (D).....\$105 ea.
 Part: MEC2205680 x2. HCPCS: E0705

Remote Stop Switch \$1,295
 Part: XRF000379

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
 Part: ACC2407046

Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
 Part: ACC2405496

16. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#).**16.1. Stealth Thoracic Laterals**

Laterals are not available with the omit back option.

1. Bracket

Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664
 Part: ST-ULPK404108B00-Q. HCPCS: E1034
 Only available with TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
 Part: ST-ULPK404104B00-Q. HCPCS: E1034
 Only available with Synergy Back and omit back configurations.

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$250
 Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250
 Part: INDPART2664. HCPCS: E0956

3. Installation

Install Laterals
 Part: PTOINDV3499

Do Not Install
 Part: PTOINDV3500

4. Lateral Extensions**REQUIRED with Synergy Back and thoracic lateral selections.**

Lateral Extension for Synergy Back (pair)..... \$58
 Part: XRF000150

16.2. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea.
 Removable Mounting Hardware
 Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Left..... \$332 ea.
 Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^

UniLink, 7" Height Adjustable, \$276 ea.
 Removable Mounting Hardware
 Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Right \$332 ea.
 Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
 Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$281 ea.
 Feature (TWBADD-L & LPHW407)
 Part: POS140428. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea.
 Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact..... \$281 ea.
 Feature (TWBADD-R & LPHW407)
 Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.
 Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^

UniLink, 4" x 12" Pad \$106 ea.
 Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^

UniLink, 4" x 8" Pad..... \$106 ea.
 Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea.
 Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$81 ea.
 Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

17. LAP TRAY

17.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q, HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
Part: ST-SSLXM010001171-Q, HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
Part: ST-SSLXL010001171-Q, HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	ML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.