



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

J4E HD Order Form with Power Positioning Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____



A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

1. BASE MODEL

REQUIRED. 400 lb. weight capacity with standard 6 mph programmable motors, 8 amp off board charger, and flat free drive wheels. Dimensions of base without legrests: 24.5"Wx32.3"L.

J4E HD 2SPHD-SS..... \$8,180

Part: J4E HD 2SPHD-SS. HCPCS: K0837

Single power option.

2. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

Red Corvette

Part: DGN152613-PRD1

Electric Blue

Part: DGN152613-PBU9

Back in Black

Part: DGN152613-PBK1

3. BATTERIES**Optional.** The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor 2 Battery \$355 ea.
 Part: BAT2004376x2. HCPCS: E2361

4. QUANTUM BACKUP CAMERA**Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
 Part: MEC1910365

1. Optional Mounting Bracket

Backup Camera Headrest Bracket.....\$21
 Part: MEC1910373
 Works with .5" x .5" vertical supports only.

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
 Part: MEC2106123

**5. TRANSIT OPTIONS****Optional.** The power base is standard without transit securement.

Unoccupied Transit No Charge
 Part: KIT130229

Allows the powerchair to be secured for **unoccupied** transit.

WC-19 Occupied Transit Option..... \$289
 Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

6. ELECTRONICS**REQUIRED.** Expandable controls are special order, please contact your inside sales representative.**6.1a. Non-Expandable Controls**

QLNES 6 Key Joystick Standard
 Part: CTL2205379
 Controls up to 2 actuators.

6.1b. Expandable Controls**1. Expandable Controls Required Components****Power Tilt, Power Lift, and Power Articulating Foot Platform combined selections require the controller module and expandable harness below.**

Q-Logic 3 Controller Module \$730
 Part: XRF000002. HCPCS: E2377

Q-Logic 3 Expandable Controller Harness..... \$458
 Part: XRF000264. HCPCS: E2313

2. Q-Logic 3 Joystick

Q-Logic 3e Joystick \$735
 Part: CTL2301467

6.1c. Specialty Controls

Specialty controls are special order, please contact your inside sales representative.

6.2. Magnetic Charger**Available with QLNES and Q-Logic 3e joystick.**

Magnetic to XLR Charger Adapter with Wall Holder \$60
 Part: ELE2103635

7. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**7.1a. Fixed, In-line Joystick Bracket**

Fixed, In-line, Bracket Standard
 Part: FRMASMB12744

7.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
 Part: ACC1908434. HCPCS: E1032

7.1c. Basic Swing Away Joystick Mount


Basic Swing Away, Q-Logic, Left Side \$245
 Part: ACC2201041. HCPCS: E1032

Basic Swing Away, Q-Logic, Right Side \$245
 Part: ACC2201042. HCPCS: E1032

7.2. Joystick Mount Location

Left Side Mount
 Part: PTOINDV1225

Right Side Mount
 Part: PTOINDV1226



8. JOYSTICK HANDLES

Optional. The joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323

Chin Cup (C) \$90
Part: ACCASMB1861. HCPCS: E2323. (E2324 when chin control)

T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323


Medium Ball Handle (E) \$69
Part: ACCASMB1859. HCPCS: E2323

Soft Ball Handle (F) \$80
Part: ACCASMB1856. HCPCS: E2323

Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323

Bodypoint 3" Goal Post Handle \$104
Part: PLSKNCP1057. HCPCS: E2323

Bodypoint 4" Goal Post Handle \$104
Part: PLSKNCP1058. HCPCS: E2323



Joystick Handles

9. POWER POSITIONING SEATING**9.1. Q4 Seat**

REQUIRED. Select one option.

Q4 Static Seat \$2,115
Part: PTOINDV4008
Requires adjustable seat lift.

Q4 Seat HD Power Tilt \$6,121
Part: PTOINDV3975. HCPCS: E1002
55° center of gravity (CG) power tilt system.

9.2. 10" Power Lift

Optional. ^ Please verify the appropriate billing code to be used for this item with all third-party payors.

Q4 Seat 10" HD Adjustable Seat Lift \$3,725
Part: PTOINDV4020. HCPCS: E2298/E2300^
10" of scissor elevate. When elevated, maximum speed will be up to 1.25 mph.

9.2. Seat Width

REQUIRED. Maximum seat width is 24". Width package range will be 20"-24" wide.

20" Wide
Part: SETWDTH1010

21" Wide
Part: SETWDTH1011

22" Wide
Part: SETWDTH1012

23" Wide
Part: SETWDTH1013

24" Wide
Part: SETWDTH1014

9.3. Seat Depth

REQUIRED. Maximum seat depth is 24". Seat depth package will be 20"-24".

20" Deep
Part: SETDPTH1010

21" Deep
Part: SETDPTH1011

22" Deep
Part: SETDPTH1012

23" Deep
Part: SETDPTH1013

24" Deep
Part: SETDPTH1014

9.4. Seat to Floor Height

REQUIRED.

Lowest Seat to Floor
Part: PTOINDV3558

Middle Seat to Floor
Part: PTOINDV3559


Highest Seat to Floor
Part: PTOINDV3560

Low	Med	High
17.5"	18.5"	19.5"

Seat to Floor Heights

9.5. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.



Power positioning seating section continued on following pages.

9. POWER POSITIONING SEATING CONTINUED**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)

Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPPReference part: XRF000118. [See cushion order form.](#)

Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPPReference part: XRF000113. [See cushion order form.](#)

Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPPReference part: XRF000112. [See cushion order form.](#)

Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPPReference part: XRF000117. [See cushion order form.](#)

Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere PReference part: XRF000122. [See cushion order form.](#)

Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SPReference part: XRF000114. [See cushion order form.](#)

Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SPReference part: XRF000115. [See cushion order form.](#)

Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity GReference part: XRF000111. [See cushion order form.](#)

Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39

Part: POP1909006-S

14"W and smaller

Medium Wedge Positioning Kit.....\$72

Part: POP1909006-M

14"W and smaller

Large Wedge Positioning Kit..... \$105

Part: POP1909006-L

14"W and smaller

Small Wedge Positioning Kit\$39

Part: POP1712340-S

15"W and larger

Medium Wedge Positioning Kit.....\$72

Part: POP1712340-M

15"W and larger

Large Wedge Positioning Kit..... \$105

Part: POP1712340-L

15"W and larger

3. Cushion Top Foam**Visco Top Foam**

Reference part: XRF000380

No Top Foam

Reference part: XRF000382

FRG Top Foam

Reference part: XRF000381

4. Cushion Liner**Stealth Urethane Liner**

Reference part: XRF000383

No Liner

Reference part: XRF000385

Stealth Polyurethane Liner

Reference part: XRF000384

5. Cushion Cover**CoolCore Cover**

Reference part: XRF000386

Stealth Tek Cover

Reference part: XRF000388

Mesh Cover

Reference part: XRF000387

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner					
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer

Reference part: XRF000389

MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner

Reference Part: XRF000390

Additional Stealth Polyurethane Liner

Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover

Reference Part: XRF000392

Additional Stealth Tek Cover

Reference Part: XRF000394

Additional Mesh Cover

Reference Part: XRF000393

9.6. Back Cane Type**REQUIRED.****Angle Adjustable No Charge**
Part: MEC140266**Flip Forward Back Canes \$184**
Part: XRF000274**2" Endomorph (1" per side) Back Cane \$74**
Part: XRF000269**4" Endomorph (2" per side) Back Cane..... \$74**
Part: XRF000271

Power positioning seating section continued on following page.

9. POWER POSITIONING SEATING CONTINUED**9.7. Back Cane Height****REQUIRED.**

16"H Cane
Part: FRMASMB15834

17"H Cane
Part: FRMASMB15835

20"H Cane
Part: FRMASMB15838

21"H Cane
Part: FRMASMB15839

24"H Cane
Part: FRMASMB15842

25"H Cane
Part: FRMASMB15843

9.8. Back Cushion Type**REQUIRED.****1a. Synergy Back Selection**

Synergy Back.....No Charge
Part: XRF000199
Select overall back height in later section. Back cushion width is determined by seat width and back cane type.

1b. TRU-Comfort 4 Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

TRU-Comfort 4 Back, Back Cane Option, 16" to 20" Width \$744
Part: TC4V10W21WBC. HCPCS: E2620

TRU-Comfort 4 Back, Back Cane Option, 22" Width..... \$855
Part: TC4V22W32WBC. HCPCS: E2621

1c. Omit Back Selection

Omit Back Pan/Cushion.....No Charge
Part: PTOINDV1287

9.9. Overall Back Cushion Height

REQUIRED with back selection. Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion.

16" Height Part: SETHGHT1029	21" Height Part: SETHGHT1034
17" Height Part: SETHGHT1030	22" Height Part: SETHGHT1035
18" Height Part: SETHGHT1031	23" Height Part: SETHGHT1036
19" Height Part: SETHGHT1032	24" Height Part: SETHGHT1037
20" Height Part: SETHGHT1033	25" Height Part: SETHGHT1038

9.10. Power Tilt Privacy Flap

Optional. A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap\$48
Part: XRF000351

9.11. Lateral Extensions

REQUIRED with thoracic lateral selections.

Lateral Extension for Synergy Back (pair)..... \$56
Part: XRF000150

10. Q4 SEAT ARMRESTS

REQUIRED. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type when applicable.

10.1a. Drop-In, Height Adjustable, Removable Armrests

Not available with toggle switch. Includes straight armpads. Do not select options from armrest size nor armpads sections.

☐ ☐ Full Length, Drop-In Armrests (9"-14" Height) No Charge L R Left part: ARM1803253. Right part: ARM156742	☐ ☐ Desk Length, Drop-In Armrests (9"-14" Height)..... No Charge L R Left part: ARM1803254. Right part: ARM156743
--	---

10.1b. Height Adjustable, Removable, Heavy Duty Armrests

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable, \$79 ea.
Heavy Duty (Sold in pairs)
Part: XRF000055. HCPCS: E0973

10"-14"H, Single Post, Height Adjustable, Removable,..... \$79 ea.
Heavy Duty (Sold in pairs)
Part: XRF000056. HCPCS: E0973

10.1c. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals

☐ ☐ Full Length, Flip Back Armrests..... \$63 L R Left part: FRMASMB7038. Right part: FRMASMB7039	☐ ☐ Desk Length, Flip Back Armrests \$63 L R Left part: FRMASMB7040. Right part: FRMASMB7041
---	--

		Back Width									
		16"	17"	18"	19"	20"	21"	22"	23"	24"	
Overall Height	16"	T	T	T	T	T	-	T	-	-	
	17"	S	T	S	T	S	T	-	S	T	S
	18"	T	T	T	T	T	-	T	-	-	
	19"	T	T	T	T	T	S	S	T	S	S
	20"	T	T	T	T	T	S	S	T	S	S
	21"	S	T	S	T	S	T	S	S	T	S
	22"	T	T	T	T	T	S	S	T	S	S
	23"	T	T	T	T	T	S	S	T	S	S
	24"	S	T	S	T	S	T	S	S	S	S
	25"	S	T	S	T	S	T	S	S	S	S

Legend: S Synergy Back T TRU-Comfort Back

Q4 Seat Back Availability Chart

Q4 Seat armrests section continued on following pages.

10. Q4 SEAT ARMRESTS CONTINUED**10.2. Armrest Height****REQUIRED with Height Adjustable, Removable, Heavy Duty Armrests.** Set left and right armrest height.**1. Left Armrest Height**

8" H PTOINDV3535	9" H PTOINDV3536	10" H PTOINDV3537	11" H PTOINDV3538	12" H PTOINDV3539	13" H PTOINDV3540	14" H PTOINDV3541
---------------------	---------------------	----------------------	----------------------	----------------------	----------------------	----------------------

2. Right Armrest Height

8" H PTOINDV3450	9" H PTOINDV3451	10" H PTOINDV3452	11" H PTOINDV3453	12" H PTOINDV3454	13" H PTOINDV3455	14" H PTOINDV3456
---------------------	---------------------	----------------------	----------------------	----------------------	----------------------	----------------------

10.3. Armrest Size**REQUIRED with Height Adjustable, Removable, Heavy Duty Armrests.** Select a left (L) and right (R) armrest size. This also affects armpad selection.

10" Length	14" Length
L R Left part: XRF000057. Right part: XRF000123	L R Left part: XRF000058. Right part: XRF000124

10.4. TRU-Balance Armpads**REQUIRED with Height Adjustable, Removable, Heavy Duty Armrests.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

No Armpad Selection..... No Charge Part: XRF000060	L R Waterfall, Armpad No Charge Left part: XRF000062. Right part: XRF000127 Not compatible with flip-back, cane mounted armrests.
L R Straight, Armpad Standard Left part: XRF000061. Right part: XRF000126	

11. HEADRESTS**Optional.** Select a headrest pad and bracket.**11.1. Headrest Pads**

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179 Part: ST-CP08CCFPM-Q. HCPCS: E0955	Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179 Part: ST-CP14CCFPM-Q. HCPCS: E0955
Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)\$179 Part: ST-CP10CCFPM-Q. HCPCS: E0955	Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179 Part: ST-CP06CCFPM-Q. HCPCS: E0955

11.2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink.....\$311 Part: ST-UHB2D33200-Q. HCPCS: E1033	Multi-Axis Removable Headrest Bracket.....\$174 Part: ST-HMO475-17M-Q. HCPCS: E1033
--	--

11.3. Q Cam Lock**Optional.**

Q Cam Lock.....\$50 Part: XRF000376
--

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport	Q-Back 2
	Standard \$50	Standard \$50
Tilt/Static/Back Canes	14"-21"H 22"-29"H	19"-24"H 25"-26"H
Power Recline Back	19"-21"H 22"-32"H	22"-25"H 26"-29"H

12. LAP BELT**REQUIRED.****12.1a. Standard Lap Belts**

60" Lap Belt Part: ACCBELT1002	70" Lap Belt Part: ACCBELT1003	80" Lap Belt Part: ACCBELT1004	90" Lap Belt Part: ACCBELT1001
-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------

12.1b. Optional Lap Belts

60" Sure-Lok Lap Belt (Non-retractable).....\$120 Part: ACCBELT1016 Required with ISO 7176-19 Compliant Transit Option.

13. LEGRESTS**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**13.1a. Base Mounted Foot Platform**

May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard Part: FRMASMB16077	Omit base mounted foot platform No Charge Part: PTOINDV3674
---	--

13.1b. Basic Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Type

Basic Center Mount..... No Charge Part: XRF000323
--

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

Legrests section continued on following pages.

13. LEGRESTS CONTINUED**2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates**

See legrest length matrix for compatibility.

Short Extension, 5"x9" Footplates No Charge
Part: MEC1906675Long Extension, 5"x9" Footplates No Charge
Part: MEC1906677Medium Extension, 5"x9" Footplates No Charge
Part: MEC1906680**3. Basic Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount \$158
Part: RIG2104369Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997
Only available with Long Length Extension.**4. Basic Center Mount Additional Accessories**Basic Center Mount Footplate Connector \$16
Part: RIG2104240**13.1c. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected. Select legrest length in following section.

1. Center Mount Foot Platform TypePower Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012Adjustable Center Mount Foot Platform No Charge
Part: PTOINDV3755**2. Extension and Angle Adjustable Footplates**

See legrest length matrix for compatibility. Do not select with Basic Center Mount.

Medium Extension, 5"x6" Footplates
Part: RIG150206Long Extension, 5"x6" Footplates
Part: RIG151505Medium Extension, 5"x8" Footplates
Part: RIG2404016Long Extension, 5"x8" Footplates
Part: RIG2404017Medium Extension, 5"x11" Footplates
Part: RIG2404018Long Extension, 5"x11" Footplates
Part: RIG2404019

Medium	Long
11"-17"	14"-17"
Power AFP Extension Range	

Medium	Long
11"-17"	14"-17"
Manually Adjustable Extension Range	

3a. Center Mount Calf Support

*Calf pads are standard with a power articulating foot platform. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019Small Calf Pads, Adjustable Center Mount Only, \$158
6"Wx6"H (pair)
Part: XRF000156Large Calf Pads, Adjustable Center Mount Only, \$158
8" taper to 6"Wx8"H (pair)
Part: XRF000157**3b. Center Mount Calf Panel**

Calf pads are standard with a power articulating foot platform.

Small Calf Panel, 18"Wx8"H \$184
Part: FRMASMB7110012Large Calf Panel, 24"Wx8"H \$237
Part: RIG123825
Only available on widths of 20" or greater.**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
Part: RIG152519
Only available with the 5x6 footplates.Outside Edge \$69
Part: RIG2404921
Available with 5"x8" and 5"x11" footplates.Width Extension for 8" \$69
Part: RIG2404919Front and Rear Cap \$16
Part: RIG2404923Width Extension for 11" \$69
Part: RIG2404920

Available with 5"x8" and 5"x11" footplates.

5. Additional AccessoriesSmall Foam Calf Pad Cover (pair) \$67
Part: ST-CSCS0775-QLarge Foam Panel Cover (pair) \$282
Part: ST-CSCP0927-QLarge Foam Calf Pad Cover (pair) \$83
Part: ST-CSCS0985-QFootplate Connector \$16
Part: RIGASMB7121023Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-Q**6. Stealth Gel Footplate Covers**5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K01085"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K01085"x11" Tapered Gel Footplate Covers \$270
Part: ST-AFPG5011E-Q. HCPCS: K01085"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-AFPG5011T-Q. HCPCS: K0108

Legrests section continued on following page.

13. LEGRESTS CONTINUED**13.1d. Heavy Duty, Drop-In, Swing-Away Legrests**

Select a configuration, lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Drop-In Legrest Selection

☐ ☐ ☐ 70° HD Drop-in for 8.5"-12" Lower Extension No Charge
 Pair part: FRMASMB4762. Left: RIG137706. Right: RIG137707

☐ ☐ ☐ 70° HD Drop-in for 13"-18" or 16"-21" Lower No Charge
 Extension
 Pair part: FRMASMB4853. Left: RIG137708. Right: RIG137709

2. Lower Extension

☐ ☐ ☐ 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103.

☐ ☐ ☐ 16"-21" HD Drop-Ins/18.5"-22.5" ELR, Ext Length
 Pair part: FRMASMB4754. L: FRMASMB4725. R: FRMASMB5393

☐ ☐ ☐ 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106

3. Angle Adjustable Footplate(s) with Heel Loops

☐ ☐ ☐ 6"Wx8"D \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.
 HCPCS: E0951

☐ ☐ ☐ 8"Wx11"D \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799.
 HCPCS: E0951

☐ ☐ ☐ 6"Wx11"D \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.
 HCPCS: E0951

13.1e. Heavy Duty, Elevating Legrests

Select a configuration, lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Heavy Duty, Elevating Legrest Selection

☐ ☐ ☐ Heavy Duty, Elevating Legrest(s) with Calf Pad(s) \$171 ea.
 Pair part: FRMASMB6702. Left: FRMASMB8444. Right: FRMASMB8445. HCPCS: E0990

2. Lower Extension

☐ ☐ ☐ 11"-14" Length
 Part: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced footplates only.

☐ ☐ ☐ 14"-17" Length
 Part: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced footplates only.

☐ ☐ ☐ 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103.
 For angle adjustable footplates only.

☐ ☐ ☐ 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable footplates only.

3a. Reinforced Footplates(s)

☐ ☐ ☐ 7.25"Wx6"D Extruded Footplate Standard
 Part pair: FRMASMB4748. Left: RIG138055.
 Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Extruded Footplate Standard
 Part pair: FRMASMB4747. Left: RIG138060.
 Right: RIG138062

3b. Angle Adjustable Footplates(s) with Heel Loops

☐ ☐ ☐ 6"Wx8"D \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.
 HCPCS: E0951

☐ ☐ ☐ 8"Wx11"D \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799.
 HCPCS: E0951
 May cause interference with seat widths less than 20"

☐ ☐ ☐ 6"Wx11"D \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.
 HCPCS: E0951

13.2. Legrest Accessories

Heel Loops (pair) \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper \$37
 Part: XRF000100

13.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	17.5" L (#3491)
9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)	18" L (#3492)

14. POWER POSITIONING ELECTRONICS**REQUIRED.****14.1. Actuator Control through Joystick**

Single Actuator Through Controller and \$2,163
Interface
Part: XRF000335. HCPCS: E2310

Multiple Actuator Through Controller and \$3,040
Interface
Part: XRF000336. HCPCS: E2311

15. ACCESSORIES**Optional.**

Accessory Bar \$132
Part: XRF000128
Large, Black, Weather Cover \$202
Part: ACCCOVR1021
Oxygen Holder \$210
Part: ACC125006 HCPCS: E2208
Includes accessory bar, a credit will be applied (PTOINDV3769).
Personal Item Pouch \$84
Part: ACC2101017

Privacy Bag \$42
Part: ACC2005723
Remote Stop Switch \$1,295
Part: XRF000379
Transfer Handle \$105 ea.
Part: MEC2205680 x2. HCPCS: E0705

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
Part: ACC2407046

Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
Part: ACC2405496

16. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#).**16.1. Stealth Thoracic Laterals**

Laterals are not available with the omit back option.

1. Thoracic Laterals Bracket

Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664
Part: ST-ULPK404108B00-Q. HCPCS: E1034
Only available with TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034
Only available with Synergy Back and omit back configurations.

2. Thoracic Laterals Pads

Curved Pads, 3"x5", pair (TWBL-3X5) \$250
Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6) \$250
Part: INDPART2664. HCPCS: E0956

3. Thoracic Laterals Installation

Install Laterals Do Not Install
Part: PTOINDV3499 Part: PTOINDV3500

16.2. Stealth Unilink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea.
Removable Mounting Hardware
Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Left \$332 ea.
Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^

UniLink, 7" Height Adjustable, \$276 ea.
Removable Mounting Hardware
Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Right \$332 ea.
Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$281 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Right \$358.87 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact \$281 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.
Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^

UniLink, 4" x 12" Pad \$106 ea.
Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^

UniLink, 4" x 8" Pad \$106 ea.
Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea.
Part: ST-UIL001150A000-Q

UniLink, 2" Link \$81 ea.
Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left Install on Right Do Not Install
Part: PTOINDV3585 Part: PTOINDV3586 Part: PTOINDV3719

17. LAP TRAY**17.1. Simple Solutions Lap Tray**

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap TraySimple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q, HCPCS: E0950Simple Solutions Large Tray with Cut Out..... \$309
Part: ST-SSLXL010001171-Q, HCPCS: E0950Simple Solutions Medium Tray with Cut Out..... \$309
Part: ST-SSLXM010001171-Q, HCPCS: E0950**2. Mounting Hardware**Swing Away Mounting Hardware..... \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	SM	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

S:Compatible

M:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

