

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

J6 Order Form with Contoured and Synergy Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering. Options noted with "XRef" have multiple possible part numbers based on system configuration. Please see [cross-reference spreadsheet](#) for the XRef part numbers.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 2 power base. 300 lb. weight capacity. Includes 5 amp off board charger. Dimensions of base without legrests: 23.19"W x 34.74"L.

J6 2S-C \$5,971

Part: J6 2S-C. HCPCS: K0823

Captains seat standard.

J6 2S-SS \$5,971

Part: J6 2S-SS. HCPCS: K0822

Solid seat pan standard.

2. BASE COLOR**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.**Red Corvette**
Part: SHDASMB4630**Back in Black**
Part: SHDASMB4632**Evolution Orange**
Part: DGN1903758-POR7**Electric Blue**
Part: DGN1903758-PBU9**Pretty in Pink**
Part: DGN1903758-PPK3**3. BASE MOUNTED FOOT PLATFORM****REQUIRED.** May interfere with seat mounted legrests. Not available with center mount foot platforms.Base Mounted Foot Platform Standard
Part: FRMASMB16077Omit Base Mounted Foot Platform No Charge
Part: PTOINDV3674Bell Shaped, Base Mounted Foot Platform \$235
Part: RIGASMB7120041**4. BATTERIES****REQUIRED.** The powerchair requires 2 batteries and will be installed if ordered.☐ U-1, 32 Amp Interceptor Batteries (Installed) \$133 ea.
Part: BAT1704611 x2 qty. HCPCS: E2365☐ No Batteries No Charge
Part: Configuration selection**5. DRIVE WHEELS****REQUIRED.** Select a type of drive wheels.10" Flat Free Black Tires Standard
Part: WHL140393 x2 qty10" Pneumatic Black Tires No Charge
Part: WHL140390 x2 qty**6. UNOCCUPIED BASE TRANSIT****Optional.**

Unoccupied Transit Securement No Charge

Limited Recline High Back part: MEC131819. Synergy part: MEC147382. Synergy Endo/Meso part: MEC147775

Only available with limited recline high back and Synergy seat.

7. ELECTRONICS**REQUIRED.** Other drive controls may be available by contacting customer service.**Non-Expandable Controls**



Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick Standard
Part: CTLDC1554NE+ 8-Key Joystick (2 Actuators & Lights) \$400
Part: CTLDC1561NE+ 6-Key Joystick (2 Actuators) \$200
Part: CTLDC1560No Joystick, Setup for NE+ Joystick No Charge
Part: Configuration selection


NE+ 6 Key Joystick

8. FULL LIGHT PACKAGE**Optional.** Available with Synergy Seating only. Requires NE+ 8-Key Joystick.☐ Seat Mounted, LED, Full Light Package \$995
Part: ELEASMB5719**9. JOYSTICK MOUNTING BRACKETS****REQUIRED.** Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket

9.1a. Fixed, In-line Joystick BracketFixed, In-line, Bracket Standard
NE/NE+ Part: FRMASMB12744Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226Swing Away Joystick
Inline Mount (Right)Swing Away Joystick
Offset Mount (Right)**9.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

Swing-away, Inline, Left Side Mount \$245
Part: See XRef. HCPCS: E1028Swing-away, 1" Offset, Left Side Mount \$245
Part: See XRef. HCPCS: E1028Swing-away, Inline, Right Side Mount \$245
Part: See XRef. HCPCS: E1028Swing-away, 1" Offset, Right Side Mount \$245
Part: See XRef. HCPCS: E1028

Retract4

9.1c. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount \$350
Part: MEC138876. HCPCS: E1028Multi-axis Mount \$100
Part: ACC142528. HCPCS: E1028

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Allows height adjustment of joystick.

Multi-Axis Joystick
Mount**Specialty Mount Side**Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226

10. JOYSTICK HANDLES

Optional. The NE/NE Plus Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B)	\$85	Soft Ball Handle (F)	\$65
Part: ACCASMB1860. HCPCS: E2323		Part: ACCASMB1856. HCPCS: E2323	
Chin Cup (C)	\$85	Extended I-Handle (G)	\$85
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)		Part: ACCASMB1857. HCPCS: E2323	
T-Bar (D)	\$120	Bodypoint 3" Goal Post Handle	\$99
Part: ACCASMB1858. HCPCS: E2323		Part: PLSKNCP1057. HCPCS: E2323	
Medium Ball Handle (E)	\$65	Bodypoint 4" Goal Post Handle	\$99
Part: ACCASMB1859. HCPCS: E2323		Part: PLSKNCP1058. HCPCS: E2323	

**11. CONTOURED SEATING**

Optional. Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.

11.1a. Captain Style Seating**1a. Depth Adjustable High Back with Limited Recline**

Full length armrests must be selected with a depth adjustable high back.

18"Wx18"-20"D with 60" Lap Belt	No Charge	20"Wx18"-20"D with 60" Lap Belt	No Charge
Part: SETASMB6895		Part: SETASMB6896	

1c. High Back with Limited Recline

16"Wx16"D	No Charge	16"Wx20"D	No Charge	20"Wx16"D	No Charge
Part: SETASMB2770		Part: SETASMB2772		Part: SETASMB2780	
16"Wx18"D	No Charge	18"Wx16"D	No Charge		
Part: SETASMB2771		Part: SETASMB2775			

1c. Medium Back

16"Wx16"	No Charge	16"Wx18"	No Charge	18"Wx16"	No Charge	18"Wx18"	No Charge
Part: SETASMB2488		Part: SETASMB2560		Part: SETASMB2159		Part: SETASMB2551	

11.1b. Contoured Seat with Solid Seat Pan**1a. High Back with Limited Recline and Solid Seat Pan**

16"Wx16"D	No Charge	18"Wx16"D	No Charge	20"Wx16"D	No Charge
Part: SETASMB3178		Part: SETASMB3181		Part: SETASMB3183	
16"Wx18"D	No Charge	18"Wx18"D	No Charge	20"Wx18"D	No Charge
Part: SETASMB3179		Part: SETASMB2848		Part: SETASMB3184	
16"Wx20"D	No Charge	18"Wx20"D	No Charge	20"Wx20"D	No Charge
Part: SETASMB3180		Part: SETASMB3182		Part: SETASMB3185	

1b. Medium Back with Solid Seat Pan

16"Wx16"	No Charge	16"Wx18"	No Charge	18"Wx16"	No Charge	18"Wx18"	No Charge
Part: SETASMB3112		Part: SETASMB3113		Part: SETASMB2492		Part: SETASMB3114	

11.1c. Special Order Contoured Seat

Complete information below. Contact inside sales specialist for pricing.

Type of Seat: _____ Width: _____ Depth: _____

11.2. Contoured Seat-to-Floor Height

REQUIRED with Contoured Seating without a Power Elevating Seat.

Low Hole Position	High Hole Position
Part: PTOINDV3558	Part: PTOINDV3560

11.3. Contoured Seat Flip-Up, Height Adjustable, Removable Armrests

REQUIRED with High Back Seat except with Depth Adjustable High Back. Select a Left and Right Armrest or to Omit them.

Full Length, Left Armrest	\$85	Full Length, Right Armrest	\$85
Part: FRMASMB8614. HCPCS: E0973		Part: FRMASMB8615. HCPCS: E0973	
Desk Length, Left Armrest	\$85	Desk Length, Right Armrest	\$85
Part: FRMASMB8612. HCPCS: E0973		Part: FRMASMB8613. HCPCS: E0973	
No Armrest	Standard		
Part: N/A			

12. SYNERGY SEATING SYSTEM**Optional.** Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 20"W x 20"D.**12.1. Synergy Frame Selection**Synergy Seating Frame \$434
Part: PTOINDV2995**12.2. Seat Dimension - Width****REQUIRED with Synergy seating.** Part number is SETWDTH plus number following # below.

10" W (#1000)	11" W (#1001)	12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.3. Seat Dimension - Width Package**REQUIRED with Synergy seating.**

10"-14" Width Part: FRMASMB18011	14"-18" Width Part: FRMASMB18012	16"-20" Width Part: FRMASMB18013
-------------------------------------	-------------------------------------	-------------------------------------

12.4. Seat Dimension - Depth**REQUIRED with Synergy seating.** Part number is SETDPTH plus number following # below.

10" D (#1000)	11" D (#1001)	12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.5. Seat Dimension - Depth Package**REQUIRED with Synergy seating.**

10"-14" Depth Package Part: FRMASMB6981	14"-18" Depth Package Part: FRMASMB6982	15"-20" Depth Package Part: FRMASMB7732
--	--	--

12.6. Synergy Seat to Floor**REQUIRED with Synergy seating.**

15.75" STF Part: PTOINDV3764	16.50" STF Part: PTOINDV1888	16.75" STF Part: PTOINDV3642	17.50" STF Part: PTOINDV1889
---------------------------------	---------------------------------	---------------------------------	---------------------------------

12.7. Synergy Back Options**REQUIRED with Synergy seating.** See the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1. Back Cane Selection**

Synergy Back Cane No Charge Part: FRMASMB15844 Select a back cane height.	2" Endomorph (1" per side) Back Cane \$195 Part: FRMASMB12115 Select a cane height.
Flip Forward Back Cane \$150 Part: FRMASMB15846 Select a back cane height. Not available with a vent tray, oxygen holder, or cane mounted and cantilever armrests.	4" Endomorph (2" per side) Back Cane \$195 Part: FRMASMB12139 Select a cane height.
2" Mesomorph (1" per side) Back Cane \$195 Part: FRMASMB12067 Select a cane height.	Manual Recline Back \$750 Part: FRMASMB12910. HCPCS: E1226 Select an overall back height. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, Shape Back, TRU-Comfort Back, rear accessories, power elevating seat, or articulating foot platform. Pivot point is 2.5" above the seat pan.
4" Mesomorph (2" per side) Back Cane \$195 Part: FRMASMB12091 Select a cane height.	

2. Synergy Back Cane Height

Omit this selection if a TRU-Comfort back is selected. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles.

14" High Canes Part: FRMASMB15832	17" High Canes Part: FRMASMB15835	20" High Canes Part: FRMASMB15838	23" High Canes Part: FRMASMB15841
15" High Canes Part: FRMASMB15833	18" High Canes Part: FRMASMB15836	21" High Canes Part: FRMASMB15839	24" High Canes Part: FRMASMB15842
16" High Canes Part: FRMASMB15834	19" High Canes Part: FRMASMB15837	22" High Canes Part: FRMASMB15840	25" High Canes Part: FRMASMB15843

3. Back Selection

Synergy Back \$275 Part: See XRef Select a back cane type and height and back height.	Sling Back \$98 Part: See XRef Select a back height. Available 18"-25" overall heights. Back cane height must match back height. Not available with flip up cantilever arms.
TRU-Comfort 4-way Stretch, 20" Width or less \$675 Part: See XRef. HCPCS: E2620 Select back height.	Omit Back No Charge Part: PTOINDV1287 Select a back cane type and height in following section.

Synergy seating section continued on following pages.

12. SYNERGY SEATING SYSTEM CONTINUED**4. Synergy Overall Back Height**

Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Parts: See XRef, configuration selection with back selection.

12" Overall Height	19" Overall Height
13" Overall Height	20" Overall Height
14" Overall Height	21" Overall Height
15" Overall Height	22" Overall Height
16" Overall Height	23" Overall Height
17" Overall Height	24" Overall Height
18" Overall Height	25" Overall Height

		Width							
		12"	14"	15"	16"	17"	18"	19"	20"
Overall Height	15"	S	S	S	-	-	-	-	-
	16"	S	S	S	-	-	-	-	-
	17"	-	S	S	S	S	S	S	S
	18"	-	S	S	S	S	S	S	S
	19"	-	S	S	S	S	S	S	S
	20"	-	S	S	S TC	S	S TC	S	S TC
	21"	-	S	S	S TC	S	S TC	S	S TC
	22"	-	S	S	S TC	S TC	S TC	S TC	S TC
	23"	-	S	S	S TC	S TC	S TC	S TC	S TC
	24"	-	S	S	S TC	S TC	S TC	S TC	S TC
	25"	-	S	S	S TC	S TC	S TC	S TC	S TC

S Synergy Back S TC Synergy or TRU-Comfort Back

Synergy and TRU-Comfort Back Availability Chart

12.8. Synergy Armrests

REQUIRED with Synergy seating. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

☐ ☐ Full Length, Flip Back Armrests..... Standard	☐ ☐ Desk Length, Flip Back Armrests Standard
☐ ☐ Left part: FRMASMB7038. Right part: FRMASMB7039	☐ ☐ Left part: FRMASMB7040. Right part: FRMASMB7041

1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐ Full Length, Two-Post Armrests.....\$150 ea.	☐ ☐ Desk Length, Two-Post Armrests\$150 ea.
☐ ☐ Left part: FRMASMB7062. Right part: FRMASMB7063. HCPCS: E0973	☐ ☐ Left part: FRMASMB7060. Right part: FRMASMB7061. HCPCS: E0973

1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H\$75 ea.	☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H.....\$75 ea.
☐ ☐ Left part: FRMASMB9982. Right part: FRMASMB9983. HCPCS: E0973	☐ ☐ Left part: FRMASMB9984. Right part: FRMASMB9985. HCPCS: E0973

1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H.....\$75 ea.	☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14"H\$75 ea.
☐ ☐ Left part: FRMASMB9986. Right part: FRMASMB9987. HCPCS: E0973	☐ ☐ Left part: FRMASMB9988. Right part: FRMASMB9989. HCPCS: E0973

1e. Pediatric, Removable, Height Adjustable Armrests

Only available with desk length armpads.

☐ ☐ Pediatric, Removable Armrests.....\$75 ea.	
☐ ☐ Left part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973	

1f. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

☐ ☐ Full Length, Cantilever Armrests.....\$97.50 ea.	☐ ☐ Desk Length, Cantilever Armrests.....\$97.50 ea.
☐ ☐ Left part: FRMASMB16188. Right part: FRMASMB16189	☐ ☐ Left part: FRMASMB16190. Right part: FRMASMB16191

1g. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests\$75 ea.	☐ ☐ Desk Length, Single Post Armrests\$75 ea.
☐ ☐ Left part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973	☐ ☐ Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973

12.9. Synergy Armpads

REQUIRED with armrest selection above. Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

Straight Armpads..... Standard	Waterfall Armpads..... No Charge
Desk/Desk part: UPSASMB1022. Full/Full part: UPSASMB1021. Desk/ Full part: UPSASMB1023	Desk/Desk part: UPSASMB1025. Full/Full part: UPSASMB1024. Desk/ Full part: UPSASMB1026
No Armpad Selection No Charge	
Part: PTOINDV1373	


12.10. Joystick Bracket Receiver

REQUIRED.

Standard Fixed Receiver Standard	Stealth Height Adjustable Receiver..... \$214
Part: FRMASMB16062	Part: ST-ARM250C-Q. HCPCS: K0108

13. HEADRESTS**Optional.** Select a headrest pad and bracket. **Not available with contoured seating.****13.1. Headrest Pads**

- ☐ Stealth 8" Comfort Plus (CP450)..... \$151.94
Part: POS148484. HCPCS: E0955
- ☐ Stealth 10" Comfort Plus (CP250) (A) \$151.94
Part: SETHEAD1057. HCPCS: E0955
- ☐ Stealth 14" Comfort Plus (CP150) \$151.94
Part: INDPART2327. HCPCS: E0955
- ☐ Pediatric Comfort Plus (6"Wx4.5"H) \$151.94
(CPN350)
Part: SETHEAD1069. HCPCS: E0955
- ☐ Stealth Small Contoured (9.75"Wx3"H) \$64.20
(STL-OTSM250)
Part: SETHEAD1052. HCPCS: E0955
- ☐ Stealth Large Contoured (11"Wx6"H) (B) \$80.25
(STL-OTLG150)
Part: SETHEAD1053. HCPCS: E0955

**13.2. Headrest Brackets**

- ☐ Stealth's Worlds Best Removable.....\$220.50
Headrest Bracket (TWB485-1) (C)
Part: FRMASMB15999. HCPCS: E1028
- ☐ Stealth Multi-Axis Removable (D).....\$160.50
Headrest Bracket (HMO475P)
Part: HDWASMB1084. HCPCS: E1028
- ☐ Pediatric Offset (HMO475P-17) (E) \$160.50
Part: HDWASMB1166. HCPCS: E1028

14. LAP BELTS**REQUIRED** except with select contoured seats with included lap belt. Select one lap belt below.**14.1a. Standard Lap Belts**

- ☐ 38" Lap Belt
Part: ACCBELT1006
- ☐ 50" Lap Belt
Part: ACCBELT1000
- ☐ 60" Lap Belt
Part: ACCBELT1002
- ☐ 70" Lap Belt
Part: ACCBELT1003
- ☐ 80" Lap Belt
Part: ACCBELT1004
- ☐ 90" Lap Belt
Part: ACCBELT1001

14.1b. Optional Lap Belts

Not available with contoured seating.

- ☐ Medium (7.5") Padded, 1.5" Belt (IFPB5213) \$86
Part: ST-IFPB5213-Q
- ☐ Extra Large Padded (11.5"), 2" Belt (IFPB5315)..... \$86
Part: ST-IFPB5315-Q
- ☐ Large Padded (9.5"), 2" Belt (IFPB5314)..... \$86
Part: ST-IFPB5314-Q

15. SYNERGY SEAT POWER AFP ELECTRONICS**REQUIRED** for Synergy Seat with Power Articulating Foot Platform.

- ☐ One Actuator Control through NE+ Joystick\$1,895
Part: ELEASMB6430. HCPCS: E2310

16. LEGRESTS**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**16.1. Center Mount Foot Platform**If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section. **Not available with contoured seating.****1. Foot Platform Assembly**

- ☐ Power Articulating Foot Platform \$2,995
Part: PTOINDV3754. HCPCS: E1012
- ☐ Manually Adjustable Center Mount\$250
Foot Platform
Part: RIG163640

**2. Lower Extension and Footplates**8"-12" Extension, 5"x6" Footplates
Part: RIG152577

Not available with power AFP.

10"-16" Extension, 5"x6" Footplates
Part: RIG15020610"-16" Extension, 5"x11" Footplates
Part: RIG1502128"-12" Extension, 5"x8" Footplates
Part: RIG152578

Not available with power AFP.

10"-16" Extension, 5"x8" Footplates
Part: RIG150203

Center Mount Foot Platform Compatibility Matrix	Synergy Seat			
	Low	Medium	High	
10"-16" Length Extension	!	✓	✓	☒ No interference issues at this seat-to-floor height. ☐ May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
13"-19" Length Extension	!	!	!	

3. Center Mount Calf Support

*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

- ☐ Small Calf Pads (pair) 6"Wx6"H \$150*
Part: RIGASMB7120018
- ☐ Small Calf Panel 18"Wx8"H \$150
Part: FRMASMB7110012
- ☐ Large Calf Pads (pair) 8" taper to 6"Wx8"H \$150*
Part: RIGASMB7120019
- ☐ Large Calf Panel 24"Wx8"H \$200
Part: RIG123824
Only available on widths of 20" or greater.

Legrests section continued on following page.

16. LEGRESTS CONTINUED**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

- | | |
|---|---|
| ☐ Width Extension for 8"D Footplates \$50
Part: RIG152517 | ☐ 4" Outside Edge (pair) \$50
Part: RIG152510
Only available with 6"D footplates and 8"D footplate extension. |
| ☐ Width Extension for 11"D Footplates \$50
Part: RIG152518 | ☐ 6" Outside Edge (pair) \$50
Part: RIG152516
Available with 8"D and 11"D footplates and 11"D footplate extension. |
| ☐ 2" Bell-Shaped Length Extension \$50
Part: RIG152519
Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates. | |

5. Additional Accessories

- | | |
|--|--|
| ☐ Reversed Dartex Calf Pad Cover (pair) \$150
Small part: RIG133717. Large part: RIG133719 | ☐ Wheel Bumpers (pair) \$12.50 ea.
Part: WHLASMB1578 |
| ☐ Reversed Dartex Panel Cover (pair) \$150
Small part: RIG134673. Large part: RIG134744 | ☐ Footplate Connector \$15
Part: RIGASMB7120023 |

16.2. 60°-90° Swing-Away Legrests

Select a legrest upper, a lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

- | | |
|--|---|
| ☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short \$81 ea.
☐ P ☐ L ☐ R
Part pair: RIG136481. Left: RIG132607. Right: RIG132609
Not available with 13"-16.5" extensions. | ☐ ☐ ☐ 60°-90° Swing-Away Legrest \$81 ea.
☐ P ☐ L ☐ R
Part pair: RIG120062. Left: FRMASMB17558. Right: FRMASMB17559
Not available with 11"-13" extensions. |
|--|---|

2. Lower Extension

- | | | |
|---|---|---|
| ☐ ☐ ☐ 11"-13" Length
☐ P ☐ L ☐ R
Part: See XRef
N/A with Multi-axis. | ☐ ☐ ☐ 11"-16.5" Length
☐ P ☐ L ☐ R
Part: See XRef
N/A with Angle Adjustable. | ☐ ☐ ☐ 13"-16.5" Length
☐ P ☐ L ☐ R
Part: See XRef
N/A with Multi-Axis. |
|---|---|---|

3a. Extruded Footplate(s)

- | | |
|--|---|
| ☐ ☐ ☐ 7.25"Wx6"D Extruded Footplate Standard
☐ P ☐ L ☐ R
Part pair: FRMASMB7847. Left: FRMASMB7861
Right: FRMASMB7862 | ☐ ☐ ☐ 8.5"Wx6"D Extruded Footplate Standard
☐ P ☐ L ☐ R
Part pair: FRMASMB7848. Left: FRMASMB7863
Right: FRMASMB7864
May cause interference on seat widths less than 18." |
|--|---|

3b. Angle Adjustable Footplate(s)¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- | | |
|--|--|
| ☐ ☐ ☐ 4"Wx6"D Angle Adjustable \$85 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560 | ☐ ☐ ☐ 6"Wx8"D ¹ Angle Adjustable \$85 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563 |
| ☐ ☐ ☐ 5"Wx6"D Angle Adjustable \$85 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561 | ☐ ☐ ☐ 6"Wx11"D ¹ Angle Adjustable \$85 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564 |
| ☐ ☐ ☐ 5"Wx8"D Angle Adjustable \$85 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562 | ☐ ☐ ☐ 8"Wx11"D ² Angle Adjustable \$85 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565 |

3c. Multi-Axis Footplate(s)¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- | | |
|---|---|
| ☐ ☐ ☐ 4"Wx6"D Multi-Axis \$247.50 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552 | ☐ ☐ ☐ 6"Wx8"D ¹ Multi-Axis \$247.50 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555 |
| ☐ ☐ ☐ 5"Wx6"D Multi-Axis \$247.50 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553 | ☐ ☐ ☐ 6"Wx11"D ¹ Multi-Axis \$247.50 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558 |
| ☐ ☐ ☐ 5"Wx8"D Multi-Axis \$247.50 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554 | ☐ ☐ ☐ 8"Wx11"D ² Multi-Axis \$247.50 ea.
☐ P ☐ L ☐ R
Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559 |

4. Swing-Away Legrest Accessories

- | | |
|---|--|
| ☐ Heel Loops (pair) \$50
Part: FRMASMB7873. HCPCS: E0951
Standard on angle adjustable and multi-axis footplates. | ☐ Gel Neoprene Footrest Sleeve (pair) \$80
Part: INDPART2243. HCPCS: K0108 |
| ☐ Wheel Bumper \$12.50 ea.
Part: See XRef | |

Legrests section continued on following page.

16. LEGRESTS CONTINUED**16.3. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ Elevating Legrest Assembly \$162.50 ea.
 P L R Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: K0195
 Select lower extension and footplates.

2. Lower Extension

☐ ☐ ☐ 11"-14" Length ☐ ☐ ☐ 12"-14.5" Length ☐ ☐ ☐ 14"-17" Length ☐ ☐ ☐ 14.5"-18" Length
 P L R Pair: FRMASMB4755 P L R Part: See XRef P L R Pair: FRMASMB4756 P L R Part: See XRef
 Left: RIG138026 Left: RIG138033 Left: RIG138033 Right: RIG138035
 Right: RIG138030 For angle adjustable only For angle adjustable only
 For reinforced only For reinforced only

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate..... Standard ☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate Standard
 P L R Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056 P L R Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

☐ ☐ ☐ 4"Wx6"D Angle Adjustable \$85 ea. ☐ ☐ ☐ 6"Wx8"D¹ Angle Adjustable with Heel Loops \$85 ea.
 P L R Pair parts: FRMASMB8812 x2. P L R Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.
☐ ☐ ☐ 5"Wx6"D Angle Adjustable \$85 ea. ☐ ☐ ☐ 6"Wx11"D¹ Angle Adjustable with Heel Loops \$85 ea.
 P L R Pair parts: FRMASMB8811 x2. P L R Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.
☐ ☐ ☐ 5"Wx8"D Angle Adjustable with Heel Loops \$85 ea. ☐ ☐ ☐ 8"Wx11"D² Angle Adjustable with Heel Loops \$85 ea.
 P L R Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792. P L R Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799.

4. Elevating Legrest Accessories

☐ Heel Loops (pair) \$50 ☐ Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$125
 Part: FRMASMB7873. HCPCS: E0951 Part: INDPART2242. HCPCS: E0953[^]
 Standard on angle adjustable footplates. ^Please verify the appropriate billing code to be used for this item
☐ Wheel Bumper \$12.50 ea. with all third-party payors.
 Part: See XRef Available on angle adjustable footplates.

16.4. Heavy Duty, Drop-In, Legrests

Available in pairs only. Select legrest uppers and footplates.

1. Heavy Duty, Drop-In Legrest Upper

☐ 70° Drop-In Legrests, 12"-14.5" Length (pair) \$420 ☐ 80° Drop-In Legrests, 12"-14.5" Length (pair) \$420
 Part: FRMASMB4762 Part: FRMASMB4762
☐ 70° Drop-In Legrests, 14.5"-20" Length (pair) \$420 ☐ 80° Drop-In Legrests, 14.5"-20" Length (pair) \$420
 Part: FRMASMB4853 Part: FRMASMB4854

2. Angle Adjustable Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ 5"Wx8"D Angle Adjustable (pair) \$170 ☐ 6"Wx11"D Angle Adjustable (pair) \$170
 Part: FRMASMB4735 Part: FRMASMB4741
☐ 6"Wx8"D Angle Adjustable (pair) \$170 ☐ 8"Wx11"D Angle Adjustable (pair) \$170
 Part: FRMASMB4738 Part: FRMASMB4744

3. Heavy Duty Drop-In Legrest Accessories

☐ Wheel Bumper \$12.50 ea.
 Part: WHLASMB1578

16.5. 90° Pediatric, Swing-Away Legrests

Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Pediatric Legrests Upper

☐ ☐ ☐ 90°, Pediatric, Swing-Away Legrest (7"-10" Length Extension) \$81 ea.
 P L R Pair part: FRMASMB8808. Left: FRMASMB8809. Right: FRMASMB8810

2. Angle Adjustable Footplate(s)

☐ ☐ ☐ 4"Wx6"D Angle Adjustable \$85 ea. ☐ ☐ ☐ 5"Wx6"D Angle Adjustable \$85 ea.
 P L R Pair part: FRMASMB8812. Left: RIG133537. Right: RIG133560 P L R Pair part: FRMASMB8811. Left: RIG133538. Right: RIG133561

3. Accessories

☐ ☐ ☐ Wheel Bumper(s) \$12.50 ea. Heel Loops (pair) \$50
 P L R Pair part: WHL169508. Left: WHL169510. Right: WHL169511 Part: FRMASMB7873. HCPCS: E0951

17. SEAT CUSHION


Optional. Specific cushion/component sizes available, part numbers, weight capacities, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

17.1. Seat Cushion Model Selection**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.

Part series: *ST-SIMGXXYY*

Simplicity G, 10"W to 20" W, E2601 \$110

**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: *ST-SPCFSPXXYY*


Spectrum Foam SPP, 16"W to 20" W, E2607 \$460

**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: *ST-SOLSPXXYY*


Solution SPP, 10"W to 20" W, E2607 \$375

**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.

Part series: *ST-GLSPXXYY*


Glacial SP, 10"W to 20"W, E2603 \$330

**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.

Part series: *ST-ZENSPXXYY*


Zen SP, 10"W to 20" W, E2603 \$315

**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Part series: *ST-TC2SPPXXYY*

TRU-Comfort 2 SPP, 16"W to 20" W, E2607 \$400

**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: *ST-SPCGSPPXXYY*

Spectrum Gel SPP, 14"W to 20" W, E2607 \$460

**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Part series: *ST-ESSSPXXYY*

Essence SPP, No Wedge, 16"-20" Wide, E2607 \$400

Essence SPP, Small Wedge, 16"-20" Wide, E2607 \$435

Essence SPP, Medium Wedge, 16"-20" Wide, E2607 . \$465

Essence SPP, Large Wedge, 16"-20" Wide, E2607 \$495

**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.

Part series: *ST-PREMPXXYY*

Premiere P, 16"W to 20" W, E2605 \$325

**17.2. Seat Cushion Dimensions**

REQUIRED with cushion model selection. Select to match or pick non-matching dimensions. See respective cushion order form for available sizes.

Match Seat Width and Depth

Cushion Width

10" W	14" W	16" W	18" W	20" W
12" W	15" W	17" W	19" W	

Cushion Depth

10" D	14" D	16" D	18" D	20" D
12" D	15" D	17" D	19" D	

17.3. Seat Cushion Substitute/Add Fluid Proof Liners

Fluid Proof Urethane Liner \$50-\$70

Part: *URL1709070-XXYY*

Not available with Simplicity, Essence, or TRU-Comfort 2.

Fluid Proof Polyurethane Liner \$55-\$75

Essence/TC2 part: *PUL1709080-XXYY*. Others part: *PUL1709030-XXYY*

Not available with Simplicity.

17.4. Seat Cushion Substitute Cover

An additional cover may be ordered or the standard cover may be substituted with a cover by checking the option in subsection 2. The "*"prefix" of the part number will be the part prefix of the cushion model selected above. For example, SPCGSPP-NBLK1816 will be an additional coolcore cover for a Spectrum Gel cushion.

Stealth Tek 2 Cover No Charge

Part: *prefix-TBLKXXYY

Not available with Simplicity, Essence, Premiere, or TRU-Comfort 2.

Mesh Cover No Charge

Part: *prefix-MBLKXXYY

Not available with Simplicity or TRU-Comfort 2.

18. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Contoured Seat Cup Holder \$20

Part: *ACCASMB2287*

Synergy Seat Cup Holder No Charge

Part: *ACCASMB2634*

Oxygen Holder for Contoured Seat \$118

Part: *ACCASMB2292*. HCPCS: E2208

Oxygen Holder for Synergy Seat \$150

Not available with Synergy Manual Recline.

Part: *FRMASMB7030*. HCPCS: E2208

Walker Holder \$103

Contoured seat part: *ACCASMB1912*. Synergy part: *FRMASMB7031*

Cane/Crutch Holder \$110

Contoured seat part: *ACCASMB2294*. HCPCS: E2207

Synergy part: *FRMASMB4822*. HCPCS: E2207

Cell Phone Holder for 12" Lift Captain Seat \$70

Part: *ACC1611300*

Rear Basket \$118

Contoured seat part: *ACCASMB1914*

Synergy part: *FRMASMB4821*

Flag Holder \$25

Contoured seat part: *ACCASMB1930*

Synergy part: *ACCASMB2415*

Accessories section continued on following page.

18. ACCESSORIES CONTINUED

Saddle Bag \$35 Part: DWR1020D010	Grip Lap Board (20005)..... \$38.95 Part: ACC152992
Stealth Cell Phone/Media Device Holder\$195 Part: ACC7120006	Grip Drink Holder (20085) \$15.50 Part: ACC152993
Stealth Large Tablet Holder \$375 Part: ACC7120005	Accessory Bar \$100 12"W-16"W part: MEC145955. 16"W-20"W part: MEC143975. 20"W-22"W part: MEC146018
Quantum Backpack \$150 Part: ACC1905853	Large, Black, Weather Cover \$202 Part: ACCCOVR1021

19. BRANCH PATCHES

Optional.

Air Force Patch \$20 Part: ACC140240	Coast Guard Patch \$20 Part: ACC140239	Navy Patch \$20 Part: ACC140244
Army Patch \$20 Part: ACC140241	Marines Patch \$20 Part: ACC140243	Purple Heart Patch \$20 Part: ACC140238

**20. POSITIONING COMPONENTS**Optional. Not available with contoured seating. More options are available on the [Positioning Components Order Form](#).**20.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, Rail Mounts, pair (TWBSTL) \$472.50 Part: INDPART2896. HCPCS: E1028
--

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5) \$210 Part: INDPART2663. HCPCS: E0956	Curved Pads, 4"x6", pair (TWBL-4X6) \$210 Part: INDPART2664. HCPCS: E0956
--	--

3. Installation

Install Laterals Part: PTOINDV3499	Do Not Install Part: PTOINDV3500
---------------------------------------	-------------------------------------

20.2. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable & Removable Hardware... \$245 ea. Part: ST-UTTP5-Q	UniLink, 7" Height Adjustable & Removable Hardware ... \$255 ea. Part: ST-UTTP7-Q
---	--

2. UniLink Pelvic/Thigh Pads

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad\$95 ea. Part: ST-UTP544056100-Q	UniLink, 4" x 12" Pad \$125 ea. Part: ST-UTP54D056100-Q
UniLink, 4" x 8" Pad..... \$105 ea. Part: ST-UTP548056100-Q	

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$75 ea. Part: ST-UIL001150A000-Q	UniLink, 2" Link..... \$75 ea. Part: ST-UIL001200A00-Q
---	---

4. UniLink Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

**20.3. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW) \$236.25 ea. Part: INDPART2804. HCPCS: E1028	Flip Down with LEFT Full Surface Contact \$236.25 ea. Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028
Removable Bracket (LP430)..... \$262.50 ea. Part: INDPART2805. HCPCS: E1028	Flip Down with RIGHT Full Surface Contact..... \$236.25 ea. Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028
Removable Bracket with 1" Offset \$400 ea. (LP430 & LHW-121) Part: INDPART2854. HCPCS: E1028	

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122).....\$26.25 Part: INDPART2821. HCPCS: K0108 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.	2 1/4" Adapter Plate (LHW-129)\$53.50 Part: POS140430. HCPCS: K0108 Only available with the flip down full surface bracket.
--	---

Positioning components section continued on following page.

20. POSITIONING COMPONENTS CONTINUED**3. Pelvic/Thigh Guide Pad**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Short Pad, 4"x4" (TWBADD-4X4) \$75 ea. Part: INDPART2801. HCPCS: E0953^/E0956	Long Pad, 4"x12" (TWBADD-4X12) \$75 ea. Part: INDPART2803. HCPCS: E0953^/E0956
Medium Pad, 4"x8" (TWBADD-4X8) \$75 ea. Part: INDPART2802. HCPCS: E0953^/E0956	

4. Pelvic/Thigh Guide InstallationInstall on Left
Part: PTOINDV3585Install on Right
Part: PTOINDV3586Do Not Install
Part: PTOINDV3719**20.4. Stealth i-Fit Chest Harness**

Upper Classic, Medium (IFPS1103) \$142 ea. Part: ST-IFPS1103-Q. HCPCS: E0960	Upper Contour, Medium (IFPS2103) \$142 ea. Part: ST-IFPS2103-Q. HCPCS: E0960
Upper Classic, Large (IFPS1104) \$142 ea. Part: ST-IFPS1104-Q. HCPCS: E0960	Upper Contour, Large (IFPS2104) \$142 ea. Part: ST-IFPS2104-Q. HCPCS: E0960

**20.5. Stealth i-Fit Chest Straps**

Chest Strap, Medium w/ Belt Anchor \$90 ea. (IFCS4113) Part: ST-IFCS4113-Q. HCPCS: E0960	Chest Strap, Large w/ Belt Anchor \$90 ea. (IFCS4114) Part: ST-IFCS4114-Q. HCPCS: E0960
--	---

20.6. Stealth Elbow Support**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket \$183.75 ea. Tube Style Arms (ARMS-ES-L) Part: POS124064. HCPCS: E1028	Right, Swing-Away, Elbow Support Bracket \$183.75 ea. Tube Style Arms (ARMS-ES-R) Part: POS124066. HCPCS: E1028
--	---

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$52.50 ea. Part: POS124065. HCPCS: K0108	4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$52.50 ea. Part: POS124067. HCPCS: K0108
--	--

20.7. Stealth Foot Boxes

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Stealth Small Foot Boxes, 5"Wx8"Dx7"H \$498.75 (SPLE107) (pair) Part: INDPART2786. HCPCS: E0954^	Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$498.75 (SPLE109) (pair) Part: INDPART2788. HCPCS: E0954^
Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$498.75 (SPLE108) (pair) Part: INDPART2787. HCPCS: E0954^	

20.8. Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D \$266 (SPLE111) (pair) Part: INDPART2783. HCPCS: E0951+E0952	Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D \$266 (SPLE114) (pair) Part: INDPART2785. HCPCS: E0951+E0952
Stealth Medium Padded Shoe Holder, 4"Wx8.25"D \$266 (SPLE112) (pair) Part: INDPART2784. HCPCS: E0951+E0952	

21. REFERENCE ORDER FORMS

Optional. Select below to note that additional order forms will accompany this order form. This will also be a link in interactive forms.

[Positioning Components Order Form](#)**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2019 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.

[youtube.com/
QuantumRehab](https://www.youtube.com/QuantumRehab)[facebook.com/
QuantumRehab](https://www.facebook.com/QuantumRehab)[twitter.com/
quantumrehab](https://twitter.com/quantumrehab)[instagram.com/
quantum_rehab](https://www.instagram.com/quantum_rehab)