

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Quantum 1450 Order Form with Contoured and Synergy Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



TRU-Balance bariatric power tilt shown.

INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering. Options noted with "XRef" have multiple possible part numbers based on system configuration. Please see [cross-reference spreadsheet](#) for the XRef part numbers.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____



A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

1. BASE MODEL

REQUIRED. Select a group 3 power base. 600 lb. weight capacity with standard 14" flat free tires, 4.5 mph motors, and 8 amp off board charger. Dimensions of base without legrests: 29.25"W x 38.25"L.

Quantum 1450 3VHD-SS\$11,095

Part: Q1450 3VHD-SS. HCPCS: K0852

Solid seat pan standard with no power option.

Quantum 1450 3VHD-C\$11,095

Part: Q1450 3VHD-C. HCPCS: K0853

Captain seat standard.

2. BASE COLOR**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.**Red Corvette**
Part: SHDASMB6684**Back in Black**
Part: SHDASMB6686**Pretty in Pink**
Part: DGN1903817**Electric Blue**
Part: SHDASMB7110018**Evolution Orange**
Part: DGN1903818**To Be Determined**
Part: N/A**3. BASE MOUNTED FOOT PLATFORM****REQUIRED.** Select a foot platform option. May interfere with seat mounted legrests.**Heavy Duty, 17"Wx13"D, Foot Platform**..... Standard
Part: FRMASMB17955**Omit Base Mounted Foot Platform** No Charge
Part: PTOINDV3674**Heavy Duty, 17"Wx16"D, Foot Platform**..... No Charge
Part: FRMASMB17956**4. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Forms](#).**4.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

VR2 6-Key Joystick Standard
Part: CTLDC1498**No Joystick, Setup for VR2 Joystick** No Charge
Part: Configuration selection**4.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers. This sections allows a selection of common configurations.

1. Q-Logic 2 Required Components

Q-Logic 2 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

Q-Logic 2 EX Controller Module \$695
Part: ELE136028. HCPCS: E2377**Q-Logic 2 Expandable Controller Harness**..... \$425
Part: DWR1111H065. HCPCS: E2313**2. Q-Logic 2 Joystick****Q-Logic 2 EX Joystick** \$700
Part: CTL123217**Stand Alone Joystick with** \$700
External Charger Port**Q-Logic 2 EX Joystick for Full Lights** \$700
Part: CTL124653

Part: ELEASMB5205

Up to 3 drive profiles and 1 seating profile when used alone. Use in conjunction with Enhanced Display for up to 5 profiles and auxiliary functions such as Bluetooth Mouse and IR. Includes charger port mounting bracket and multiplier box harness.

No Joystick, Setup for Q-Logic
Part: Configuration selection.**3. Q-Logic 2 Joystick Armor****REQUIRED with Q-Logic 2 EX Joystick selection above.** Reference colors shown in section 2.**Black, Soft Touch** Standard
Part: DGN124307**Match Joystick Armor to Base Color** \$50
Part: PTOINDV3421**Non-matching Joystick Armor** \$50
Part: PTOINDV3422

Select color below

Red Corvette
Part: DGN124336-PRD1**Back in Black**
Part: DGN124336-PBK1**Pretty in Pink**
Part: DGN124336-PPK3**Electric Blue**
Part: DGN124336-PBU9**Evolution Orange**
Part: DGN124336-POR7**To Be Determined**
Part: N/AQ-Logic 2 EX
Controller Module

Stand Alone Joystick



Q-Logic 2 EX Joystick

5. JOYSTICK HANDLES**REQUIRED.** The Q-Logic 2 EX joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.**Mushroom Handle (B)** \$85
Part: ACCASMB1860. HCPCS: E2323**Soft Ball Handle (F)** \$65
Part: ACCASMB1856. HCPCS: E2323**Chin Cup (C)** \$85
Part: ACCASMB1861.
HCPCS: E2323. (E2324 when chin control)**Extended I-Handle (G)** \$85
Part: ACCASMB1857. HCPCS: E2323**T-Bar (D)** \$120
Part: ACCASMB1858. HCPCS: E2323**Bodypoint 3" Goal Post Handle** \$99
Part: PLSKNCP1057. HCPCS: E2323**Medium Ball Handle (E)** \$65
Part: ACCASMB1859. HCPCS: E2323**Bodypoint 4" Goal Post Handle** \$99
Part: PLSKNCP1058. HCPCS: E2323

Joystick Handles


6. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket


6.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket Standard


Q-Logic Part: FRMASMB10617. VR2 Part: FRMASMB9532

Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**6.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" is used with TRU-Balance 3 arms when the joystick is inline with armpad. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

Swing-away, Inline, Left Side Mount \$245
Part: See XRef. HCPCS: E1028Swing-away, 1" Offset, Left Side Mount \$245
Part: See XRef. HCPCS: E1028Swing-away, Inline, Right Side Mount \$245
Part: See XRef. HCPCS: E1028Swing-away, 1" Offset, Right Side Mount \$245
Part: See XRef. HCPCS: E1028**6.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$350
Q-Logic Part: MEC138860. HCPCS: E1028
Not available with VR2 joystick. Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.Stealth Swing-Away/Flip Down Bracket \$294.25
Q-Logic Part: ST-TWBM480QL-Q. HCPCS: E1028
SAJ Part: ST-TWBM480CJ-Q. HCPCS: E1028
Not available with VR2 joystick. Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.Multi-axis Mount \$100
VR2 Part: ACC142531. HCPCS: E1028
Q-Logic Part: ACC131216. HCPCS: E1028
Allows height adjustment of joystick.**Specialty Mount Side**☐ Left Side Mount
Part: PTOINDV1225☐ Right Side Mount
Part: PTOINDV1226**7. ELECTRONICS ACCESSORIES****Optional.**Q-Logic EX/NE Programmer \$545
Part: ELEASMB6976Q-Logic PC Programming Disk with Cable \$750
Part: ELEASMB5216**8. CONTOURED SEATING****Optional.** Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.**8.1. Seat to Floor Height**Lowest Seat to Floor
Part: PTOINDV3558Middle Seat to Floor
Part: PTOINDV3559Highest Seat to Floor
Part: PTOINDV3560**8.2a. Captain Seating****1a. Depth Adjustable, Fold Flat, High Back Seat with Limited Recline**

Full length armrests must be selected with a depth adjustable high back.

20"Wx18"-20"D with 60" Lap Belt No Charge
Part: SETASMB6896
Max weight capacity is 450 lbs. when seat is set at 18"D.22"Wx20"-22"D with 70" Lap Belt No Charge
Part: SETASMB6897**1b. High Back Seat with Limited Recline**20"Wx22"D Standard
Part: SETASMB4889**1c. Pinchless Hinge, High Back Seat**20"Wx24"D Standard
Part: SETASMB238022"Wx24"D Standard
Part: SETASMB238424"Wx24"D Standard
Part: SETASMB238922"Wx20"D Standard
Part: SETASMB238224"Wx20"D Standard
Part: SETASMB238722"Wx22"D Standard
Part: SETASMB238324"Wx22"D Standard
Part: SETASMB2388

Contoured seating section continued on following page.

8. CONTOURED SEATING CONTINUED**8.2a. Captain Seating continued****1c. Heavy Duty, High Back Seat**

Seats include a lumbar cutout. Select armrests and armpads in Armrest for Heavy Duty Contoured Seating section.

26"Wx20"D.....\$1,000 Part: SETASMB4566	28"Wx24"D\$1,000 Part: SETASMB4574	30"Wx28"D\$1,000 Part: SETASMB4582
26"Wx22"D.....\$1,000 Part: SETASMB4567	28"Wx26"D.....\$1,000 Part: SETASMB4575	32"Wx20"D\$1,000 Part: SETASMB4584
26"Wx24"D.....\$1,000 Part: SETASMB4568	28"Wx28"D\$1,000 Part: SETASMB4576	32"Wx22"D.....\$1,000 Part: SETASMB4585
26"Wx26"D.....\$1,000 Part: SETASMB4569	30"Wx20"D\$1,000 Part: SETASMB4578	32"Wx24"D.....\$1,000 Part: SETASMB4586
26"Wx28"D.....\$1,000 Part: SETASMB4570	30"Wx22"D\$1,000 Part: SETASMB4579	32"Wx26"D.....\$1,000 Part: SETASMB4587
28"Wx20"D\$1,000 Part: SETASMB4572	30"Wx24"D\$1,000 Part: SETASMB4580	32"Wx28"D.....\$1,000 Part: SETASMB4588
28"Wx22"D\$1,000 Part: SETASMB4573	30"Wx26"D.....\$1,000 Part: SETASMB4581	

8.2b. Contoured Seat with Solid Seat Pan**1c. Limited Recline, High Back Seat with Solid Seat Pan**

20"Wx20"D Standard Part: SETASMB3185	20"Wx22"D Standard Part: SETASMB5008
---	---

1c. Pinchless Hinge, High Back Seat with Solid Seat Pan

20"Wx24"D Standard Part: SETASMB2982	22"Wx24"D Standard Part: SETASMB2987	24"Wx24"D Standard Part: SETASMB2992
22"Wx20"D Standard Part: SETASMB2985	24"Wx20"D Standard Part: SETASMB2990	
22"Wx22"D Standard Part: SETASMB2986	24"Wx22"D Standard Part: SETASMB2991	

1c. Heavy Duty, High Back Seat with Solid Seat Pan

Seats include a lumbar cutout. Select armrests and armpads in Armrest for Heavy Duty Contoured Seating section.

26"Wx20"D.....\$1,000 Part: SETASMB5479	28"Wx24"D\$1,000 Part: SETASMB5487	30"Wx28"D\$1,000 Part: SETASMB5495
26"Wx22"D.....\$1,000 Part: SETASMB5480	28"Wx26"D.....\$1,000 Part: SETASMB5488	32"Wx20"D\$1,000 Part: SETASMB5497
26"Wx24"D.....\$1,000 Part: SETASMB5481	28"Wx28"D\$1,000 Part: SETASMB5489	32"Wx22"D.....\$1,000 Part: SETASMB5498
26"Wx26"D.....\$1,000 Part: SETASMB5482	30"Wx20"D\$1,000 Part: SETASMB5491	32"Wx24"D.....\$1,000 Part: SETASMB5499
26"Wx28"D.....\$1,000 Part: SETASMB5483	30"Wx22"D\$1,000 Part: SETASMB5492	32"Wx26"D.....\$1,000 Part: SETASMB5500
28"Wx20"D\$1,000 Part: SETASMB5485	30"Wx24"D\$1,000 Part: SETASMB5493	32"Wx28"D.....\$1,000 Part: SETASMB5501
28"Wx22"D\$1,000 Part: SETASMB5486	30"Wx26"D.....\$1,000 Part: SETASMB5494	

8.2c. Special Order Contoured Seat

Complete information below. Contact inside sales specialist for pricing.

Type of Seat: _____ Width: _____ Depth: _____

8.3a. Armrests for Contoured Seating

Not available with heavy duty contoured seating. Includes arm pads.

L R Full Length, Flip-up, Height Adjustable, Removable \$85 ea. Left part: FRMASMB8614. Right part: FRMASMB8615 HCPCS: E0973	L R Desk Length, Flip-up, Height Adjustable, Removable \$85 ea. Left part: FRMASMB8612. Right part: FRMASMB8613 HCPCS: E0973 Not available with depth-adjustable seats.
---	--

Contoured seating section continued on following page.

8. CONTOURED SEATING CONTINUED**8.3b. Armrest for Heavy Duty Contoured Seating****1a. Bariatric, Removable, Height Adjustable Armrests**

Only available with bariatric armpad.

L R 8"-12", Full Length, Bariatric, Height Adjustable, No Charge
Removable
Left part: FRMASMB6728. Right part: FRMASMB6729
HCPCS: E0973

L R 10"-14" Full Length, Bariatric, Height Adjustable..... No Charge
Removable
Left part: FRMASMB6730. Right part: FRMASMB6731
HCPCS: E0973

1b. Single Post, Removable, Height Adjustable, Heavy Duty Armrests

L R 8"-12", Full Length, Heavy Duty, Single Post, \$85 ea.
Height Adjustable and Removable
Left part: FRMASMB9982. Right part: FRMASMB9983
HCPCS: E0973

L R 10"-14", Full Length, Heavy Duty, Single Post, \$85 ea.
Height Adjustable and Removable
Left part: FRMASMB9986. Right part: FRMASMB9987
HCPCS: E0973

L R 8"-12", Desk Length, Heavy Duty, Single Post, \$85 ea.
Height Adjustable and Removable
Left part: FRMASMB9984. Right part: FRMASMB9985
HCPCS: E0973

L R 10"-14", Desk Length, Heavy Duty, Single Post, \$85 ea.
Height Adjustable and Removable
Left part: FRMASMB9988. Right part: FRMASMB9989
HCPCS: E0973

2. Armrests for Single Post, Heavy Duty Armrests

Select a style of armrest for heavy duty, contoured seat armrests. Armrests are configured as pairs.

Bariatric Armrests No Charge
Part: UPSASMB1025
Only available with bariatric armrests.

Straight Armrests No Charge
Desk/Desk part: UPSASMB1022. Full/Full part: UPSASMB1021.
Desk/Full part: UPSASMB1023
Not available with bariatric armrests.

Waterfall Armrests No Charge
Desk/Desk part: UPSASMB1025. Full/Full part: UPSASMB1024.
Desk/Full part: UPSASMB1026
May interfere with seat back. Not available with bariatric armrests.

Omit Armrests No Charge
Part: PTOINDV1373

9. SYNERGY SEATING SYSTEM**Optional.** Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 32"W x 28"D.**9.1. Synergy Frame Selection**

Synergy Seating System \$800
Part: PTOINDV2995

9.2. Seat Dimension - Width**REQUIRED with Synergy seating.** Maximum seat width is 32". Part number is SETWDTH plus number following # below.

20" W (#1010)	22" W (#1012)	24" W (#1014)	26" W (#1054)	28" W (#1055)	30" W (#1059)	32" W (#1061)
21" W (#1011)	23" W (#1013)	25" W (#1056)	27" W (#1057)	29" W (#1058)	31" W (#1060)	

9.3. Seat Dimension - Width Package**REQUIRED with Synergy seating.**

20"-24" Width Package Part: FRMASMB7120014	26"-28" Width Package Part: FRMASMB7120019	30"-32" Width Package Part: FRMASMB7120021
24"-26" Width Package Part: FRMASMB7120015	28"-30" Width Package Part: FRMASMB7120020	

9.4. Seat Dimension - Depth**REQUIRED with Synergy seating.** Maximum seat depth is 28". Part number is SETDPTH plus number following # below.

20" D (#1010)	22" D (#1012)	24" D (#1014)	26" D (#1061)	28" D (#1063)
21" D (#1011)	23" D (#1013)	25" D (#1060)	27" D (#1062)	

9.5. Seat Dimension - Depth Package**REQUIRED with Synergy seating.**

15"-20" Depth Package Part: FRMASMB7732	20"-24" Depth Package Part: FRMASMB7733	24"-28" Depth Package Part: FRMASMB6605
--	--	--

9.6. Synergy Back Options**REQUIRED with Synergy seating.** See the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1. Back Selection**

Synergy Back \$275
Part: See XRef
Select a back cane type and height and back height.

TRU-Comfort 4-way Stretch, 20" Width or less \$675
Part: See XRef. HCPCS: E2620
Select back height.

TRU-Comfort 4-way Stretch, 22" Width or greater \$775
Part: See XRef. HCPCS: E2621
Select back height.

HD Sling Back \$98
Part: See XRef
Select a back height. Available with 24"-28" Width and 18"-25" overall heights. Back cane height must match back height. Not available with flip up cantilever arms.

Omit Back No Charge
Part: PTOINDV1287
Select a back cane type and height in following section.

Synergy seating continued on following page.

9. SYNERGY SEATING SYSTEM CONTINUED

		Back Width																
		16"	17"	18"	19"	20"	21"	22"	23"	24"	25"	26"	27"	28"	29"	30"	31"	32"
Overall Height	17"	S	S	S	S	S	-	S	-	S	-	-	-	-	-	-	-	-
	18"	S	S	S	S	S	-	S	-	S	-	-	-	-	-	-	-	-
	19"	S	-	S	S	S	S	S	S	S	-	S	S	S	S	-	S	-
	20"	STC	-	STC	STC	STC	S	S	S	S	-	S	S	S	S	-	S	-
	21"	STC	S	STC	STC	STC	S	S	S	S	S	-	S	S	-	S	-	-
	22"	STC	STC	STC	STC	STC	S	S	STC	S	S	S	-	S	S	-	S	-
	23"	STC	STC	STC	STC	STC	S	S	STC	S	S	S	S	-	S	-	S	S
	24"	STC	STC	STC	STC	STC	S	STC	STC	STC	S	S	S	S	-	-	S	S
	25"	STC	STC	STC	STC	STC	S	STC	STC	STC	S	S	S	S	-	-	S	S
	27"	TC	-	TC	TC	TC	-	TC	-	TC	-	-	-	-	-	-	-	-
28"	TC	-	TC	TC	TC	-	TC	-	TC	-	-	-	-	-	-	-	-	

S

 Synergy Back

TC

 TRU-Comfort Back

STC

 Synergy/TRU-Comfort Back

Synergy and TRU-Comfort Back Availability Chart

2. Back Cane Selection

Synergy Back Cane No Charge
 Part: FRMASMB15844
 Select a back cane height.

Flip Forward Back Cane \$150
 Part: FRMASMB15846
 Select a back cane height. Not available with oxygen holder, or cane mounted and cantilever armrests.

2" Mesomorph (1" per side) Back Cane \$195
 Part: FRMASMB12067
 Select a cane height.

4" Mesomorph (2" per side) Back Cane \$195
 Part: FRMASMB12091
 Select a cane height.

2" Endomorph (1" per side) Back Cane \$195
 Part: FRMASMB12115
 Select a cane height.

4" Endomorph (2" per side) Back Cane \$195
 Part: FRMASMB12139
 Select a cane height.

Manual Recline Back \$750
 Part: FRMASMB12910. HCPCS: E1226
 Select an overall back height. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, rear accessories, or articulating foot platform.

3. Synergy Push Handles

Push Handles Standard
 Part: Configuration selection with back cane height

Omit Push Handles \$125
 Part: Configuration selection with back cane height
 Not available with Sling Back. Minimum cane height is 16".

4. Synergy Accessory Bar

Accessory Bar \$100
 12"W-16"W Part: MEC145955. 16"W-20"W Part: MEC143975.

5. Synergy Overall Back Height

Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Parts: See XRef, configuration selection with back selection.

17" High Back	19" High Back	21" High Back	23" High Back	25" High Back
18" High Back	20" High Back	22" High Back	24" High Back	

6. Synergy Back Cane Height

Omit this selection if a TRU-Comfort back is selected. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles. Parts: See XRef, configuration selection with push handles selection.

14" High Canes	17" High Canes	20" High Canes	23" High Canes
15" High Canes	18" High Canes	21" High Canes	24" High Canes
16" High Canes	19" High Canes	22" High Canes	25" High Canes

9.7. Synergy Armrests

REQUIRED with Synergy seating. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Bariatric, Single Post, Heavy Duty, Height Adjustable, Removable Armrests

May interfere with swing-away laterals.

☐ ☐ Full Length, Bariatric Armrests, 8"-12" H No Charge
L R Left part: FRMASMB6728. Right part: FRMASMB6729

☐ ☐ Full Length, Bariatric Armrests, 10"-14"H No Charge
L R Left part: FRMASMB6730. Right part: FRMASMB6731.

1b. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H \$75 ea.
L R Left part: FRMASMB9982. Right part: FRMASMB9983. HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H \$75 ea.
L R Left part: FRMASMB9984. Right part: FRMASMB9985. HCPCS: E0973

1c. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14" H \$75 ea.
L R Left part: FRMASMB9986. Right part: FRMASMB9987. HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14"H \$75 ea.
L R Left part: FRMASMB9988. Right part: FRMASMB9989. HCPCS: E0973

Synergy seating continued on following page.

9. SYNERGY SEATING SYSTEM CONTINUED**9.7. Synergy Armrests continued****1d. Height Adjustable, Cantilever Armrests**

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

Full Length, Cantilever Armrests \$97.50 ea.
 Left Part: FRMASMB16188. Right part: FRMASMB16189

Desk Length, Cantilever Armrests \$97.50 ea.
 Left part: FRMASMB16190. Right part: FRMASMB16191

9.8. Synergy Armpads**REQUIRED with armrest selection above.** Armpads are configured as pairs.

Bariatric Armpads No Charge
 Part: UPSASMB1025.
 Only available with Bariatric Armrests.

Straight Armpads No Charge
 Desk/Desk part: UPSASMB1022. Full/Full part: UPSASMB1021.
 Desk/Full part: UPSASMB1023
 Not available with bariatric armrests.

Waterfall Armpads No Charge
 Desk/Desk part: UPSASMB1025. Full/Full part: UPSASMB1024.
 Desk/Full part: UPSASMB1026
 May cause interference with seat back. Not available with bariatric armrests.

Omit Synergy Armpads No Charge
 Part: PTOINDV1373

9.9. Synergy Seat to Floor**REQUIRED with Synergy seating.**

16.50" STF
 Part: PTOINDV1888

17.50" STF
 Part: PTOINDV1889

18.50" STF
 Part: PTOINDV1890

10. JOYSTICK BRACKET RECEIVER**REQUIRED.**

Standard Fixed Receiver Standard
 Part: ACC123525

Stealth Height Adjustable Receiver \$214
 Part: ST-ARM250C-Q. HCPCS: K0108

**11. HEADRESTS****Optional.** Select a headrest pad and bracket. **Not available with contoured seating.****11.1. Headrest Pads**

Stealth 8" Comfort Plus (CP450) \$151.94
 Part: POS148484. HCPCS: E0955

Stealth 10" Comfort Plus (CP250)(A) \$151.94
 Part: SETHEAD1057. HCPCS: E0955


Stealth 14" Comfort Plus (CP150) \$151.94
 Part: INDPART2327. HCPCS: E0955

Pediatric Comfort Plus (6"Wx4.5"H) \$151.94
 (CPN350)
 Part: SETHEAD1069. HCPCS: E0955

Stealth Small Contoured (9.75"Wx3"H) \$64.20
 (STL-OTSM250)

Part: SETHEAD1052. HCPCS: E0955

Stealth Large Contoured (11"Wx6"H) (B) \$80.25
 (STL-OTLG150)
 Part: SETHEAD1053. HCPCS: E0955


**11.2. Stealth Headrest Brackets**

TWB Removable Headrest Bracket (TWB485-1) (C) \$220.50
 Part: FRMASMB15999. HCPCS: E1028

Multi-Axis Removable Headrest Bracket (HMO475P) (D) \$160.50
 Part: HDWASMB1084. HCPCS: E1028

12. FULL LIGHT PACKAGE**Optional.** Requires Q-Logic 2 Joystick for Full Light Package.

Seat Mounted, LED, Full Light Package \$995
 Part: ELEASMB5719

**13. LAP BELTS****REQUIRED.** Select a lap belt below.**13.1a. Standard Lap Belts**

38" Lap Belt
 Part: ACCBELT1006

50" Lap Belt
 Part: ACCBELT1000

60" Lap Belt
 Part: ACCBELT1002

70" Lap Belt
 Part: ACCBELT1003

80" Lap Belt
 Part: ACCBELT1004

90" Lap Belt
 Part: ACCBELT1001

13.1b. Optional Lap Belts

i-Fit 60" Single Pull, Push Button, Medium (7.5") Padded, 1.5" Belt (IFPB5213) \$86
 Part: ST-IFPB5213-Q

i-Fit 60" Single Pull, Push Button, Large Padded (9.5"), 2" Belt (IFPB5314) \$86
 Part: ST-IFPB5314-Q

i-Fit 60" Single Pull, Push Button, Extra Large Padded (11.5"), 2" Belt (IFPB5315) \$86
 Part: ST-IFPB5315-Q

14. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

14.1. Manually Adjustable Center Mount Foot Platform

Not available with contoured seating.

1. Foot Platform Assembly

Center Mount Foot Platform..... \$250
20"-24"W part: RIG155637. 24"-28"W part: RIG155527. 28"-33"W part: RIG155528

2. Footplates

17"Wx13"D One Piece Foot Platform
Part: RIG123289

17"Wx16"D One Piece Foot Platform
Part: RIGASMB7120011

3. Center Mount Calf Support

Calf pad and panel covers include 1/2" of soft visco foam.

Large Calf Pads (pair) 8" taper to 6"Wx8"H \$150
Part: RIGASMB7120019

Large Calf Panel 24"Wx8"H \$200
Part: RIG123824

4. Additional Accessories

Reversed Dartex Calf Pad Cover (pair) \$150
Small part: RIG133717. Large part: RIG133719

Wheel Bumpers (pair)..... \$12.50 ea.
Part: WHLASMB1578 x2

Reversed Dartex Panel Cover (pair) \$150
Small part: RIG134673. Large part: RIG134744

Footplate Connector \$15
Part: RIGASMB7120023

5. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers..... \$205
Part: ST-FPG5008T-Q

6.5"x8" Tapered Gel Footplate Covers for Width Extension..... \$225
Part: ST-FPG6508T-Q

5"x11" Tapered Gel Footplate Covers..... \$225
Part: ST-FPG5011T-Q

7"x11" Tapered Gel Footplate Cover for Width Extension \$250
Part: ST-FPG7011T-Q

14.2. Heavy Duty, Drop-In, Swing-Away Legrests

Select a configuration, a legrest with extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Drop-In Legrest Selection

☐ ☐ ☐ 70° HD Drop-in with 12"-14.5" Length No Charge
Lower Extension
Pair part: FRMASMB4762. Left: RIG137706. Right: RIG137707
Only available with heavy duty high back and Synergy seat.

☐ ☐ ☐ 80° HD Drop-in with 12"-14.5" Length No Charge
Lower Extension
Pair part: FRMASMB4763. Left: RIG137712. Right: RIG137713
Only available with Synergy seat.

☐ ☐ ☐ 70° HD Drop-in with 14.5-18" Length No Charge
Lower Extension
Pair part: FRMASMB4853. Left: RIG137708. Right: RIG137709
Only available with heavy duty high back and Synergy seat.

☐ ☐ ☐ 80° HD Drop-in with 14.5-18" Length No Charge
Lower Extension
Pair part: FRMASMB4854. Left: RIG137710. Right: RIG137711
Only available with Synergy seat.

2. Angle Adjustable Footplate(s) with Heel Loops

☐ ☐ ☐ 6"Wx8"D \$85 ea.
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.
HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D \$85 ea.
Pair part: FRMASMB4744. Left: RIG137708. Right: RIG137709.
HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx11"D \$85 ea.
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.
HCPCS: K0040+E0951

3. Accessories

☐ ☐ ☐ Wheel Bumper(s) \$12.50 ea.
Pair part: WHLASMB1578. Left part: WHL169510. Right part: WHL169511

14.3. Heavy Duty, Elevating Legrests

Not available with contoured seating. Select a configuration, a legrest with extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Heavy Duty, Elevating Legrest Selection

☐ ☐ ☐ Heavy Duty, Elevating Legrest(s) with Calf Pad(s)..... \$162.50 ea.
Pair part: FRMASMB6702. Left: FRMASMB8444. Right: FRMASMB8445. HCPCS: E0990

2. Lower Extension

☐ ☐ ☐ 11"-14" Length
Part: FRMASMB4755. Left: RIG138026. Right: RIG138030
For reinforced footplates only.

☐ ☐ ☐ 14"-17" Length
Part: FRMASMB4756. Left: RIG138033. Right: RIG138035
For reinforced footplates only.

☐ ☐ ☐ 12"-14.5" Length
Part: FRMASMB4752. Left: FRMASMB5391. Right: FRMASMB4723
For angle adjustable footplates only.

☐ ☐ ☐ 14.5"-18" Length
Part: FRMASMB4753. Left: FRMASMB4724. Right: FRMASMB5392
For angle adjustable footplates only.

Legrests section continued on following page.

14. LEGRESTS CONTINUED**3. Reinforced Footplates(s)**

7.25"Wx6"D Extruded Footplate..... Standard
 P L R Part pair: FRMASMB4748. Left: RIG138055.
 Right: RIG138056

8.5"Wx6"D Extruded Footplate Standard
 P L R Part pair: FRMASMB4747. Left: RIG138060.
 Right: RIG138062

4. Angle Adjustable Footplates(s) with Heel Loops

6"Wx8"D \$85 ea.
 P L R Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.
 HCPCS: K0040+E0951

8"Wx11"D \$85 ea.
 P L R Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799.
 HCPCS: K0040+E0951

6"Wx11"D \$85 ea.
 P L R Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.
 HCPCS: K0040+E0951

May cause interference with seat widths less than 20"

15. SET LEGREST LENGTH

REQUIRED with matching left and right legrest configuration. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	


16. SEAT CUSHION

Optional. Specific cushion/component sizes available, part numbers, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

16.1. Seat Cushion Model Selection**Stealth Simplicity G Cushion**


General use cushion with Stealth Tek 2, 2-way stretch, breathable cover. **Weight capacity for 16" to 20" W is 350 lbs, over 20"W is 400 lbs.**
 Part series: ST-SIMGXXYY

Simplicity G, E2601, 20"W to 21"W \$110
 Simplicity G, E2602, 22"W to 24"W \$165

**Stealth Spectrum Foam SPP Cushion**


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner. **Weight capacity for 20" W is 350 lbs, over 20"W is 650 lbs.**
 Part series: ST-SPCFSPXXYY

Spectrum Foam SPP, E2607, 20"W \$460
 Spectrum Foam SPP, E2608, 22"W to 32"W \$530

**Stealth Solution SPP Cushion**


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner. **Weight capacity for 20" W is 350 lbs, over 20"W is 550 lbs.**
 Part series: ST-SOLSPXXYY

Solution SPP, E2607, 20"W \$375
 Solution SPP, E2608, 22"W to 24"W \$430

**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover. **Weight capacity is for 20"W is 350 lbs, over 20"W is 450 lbs.**
 Part series: ST-GLSPXXYY

Glacial SP, E2603, 20"W to 21"W \$330
 Glacial SP, E2604, 22"W to 24"W \$440

**Stealth Zen SP Cushion**


Skin protection with Coolcore, 4-way stretch, breathable cover. **Weight capacity for 20"W-21"W is 350 lbs, 22"W and over is 400 lbs.**
 Part series: ST-ZENSPXXYY

Zen SP, E2603, 20"W to 21" W \$315
 Zen SP, E2604, 22"W to 24" W \$425

**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. **Weight capacity for 16" to 20" W is 350 lbs, over 20"W is 450 lbs.**
 Part series: ST-TC2SPPXXYY

TRU-Comfort 2 SPP, E2607, 20"W \$400
 TRU-Comfort 2 SPP, E2608, 22"W \$500

**Stealth Spectrum Gel SPP Cushion**


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner. **Weight capacity for 20" W is 350 lbs, over 20"W is 550 lbs.**
 Part series: ST-SPCGSPPXXYY

Spectrum Gel SPP, E2607, 20"W \$460
 Spectrum Gel SPP, E2608, 22"W to 24"W \$530

**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. **Weight capacity for 20" W is 350 lbs, over 20"W is 550 lbs.**
 Part series: ST-ESSSPXXYY

Essence SPP, No Wedge, 20" Wide, E2607 \$400
 Essence SPP, No Wedge, 22" Wide, E2608 \$500
 Essence SPP, Small Wedge, 20" Wide, E2607 \$435
 Essence SPP, Small Wedge, 22" Wide, E2608 \$535
 Essence SPP, Medium Wedge, 20" Wide, E2607 \$465
 Essence SPP, Medium Wedge, 22" Wide, E2608 \$565
 Essence SPP, Large Wedge, 20" Wide, E2607 \$495
 Essence SPP, Large Wedge, 22" Wide, E2608 \$595

**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover. **Weight capacity for 20"W is 350 lbs, 22"W is 500 lbs.**
 Part series: ST-PREMPXXYY

Premiere P, 20"W, E2605 \$325
 Premiere P, 22"W, E2606 \$475



Seat cushion section continued on following page.

16. SEAT CUSHION CONTINUED**16.2. Seat Cushion Dimensions****REQUIRED with cushion model selection.** Select to match or pick non-matching dimensions. See respective cushion order form for available sizes.

Match Seat Width and Depth

Cushion Width

20" W 22" W 26" W 30" W
 21" W 24" W 28" W 32" W

Cushion Depth

20" D 21" D 22" D 23" D 24" D 26" D 28" D

16.3. Seat Cushion Substitute/Add Fluid Proof Liners

Fluid Proof Urethane Liner \$50-\$70
 Part: URL1709070-XXYY
 Not available with Simplicity, Essence, or TRU-Comfort 2.

Fluid Proof Polyurethane Liner \$55-\$75
 Essence/TC2 part: PUL1709080-XXYY. Others part: PUL1709030-XXYY
 Not available with Simplicity.

16.4. Seat Cushion Substitute Cover

An additional cover may be ordered or the standard cover may be substituted with a cover by checking the option in subsection 2. The "prefix" of the part number will be the part prefix of the cushion model selected above. For example, SPCGSPP-NBLK1816 will be an additional coolcore cover for a Spectrum Gel cushion.

Stealth Tek 2 Cover No Charge
 Part: *prefix-TBLKXXYY
 Not available with Simplicity, TRU-Comfort 2, Essence, and Premiere.

Mesh Cover No Charge
 Part: *prefix-MBLKXXYY
 Not available with Simplicity or TRU-Comfort 2.

17. PROVIDER INSTALLED ACCESSORIES**Optional.**

Medium, Black, Weather Cover \$168
 Part: ACCCOVR1022

Large, Black, Weather Cover \$202
 Part: ACCCOVR1021

Oxygen Holder \$118
 Contoured seat part: ACCASMB2292. Synergy part: FRMASMB7030.
 HCPCS: E2208
 Not available with Synergy Manual Recline.

Contoured Seat Cup Holder \$20
 Part: ACCASMB2287

Synergy Seat Cup Holder No Charge
 Part: ACCASMB2634

Rear Basket \$118
 Contoured seat part: ACCASMB1914
 Synergy part: FRMASMB4821

Walker Holder \$103
 Contoured seat part: ACCASMB1912. Synergy part: FRMASMB7031.

Saddle Bag \$35
 Part: DWR1020D010

Flag Holder \$25
 Contoured seat part: ACCASMB1930
 Synergy part: ACCASMB2415

Cane/Crutch Holder \$110
 Contoured seat part: ACCASMB2294. HCPCS: E2207
 Synergy part: FRMASMB4822. HCPCS: E2207

Stealth Large Tablet Holder \$375
 Part: ACC7120005

Stealth Cell Phone/Media Device Holder \$195
 Part: ACCINDV1102

18. BATTERIES**REQUIRED.** The powerchair requires 2 batteries and will be installed if ordered.

Group 24, Gel Cell, Interceptor Battery \$365 ea.
 Part: BAT1704614. HCPCS: E2363

Group 24, Gel Cell, MK Battery \$495 ea.
 Part: BATGEL1009. HCPCS: E2363

No batteries. Provider will assume responsibility of batteries.
 Part: Configuration selection

19. BRANCH PATCHES**Optional.**

Air Force Patch \$20
 Part: ACC140240

Coast Guard Patch \$20
 Part: ACC140239

Navy Patch \$20
 Part: ACC140244

Army Patch \$20
 Part: ACC140241

Marines Patch \$20
 Part: ACC140243

Purple Heart Patch \$20
 Part: ACC140238



Military Patches

20. MANUFACTURER INSTALLED POSITIONING COMPONENTS**Optional. Only available with the Synergy Seating system.** More options are available on the [Positioning Components Order Form](#).**20.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, Rail Mounts, pair (TWBTLTC) \$472.50
 Part: INDPART2677. HCPCS: E1028

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5) \$210
 Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6) \$210
 Part: INDPART2664. HCPCS: E0956

3. Installation

Install Laterals
 Part: PTOINDV3499

Do Not Install
 Part: PTOINDV3500

Manufacturer installed positioning components section continued on following page.

20. MANUFACTURER INSTALLED POSITIONING COMPONENTS CONTINUED**20.2. Stealth UniLink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable & Removable Hardware... \$245 ea.
 QTY Part: ST-UTTP5-Q

UniLink, 7" Height Adjustable & Removable Hardware... \$255 ea.
 QTY Part: ST-UTTP7-Q

2. UniLink Pelvic/Thigh Pads

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad\$95 ea.
 QTY Part: ST-UTP544056100-Q

UniLink, 4" x 12" Pad \$125 ea.
 QTY Part: ST-UTP54D056100-Q

UniLink, 4" x 8" Pad..... \$105 ea.
 QTY Part: ST-UTP548056100-Q

**3. UniLink Pelvic/Thigh Additional Links**

UniLink, 1.5" Link \$75 ea.
 QTY Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$75 ea.
 QTY Part: ST-UIL001200A00-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

20.3. Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW) \$236.25 ea.
 QTY Part: INDPART2804. HCPCS: E1028

Removable Bracket (LP430)..... \$262.50 ea.
 QTY Part: INDPART2805. HCPCS: E1028

Removable Bracket with 1" Offset \$400 ea.
 QTY (LP430 & LHW-121)
 Part: INDPART2854. HCPCS: E1028

Flip Down with LEFT Full Surface Contact \$236.25 ea.
 Feature (TWBADD-L & LPHW407)
 Part: POS140428. HCPCS: E1028

Flip Down with RIGHT Full Surface Contact..... \$236.25 ea.
 Feature (TWBADD-R & LPHW407)
 Part: POS140429. HCPCS: E1028

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122).....\$26.25
 Part: INDPART2821. HCPCS: K0108
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129)\$53.50
 Part: POS140430. HCPCS: K0108
 Only available with the flip down full surface bracket.

3. Pelvic/Thigh Guide Pad

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Short Pad, 4"x4" (TWBADD-4X4)..... \$75 ea.
 QTY Part: INDPART2801. HCPCS: E0953^/E0956

Long Pad, 4"x12" (TWBADD-4X12)..... \$75 ea.
 QTY Part: INDPART2803. HCPCS: E0953^/E0956

Medium Pad, 4"x8" (TWBADD-4X8)..... \$75 ea.
 QTY Part: INDPART2802. HCPCS: E0953^/E0956

4. Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

21. PROVIDER INSTALLED POSITIONING COMPONENTS

Optional. Only available with the Synergy Seating system. More options are available on the [Positioning Components Order Form](#).

21.1. Stealth Foot Boxes

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$498.75
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954^

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$498.75
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954^

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$498.75
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954^

21.2. Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D \$266
 (SPLE111) (pair)
 Part: INDPART2783. HCPCS: E0951+E0952

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D \$266
 (SPLE114) (pair)
 Part: INDPART2785. HCPCS: E0951+E0952

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D \$266
 (SPLE112) (pair)
 Part: INDPART2784. HCPCS: E0951+E0952

Provider installed positioning components section continued on following page.

21. PROVIDER INSTALLED POSITIONING COMPONENTS CONTINUED**21.3. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket\$183.75 ea.
 QTY Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket\$183.75 ea.
 QTY Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

Left, Fixed Elbow Support Bracket.....\$183.75 ea.
 QTY TRU-Balance 3 Arms (ARMS-ESFTB3-L)
 Part: POS133465

Right, Fixed Elbow Support Bracket\$183.75 ea.
 QTY TRU-Balance 3 Arms (ARMS-ESFTB3-R)
 Part: POS133466

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$52.50 ea.
 QTY Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)\$52.50 ea.
 QTY Part: POS124067. HCPCS: K0108

21.4. Stealth i-Fit Chest Harness

Upper Classic, Medium (IFPS1103) \$142 ea.
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$142 ea.
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$142 ea.
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$142 ea.
 Part: ST-IFPS2104-Q. HCPCS: E0960

21.5. Stealth i-Fit Chest Straps

Chest Strap, Medium w/ Belt Anchor \$90 ea.
 (IFCS4113)
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor \$90 ea.
 (IFCS4114)
 Part: ST-IFCS4114-Q. HCPCS: E0960

21.6. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. The 1" ball kit must be selected to mount the facial lateral to the headrest bracket. Also select a facial pad(s) with the facial lateral(s).

1. Facial Lateral Hardware

Facial Lateral for Comfort Series.....\$294.25 ea.
 Right Side (CPS9-R-HW)
 Part: INDPART2647. HCPCS: E1028

Facial Laterals for Comfort Series \$347.75 ea.
 Left and Right Side (CPS9-B-HW)
 Part: INDPART2649. HCPCS: E1028

Facial Lateral for Comfort Series.....\$294.25 ea.
 Left Side (CPS9-L-HW)
 Part: INDPART2648. HCPCS: E1028

2. Facial Lateral Pads and Accessories

1" Ball Kit for Facial Laterals (TWB-1BK) (D) \$80.25
 Part: POS150532
 Required with facial laterals to mount to comfort headrest brackets.

Facial Pad, 4"x2" (900-Q)..... \$53.50 ea.
 QTY Part: ST-900-Q. HCPCS: K0108

Facial Pad, 6"x2" (920-Q)..... \$53.50 ea.
 QTY Part: ST-920-Q. HCPCS: K0108

Egg Switch Mount (ESM) (A) \$52.50 ea.
 QTY Part: INDPART2815

Swing-Away Joystick Mount,.....\$250 ea.
 QTY 18" Rod (C)

Figure 8 Bracket with 13" Rod..... \$267.50 ea.
 QTY and Egg Switch Mount (B)
 Part: ST-F8675-13-Q

Allows mounting of facial pads and attaches to the headrest




Contoured and Classic Chest Harness



Chest Strap



Left and Right Facial Lateral with Pads



Headrest Accessories

22. SPECIALTY CONTROLS

Optional. Only available with Synergy Seat. Select a type of specialty controls below or to refer to the specialty controls order form.

22.1. SCIM Module

Optional.

SCIM Module Kit for Q-Logic 2.....\$525
 Part: CTL168055

Available with Q-Logic 2 EX Joystick selection above. Allows use of joystick display for specialty devices and offers up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

22.2. Q-Logic 2 Stand Alone Joystick

Q-Logic 2 EX Stand Alone Joystick.....\$700
 Part: ELEASMB5205

Up to 3 drive profiles and 1 seating profile when used alone. Use in conjunction with Enhanced Display for up to 5 profiles and auxiliary functions such as Bluetooth Mouse and IR. Includes charger port mounting bracket and multiplier box harness.



Stand Alone Joystick

Specialty controls section continued on following page.

22. SPECIALTY CONTROLS CONTINUED**22.3. Enhanced Display**

Q-Logic 2 Enhanced Display Kit.....\$1,720
 Part: KIT140039

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 2 Enhanced Display Shroud

Required with Q-Logic 2 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft TouchStandard **Match Joystick Shroud to Base Color** \$50
 Part: DGN137728 Part: See color below

Non-matching Joystick Shroud..... \$50
 Select color below

Red Corvette
 Part: DGN137729-PRD1

Electric Blue
 Part: DGN137729-PBU9

Back in Black
 Part: DGN137729-PBK1

Evolution Orange
 Part: DGN137729-POR7

Pretty in Pink
 Part: DGN137729-PPK3



Enhanced Display

22.4. Head Array

Adult, i-Drive Tri-Array Series, 3-Switch \$4,414
Head Array (IDHBT200-1)
 Part: ST-IDHBT200-1-Q. HCPCS: E2330
 Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.

Adult, i-Drive Pro Series Head Array \$4,461
(IDHBT300-1)
 Part: ST-IDHBT300-1-Q. HCPCS: E2330
 Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

22.5. Sip and Puff Drive**1. Sip and Puff Module**

Sip and Puff Module for Q-Logic 2 (module and bus cable).....\$1,895
 Part: KITASMB1697. HCPCS: E2325

2. Sip and Puff Straw and Mount

Stealth Sip and Puff for Comfort Series, Left Mount (SP625-CL)\$294.25
 Part: ST-SP625-CL-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, Head Array (SP625-U)\$294.25
 Part: ST-SP625-U-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series, Right Mount (SP622-CR)\$294.25
 Part: ST-SP625-CR-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Sip & Puff Module



Stealth Sip & Puff

22.6. Egg Switches

Red Egg Switch (ESRED)..... \$80.85 ea.
 Part: SWTMCRO1077
 QTY _____
Blue Egg Switch (ESBLU)..... \$80.85 ea.
 Part: SWTMCRO1079
 QTY _____
Black Egg Switch (ESBLK)..... \$80.85 ea.
 Part: SWTMCRO1081
 QTY _____

Green Egg Switch (ESGRN) \$80.85 ea.
 Part: SWTMCRO1078
 QTY _____
Yellow Egg Switch (ESYEL) \$80.85 ea.
 Part: SWTMCRO1080
 QTY _____



Egg Switch

22.7. Switch Adapter

2 to 1 Switch Adapter..... \$13.39
 Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

23. REFERENCE ORDER FORMS

Optional. Select below to note that additional order forms will accompany this order form. This will also be a link in interactive forms.

[Alternative Controls & Electronics Order Form](#)

[Positioning Components Order Form](#)

Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2019 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.



youtube.com/
QuantumRehab



facebook.com/
QuantumRehab



twitter.com/
quantumrehab



instagram.com/
quantum_rehab