




Stealth Solution SPP Cushion Order Form

401 York Ave., Duryea, PA 18642 | Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Account Number: _____ Date: _____
 Provider Name: _____
 Contact: _____
 Phone: _____ Fax: _____
 Email: _____
 PO Number: _____
 Marked for: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____



This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. All prices are on-chair MSRP. Prices, specifications, part numbers, and availability are subject to change without notice.

SOLUTION SPP CUSHION

This cushion offers skin protection and positioning with a Coolcore, 4-way stretch, breathable cover and an urethane liner. Weight capacities are 175 lbs for 10"-14" width, 350 lbs for 15"-21" width, and 550 lbs for 22" width and greater. The XXYY portion of all part numbers on this page directly relate to the dimensions of the cushion, where XX represents the width and YY represents the depth. For example, a 16" wide and 18" deep cushion would be ST-SOLSPP1618. Available sizes are shown on the chart to the right. Part numbers and prices for options with less predictable part numbers or prices can be found on the [parts list attachment](#). Multiple cushion orders can be configured on the second page.

1. Cushion Selection

- Solution SPP Cushion, 10"-20" Wide, E2607..... \$400
- Solution SPP Cushion, 22"-24" Wide, E2608..... \$500

Cushion Part: ST-SOLSPPXXYY

Cushion Width: _____
 Cushion Depth: _____

| | | Width | | | | | | | | | |
|-------|-----|-------|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| | | 10" | 12" | 14" | 16" | 17" | 18" | 19" | 20" | 22" | 24" |
| Depth | 10" | ✓ | ✓ | - | - | - | - | - | - | - | - |
| | 12" | ✓ | ✓ | ✓ | - | - | - | - | - | - | - |
| | 14" | ✓ | - | ✓ | - | - | - | - | - | - | - |
| | 16" | - | - | ✓ | ✓ | ✓ | ✓ | - | ✓ | - | - |
| | 17" | - | - | - | ✓ | ✓ | - | - | - | - | - |
| | 18" | - | - | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | - |
| | 19" | - | - | - | ✓ | ✓ | ✓ | ✓ | ✓ | - | - |
| | 20" | - | - | - | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| | 22" | - | - | - | - | - | ✓ | - | ✓ | ✓ | ✓ |

2. Substitute Fluid Proof Liner

The Solution SPP cushion is standard with an urethane liner but can be substituted with selection below. See [parts list attachment](#) for specific prices.

- Fluid Proof Polyurethane Liner \$55-\$75
 Part: PUL1709030-XXYY

3. Add Solid Base

- ABS Flat Rigidizer \$45-\$65
 Part: See [parts list attachment](#).
- ACM Rigidizer \$95
 Part: See [parts list attachment](#).
 For use with manual wheelchair only.

4. Substitute Cushion Cover

The Solution SPP cushion is standard with a Coolcore cover (TC2SPP-NBLKXXYY). The standard can be substituted with a Tek 2 or Mesh Cover with the selection below.

Substitute with Tek 2 Cover
 Part: TC2SSP-TBLKXXYY

Substitute with Mesh Cover
 Part: TC2SPP-MBLKXXYY

5. Notes

MULTIPLE ORDER SOLUTION SPP CUSHION

This section allows multiple cushions of various sizes to be ordered. Reference parts, prices, and available sizes on the previous page.


| | Quantity | Size | Width | Depth | Liner | Urethane | Polyur. | Rigidizer | ABS | ACM | Cover | Coolcore | Tek 2 | Mesh |
|---|----------|------|-------|-------|-------|----------|---------|-----------|-----|-----|-------|----------|-------|------|
| A | | | | | | | | | | | | | | |
| B | | | | | | | | | | | | | | |
| C | | | | | | | | | | | | | | |
| D | | | | | | | | | | | | | | |
| E | | | | | | | | | | | | | | |
| F | | | | | | | | | | | | | | |

SPECIAL ORDER SOLUTION SPP CUSHION


Please detail your special order cushion request. All special order cushion requests must be approved with the Quantum Custom Department. Please contact Quantum Inside Sales or Custom Department with questions. Extended lead time will apply. Cushion part numbers involved will be the cushion part number listed in previous sections modified to meet special order request stated in details below. Drawings are available below for sketch details.

Special Order Solution SPP Cushion. E2607/E2608..... \$525-\$585

Details



Top View



Rear View



Right Side View