

Tru-Balance® 2 Power Positioning Order Form



Discount: _____
Terms: _____
Freight: _____
Taken By: _____

Account Number: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
E-Mail: _____
Date: _____
PO#: _____
Marked For: _____
Ship To Address: _____
City: _____ State: _____ Zip: _____

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

BASES: SEATING SPECIFICATIONS:

- 250lb Weight Capacity. Max Dimensions 20"W x 20"D**
 J6
- 300lb Weight Capacity. Max Dimensions 20"W x 20"D**
 Q610
- 300lb Weight Capacity. Max Dimensions 22"W x 22"D**
 Rival
- 300lb Weight Capacity. Max Dimensions 22"W x 22"D**
 Q6 Edge

- Seat Width:**
 14" 15" 16" 17" 18"
 19" 20" 21" 22" 23" 24"
- Seat Width Package:**
 14"W-18"W 16"W-20"W 20"W-24"W
- Seat Depth:**
 14" 15" 16" 17" 18"
 19" 20" 21" 22" 23" 24"
- Seat Depth Package:**
 14"D-18"D 15"D-20"D 20"D-24"D
- Back Cane Height:**
 14" 15" 16"¹ 17"¹
 18" 19" 20" 21"
 22" 23" 24" 25"

****For more detailed information, please refer to the base order forms****

For Discontinued Bases, Contact Inside Sales Specialist

1. Selection of a Synergy® Back & Flip Up Armrests disqualified 16" or 17"H back canes.

TRU-BALANCE® 2 SEATING SYSTEMS BACK OPTIONS

Description	Price	HCPCS Code
List base serial # for off chair seating		
S/N: _____		
<input type="checkbox"/> TRU-Balance® 2 Power Tilt (55°) <i>55 degree center of gravity (CG) power tilt seating system.</i>	\$4,995	E1002
<input type="checkbox"/> TRU-Balance® 2 Power Recline (168°) <i>168 degrees of power recline with full drive lock-out set at 20 degrees</i>	\$5,595	E1004
<input type="checkbox"/> TRU-Balance® 2 Power Tilt (55°) & Recline (168°) <i>55 degree center of gravity (CG) power tilt seating system with 168 degrees of power recline Full drive lock-out at 20 degrees</i>	\$9,595	E1007
<input type="checkbox"/> TRU-Balance® Lift (10") ¹ <i>Must be selected in conjunction with TB Power Tilt or TB Power Tilt & Recline</i>	\$2,895	E2300
<input type="checkbox"/> TRU-Balance® 6" Power Elevate ² <i>Includes 1/4 Speed Inhibit.</i>	\$1895	
<input type="checkbox"/> TRU-Balance® 7" Power Elevate ³ <i>Includes 1/4 Speed Inhibit.</i>	\$1895	
<input type="checkbox"/> TRU-Balance® 400lb Power Tilt (55°) ⁴ <i>55 degree center of gravity (CG) power tilt seating system. 400lb wt cap.</i>	\$4,995	E1002

Description	Price	HCPCS Code
Tilt Back:		
<input type="checkbox"/> Omit Back	No Charge	
<input type="checkbox"/> Synergy® Back with 1.5" Thick Back Cushion	No Charge	
Overall Back	<input type="checkbox"/> 18" H <input type="checkbox"/> 19" H <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H	
Height:	<input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H	
Choose Lateral Extensions ¹ :		
<input type="checkbox"/> 6" for 13"-16"H Back	\$26.25	
<input type="checkbox"/> 10" for 17"-18"H Back	\$26.25	
<input type="checkbox"/> 12" For 19"-22"H Back	\$26.25	
<input type="checkbox"/> 16" For 23"-25"H Back	\$26.25	
<input type="checkbox"/> TRU-Comfort™ Back 4-way Stretch Material Cover		
Width:	<input type="checkbox"/> 16" W <input type="checkbox"/> 18" W <input type="checkbox"/> 20" W	\$675 E2620
Overall Back	<input type="checkbox"/> 20" H <input type="checkbox"/> 21" H <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H	
Height ² :	<input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H	
Width:	<input type="checkbox"/> 17" W <input type="checkbox"/> 19" W	\$675 E2620
Overall Back	<input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H	
Height ² :	<input type="checkbox"/> 27" H <input type="checkbox"/> 28" H	
Width:	<input type="checkbox"/> 22" W <input type="checkbox"/> 24" W	\$775 E2621
Overall Back	<input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H	
Height ² :		

1. 10" of scissor Lift, must select TB Tilt or TB Tilt & Recline in conjunction with lift, 1/4 Speed Inhibit, 300 lbs. max. wt., Allowable dimensions 16" x 16" - 20" x 20" on most bases, Requires Q-Logic 2 EX, Not available with Manual Recline or Power Recline Only.
2. Only available on Q6 Edge. Maximum seat size 20"W x 20"D. Available with all power positioning except 10" Lift. Increases STF Height. Includes 1/4 speed inhibit.
3. Only available on Q6000Z and Q6 Edge. Maximum seat size 20"W x 20"D. Available with all power positioning except 10" Lift. Increases STF Height. Includes 1/4 speed inhibit.
4. Q6000Z HD only.

1. Lateral extensions are required to mount lateral thoracic supports.
2. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".

BACK OPTIONS (CONT'D)

Description	Price	HCPCS Code
Tilt Back Options:		
<input type="checkbox"/> Standard Back Canes	Standard	
<input type="checkbox"/> Omit Push Handles	\$125	
<input type="checkbox"/> Flip Forward Back Canes ¹	\$150	
<input type="checkbox"/> Omit Push Handles	\$125	
<input type="checkbox"/> Offset (Meso) ²	\$195	
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)		
<input type="checkbox"/> Inset (Endo) ³	\$195	
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)		
<input type="checkbox"/> Attendant Manual Recline ³	\$750	E1226
Power Recline Back:		
<input type="checkbox"/> Synergy® Reclining Back with 1.5" Thick Back Cushion	Standard	
<input type="checkbox"/> TRU-Comfort™ Back 4-way Stretch Material Cover		
Width: <input type="checkbox"/> 16"W <input type="checkbox"/> 17" W ⁴ <input type="checkbox"/> 18"W	\$675/\$775	E2620
<input type="checkbox"/> 19"W ⁵ <input type="checkbox"/> 20"W		
<input type="checkbox"/> 22"W	\$675/\$775	E2621
Overall Back Height ⁵ : <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H <input type="checkbox"/> 22"H <input type="checkbox"/> 23"H		
<input type="checkbox"/> 24"H <input type="checkbox"/> 25"H <input type="checkbox"/> 27"H <input type="checkbox"/> 28"H		
<input type="checkbox"/> After-Market Back Interface ⁶	\$200	
Recline Back Options:		
<input type="checkbox"/> Offset (Meso) ³	\$195	
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)		
<input type="checkbox"/> Inset (Endo) ³	\$195	
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)		
Recline Color Shroud Option:		
<input type="checkbox"/> Black (Standard color option)		
Metallic Colors:N/C		
<input type="checkbox"/> Candy Apple Red	<input type="checkbox"/> Emerald Green	
<input type="checkbox"/> Viper Blue	<input type="checkbox"/> Sapphire Blue	
<input type="checkbox"/> Sunburst Orange	<input type="checkbox"/> Tru Blue ¹	
<input type="checkbox"/> Pearl Pink	<input type="checkbox"/> Orion Yellow ¹	
<input type="checkbox"/> Silver		

1. Not available with Attendant Manual Recline, Vent Tray, or O2 Holder or 400LB tilt. May interfere with Lateral Supports and Waterfall Armpads.
 2. Dual Post Flip-back Armrests interfere with offset canes. Not available with Manual Recline.
 3. Not available with Flip Back Arms, Cantilever Arms, Shape Back, Tru Comfort back, AFP, power elevating seat or 10" Lift. Requires 22"H canes. Pivot point 2.5" above seat pan.
 4. SOF is required with no extended lead time
 5. Overall Back height measurement is measured from the seat pan to the top of the back. All back heights except 25" and 28" are SOF with no extended lead time.
 6. Overall Seat depth may need to be adjusted to accommodate certain backs. Contact your Quantum Specialist for details. Does not include brackets to mount after market back.

SYNERGY® ARMRESTS

Description	Price	HCPCS Code
Tilt Armrests:		
<input type="checkbox"/> Flip Back, Cane Mounted Armrests ¹	No Charge	
<input type="checkbox"/> 2-Post, Flip-back, Height Adjustable, Removable ²	No Charge	
<input type="checkbox"/> Single Post Height Adjustable Removable ³	No Charge	
<input type="checkbox"/> Quad Release Lever		
<input type="checkbox"/> Removable, Single Post, Height Adjustable, Heavy Duty	No Charge	
Height: <input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H - 14"H		
<input type="checkbox"/> Pediatric, Removable ⁴	No Charge	
<input type="checkbox"/> Quad Release Lever		
<input type="checkbox"/> Height Adjustable Cantilever Armrests ⁵	\$195	
<i>Angle Adjustable: -30° through +30°</i>		
Power Recline Armrests:		
<input type="checkbox"/> Reclining, Flip Back	Standard	
Height: <input type="checkbox"/> 8"H - 9"H <input type="checkbox"/> 10"H - 11"H		
<input type="checkbox"/> Height Adjustable Cantilever Armrests ¹	\$195	
<i>Angle Adjustable: -30° through +30°</i>		

1. Interferes with Swing-Away laterals. Not available w/ manual recline or flip forward back canes.
 2. Height adjustment from 9"-13" in 1" increments. Interferes with swing-away laterals. Not available with a depth less than 16".
 3. Height adjustable from 11" - 17" in 1/2" increments.
 4. Height adjustment from 7.75" - 9.75" in 1/2" increments. Not available with full length armpads.
 5. Height adjustment from 11"-17" depends on cane height. Not available with Manual Recline, Flip Forward Back Canes, Cane Clamp Pushbar, or Shape Back.

SYNERGY® ARMPADS

Description	Price	HCPCS Code
<input type="checkbox"/> Straight, Full Length (2"x14")	Standard	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Straight, Desk Length (2"x10")	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Waterfall, Full Length (2"x14") ¹	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Waterfall, Desk Length (2"x10") ¹	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Falcon Gel, Full Length, (3.5"x14")*	\$92.50 ea.	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Full Length Gel Armpad Mounting Plate (Required)*	\$56.69 ea.	
<input type="checkbox"/> Falcon Gel, Desk Length, (3.5"x10")*	\$70.00 ea.	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Desk Length Gel Armpad Mounting Plate (Required)*	\$53.52 ea.	

1. Not available w/ Cane Mounted, Flip Back Armrests combined with a Synergy Back.

JOYSTICK RECEIVER MOUNT

Description	Price	HCPCS Code
<input type="checkbox"/> Universal Joystick Receiver Mount	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Omit Universal Joystick Receiver Mount	No Charge	

POWER POSITIONING ELECTRONICS

Choose Power Positioning Electronics From The Section That Matches The Electronics On The Base.

Q-Logic 2 Electronics

Description	Price	HCPCS Code
Single Actuator System:		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through Q-Logic 2 EX Controller		
<input type="checkbox"/> Single Actuator Module (AM1)	\$1,175	E2310 ¹
<input type="checkbox"/> Multiple Actuator Module (AAM) ²	\$2,375	E2311 ¹
Multiple Actuator System:		
<input type="checkbox"/> Through Q-Logic 2 EX Controller	\$2,375	E2311 ¹
<input type="checkbox"/> Through 4-Way Toggle	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through 4-Way Button Toggle	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through Quad Pushbutton	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		

1. Applies to covered actuators only.
 2. For adding up to 4 additional actuator functions in the field.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

POWER POSITIONING ELECTRONICS (CONT'D) NON-POWER LEGRESTS

Description	Price	HCPCS Code
<u>Remote Plus Electronics</u>		
Single Actuator System:		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through Remote+ Controller	\$1,175	E2310 ¹
Multiple Actuator System:		
<input type="checkbox"/> Through Remote+ Controller	\$2,375	E2311 ¹
<input type="checkbox"/> Through 4-Way Toggle	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through 4-Way Button Toggle	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through Quad Pushbutton	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<u>NE+ Electronics</u>		
Single Actuator System:		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through NE+ 6-Key Joystick	\$1,175	E2310 ¹
Two Actuator System:		
<input type="checkbox"/> Through NE+ 6-Key Controller	\$2,375	E2311 ¹
<u>NE Electronics</u>		
Tilt Only:		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<u>VR2 Electronics</u>		
Tilt Only:		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through VR2 6-Key Joystick	No Charge	
<u>VSI Electronics</u>		
Tilt Only:		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through VSI 6-Key Joystick		

1. Applies to covered actuators only.

HEADRESTS




Description	Price	HCPCS Code
<input type="checkbox"/> No Headrest Required		
Stealth:		
<input type="checkbox"/> Small (9.75"W x 3"H) (STL-OTSM250)	\$64.20	E0955
<input type="checkbox"/> Large (11"W x 6"H) (STL-OTLG150)	\$80.25	E0955
<input type="checkbox"/> Pediatric Comfort Plus (6"W) (CPN350)	\$151.94	E0955
<input type="checkbox"/> Small Comfort Plus (10"W) (CP250)	\$151.94	E0955
<input type="checkbox"/> Large Comfort Plus (14"W) (CP150)	\$151.94	E0955
<input type="checkbox"/> Removable Headrest Bracket ² (HMO475P)	\$160.50	E1028
Aftermarket Bracket:		
<input type="checkbox"/> Stealth's World's Best Removable Headrest Hardware (TWB485)	\$220.50	E1028

Description	Price	HCPCS Code
<input type="checkbox"/> 60-90° Swing-away Legrests¹	Standard	
(Adjustable in 5° increments)		
<input type="checkbox"/> Standard Extruded Footplates	No Charge	
Standard Extruded Extensions:		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
Standard Extruded Footplates:		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle-Adjustable Footplates ²	\$85 ea.	K0040 ³
Angle Adjustable Extensions:		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
Angle-Adjustable Footplates:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ⁴ <input type="checkbox"/> 6" x 11" ⁴ <input type="checkbox"/> 8" x 11" ⁵		
<input type="checkbox"/> Multi-Axis Footplates	\$495	K0040 ³
Multi-Axis Extensions:		
<input type="checkbox"/> 15.5"-19.5"		
Multi-Axis Footplates:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ⁴ <input type="checkbox"/> 6" x 11" ⁴ <input type="checkbox"/> 8" x 11" ⁵		
<input type="checkbox"/> 90° Pediatric Swing-Away Legrests⁶	\$162	
<input type="checkbox"/> Angle Adjustable Footplates ²		
Angle Adjustable Extension Length:		
<input type="checkbox"/> 7"-10"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6"		
<input type="checkbox"/> Heavy Duty Drop-In Legrests⁶	\$420	
<input type="checkbox"/> Angle Adjustable Footplates ²	\$85 ea.	K0040 ³
Angle Adjustable Extension Length:		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ⁴ <input type="checkbox"/> 6" x 11" ⁴ <input type="checkbox"/> 8" x 11" ⁵		
<input type="checkbox"/> Elevating Legrests⁶	\$325	E0990
<input type="checkbox"/> Standard Footplates		
Standard Extruded Extension Length:		
<input type="checkbox"/> 11"-14" <input type="checkbox"/> 14"-17"		
Standard Extruded Footplate Size:		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle Adjustable Footplates ²	\$85 ea.	K0040 ³
Angle Adjustable Extension Length:		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ⁴ <input type="checkbox"/> 6" x 11" ⁴ <input type="checkbox"/> 8" x 11" ⁵		
<input type="checkbox"/> High Mount, Clamp-On Footplates*		
<input type="checkbox"/> Angle Adjustable		K0040 ³
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<input type="checkbox"/> Multi-Axis		K0040 ³
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	

1. Possible caster wheel interference when adjusted to less than 70 degrees.
 2. Includes Heel Loops
 3. K0040 is not separately reimbursable on Group 2 units.
 4. Causes interference on seat widths less than 18".
 5. Causes interference on seat widths less than 20".
 6. Extensions are measured from the seatpan to footplate.




Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

ADJUSTABLE CENTER MOUNT FOOT PLATFORM

Description	Price	
<input type="checkbox"/> Adjustable Center Mount Foot Platform ¹	\$250	
Footplates :		
<input type="checkbox"/> 10"W x 11"D Split Footplates		
<input type="checkbox"/> 10"W x 8"D Split Footplates		
Extensions: ²		
<input type="checkbox"/> 10"-16" Ext.		<input type="checkbox"/> 13"-19" Ext
Calf Supports ³ :		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small (6"H x 6"W)	\$150	
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)	\$150	
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) ⁴	\$200	
Optional Modular Footplate Extension Components/Kits ^{3,5}		
<input type="checkbox"/> 4" Width Extension Kit	\$50	
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/> 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75	
Additional Optional Accessories ⁵ :		
<input type="checkbox"/> Footplate Connector	\$15	
<input type="checkbox"/> Wheel Bumpers, pair	\$12.50 ea.	

1. Not available on 614 or Select Series.
 2. Extension length is measured from top of the seat pan to footplate.
 3. Choose in Rigging Accessory Section in iQuantum.
 4. Only available on 20"W or greater.
 5. Assembly required. May require seating configuration adjustments to avoid front caster interference. N/A with the 10"W x 8"D footplates.

CENTER MOUNT POWER ARTICULATING FOOT PLATFORM

Description	Price	HCPCS Code
<input type="checkbox"/> Center Mount Power Articulating Foot Platform	\$2895	K0108
Footplates :		
<input type="checkbox"/> 10"W x 11"D Split Footplates		
<input type="checkbox"/> 10"W x 8"D Split Footplates		
Extensions: ¹		
<input type="checkbox"/> 10"-16" Ext.		<input type="checkbox"/> 13"-19" Ext
Calf Supports ² :		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small (6"H x 6"W)		
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)		
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) ³	\$200	
Optional Modular Footplate Extension Components/Kits ^{2,4}		
<input type="checkbox"/> 4" Width Extension Kit	\$50	
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/> 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75	
Additional Optional Accessories ⁴ :		
<input type="checkbox"/> Footplate Connector	\$15	
<input type="checkbox"/> Wheel Bumpers, pair	\$12.50 ea.	

1. Extension length is measured from top of the seat pan to footplate.
 2. Choose in Rigging Accessory Section in iQuantum.
 3. Only available on 20"W or greater.
 4. Assembly required. May require seating configuration adjustments to avoid front caster interference.

POWER ELEVATING LEGRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> 70° Articulating Swing-Away Legrests, 4" Articulation Range ¹	\$2,395	E1010
<input type="checkbox"/> Quad Release Lever		
Calf Pads ² :		
<input type="checkbox"/> 6" x 5.35"		<input type="checkbox"/> 8" x 5.86"
Angle Adjustable Extension Length ³ :		
<input type="checkbox"/> 11"-13"		<input type="checkbox"/> 14.5"-18"
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6"	<input type="checkbox"/> 5" x 6"	<input type="checkbox"/> 5" x 8"
<input type="checkbox"/> 6" x 8" ⁴	<input type="checkbox"/> 6" x 11" ⁴	<input type="checkbox"/> 8" x 11" ⁵
Operational Method:		
<input type="checkbox"/> Independent		<input type="checkbox"/> Combined

1. Not available on J6 or Q610 Series
 2. Choose in Rigging Accessory Section in iQuantum.
 3. Extensions are measured from top of seat pan to footplate.
 4. Only available with seat width 18" or greater.
 5. Only available with seat width 20" or greater.

FOOT RIGGING ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Wheel Bumpers	\$12.50 ea.	
<input type="checkbox"/> Heel Loops (pair) ¹	\$50	E0951
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) ²	\$80	K0108
<input type="checkbox"/> ELR/ALR Adductor Buttons (pair)	\$125	K0108
<input type="checkbox"/> Residual Limb Support	\$230 ea.	E1020
<input type="checkbox"/> Left		<input type="checkbox"/> Right

1. Standard with angle adjustable footplates.
 2. Available only with swing-away legrests.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

SYNERGY® SEAT CUSHIONS

Description	Price	HCPCS Code
<input type="checkbox"/> No Cushion Required <i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
Simplicity Cushion <i>General Use</i>		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$98	E2601
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$98	E2602
Solution® 1 Cushion <i>Skin Protection:</i>		
<input type="checkbox"/> 14" x 14" - 20" x 24"	\$430	E2603
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$480	E2604
Solution® Cushion <i>Skin Protection & Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 24"	\$430	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$480	E2608
Spectrum® Gel Cushion <i>Positioning</i>		
<input type="checkbox"/> 12" x 10" - 20" x 24"	\$530	E2605
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$580	E2606
Spectrum® Visco Cushion <i>Skin Protection & Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 24"	\$530	E2607
<input type="checkbox"/> 22" x 16" - 32" x 28"	\$580	E2608
TRU-Comfort™ 2 4-Way Stretch <i>Skin Protection</i>		
<input type="checkbox"/> 16" x 14" - 20" x 20"	\$400	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$500	E2604

MISCELLANEOUS POSITIONING COMPONENTS

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Foot Boxes (pair)* ○ Small ○ Medium <input type="checkbox"/> Large	\$498.75	K0108
<input type="checkbox"/> Stealth Padded Shoe Holders (pair)* <input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Large	\$266	K0108

Please see positioning components order form for additional types, sizes and details.

ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Quantum Backpack ¹	\$75	
<input type="checkbox"/> Oxygen Holder ²	\$150	E2208
<input type="checkbox"/> Recline Accessory Mounting Bracket ³	\$50	
<input type="checkbox"/> Weather Cover, Medium, Black	\$168	
<input type="checkbox"/> Cup Holder <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge	
<input type="checkbox"/> Tablet Holder	\$315	
<input type="checkbox"/> Deluxe Cell Phone/Media Holder	\$140	

1. Backpack may interfere w/ other rear accessories.
2. Not available with a Power Elevating Seat or a Lift/Tilt. Oxygen holder on Tilt with Tru Comfort back requires SOF. N/A on Power Recline.
3. 1" adaptor tube allows cane mounted accessories to mount to power recline back or siderail.

LATERALS

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Thoracic Laterals (pair)* <i>Bracket:</i>		
<input type="checkbox"/> 125° Swing-Away Bracket	\$472.50	E1028
<i>For Use on TRU-Comfort Backs w/ Recline</i>		
<input type="checkbox"/> Swing-Away Rail Mounted	\$472.50	E1028
<i>For Use on Synergy Back and Tru-Comfort</i>		
<i>Pad:</i>		
<input type="checkbox"/> Curved Pad, 3"x5", 4-Way Stretch Cover	\$210	E0956
<input type="checkbox"/> Curved Pad, 4"x6", 4-Way Stretch Cover	\$210	E0956
<input type="checkbox"/> Stealth Pelvic/Thigh Laterals* <i>Side:</i>		
<input type="checkbox"/> Left Side <input type="checkbox"/> Right Side		
<i>Bracket:</i>		
<input type="checkbox"/> Removable Bracket (ea.)	\$262.50	E1028
<input type="checkbox"/> Removable Bracket w/ 1" Offset (ea.)	\$400	E1028
<input type="checkbox"/> Flip Down Bracket (ea.)	\$236.25	E1028
<i>Pad:</i>		
<input type="checkbox"/> Short Pad, 4"x4" (ea.)	\$75	E0956
<input type="checkbox"/> Medium Pad, 4"x8" (ea.)	\$75	E0956
<input type="checkbox"/> Long Pad, 4"x12" (ea.)	\$75	E0956
<i>Adapter:</i>		
<input type="checkbox"/> 2" Depth/Height Adapter Plate (ea.) ²	\$26.25	K0108

1. Swing-away laterals interfere with flip-up armrests.
2. Height adapter plate is needed w/ Tru Comfort seat cushion. Not available on the flip down bracket.

LAP BELTS

Description	Price
<input type="checkbox"/> 38" Lap Belt	Standard
<input type="checkbox"/> 50" Lap Belt	Standard
<input type="checkbox"/> 60" Lap Belt	Standard
<input type="checkbox"/> 60" Padded Lap Belt	\$125
<input type="checkbox"/> 70" Lap Belt	Standard
<input type="checkbox"/> 80" Lap Belt	Standard
<input type="checkbox"/> 90" Lap Belt	Standard