

THE POWER OF CHOICE

KNOW YOUR OPTIONS

The Q-Logic 2 Advantage

- ✓ Comes standard with Bluetooth to connect to Apple or Android devices
- ✓ Infrared signal output is built into the Enhanced Display (**no additional cost**)
- ✓ Mouse emulation is built into the hand control and Enhanced Display (**no additional cost**)
- ✓ Ability to adapt to technology advancements with firmware updates
- ✓ Ability to customize auxiliary function templates with an ECU Module
- ✓ Turn sensitivity parameter available to increase or decrease joystick sensitivity
- ✓ Up to four alternate input devices can be connected (joystick, stand alone joystick, attendant control and specialty control)
- ✓ Single-switch scanning (12 scanner options: single or two stage option)
- ✓ Preset maintenance reminders



Q-LOGIC 2
DRIVE CONTROL SYSTEM



THE POWER OF CHOICE

KNOW YOUR OPTIONS

Q-Logic 2 Features

- ✓ Joystick & Enhanced Display have two switch jacks: on/off and mode
- ✓ Shortcut buttons available on the joystick (I and II)
- ✓ Profiles can be assigned as pre-set text, icons, or colors and shapes
- ✓ Ability to download digital photos to hand control or Enhanced Display (Photo management software and PC programming cable required)
- ✓ Ability to program reminders for pressure relief with a set time or as an interval; optional audible alerts
- ✓ Joystick can be programmed to be a proportional switch
- ✓ Consumer can program the 24-hour time clock, language, trip odometer reset, reminders & backlight
- ✓ Decibel level of horn can be increased in menu settings
- ✓ Caregivers can access buttons on the Enhanced Display that control power, mode, directional arrows and operate seat functions
- ✓ Enhanced Display is not required for alternative drive control access (SCIM)



Enhanced
LCD Display



Q-Logic 2
Controller



Q-LOGIC 2
DRIVE CONTROL SYSTEM