








Throw Force and Distance Chart

		IDM-HD 650 grams/1.43lbf 40mm / 1.57in
IDM-ALLAROUND-9 249 grams/.55lbf 19.5mm / .74in		
		MJM250 227 grams/.5lbf 19mm / .74in
IDM-ARLITE-9 120 grams/.26lbf 7mm / .27in		
IDPMP-9 / IDPMP-1 43 grams/.094lbf 9.6mm / .37in (castle) 10.3mm / .40in (cone) 14.4mm / .56in (tip)		
		IDM-MULTI-9 48.89 grams/.11lbf 8mm / .31in
		IDM-MICRO-9 8.5 grams/.018lbf 3.3mm / .13in



Stealth Products Joystick Matrix

Disclaimer: HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage, and bundling guidelines may apply. Stealth Products, LLC's review of Medicare guidance suggests that accessories acquired at zero-dollar net cost may not be eligible for separate Medicare reimbursement despite being separately coded and eligible to be unbundled at the time of initial purchase. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Stealth Products, LLC recommends that providers seek competent legal advice to assure a compliant Medicare billing protocol.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

stealthproducts.com

Stealth Products Joystick Matrix



Name	Micro	PMPJ	Multi	All-Around Lite	Mushroom	All-Around	All-Around HD
Part #	IDM-MICRO-9	IDPMP-9 / IDPMP-1	IDM-MULTI-9	IDM-ARLITE-9	MJM250	IDM-ALLAROUND-9	IDM-HD
Manufacturer	Mo-Vis	Stealth Products	Mo-Vis	Mo-Vis	Stealth Products	Mo-Vis	Mo-Vis
MSRP	\$3229.00	IDPMP-9 - \$3229.00 IDPMP-1 - \$3964.00	\$3229.00	\$1800.00	\$1995.00- \$2495.00	\$1575.00	\$2500.00
HCPC	E2312	E2312	E2312	E2373	K0108	E2373	E2373
Throw Force	8.5 grams / .018lbf	43 grams / .094lbf 173 grams / .38lbf	48.89 grams / .11lbf	120 grams / .26lbf	227 grams / .5lbf	249 grams / .55lbf	650 grams / 1.43lbf
Throw Distance	3.3mm / .13in	9.6mm / .37in (castle) 10.3mm / .40in (cone) 14.4mm / .56in (tip)	8mm / .31in	7mm / .27in	19mm / .74in	19.5mm / .77in	40mm / 1.57in
Throw Angle(Center to End)	10°	25°	15°	8°	18°	18°	18°
Length of stick from mechanical pivot	19mm / .74in	22.17mm / .87in (castle) 23.88mm / .94in (cone) 33.19mm / 1.3in (tip)	30mm / 1.18in	50mm / 1.97in	60mm / 2.36in	60mm / 2.36in	120mm / 4.72in
Expected Life	2,800,000 cycles	1,000,000 cycles	2,800,000	>2,800,000	N/A	15,000,000 cycles	>2,800,000 cycles
Max. Horizontal Load	4,078.86 grams / 8.99lbf	15,286 grams / 33.7lbf	5,098 grams / 11.2lbf	15,286 grams / 33.7lbf	N/A	30,591 grams / 67lbf	76,478 grams / 168lbf
Max. Vertical Load	7,647.87 grams / 16.8lbf	20,411.6 grams / 45lbf	10,197 grams / 22.48lbf	30,591 grams / 67lbf	N/A	30,591 grams / 67lbf	76,478 grams / 168lbf
Mass	225 grams / 7.93 ounces	.64oz/18.25g (.18oz / 5g)	250 grams / 8.81 ounces	315 grams / 11.1 ounces	930 grams / 32.8 ounces	315 grams / 11.1 ounces	800 grams / 28.2 ounces
Interface Cord Length	1800mm / 70.8in	2m / 6ft	1800mm / 70.8in	1800mm / 70.8in	=1.5m / 5ft	1800mm / 70.8in	1800mm / 70.8in
Operating Temp	-25°C to 50°C / -13°F to 122°F	-40°C to 85°C / -40°F to 185°F	-25°C to 50°C / -13°F to 122°F	-25°C to 50°C / -13°F to 122°F	-40°C to 65°C / -40°F to 149°F	-25°C to 50°C / -13°F to 122°F	-25°C to 50°C / -13°F to 122°F
Storage Temp	-40°C to 65°C / -40°F to 149°F	-40°C to 85°C / -40°F to 185°F	-40°C to 65°C / -40°F to 149°F	-40°C to 50°C / -40°F to 122°F	-40°C to 65°C / -40°F to 149°F	-40°C to 50°C / -40°F to 122°F	-40°C to 50°C / -40°F to 122°F
Compatible Electronics	Q-Logic, R-Net, e-Fix Alber, Ottobock Curtis, DX	Q-Logic, R-Net, Curtis, DX, MK6i	Q-Logic, R-Net, e-Fix Alber, Ottobock Curtis, DX	Q-Logic, R-Net, e-Fix Alber, Ottobock Curtis, DX	Q-Logic, R-Net, e-Fix, DX, MK6i, Mark V, Mark IV	Q-Logic, R-Net, e-Fix Alber, Ottobock Curtis, DX	Q-Logic, R-Net, e-Fix Alber, Ottobock Curtis, DX
Mounts	GAT451M ARMS260-MICRO SM651 GAT441M IDCH IDM-PCB	GAT486 GAT497 SM670 ARMS260-PMPJ IDCH IDM-PCB	GAT461M ARMS260 SM661M IDCH IDM-PCB	GAT472P-M SM660 ARMS260-AR IDM-PCB	MJM250 MJMS250 MJM260	GAT472P-M SM660 ARMS260-AR IDM-PCB	IDM-HDMT
Typical Site	Finger/Chin/Lip/Tongue	Finger/Chin	Finger/Chin/Lip	Chin/Hand	Hand/Foot	Chin/Hand/Attendant	Hand/Foot