


**Quantum® 6000 XL
Series (HD)
Group 2**
Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.
Pride Mobility Products Corporation®
182 Susquehanna Ave., Exeter, PA 18643
Sales: Phone: (866) 800-2002 • Fax: (866) 707-3422



This product is available for purchase online at www.mypride.com

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:	P.O. #:	Discount:	
Address:	Address:	Terms:		
Phone Contact:	Marked For:	Freight:		
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.		Taken By:		

All pricing is MSRP.

☐ QUANTUM 6000 2SPHD-SS REFERENCE GUIDE #1-2

- Quantum 6000 XL 2SPHD-SS \$9,386 -Flip-back, cane mounted, full length arms
- Group 2 / Single power/ Heavy duty -100 amp, Remote Plus controller
- K0837** -14" Solid tires
- 400 lbs. weight capacity -8 amp off-board charger
- 60" lapbelt
- Swing-away legrests
- Reinforced footplates
- Seat Options: Omit Solid Tires (-\$213)
- TB Tilt E1002 Omit contoured seat (-\$275)

Select seating , front rigging options in appropriate section & base color.

☐ QUANTUM 6000 2SPHD-C REFERENCE GUIDE #2-2

- Quantum 6000 2SPHD-C \$9,386 -20" W x 22" D Ltd. Recline Gray vinyl high-back
- Group 2 / Single Power / Heavy Duty -Full length armrests
- K0838** -100 amp, Remote Plus controller
- 400 lbs. weight capacity -14" Solid tires
- 8 amp off-board charger
- 60" lapbelt
- Foot Platform
- Omit Solid Tires (-\$213)
- Omit contoured seat (-\$275)

Select seating , front rigging options in appropriate section & base color.

When configuring power chair use reference guides to determine accessories included in the MSRP pricing package

Model number with suffix must be used for all chairs billed through Medicare.

SPECIFIC CONFIGURATION

- Quantum 6000 XL \$9,173
- Quantum 6000 2SPHD-SS \$9,386
- Quantum 6000 2SPHD-C \$9,386
- Omit contoured seat²⁰ (-\$275)

FRAME COLOR

These colors are on a painted metal frame.

- Metallic Red NC
- Metallic Blue NC
- Sunburst Orange NC
- Metallic Black NC
- Metallic Silver NC

Due to the nature of the painted frame design, color substitutions may not be made in the field.

BATTERIES (2 required)

- Group 24, 70 amp, Gel E2363 \$315
- NF-22 Battery Required w/ Vent Pkg.⁷ (Select below)
- NF-22, 50-amp Batteries E2361
 - Maintenance-free AGM \$225
 - Maintenance-free Gel Cell \$265
- Ship unit with batteries installed \$200
- Ship unit with batteries, but do not install NC

BATTERY CHARGER

- 8 amp Off-board Charger Std.

CONTROLLER

Expandable Controllers:

- (Applicable for use w/ Tilt)
- 3-Key Remote Plus 100 amp^{3,5} Std.
(Select Box for Profile)
- 1 Profile Multiple Profile
- 7-Key Access. Remote Plus 100 amp^{1,3,4} \$275
- Remote Plus 100 amp with Speed Dial^{3,5} \$275
- Remote Plus 100 amp with Toggles^{3,5} \$275
- Remote Plus 100 amp with Switch Jacks^{3,5} \$275
- Microdrive, 100 amp^{3,5,6} \$275
- Omni Plus²⁵ \$1721.63

PROGRAMMER

- PG Programmer \$423

Quantum® 6000 XL Series (HD) Group 2

Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

CONTROLLER MOUNT

<input type="checkbox"/> Fixed, In-Line Joystick Mount	Std.
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Quick Disconnect, Joystick Mount ²²	NC
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Swing-away Joystick Mount ²³ E1028	\$196
<input type="checkbox"/> 4" L ²⁶ <input type="checkbox"/> 7" L	
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Multi-axis Joystick Mount K0108	\$100
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<i>(Fixed, does not swing-away)</i>	
<input type="checkbox"/> Vertical Adj. Clamp ⁸ K0108	\$95
Provides vertical joystick adjustment w/o adjusting armrest height.	

JOYSTICK HANDLE OPTIONS

PG Controls

<input type="checkbox"/> Std. Knob	Std.
<input type="checkbox"/> Mushroom Handle E2323	\$71
<input type="checkbox"/> T-Bar E2323	\$107
<input type="checkbox"/> Chin Cup E2323	\$71
<input type="checkbox"/> Extended I-Handle E2323	\$71
<input type="checkbox"/> Soft Ball Handle E2323	\$50
<input type="checkbox"/> Medium Ball Handle E2323	\$50
<input type="checkbox"/> Goal Post Handle E2323	\$66
<input type="checkbox"/> 3" W <input type="checkbox"/> 4" W	

DRIVE WHEELS

<input type="checkbox"/> 14" Pneumatic (pair) ²⁴	NC
<input type="checkbox"/> 14" Flat-Free (pair)	\$213

FRONT RIGGINGS

<input type="checkbox"/> Swing-away footrests ²¹	Std.
(\$162 w/ captains seating)	
<input type="checkbox"/> Elevating Leg Rests E0990	\$325
Adj. Calf Pads w/ ELRS (1 hole std.)	
<input type="checkbox"/> 1 hole <input type="checkbox"/> 3 holes* <input type="checkbox"/> 6 holes*	
*Potential for interference w/ shrouds	
<input type="checkbox"/> HD Drop-Ins, 70 ^{o2,8,15} K0108	\$420
(HD Drop-ins N/A w/ contoured or Specialty seating)	
<input type="checkbox"/> Foot Platform (Substitute for S/A)	NC
<input type="checkbox"/> No Front Rigging	NC
<input type="checkbox"/> Power-elevating, high-pivot ^{16,17} E1010	\$1695
Legrests, (calf pads included) ¹⁸	
<input type="checkbox"/> Independent <input type="checkbox"/> Combined ¹⁶	
<input type="checkbox"/> Power-elevating, Articulating ^{15,16,17} E1010	\$2395
Articulating Range: <input type="checkbox"/> 4" L Std. <input type="checkbox"/> 6" L NC	
Set at: <input type="checkbox"/> 70°	
<input type="checkbox"/> Independent <input type="checkbox"/> Combined ¹⁶	
<input type="checkbox"/> Swing-away Power-elevating, E1010	\$2395
Articulating ^{15,16,17}	Std.
Articulating Range: <input type="checkbox"/> 4" L	
Set at: <input type="checkbox"/> 70°	
<input type="checkbox"/> Quad Release Lever	NC
<input type="checkbox"/> Independent <input type="checkbox"/> Combined ¹⁶	
<input type="checkbox"/> Calf Pads Left & Right (6 x 5.35)	Std.
<input type="checkbox"/> Calf Pads Left & Right (8.00 x 5.86)	NC
Footplates (pair)	
<input type="checkbox"/> Reinforced Footplate	Std.
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5 x 6"	
<input type="checkbox"/> 5"x8" Angle-Adj. Footplate ¹⁰	\$150
<input type="checkbox"/> 6"x8" Angle-Adj. Footplate ⁹	\$150
<input type="checkbox"/> 6"x11" Angle-Adj. Footplate ⁹	\$150
<input type="checkbox"/> 8"x11" Angle-Adj. Footplate ¹¹	\$150
<input type="checkbox"/> Wheel Bumpers (Ang-adj. only)	\$25
Extensions:²²	
<input type="checkbox"/> 11"-14" Ext. ¹⁹	NC
<input type="checkbox"/> 14"-17" Ext	Std.
<input type="checkbox"/> High mount Clamp-on; Ang. Adj.	\$326
<input type="checkbox"/> Small <input type="checkbox"/> Large	
<input type="checkbox"/> High mount Clamp-on; Multi-axis	\$412
<input type="checkbox"/> Small <input type="checkbox"/> Large	
<input type="checkbox"/> High mount footplatform ⁸ K0108	\$450
<input type="checkbox"/> Mid-wheel 6 foot platform	
Lower Extensions	
<input type="checkbox"/> 11"-14" L	
<input type="checkbox"/> 14"-17" L	
Rigging Options	
<input type="checkbox"/> Heel Loops (pair) ¹³ E0951	\$20
<input type="checkbox"/> Gel Neoprene Footrest K0108	
Sleeve (pair) ¹⁴	\$65
<input type="checkbox"/> ELR / ALR Adductor Buttons (pair) K0108	\$100
<input type="checkbox"/> POSAlinc Knee Adductor (pair)	
for S/A legrests K0108	\$464
<input type="checkbox"/> Residual limb support left E1020	\$300
<input type="checkbox"/> Residual limb support right E1020	\$300

Footnotes for pages 1 & 2

- Std. w/ one programmable drive profile.
- May cause interference w/ front casters w/ certain seat widths and lower leg lengths. Call Customer Service.
- Needed for controlling power positioning systems.
- Needed for lighting packages.
- Avail. w/ five individual drive profiles.
- Requires multi-axis joystick mount; N/C for mount.
- Vent configuration consists of (2) Group 24 batteries and (1) NF-22.
- N/A on Specialty, contoured seating or waterfall armpads.
- Only avail. on 16" W or greater seats.
- Only avail. on 14" W or greater seats.
- Only avail. on 20" W or greater seats.
- Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured or rehab seating w/ cushions is selected. The overall length may vary approx. 3" depending on choice of seat and /or cushion.
- Std. w/ ang-adj. footplates.
- Avail. w/ S/A footrests & 70 degree HD drop-ins.
- Standard with heel loops and ang-adj. footplates, check sizes below.
- Requires the addition of an ALM.
- Power legrests are only available on TRU-Balance and Synergy Seating Systems. Only the combined & independent motion can be selected on the Synergy Seat.
- Std. calf pad is not depth adjustable.
- Special order on ALR's.
- Applicable when ordering no seat, Synergy seat or any power positioning system.
- Please note swing away footrests are standard when ordered with Power Positioning.
- N/A w/ waterfall armpads; only available on Quantum Seating Systems and only one option maybe mounted.
- Can only mount one armrest accessory in combination w/ S/A joystick.
- If ordered w/ reference guides or configuration packages discounts may apply.
- Please consult specialty controls order form for a variety of secondary controls and mounting options.
- N/A with waterfall armpads and dual post armrests.

**Quantum® 6000 XL
Series (HD)
Group 2**

**Order Form / Price List
Revised: February 2007**

EFFECTIVE NOVEMBER 15TH 2006.

STANDARD SEATING

Medium-Back NC

**Medium-Back
w/ Solid Seat Pan^{1,2}** NC

- Brown Vinyl
 Beige Vinyl
 Gray Vinyl
 Black Vinyl
 Gray Fabric
 Brown Fabric

Medium-Back Seating Dimensions

(Dimensions 18"W x 18"D)

Width: 16" NC 18" Std.
Depth: 16" NC 18" Std.

UPGRADED SEATING

Ltd. Recline High-Back w/ Headrest NC

**Ltd. Recline High-Back with Headrest^{1,2}
w/ Solid Seat Pan** NC

- Black Vinyl
 Gray Vinyl
 Beige Vinyl
 Gray Fabric
 Brown Fabric

High-Back Seating Dimensions

(Dimensions 18"W x 18"D)

Width: 16" NC 18" Std. 20" NC 22" NC 24" NC
Depth: 16" NC 18" Std. 20" NC 22" NC 24" NC

22" W AND ABOVE WILL NOT BE LTD RECLINE.

ARMRESTS

For Medium and High-Back Seats

- Flip-up, removable Std.
 Flip-up, height-adjustable \$170
removable⁴ E0973

Armpads

- Desk Length, 10" NC
 Left NC Right NC
 Full Length, 14" NC
 Left Std. Right Std.

MANUAL RECLINER^{2,5}

Manual Recline High Back E1226 \$600

**Manual Recline High Back¹
w/ Solid Seat Pan** E1226 \$600
(180 degrees of Recline)
(Dimensions 18"W x 18"D)

- Gray Vinyl Beige Vinyl
Width: 16" NC 18" Std. 20" NC
Depth: 16" NC 18" Std. 20" NC

**When ordering a lap belt with this seat, an
additional 10" is added to accommodate**

CUSHIONS³

- Simplicity Cushion**
16 x 16 - 20 x 20 E2601 \$98
 Solution Cushion
10 x 10 - 20 x 18 E2607 \$380
20 x 20 - 20 x 20 E2607 \$430
22 x 16 - 22 x 20 E2608 \$430
 Spectrum Cushion
 Gel
10 x 10 - 20 x 18 E2605 \$480
20 x 20 - 20 x 20 E2605 \$530
22 x 16 - 22 x 20 E2606 \$530
 Visco
10 x 10 - 20 x 18 E2607 \$480
20 x 20 - 20 x 20 E2607 \$530
22 x 16 - 22 x 20 E2608 \$530
 Structure Cushion
 Gel Visco
10 x 10 - 20 x 18 E2607 \$560
20 x 20 - 20 x 20 E2607 \$610
22 x 16 - 22 x 20 E2608 \$610
 TRU-Comfort Cushion
16 x 16 - 20 x 18 \$380
20 x 20 - 20 x 22 \$380
22 x 16 - 22 x 24 \$480

Footnotes for pages 3:

1. Includes black vinyl seat pan, back available in all listed colors.
2. 18" x 18" std. Special measurements are not available.
3. See Synergy Cushion order forms for detailed pricing, sizes & codes.
4. N/A on manual recliner.
5. Rear mount accessories are N/A.

Quantum® 6000 XL Series (HD) Group 2

Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

Synergy Seating Systems¹⁰

Synergy Seating™ System \$434
(Price with Base)

WIDTH

Set Seat Width At:

<input type="checkbox"/> 10" NC	<input type="checkbox"/> 11" NC
<input type="checkbox"/> 12" NC	<input type="checkbox"/> 13" NC
<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	

Seat Width Package

<input type="checkbox"/> 10"W-14"W	<input type="checkbox"/> 14"W-18"W
<input type="checkbox"/> 16"W-20"W	<input type="checkbox"/> 20"W-24"W

DEPTH

Set Seat Depth At:

<input type="checkbox"/> 10" NC	<input type="checkbox"/> 11" NC
<input type="checkbox"/> 12" NC	<input type="checkbox"/> 13" NC
<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	

Seat Depth Package

<input type="checkbox"/> 10"D-14"D	<input type="checkbox"/> 14"D-18"D
<input type="checkbox"/> 16"D-20"D	<input type="checkbox"/> 20"W-24"W

BACK CANES

Back Cane Height

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

BACK OPTIONS

<input type="checkbox"/> Solid Plate Back up to 20" W	\$275
<input type="checkbox"/> Solid Plate Back 21-24" W	\$325

Back Height:¹

<input type="checkbox"/> 18"H	<input type="checkbox"/> 19"H	<input type="checkbox"/> 20"H	<input type="checkbox"/> 21"H
<input type="checkbox"/> 22"H	<input type="checkbox"/> 23"H	<input type="checkbox"/> 24"H	<input type="checkbox"/> 25"H
<input type="checkbox"/> Sling Back	\$98		

Back height will be the same as back cane height.
Includes gap for cushion.

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

TRU-Comfort Back⁹

Back Height:	<input type="checkbox"/> 22"	<input type="checkbox"/> 25"
<input type="checkbox"/> Back up to 20" W	\$595	
<input type="checkbox"/> 22" W	\$695	

ARMRESTS

<input type="checkbox"/> Flip-back, cane mounted ²	Std.
<input type="checkbox"/> 2-Post, flip-back, height ²	E0973 \$150ea
adjustable, removable armrests	
<input type="checkbox"/> Removable, adjustable ³	E0973 \$75ea
quick height	
<input type="checkbox"/> Quad Release Lever	NC
(3" of height adjustment in 3/8" increments)	
<input type="checkbox"/> Removable, height adjustable	E0973 \$75ea
heavy duty armrests	
<input type="checkbox"/> 8"H-12"H	<input type="checkbox"/> 10"H-14"H

ARMPADS

<input type="checkbox"/> Full length, Straight, 2" X 14"	Std.
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Desk Length, Straight, 2" X 10"	Std.
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Full Length, Waterfall, 2" X 14" ⁴	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Desk Length, Waterfall, 2" X 10" ^{4,5}	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right

HEAD REST⁶

<input type="checkbox"/> Stealth ^{7,8}	E0955
Large (11" wide x 6" high)	
	\$210
<input type="checkbox"/> Stealth ^{7,8}	E0955
Small (9 3/4" wide x 3" high)	
	\$210
<input type="checkbox"/> Stealth, swing-away (Flip-Down)	
Head Rest Bracket	
	E1028 \$150
Interfaces w/ any headrest that uses 1/2" square down-tube.	
<input type="checkbox"/> Medical Applications ^{7,8}	E0955
Large (11" wide x 6" high)	
	\$275
<input type="checkbox"/> Medical Applications ^{7,8}	E0955
Small (9 3/4" wide x 3" high)	
	\$275
<input type="checkbox"/> Fixed, Offset Mounting Bracket	NC
<input type="checkbox"/> TRU-Comfort ^{7,8}	E0955
<input type="checkbox"/> 8" Wide w/ Cobra hardware	
	\$417
<input type="checkbox"/> 10" Wide w/ Cobra hardware	
	\$417
<input type="checkbox"/> Bracket for after market headrest	\$180

Footnotes for pages 4:

1. Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion
2. S/A laterals interfere w/ flip-back arms.
3. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
4. Only avail. on removable, adjustable, quick-height arms.
5. N/A w/ S/A joystick mount.
6. Available w/ solid backs only.
7. Comes with swing-away / removable hardware; HCPCS code is E1028.
8. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
9. Cane mounted, flip-back arms are N/A.
10. Manual systems do not include headrest.

Quantum® 6000 XL Series (HD) Group 2

Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

TRU-BALANCE TILT

TRU-Balance Power Tilt E1002 \$4995
55 degree center of gravity (CG) power tilt seating system. 400 lbs. weight capacity

Tower Selection

- Adjustable-angle Std.
Will raise minimum STF by 7/8".
 Fixed NC

SYSTEM SPECIFICATIONS

Set Seat Width At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	

Seat Width Package

<input type="checkbox"/> 14"W-18"W	<input type="checkbox"/> 16"W-20"W
<input type="checkbox"/> 20"W-24"W	

Set Seat Depth At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	

Seat Depth Package

<input type="checkbox"/> 14"D-18"D	<input type="checkbox"/> 16"D-20"D
<input type="checkbox"/> 20"D-24"D	

Push Bar selection must match back width;
center tube standard when selecting flip-up arms.

Set Back Cane Height At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

Seat Overall Back Height¹

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

Back Type

- Solid Plate Back Std.
 TRU-Comfort Back⁴
Back Height: 22" 25"
 Back up to 20" W \$595
 22" W \$695
 Omit Back -(\$100)

Quantum_6000_XL_Series_ORD-01-25-07_Group_2

ARMRESTS

- Flip-back, cane mounted⁹ Std.
Mounted on back canes with full width spreader bar.
 2-Post, flip-back, hgt. adj, removable^{9,10} NC
 Removable, adj. quick height^{11,14} NC
 Quad Release Lever NC
(3" of height adjustment in 3/8" increments)
 Removable, single-post, height adjustable, heavy duty¹⁴ NC
 8"-12" Adjustment
 10"-14" Adjustment

Armrest Pads

- Full Length, straight, 2"x14" NC
 Left Right
 Desk Length, straight, 2"x10" NC
 Left Right
 Full Length, Waterfall, 2"x14"¹² NC
 Left Right
 Desk Length, Waterfall, 2"x10"^{12,13} NC
 Left Right

HEAD REST

- Stealth^{7,8} E0955
Large (11" wide x 6" high) \$210
 Stealth^{7,8} E0955
Small (9 3/4" wide x 3" high) \$210
 Stealth, swing-away (Flip-Down)
Head Rest Bracket E1028 \$150
Interfaces w/ any headrest that uses 1/2" square down-tube.
 Medical Applications^{7,8}
Large (11" wide x 6" high) E0955 \$275
 Medical Applications^{7,8}
Small (9 3/4" wide x 3" high) E0955 \$275
 Omit Head Rest - To be supplied by provider
 Fixed, Offset Mounting Bracket NC
 TRU-Comfort^{7,8} E0955
 8" Wide w/ Cobra hardware \$417
 10" Wide w/ Cobra hardware \$417
 Bracket for after market headrest \$180

POWER POSITIONING ELECTRONICS

Select Electronics to be used

- PG Remote Plus

Electronics for Tilt Only

- Through Toggle Box¹⁵ NC
 Left Right (mounts to side of armrest)
 Through Single Switch Tilt Module^{16,17} E2310 \$500

Single Motor Controller:

- Through Remote Plus Controller E2310* \$1175
 Through single or multiple motor controller w/ 2 Jacks^{16,17} E2310* \$1175

*E2310 coding may vary due to configuration

ACCESSORIES

- Rear Basket^{2,3,5} \$103
 Cane/Crutch Holder^{2,3,5} \$98
 Oxygen Holder^{2,3,5} E2208 \$103
 Oxygen Holder^{2,3,4} E2208 \$150
(TRU-Balance Tilt & Synergy Seat only)
 Walker Holder^{2,3,5} \$92
 Flag Holder (SOF on rehab seating) \$17
 Weather Cover, Large, Black \$202

ACCESSORIES CONT.

- Saddle bag (contoured seating only) \$28
 Cup holder \$17
 Left Right
N/A on desk length arm with S/A mount or toggle box.
 Lap Belt Std.
 Pediatric 38" N/A w/ limited recline seating
 50"³ 60"³ 70" 80" 90"
 Vent battery storage kit⁶ \$250
(Includes extended battery cover, strap & bracket, for use with third battery in base)
Can accommodate up to NF-22 battery-not included
N/A w/ footplatform or high mount footplatform

Footnotes for page 5:

- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- Includes mounting bracket and hardware.
- N/A on manual recliner high-back.
- N/A w/ cane mounted, flip-back arms.
- N/A on TB Tilt.
- Vent Tray is SOF due to vent and seating type.
- Comes with swing-away / removable hardware; HCPCS code is E1028.
- These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
- S/A laterals interfere w/ flip-back arms.
- May cause interference with certain dimension configurations. Call Quantum Customer Service.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- N/A w/ S/A joystick mount.
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- N/A w/ waterfall armpads.
- Consult specialty control order form for a wide variety of switch types and sizes--2 required.
- Does not include switch.

Lowest seat-to-floor height

Synergy Seat:	16.625"
TB Tilt:	19.5"
(w/ fixed towers towers)	
TB Tilt:	17.375"
(w/ ang-adj. towers towers)	
High-back Seat:	21.5"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers (please see Tilt information above), 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Increased widths may increase seat-to-floor height. Contact customer service for details. Listed measurements are in the lowest position. It does not include elevating systems.

Quantum Custom Delivery Service

Includes: liftgate service and delivery of your product(s) inside your facility
Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice