



**Quantum Rehab®**  
 A Division of Pride Mobility Products Corporation®  
 182 Susquehanna Ave., Exeter, PA 18643  
 Phone: (866)800-2002 • Fax: (866) 707-3422

# Order Form

Discount: \_\_\_\_\_  
 Terms: \_\_\_\_\_  
 Freight: \_\_\_\_\_  
 Taken By: \_\_\_\_\_



Account Number: \_\_\_\_\_  
 Provider Name: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 PO#: \_\_\_\_\_  
 Marked For: \_\_\_\_\_  
 Ship To Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

*DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.*

## QUICK SHIP OPTION 1

Description	Price	HCPCS Code
<input type="checkbox"/> Stylus LS Option 1	\$965	K0004
250 lbs. Weight Capacity 17.5" Seat-to-Floor Height Includes:		
- Tension Adjustable, Fluid Resistant Seat		
- Non-Tension Adjustable Back		
- 17"-19" Back Upholstery Adjustment		
- 70° Swing-Away Foot Riggings w/ 14-18" Toolless Ext.		
- Reinforced Footplates w/ Heel Loops		
- Quick Release Axle		
- 24" Composite Wheels w/ Smooth Poly Tire & Aluminum Hand Rim		
- 8" Flat Free Casters		
- Push-to-Lock Wheel Locks		
Adjustable Rear Anti-Tip Wheels <sup>1</sup>	\$80	E0971
Full Length Height Adj. Armrests, 9"-14" <sup>2</sup>	\$245	E0973
48" Lap Belt	\$50	E0978
<b>Dimensions:</b>		
<b>Qty:</b>		
___ 16x16 - STYLUSLS1007	No Charge	
___ 16x18 - STYLUSLS1008	No Charge	
___ 18x16 - STYLUSLS1009	No Charge	
___ 18x18 - STYLUSLS1010	No Charge	
___ 20x16 - STYLUSLS1011	\$415	E2201
___ 20x18 - STYLUSLS1012	\$415	E2201

1. Anti-tip wheels must be used in conjunction with O2 holder. Non-billable.  
 2. Desk Length Armrests Included.

## QUICK SHIP OPTION 2

Description	Price	HCPCS Code
<input type="checkbox"/> Stylus LS Option 2	\$965	K0004
250 lbs. Weight Capacity 16.5" Seat-to-Floor Height Includes:		
- Tension Adjustable, Fluid Resistant Seat		
- Non-Tension Adjustable Back		
- 17"-19" Back Upholstery Adjustment		
- 70° Swing-Away Foot Riggings w/ 14-18" Toolless Ext.		
- Reinforced Footplates w/ Heel Loops		
- Quick Release Axle		
- 22" Composite Wheels w/ Smooth Poly Tire & Aluminum Hand Rim		
- 6" Flat Free Casters		
- Push-to-Lock Wheel Locks		
Adjustable Rear Anti-Tip Wheels <sup>1</sup>	\$80	E0971
Full Length Height Adj. Armrests, 9"-14" <sup>2</sup>	\$245	E0973
48" Lap Belt	\$50	E0978
<b>Dimensions:</b>		
<b>Qty:</b>		
___ 16x16 - STYLUSLS1001	No Charge	
___ 16x18 - STYLUSLS1002	No Charge	
___ 18x16 - STYLUSLS1003	No Charge	
___ 18x18 - STYLUSLS1004	No Charge	
___ 20x16 - STYLUSLS1005	\$415	E2201
___ 20x18 - STYLUSLS1006	\$415	E2201

1. Anti-tip wheels must be used in conjunction with O2 holder. Non-billable.  
 2. Desk Length Armrests Included.

## ADDITIONAL OPTIONS

Description	Price	HCPCS Code
<b>Armrests:</b>		
<input type="checkbox"/> Desk Length, Flip-Up <sup>1</sup>	No Charge	
<b>Back Canes:</b>		
<input type="checkbox"/> Angle Adjustable Back Canes <sup>2</sup>	\$100	
<b>Back Options:</b>		
<input type="checkbox"/> Tension Adjustable Seat Back	\$350	E2611
<input type="checkbox"/> Shape GP Back	\$465	E2611
<input type="checkbox"/> Comfort Company Back	Various	
<b>See Comfort Company Order Form</b>		
<b>Foot Riggings:</b>		
<input type="checkbox"/> Elevating Legrests, 13-17" Extensions w/ Calf Pads <sup>3,4</sup>	\$275	K0195
<b>Accessories:</b>		
<input type="checkbox"/> O2 Holder	\$150	E2208
<input type="checkbox"/> Wheel Lock Extensions	\$18.50 ea.	E0961
<input type="checkbox"/> Pull-To-Lock Wheel Locks <sup>3</sup>	\$30	
<b>Seat Extension:</b>		
<input type="checkbox"/> 2" Seat Depth Extension Kit <sup>4,5</sup>	\$79	
<b>Ultra Hemi:</b>		
<input type="checkbox"/> Ultra Hemi Kit	\$249	K0056 <sup>6</sup>
<i>Includes 20" Tires (WHLASMB2175) and 5" Casters (FRMASMB16919) w/ Short Fork To Achieve 15" Seat-to-Floor Height</i>		

1. Additional lead times may apply. Please consult your Quantum Specialist.  
 2. Adjustable from 85° to 105° in 10° increments. Not compatible with standard Non Tension Adj back.  
 3. Pull-to-Lock wheel locks are recommended when choosing ELR's to avoid interference with locking mechanism.  
 4. Seat extensions are not available w/ ELR's  
 5. Use HCPCS code E2203 when extending from 18"D to 20"D.  
 6. HCPCS Code applies only to 20" tires..

## SYNERGY® CUSHIONS

Description	Price	HCPCS Code
<input type="checkbox"/> No Cushion Required		
<i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
<b>Simplicity Cushion</b> <i>General Use</i>		
<input type="checkbox"/> 16" x 16" - 20" x 18"	\$98	E2601
<b>Simplicity GP Cushion</b> <i>General Use</i>		
<input type="checkbox"/> 16" x 16" - 20" x 18"	\$98	E2601
<b>Solution® 1 Cushion</b> <i>Skin Protection</i>		
<input type="checkbox"/> 14" x 14" - 20" x 18"	\$380	E2603
<b>Solution® Cushion</b> <i>Skin Protection &amp; Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 18"	\$380	E2607
<b>Spectrum® Gel Cushion</b> <i>Positioning</i>		
<input type="checkbox"/> 12" x 10" - 20" x 19"	\$480	E2605
<b>Spectrum® Visco Cushion</b> <i>Skin Protection &amp; Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 18"	\$480	E2607

## SEAT-TO-FLOOR HEIGHTS

STF Height	Wheel Position				Ultra-Hemi STF Height	Wheel Position	
	24"	8" Tall Fork	22"	6" Tall Fork		20"	5" Short Fork
19.5"	BOT	BOT			16.0"	MID	BOT
18.5"	MID	MID	BOT	BOT	15.5"	MID/TOP <sup>5</sup>	MID
17.5"	TOP	TOP	MID	MID	15.0"	TOP	TOP
16.5"			TOP	TOP			